

Dynabeads[®] M-270 Epoxy

For research use only

INDEX

1. PRODUCT DESCRIPTION
2. PRINCIPLE
3. INSTRUCTIONS FOR USE
 - 3.1. Weight of Dynabeads M-270 Epoxy
 - 3.2. Washing Procedure
 - 3.3. Resuspension in Organic Solvent
 - 3.4. Coating Procedure
 - 3.5. Isolation of Target Molecule
 - 3.6. Target Elution Procedure
4. GENERAL INFORMATION
 - 4.1. Product Characteristics
 - 4.2. References
 - 4.3. Additional Material Needed
 - 4.4. Recommended Buffers/Solutions
 - 4.5. Related Dynabeads Products
 - 4.6. Storage & Stability
 - 4.7. Warnings & Limitations
 - 4.8. Trademarks & Patents
 - 4.9. Intellectual Property Disclaimer
 - 4.10. Limited Use Label License
 - 4.11. Warranty

1. PRODUCT DESCRIPTION

Dynabeads[®] M-270 Epoxy are uniform, superparamagnetic beads composed of highly cross-linked polystyrene with magnetic material precipitated in pores evenly distributed throughout the particles. The beads are further coated with a hydrophilic layer of glycidyl ether (epoxy) functional groups which seals the iron oxide inside the beads. These surface reactive epoxy groups allow for binding of proteins, peptides or other ligands, with covalent bond formation at neutral pH. Binding of ligands through amine or thiol groups occurs with no further activation of the surface.

2. PRINCIPLE

Dynabeads M-270 Epoxy act as a solid support in a wide variety of biomagnetic separations and manipulations. The hydrophilic surface ensures low non-specific binding, excellent dispersion abilities and easy handling of the beads in a wide variety of buffers.

Their size makes them particularly suitable for protein isolation for sample preparation, bioassays, selection of affinity binders etc. Due to the gentle pull of the beads to the magnet, they can also be used for selection of fragile cells. For cell separation in general, the larger Dynabeads (4.5 µm) are recommended.

Via primary amino or sulf-hydryl groups in the ligand, Dynabeads can be coated directly with proteins, peptides, antibodies, enzymes or other target specific molecules for the isolation of targets such as hormones, receptors, disease markers,

bacteriophages etc. Alternatively, the bead-surface can be further modified to bind additional functional groups.

Once coupled with ligand, the beads are added to the sample containing your target molecule. After a short incubation allowing affinity capture of the target, the beads are pulled to the side of the test tube by the use of a magnet (Dynal MPC[™]), allowing aspiration of unbound material. The magnet facilitates washing and target-concentration. Bead-bound targets can be used directly in bioassays, boiled in application buffer and analysed on SDS-PAGE. Alternatively, the target molecule can be eluted off the beads with conventional elution methods (high salt, low pH etc.).

3. INSTRUCTIONS FOR USE

The freeze dried Dynabeads M-270 Epoxy must be equilibrated in an appropriate buffer before coating. A sample of beads can be weighed out directly from the vial, using a microscale balance. Due to the short stability of hydrophilic epoxy groups in aqueous buffers we recommend to only prepare the wanted amount of beads for each experiment.

The beads can also be resuspended in organic solvents like diglyme or DMF, in which the beads are stable for at least one year at 2-8°C (see section 3.3. below).

3.1. Weight of Dynabeads M-270 Epoxy

Use the table below to be able to weigh out the wanted number of beads.

No. of beads	Weight (mg)
10 ⁹	15
2 x 10 ⁹	30
3 x 10 ⁹	45
4 x 10 ⁹	60 (content of 143.01)
5 x 10 ⁹	75
10 ¹⁰	150
1,5 x 10 ¹⁰	225
2 x 10 ¹⁰	300 (content of 143.02D)

3.2. Washing Procedure

A washing step is necessary to equilibrate the beads in appropriate buffer.

1. Add 0.1 M sodium phosphate buffer, pH 7.4 to the dry beads, to give a concentration of approximately 10⁹ beads per ml.
2. Vortex for 30 seconds and incubate with mixing for 10 minutes.
3. Place the tube on a magnet (Dynal MPC) for 2 minutes and carefully pipette off the supernatant, leaving beads undisturbed.
4. Remove the test tube from the magnet and resuspend the beads carefully in the same volume with 0.1 M sodium phosphate buffer, pH 7.4 again. Vortex to mix properly.
5. Place the tube on the magnet for 2 minutes and pipette off the supernatant.
6. Resuspend the washed beads in the same buffer to a volume that will give the recommended bead concentration and ammonium sulfate concentration after addition of ligand solution and stock solution of ammonium sulfate (see below for recommended values). The beads are now washed and ready for coating.

3.3. Resuspension in Organic Solvent

Resuspension in an organic solvent like DMF (dimethyl formamide) makes it easy to withdraw samples from the vial.

1. Add 2 ml diglyme to the 143.01-vial (contains 4 x 10⁹ beads) or 10 ml to the 143.02-vial (contains 2 x 10¹⁰ beads) to give a final concentration of 2 x 10⁹ beads/ml.
2. When samples are to be used: Resuspend the beads well in diglyme by vortexing for 1-2 minutes and transfer the wanted amount of beads to a tube with cap.
3. Place the tube on the magnet for 4 minutes. Pipette off the supernatant carefully, leaving beads undisturbed.
4. Equilibrate the beads (see procedure 3.2. above).

3.4. Coating Procedure

For coating with proteins or other ligands, it is recommended to use > 3 µg pure ligand per 10⁷ Dynabeads M-270 Epoxy for coupling of small ligands (such as peptides), and slightly higher concentrations for coupling of larger ligands. An overall final concentration of 1-2 x 10⁹ beads per ml during incubation is recommended.

To enhance binding, a final concentration of ammonium sulfate of 1-3 M should be used during coating. The optimal ammonium sulfate concentration will depend on the nature of the ligand. A hydrophilic or small ligand (peptide) requires high ammonium sulfate concentration (up to 3 M), whereas most proteins will be sufficiently coated with 1 M ammonium sulfate. Some biomolecules may even lose their function (some Igs) at molarities > 1.5 M.

The ligand should be dissolved in PBS or a similar buffer, without other proteins or stabilizers. Avoid buffers with amino or sulfo groups (e.g. Tris). Buffers with slightly higher pH (e.g. borate buffer pH 9.5) may speed up the formation of covalent bonds between amino-groups in the ligand and the bead-surface. Such buffers can also be used to adjust the volume to give the wanted concentrations of beads, ligand and ammonium sulfate if necessary.

Optimization of the protocol for each ligand is always recommended.

NOTE: Sugars or stabilizers may disturb the binding and should be removed from the ligand preparation by dialysis, spin columns, gel chromatography etc. prior to coupling.

Calculation example: 3 mg beads = 2 x 10⁸ beads. The ligand requirement for 2 x 10⁸ beads using 3 µg/10⁷ beads (20 µg/mg beads) is 60 µg ligand. At a ligand concentration of 1 mg/ml, aliquote a volume of 60 µl. Wash and resuspend the beads in 60 µl 0.1 M phosphate buffer to adjust the bead-concentration and mix ligand and beads thoroughly. The needed ammonium sulfate concentration is 1 M. Add 60 µl (1/3 the final volume) of a 3 M stock-solution. The final volume is 180 µl, which gives a final bead concentration of 1.2 x 10⁹ beads/ml (the ligand volume is then 1/3 of the final volume).

1. Make a homogeneous suspension of the washed beads as described in 3.2 above.
2. Add the calculated amount of ligand solution to the bead suspension. Vortex to ensure good mixing before adding the calculated ammonium sulfate stock solution.
3. Incubate for 16-24 hours at 37°C with slow tilt rotation. Incubation with temperatures down to 4°C may be used for temperature sensitive ligands, but be aware that the covalent bond formation is slower and less efficient at low temperatures, and an additional 24 hours should be used to ensure covalent coupling. Do not let the beads settle during the incubation period.
4. Place the tube on the magnet for 4 minutes for magnetic separation. Carefully turn the magnet (with the tube in place) upside-down twice, to ensure collection of any beads that might remain in the cap. Remove the supernatant.
5. Wash the coated beads a total of four times with PBS or PBS with blocking protein. Blocking protein like BSA or skimmed milk powder should be added to 0.1-0.5% when this does not interfere with your downstream application.

6. If the downstream application involves elution steps, physically adsorbed ligand can be removed by washing for 10 minutes in 0.5-1% (w/v) Tween 20/Triton X-100 or similar non-ionic detergent.
7. Resuspend the coated beads to the desired concentration in PBS or PBS with blocking protein. The beads are now ready for use.

3.5. Isolation of Target Molecule

Efficient isolation of target molecules is dependent on bead concentration, target molecule concentration, the ligand's affinity for the target molecule and time. Equilibrium target-ligand binding is reached after approximately 1 hour. Binding is performed at 2-8°C, at a recommended bead concentration of 1-10 x 10⁹ beads/ml.

1. Add sample containing target molecule to the coated beads (2 x 10⁸ beads). For 100 kD protein use volume containing approximately 25 µg target molecule to assure an excess of this molecule.
2. Incubate with tilting and rotation for one hour to capture the target (incubation times as low as 10 minutes can be used with concentrated protein samples).
3. Place the tube on the magnet for 4 minutes to collect the beads at the tube wall (for viscous samples, increase the time on the magnet). Pipette off the supernatant.
4. Wash beads 3 x using 1 ml PBS each time (see procedure 3.2. above).

3.6. Target Elution Procedure

Conventional elution methods can be applied for the elution of target protein from the beads. Low pH (2.8-3.5), change in ionic strength, affinity elution, electrophoresis, polarity reducing agents, deforming eluants can be applied, or even boiling the beads in SDS-PAGE application buffer for direct characterization of protein on SDS-PAGE. The method of choice depends on affinity of target molecule to the ligand coated onto the beads, stability of target molecule and downstream application and detection methods. Most proteins will be eluted off at pH 3.1 following the procedure described below, but some protein functionality might be lost under such harsh conditions. If maintaining functionality of the target molecule is essential, try milder elution conditions first (such as high salt, e.g. 2M NaI, or stepwise elution reducing pH from 6 to 3). Use the same if the bead-bound ligand must remain functional to allow reuse of the beads.

1. Add an appropriate amount (e.g. 60 µl) of 0.1 M citrate (pH 3.1) to beads with immobilized target.
2. Mix well by tilting and rotation for 2 minutes.
3. Place the test tube on the magnet and transfer the supernatant, containing purified target, to a clean tube.

To ensure reuse of the beads and functionality of the isolated target molecule, bring both beads and target back to physiological pH (7.4) immediately after elution.

4. GENERAL INFORMATION


Invitrogen Dynal AS complies with the Quality System Standards ISO 9001:2000 and ISO 13485:2003.

4.1. Product Characteristics

Dynabeads M-270 Epoxy are sold as freeze dried beads. Each mg contains 6-7 x 10⁷ beads. The product is available in two formats, 60 mg (Cat. no. 143.01) and 300 mg (Cat. no. 143.02D).

Typical bead-characteristics for any given lot of this product:

Diameter: 2.8 µm
 Specific surface area: 2-5 m²/g
 Active chemical functionality: 0.1 - 0.2 mmol/g



4.2. References:

1. Marfatia KA *et al.* Identification and characterization of the human MOG1 gene. *Gene* 2001, 266(1-2):45-56.
2. Hamm J, Alessi DR and Biondi RM. Bi-functional, substrate mimicking RNA inhibits MSK1-mediated cAMP-response element-binding protein phosphorylation and reveals magnesium independent conformational changes of the kinase. *J.Biol.Chem.* 2002, 277(48):45793-45802.
3. Laine S *et al.* In vitro and in vivo interactions between the hepatitis B virus protein P22 and the cellular protein gC1qR J. *Virology*. 2003 Dec; 77(23):12875-80.

4.3. Additional material needed

- Magnetic device (Dynal MPC™, Magnetic Particle Concentrator)
- Mixing/rotation device (Dynal® Sample Mixer or similar)
- Test tubes, glassware and pipettes
- Antibody/other selecting molecule
- Buffers/solutions (see below)

4.4. Recommended buffers/solutions

Invitrogen Dynal recommends the use of 0.1M Phosphate buffer pH 7.4, PBS and PBS Tween. Other buffers might be used, but buffers containing amino groups (e.g. Tris) should not be used for coating of ligand.

High pH and high temperature during coating procedures are optimal conditions for quick formation of chemical bonds. (The upper pH and temperature limit is determined by the ligand). The molarity of salt in the final coating solution should never be less than 0.05M.

0.1M sodium phosphate buffer (pH 7.4):

- 2.62 g NaH₂PO₄ x H₂O (MW 137.99)
- 14.42 g Na₂HPO₄ x 2H₂O (MW 177.99)

Dissolve in distilled water, adjust pH if necessary and adjust to 1 litre. This buffer is used for pre-washing of the beads and should not be added any protein, sugar etc.

3M ammonium sulfate (stock solution):

- 39.6 g (NH₄)₂SO₄ (MW 132.1)
- Dissolve in 0.1M sodium phosphate buffer (pH 7.4) and adjust to 100 ml.

0.1 M citrate pH 3.1:

- 2.10 g citric acid (C₆H₈O₇ x H₂O, MW 210.14)
- Dissolve in 90 ml distilled water, adjust to pH 3.1 and adjust to 100 ml.

2 M NaI:

- 3 g NaI (MW 149.9) to 10 ml distilled water.

PBS pH 7.4 (phosphate buffered saline):

- 0.26 g NaH₂PO₄ x H₂O (MW 137.99)
- 1.44 g Na₂HPO₄ x 2H₂O (MW 177.99)
- 8.78 g NaCl (MW 58.5)

Dissolve in 900 ml distilled water, adjust pH if necessary and adjust to 1 litre

PBS with 0.1% BSA/HSA/skimmed milk:

- Include 0.1% (w/v) BSA/HSA/skimmed milk (0.1g) in 100 ml PBS (above).

PBS/Tween 20/Triton X:

- Include 0.5-1.0 % (w/v) Tween 20/Triton X (50-100 mg) in 100 ml PBS (above).

4.5. Related Dynabeads Products

Cat. no.	Product name
142.03/04 & 301.01	Dynabeads® M-280 Tosylactivated
143.05D/06D	Dynabeads® M-270 Carboxylic Acid
143.07D/08D 302.12	Dynabeads® M-270 Amine Dynabeads® M-450 Tosylactivated
302.02D/03D 655.01/02/03	Dynabeads® M-450 Epoxy Dynabeads® MyOne™ Tosylactivated
650.11/12/13	Dynabeads® MyOne™ Carboxylic Acid

4.6. Storage & Stability

Precautions should be taken to prevent bacterial contamination of the beads.

When stored in unopened vials at 2-8°C,

Dynabeads M-270 Epoxy are stable until the expiration date printed on the label.

The beads should be washed once before use. Beads should not be autoclaved, but can be incubated with ethanol (70%, 1 hour) or gamma irradiated.

Coated beads may be stored at 2-8°C for several weeks or even months, depending on the stability of the immobilized ligand. Coated beads should be washed once for 5 min in PBS/BSA before use. Use the magnet to collect the beads according to the washing procedure.

If a preservative is needed for storage of coated beads, a final concentration of 0.02% (w/v) sodium azide (NaN₃) may be added to the storage buffer. Carefully remove before use by washing (see 3.2 above). Required safety precautions must be followed when handling this cytotoxic material.

4.7. Warnings & Limitations

This product is for research use only. Not intended for any animal or human therapeutic or diagnostic use unless otherwise stated.

Sodium azide is toxic if ingested. Avoid pipetting by mouth. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. When disposing through plumbing drains, flush with large volumes of water to prevent azide buildup.

Certificate of Analysis (CoA) is available upon request.

Material Safety Data Sheet (MSDS) is available at <http://www.invitrogen.com>.

4.8. Trademarks & Patents

Dynal®, Dynabeads® and Dynal MPC™ are either registered trademarks or trademarks of Invitrogen Dynal AS, Oslo, Norway. Any registration or trademark symbols used herein denote the registration status of trademarks in the United States. Trademarks may or may not be registered in other countries.

4.9. Intellectual Property Disclaimer

Invitrogen Dynal will not be responsible for violations or patent infringements that may occur with the use of our products.

4.10. Limited Use Label License

No. 5: Invitrogen Technology – The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of patents owned or controlled by Invitrogen Corporation which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser

is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact

Licensing Department,
Invitrogen Corporation,
1600 Faraday Avenue, Carlsbad,
California 92008.
Phone (760) 603-7200.
Fax (760) 602-6500.
Email: outlicensing@invitrogen.com

4.11. Warranty

The products are warranted to the original purchaser only to conform to the quantity and contents stated on the vial and outer labels for the duration of the stated shelf life. Invitrogen Dynal's obligation and the purchaser's exclusive remedy under this warranty is limited either to replacement, at Invitrogen Dynal's expense, of any products which shall be defective in manufacture, and which shall be returned to Invitrogen Dynal, transportation prepaid, or at Invitrogen Dynal's option, refund of the purchase price.

Claims for merchandise damaged in transit must be submitted to the carrier.

This warranty shall not apply to any products which shall have been altered outside Invitrogen Dynal, nor shall it apply to any products which have been subjected to misuse or mishandling. ALL OTHER WARRANTIES, EXPRESSED, IMPLIED OR STATUTORY, ARE HEREBY SPECIFICALLY EXCLUDED, INCLUDING BUT NOT LIMITED TO WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. Invitrogen Dynal's maximum liability is limited in all events to the price of the products sold by Invitrogen Dynal. IN NO EVENT SHALL INVITROGEN DYNAL BE LIABLE FOR ANY SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES. Some states do not allow limits on warranties, or on remedies for breach in certain transactions. In such states, the limits set forth above may not apply.

Invitrogen Dynal is a part of the Invitrogen Group.

Contact details for your local Invitrogen sales office/technical support can be found at <http://www.invitrogen.com/contact>

© Copyright 2008 Invitrogen Dynal AS, Oslo, Norway. All rights reserved.

SPEC-05660