

# One Shot® Mach1™-T1<sup>R</sup> Chemically Competent *E. coli*

Cat. No. C8620-03

Size: 20 reactions

## Shipping and Storage

The One Shot® Mach1™-T1<sup>R</sup> Chemically Competent *E. coli* kit is shipped on dry ice. Upon receipt, store at -80°C.

## Kit Contents

Each kit contains the reagents listed below. Transformation efficiency is greater than  $1 \times 10^9$  cfu/μg DNA.

Item	Composition	Amount
S.O.C. Medium (store at room temperature or +4°C)	2% Tryptone, 0.5% Yeast Extract, 10 mM NaCl, 2.5 mM KCl, 10 mM MgCl <sub>2</sub> , 10 mM MgSO <sub>4</sub> , 20 mM glucose	6 ml
Mach1™-T1 <sup>R</sup> Cells	---	21 x 50 μl
pUC19 Control DNA	10 pg/μl in 5 mM Tris-HCl, 0.5 mM EDTA, pH 8	50 μl

## Genotype

F: φ80(*lacZ*)Δ*M15* Δ*lacX74* *hsdR*(r<sub>K</sub>m<sub>K</sub><sup>+</sup>) Δ*recA1398* *endA1* *tonA*

## Information for Non-U.S. Customers

*For European Customers*

The Mach1™-T1<sup>R</sup> *E. coli* strain is genetically modified to carry the *lacZ*Δ*M15* *hsdR* *lacX74* *recA* *endA* *tonA* genotype. As a condition of sale, this product must be in accordance with all applicable local legislation and guidelines including EC Directive 90/219/EEC on the contained use of genetically modified organisms.

*For All Non-U.S. Customers*

The parental strain of Mach1™-T1<sup>R</sup> *E. coli* is the non-K-12, wild-type W strain (ATCC #9637, S. A. Waksman). Although the parental strain is generally classified as Biosafety Level 1 (BL-1), we recommend that you consult the safety department of your institution to verify the Biosafety Level.

## Product Qualification

Competent cells are tested for transformation efficiency using the control plasmid included in the kit. Transformed cultures are plated on LB plates containing 100 μg/ml ampicillin and the transformation efficiency is calculated. Test transformations are performed in duplicate. Transformation efficiency should be greater than  $1 \times 10^9$  cfu/μg plasmid DNA. In addition, untransformed cells are tested for appropriate antibiotic sensitivity, the absence of phage contamination, and resistance to phage T5 (a standard test that demonstrates resistance to phage T1).

## Features of the Strain

The Mach1™-T1<sup>R</sup> *E. coli* strain is modified from the wild-type W strain (ATCC #9637, S. A. Waksman) and has a faster doubling time compared to other standard cloning strains. With Mach1™-T1<sup>R</sup> cells, you can visualize colonies 8 hours after plating on ampicillin selective plates. You can also prepare plasmid DNA 4 hours after inoculating a single, overnight-grown colony in the selective media of choice. Note that this feature is not limited to ampicillin selection. Additional key features of the Mach1™-T1<sup>R</sup> *E. coli* strain include:

- *lacZ*Δ*M15* for blue/white color screening of recombinants
- *hsdR* mutation for efficient transformation of unmethylated DNA from PCR applications
- Δ*recA1398* mutation for reduced occurrence of homologous recombination in cloned DNA
- *endA1* mutation for increased plasmid yield and quality
- *tonA* mutation to confer resistance to T1 and T5 phage

## General Guidelines

Follow these guidelines when using One Shot® Mach1™-T1<sup>R</sup> Chemically Competent *E. coli*.

1. Handle competent cells gently as they are highly sensitive to changes in temperature or mechanical lysis caused by pipetting. Thaw One Shot® competent cells on ice, and transform cells immediately following thawing. After adding DNA, mix by swirling or tapping the tube gently. **Do not mix cells by pipetting.**
2. One Shot® Mach1™-T1<sup>R</sup> cells do not require IPTG to induce expression from the *lac* promoter. If blue/white screening is required to select for transformants, spread 40 μl of 40 mg/ml X-Gal in dimethylformamide on top of the agar. Let the X-Gal diffuse into the agar for approximately 1 hour.

## Transforming Competent Cells

Perform the following before starting the transformation procedure:

- Equilibrate a water bath to 42°C.
- Warm the vial of S.O.C. Medium (supplied with the kit) to room temperature.
- Spread X-Gal onto LB agar plates containing antibiotic, if desired.
- Warm the selective plates in a 37°C incubator for 30 minutes (use one plate for each transformation). If you are including the pUC19 control, make sure that you have one LB agar plate containing 100 µg/ml ampicillin. **Note:** For optimal growth of Mach1™-T1<sup>R</sup> *E. coli* cells, it is essential that selective plates are prewarmed to 37°C prior to spreading.

### Transformation Procedure

We recommend including the pUC19 control plasmid DNA supplied with the kit in your transformation experiment to verify the efficiency of the competent cells. **Do not** use these cells for electroporation.

1. Thaw, on ice, one vial of One Shot® Mach1™-T1<sup>R</sup> Chemically Competent *E. coli* for each transformation.
2. Add 1 to 5 µl of the DNA (10 pg to 100 ng) into a vial of One Shot® cells and mix gently. **Do not mix by pipetting up and down.** If you are transforming the pUC19 control, add 1 µl (10 pg) into a separate vial of One Shot® cells and mix gently.
3. Incubate the vial(s) on ice for 30 minutes.
4. Heat-shock the cells for 30 seconds at 42°C without shaking.
5. Remove the vial(s) from the 42°C bath and place them on ice for 2 minutes.
6. Add 250 µl of room temperature S.O.C. Medium to each vial.
7. Cap the vial(s) tightly and shake horizontally at 37°C for 1 hour at 225 rpm in a shaking incubator.
8. Spread 25-100 µl of the transformation mix on a **prewarmed** selective plate. Store the remaining transformation mix at +4°C. Additional cells may be plated out the next day, if desired.
9. Invert the plate(s) and incubate at 37°C. If you are using ampicillin selection, visible colonies should appear within 8 hours, and blue/white screening can be performed after 12 hours. If you are selecting transformants with an antibiotic other than ampicillin, incubate plates overnight.
10. Select overnight-grown colonies and analyze by plasmid isolation, PCR, or sequencing. For plasmid isolation, inoculate a single, overnight-grown colony in 2 ml of **prewarmed** selective media (e.g. LB + ampicillin, LB + kanamycin, LB + Zeocin™, etc.). For optimal results, we recommend inoculating as much of the single colony as possible. Shake at 37°C for 4 hours before isolating the plasmid.

### Calculating Transformation Efficiency

Use the following formula to calculate the transformation efficiency as transformants (in cfu) per µg of plasmid DNA.

$$\frac{\# \text{ of colonies}}{10 \text{ pg transformed DNA}} \times \frac{10^6 \text{ pg}}{\mu\text{g}} \times \frac{300 \mu\text{l total transformation volume}}{X \mu\text{l plated}} = \frac{\# \text{ transformants}}{\mu\text{g plasmid DNA}}$$

This product is the subject of one or more of U.S. Patent No. 4,981,797 and foreign equivalents owned by Invitrogen Corporation. The purchase of this product conveys to

### Limited Use Label License No. 5: Invitrogen Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of patents owned or controlled by Invitrogen Corporation which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Invitrogen Corporation, 1600 Faraday Avenue, Carlsbad, California 92008. Phone (760) 603-7200. Fax (760) 602-6500. Email: [outlicensing@invitrogen.com](mailto:outlicensing@invitrogen.com)

©2003-2006 Invitrogen Corporation. All rights reserved. For research use only. Not intended for any animal or human therapeutic or diagnostic use.