



One Shot[®] *ccdB* Survival[™] 2 T1^R Chemically Competent Cells

Cat. No. A10460

Size: 10 reactions
Store at -80°C

Description

The One Shot[®] *ccdB* Survival[™] 2 T1^R strain is resistant to the toxic effects of the *ccdB* gene (Bernard and Couturier, 1992; Bernard *et al.*, 1993) and may be used to propagate and maintain vectors containing the *ccdB* gene (e.g. Gateway[®] Technology vectors). The transformation efficiency of One Shot[®] *ccdB* Survival[™] 2 T1^R chemically competent cells is greater than 1×10^9 cfu/ μ g pUC19 Control DNA.

Components Supplied

Components Supplied	Amount
<i>ccdB</i> Survival [™] 2 T1 ^R Competent Cells	11 \times 50 μ l
pUC19 Control DNA (10 pg/ μ l)	50 μ l
S.O.C. Medium	6 ml

Genotype

F⁻ *mcrA* Δ (*mrr-hsdRMS-mcrBC*) Φ 80*lacZ* Δ M15 Δ *lacX74* *recA1* *ara* Δ 139 Δ (*ara-leu*)7697 *galU galK rpsL* (Str^R) *endA1 nupG fhuA::IS2*

Product Qualification

The Certificate of Analysis (CofA) provides detailed quality control information for each product. The CofA is available on our website at www.invitrogen.com/cofa, and is searchable by product lot number, which is printed on each box.

Part No. 100003536

Rev. Date: 28 May 2008

General Guidelines

- Handle competent cells gently as they are highly sensitive to changes in temperature and mechanical lysis caused by pipetting.
- Thaw competent cells on ice, and transform cells immediately. Add DNA to tube and mix contents by gentle swirling or tapping. **Do not mix cells by pipetting.**

Cloning Unstable Inserts

When cloning unstable inserts such as lentiviral DNA containing direct repeats (e.g. Invitrogen's ViraPower™ Lentiviral Expression Kits), we recommend using the following modifications to reduce the chance of recombination between direct repeats:

- Select and culture transformants at 25°C–30°C.
- **Do not use** “rich” bacterial medias as they tend to give rise to a greater number of unwanted recombinants.
- If your plasmid confers chloramphenicol resistance, select and culture transformants using LB medium containing 15–30 µg/ml chloramphenicol **in addition to** the antibiotic appropriate for selection of your plasmid.

Transforming Competent Cells

Perform the following before starting the transformation procedure:

- Equilibrate a water bath to 42°C.
- Warm the vial of S.O.C. Medium (supplied with the kit) and LB Medium (if needed) to room temperature.
- Warm selective plates in a 37°C incubator for 30 minutes (use 1 or 2 plates for each transformation). If you are including the pUC19 control, make sure that you have one LB agar plate containing 100 µg/ml ampicillin.

Transformation Procedure

Use this procedure to transform One Shot[®] *ccdB* Survival[™] 2 T1^R chemically competent cells. We recommend including the pUC19 control plasmid DNA to verify the transformation efficiency. **Do not** use these cells for electroporation.

1. Thaw one vial of One Shot[®] cells on ice for each transformation.
2. Add 1–5 μ l of DNA (10 pg to 100 ng) into a vial of One Shot[®] cells and mix gently. **Do not mix by pipetting up and down.** For the pUC19 control, add 1 μ l (10 pg) of DNA into a separate vial of One Shot[®] cells and mix gently.
3. Incubate the vials on ice for 30 minutes.
4. Heat-shock the cells for 30 seconds at 42°C without shaking.
5. Remove the vials from the 42°C bath and place them on ice for 2 minutes.
6. Add 250 μ l of pre-warmed S.O.C. Medium to each vial.
7. Cap the vials tightly and shake horizontally at 37°C for 1 hour at 225 rpm in a shaking incubator. *
8. Spread 25–100 μ l from each transformation on a pre-warmed selective plate. We recommend that you plate two different volumes to ensure that at least one plate will have well-spaced colonies. For the pUC19 control, dilute the transformation mix 1:10 into LB Medium and plate 25–100 μ l.
9. Store the remaining transformation mix at +4°C. Additional cells may be plated out the next day, if desired.
10. Incubate plates overnight at 37°C. *

* **Note:** When cloning unstable inserts, incubate cells at 25°C–30°C.

Calculating Transformation Efficiency

Use the following formula to calculate the transformation efficiency as the number of transformants (in cfu) per μg of plasmid DNA.

$$\frac{\text{\# of colonies}}{\text{pg pUC19 DNA}} \times \frac{10^6 \text{ pg}}{\mu\text{g}} \times \frac{\text{volume of transformants}}{X \mu\text{l plated}} \times \text{dilution factor}$$

For example, if transformation of 10 pg of pUC19 DNA yields 50 colonies when 25 μl of a 1:10 dilution is plated, then the transformation efficiency is:

$$\frac{50 \text{ colonies}}{10 \text{ pg DNA}} \times \frac{10^6 \text{ pg}}{\mu\text{g}} \times \frac{300 \mu\text{l}}{25 \mu\text{l plated}} \times 10 = 6 \times 10^8$$

References

Bernard, P., Couturier, M., "Cell killing by the F plasmid CcdB protein involves poisoning of DNA-topoisomerase II complexes". *J. Mol. Biol.* (1992) Aug 5; 226(3):735-45.

Bernard, P., Kezdy, K.E., Van Melderen, L., Steyaert, J., Wyns, L., Pato, M.L., *et al.* (1993) The F plasmid CcdB protein induces efficient ATP-dependent DNA cleavage by gyrase. *J Mol Biol* 234: 534-541.

©2008 Invitrogen Corporation. All rights reserved.

For research use only. Not intended for any animal or human therapeutic or diagnostic use.