Mouse anti-TTF-1, Clone 8G7G3/1

For In Vitro Diagnostic Use

Lot No.

[] 08-1221 6.0 mL 2nd Gen Predilute Antibody, Ready-To-Use
 [X] 18-0221 1.0 mL Concentrate Antibody

INTENDED USE

For In Vitro Diagnostic Use

Invitrogen's monoclonal Mouse anti-TTF-1 antibody is intended to qualitatively detect TTF-1 in formalin-fixed paraffin-embedded (FFPE) or frozen tissue sections. Interpretation must be made within the context of the patient's clinical history and other diagnostic test by a qualified pathologist.

BACKGROUND

TTF-1 (38 kD) is a member of the *Nkx2* family of homeodomain-containing transcription factors. TTF-1 has been identified as a mediator of thyroid-specific gene transcription¹, and is involved in the activation of pulmonary-specific genes.²

TTF-1 is expressed in the epithelial cells of the lung and thyroid gland. Clone 8G7G3/1 detects a single ~40 kD band on Western blots of cell lines known to express TTF-1 mRNA and was negative on cell lines known not to contain TTF-1. Immunohistochemical staining of tissue strongly correlated with staining of an established TTF-1 antibody and staining by in situ hybridization.³

TTF-1 is useful as part of a panel of antibodies used to differentiate lung and thyroid primary carcinomas from other carcinomas that metastasize to the lung such as breast and colon carcinomas³. It may also be useful as part of a panel of antibodies, including anti-calretinin, used to differentiate primary lung adenocarcinomas from mesotheliomas.⁴

REAGENT PROVIDED

Mouse anti-TTF-1 Antibody is diluted in phosphate buffered saline (PBS), pH 7.4, and 1% bovine serum albumin (BSA) with 0.1% sodium azide (NaN₃) as a preservative.

Immunogen: Recombinant protein derived from the rat TTF-1 sequence Total protein concentration: g/L

Clone: 8G7G3/1³ Isotype: IgG₁ Antibody concentration: mg/L

STORAGE: 2-8°C

POSITIVE CONTROL TISSUE: Lung carcinoma or thyroid gland

EXPECTED STAINING PATTERN: Nucleus

INSTRUCTIONS FOR USE

PRETREATMENT REQUIREMENTS:

Epitope Retrieval: Required (EDTA pH 8.0) (See page 2 for protocol)

Enzyme Digestion: Not Required

Mouse anti-TTF-1 may be diluted according to Table 1 when using the Invitrogen detection systems below.

Table 1. Dilution Table

Invitrogen Kit	Predilute Antibody	Dilution for Concentrate	Incubation Time
Histostain-SP or SAP kits*	Ready-To-Use	1: 50 - 1: 100	60 min.
Histostain [®] -Plus Kits	Ready-To-Use	1: 100 - 1: 200	30-60 min.
SuperPicTure [™] Polymer Kits	Ready-To-Use	1: 100 - 1: 200	30-60 min.
Cap-Plus [™] Kits	Ready-To-Use	1: 100 - 1: 200	30-60 min.

* Use Histostain-SP or SAP kits only for Cat. No. 08-0XXX and 18-X001 to 18-X200 primary antibodies.

This is a guideline only. Optimal antibody concentrations may vary based on specimen and preparation method used, and should be determined by each individual laboratory.

SPECIMEN PREPARATION

1. Use tissue fixed in 10% neutral buffered formalin or other fixative on regular basis, or frozen tissue sections.
2. Cut 3-4 µm sections and place on positively charged slides.
3. Dry overnight at 37° C or for 2-4 hrs at 58°C.

PRETREATMENT

Heat Induced Epitope Retrieval (HIER), if required

1. Deparaffinize slides. Tissue sections should be mounted on silane, poly-L-Lysine, or HistoGrip (Cat. No. 00-8050) coated slides.
2. Wash slides with distilled water 3 times for 2 min each.
3. Place a 1L glass (Pyrex) beaker containing 500 mL of 0.01 M citrate buffer (Cat. No. 00-5000) or EDTA solution (Cat. No. 00-5500) on a hot plate. Heat the buffer solution until it boils. (*This step may be prepared before slide deparaffinization, as the buffer may take several minutes to boil*).
4. Put the slides in a slide rack and place in the beaker with boiling buffer. Keep it boiling for 15 minutes.
5. After heating, remove beaker from the hot plate and allow it to cool down for at least 15-20 minutes at room temperature.
6. Rinse slides with PBS (Cat. No. 00-3000) and begin the immunostaining protocol.

Enzyme Digestion, if required

1. Prewarm the enzyme of choice at 37°C for 10 min.
2. Add the prewarmed enzyme to a tissue section and incubate at 37°C for 10 min.
3. Wash in several changes of PBS (Cat. No. 00-3000) and begin the immunostaining protocol.

RECOMMENDED MANUAL STAINING PROCEDURE

1. Submerge slides in peroxidase quenching solution and rinse with PBS.
2. Apply serum blocking solution.
3. Apply primary antibody and incubate for 30-60 min at room temperature; rinse with PBS.
4. Apply secondary antibody and incubate for 10 min at room temperature; rinse with PBS.
5. Apply enzyme conjugate and incubate for 10 min at room temperature; rinse with PBS.
6. Apply chromogen and incubate for 5-10 min at room temperature; rinse with PBS.

MATERIALS REQUIRED BUT NOT PROVIDED

Reagent	Catalog No.
1. HistoGrip™	00-8050
2. PAP Pen	00-8888
3. Isotype Control for Rabbit or Mouse Primary Antibody	08-6199 or 08-6599
4. Antibody Diluent	00-3118
5. PBS (0.01 M PBS)	00-3000
6. Mayer's Hematoxylin	00-8011
7. Citrate Buffer pH 6.0 (if required for HIER)	00-5000
8. EDTA Solution (if required for HIER)	00-5500
9. Digest-All™ 1, Digest-All™ 2, or Digest-All™ 3 (if required for Enzyme Digestion)	00-3007 or 00-3008 or 00-3009
10. SuperPicTure™ polymer kit, or LAB-SA kit (Histostain®-Plus and Cap-Plus™) (see Table 2).	
11. Chromogen/substrate (if not included in detection kit): Single Solution AEC (Cat. No. 00-1111), or DAB (Cat. No. 00-2014), or Fast-Red (Cat. No. 00-2234).	
12. Mounting solution: Histomount™ (for DAB: Cat. No.00-8030), GVA (for AEC, or Fast-Red: Cat. No. 00-8000), or Clearmount™ (for AEC, DAB, or Fast-Red: Cat. No. 00-8010).	

REFERENCES

1. Ikeda, K, et al. *J. Biol. Chem.* 270: 8108-8114, 1995.
2. Bohinski RJ, et al. *Mol Cell Bio.* 14(9):5671-5681, 1994.
3. Holzinger A, et al. *Hybridoma* 15(1):49-53, 1996.
4. Khor A, et al. *Hum Pathol* 30(6):695-700, 1999.

TRADEMARKS

Cap-Plus™, Clearmount™, Digest-All™, HistoGrip™, Histomount™, and Histostain®, HistoST5050™, NBA™, PicTure™, and Zymed® are trademarks of Zymed Laboratories, Inc.

Authorized Representative for IVDD 98/79/EC: Invitrogen Ltd.
Inchinnan Business Park
3 Fountain Drive
Paisley
PA4 9RF
UK

www.invitrogen.com

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288 • E-mail: techsupport@invitrogen.com

PIN: 30604

(Rev 08/08) DCC-08-1089

© 2007 Invitrogen Corporation. All right reserved. These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or www.invitrogen.com). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses.

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