

Dynabeads® mRNA DIRECT™ Micro Kit

Isolate mRNA for reverse transcription and PCR

Cat. no. 610.21

Initial Check-List

- All buffers, except the 10 mM Tris-HCl, should be brought to room temperature prior to use. The 10 mM Tris-HCl should be stored on ice or 2-8°C prior to use.
- Ensure that the Dynabeads® Oligo (dT)₂₅ have been fully resuspended before use. Resuspend by brief vortexing or careful pipetting.
- Check that your Lysis/Binding Buffer has not precipitated. If any precipitation is observed, warm to room temperature and shake to full resuspension.
- A rapid lysis in Lysis/Binding Buffer is critical for obtaining un-degraded mRNA. Thawing of frozen material prior to lysis must be avoided.
- When working with cells isolated with cell-specific Dynabeads®, make sure all the cell-specific Dynabeads® are removed from the lysate before adding Dynabeads® Oligo (dT)₂₅.
- Prepare your reverse transcription reaction mix prior to mRNA isolation and keep on ice.
- We recommend to immediately use the Dynabeads®-mRNA complex for reverse transcription and PCR.
- Use sterile, RNase-free microtubes and pipette tips.
- Wear disposable gloves and change them frequently.
- Refer to the package insert for detailed protocols.



Add 100 μ l Lysis/Binding Buffer to your starting sample and homogenize (e.g. by pipetting) until complete lysis is obtained.




Transfer the lysate to a tube containing 20 μ l pre-washed Dynabeads[®] Oligo (dT)₂₅



Mix 5 minutes on a roller at room-temperature for mRNA to hybridize to the oligo-dT on the beads.



Place on magnet and remove supernatant. Add 100 μ l Washing Buffer A and resuspend by careful pipetting.



Wash once in 100 μ l Washing Buffer A and twice in 100 μ l Washing Buffer B. The magnet allows for easy washing and change of buffers.



Resuspend the Dynabeads[®]-mRNA complex in 100 μ l ice-cold Tris-HCl and store on ice for reverse transcription and downstream applications. (Elution is optional at this step.)



Immediately prior to adding the reverse transcription PCR-mix to the Dynabeads[®]-mRNA complex, apply the magnet and remove the supernatant.



The Dynabeads[®]-mRNA complex is now ready for reverse transcription and PCR.