



Cell Therapy Systems
Recombinant Human
Tumor Necrosis Factor- α CTS™
(TNF- α)

PRODUCT ANALYSIS SHEET

Catalog Number:	CTP3011	CTP3013
Quantity:	100 μ g	1 mg
Lot Number:	See product label	
Molecular Weight:	17.5 kDa	
Purity:	>95% as determined by SDS-PAGE analysis.	
Amino Acid Sequence:	VRSSSRTPSD KPAHVVANP QAEGQLQWLN RRANALLANG VELRDNQLVV PSEGLYLIYS QVLFKGGQCP STHVLLTHTI SRIAVSYQTK VNLLSAIKSP CQRETPEGAE AKPWYEPIYL GGVFQLEKGD RLSAEINRPD YLDFAESGQV YFGIIL	
Biological Activity:	ED ₅₀ range = 0.02 to 0.05 ng/mL (Specific Activity: 2×10^7 to 5×10^7 units/mg), determined by the dose dependent cytotoxic effect on L929 cells in the presence of actinomycin D.	
Formulation:	Lyophilized, carrier free.	
Sterility:	Filtered prior to lyophilization through a 0.22 micron sterile filter.	
Endotoxin:	<0.1 ng/ μ g	
Production:	Recombinant human TNF- α is produced in <i>E. coli</i> and purified via sequential chromatography.	
Reconstitution Recommendation:	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute lyophilized recombinant human TNF- α in sterile, distilled water to a concentration of 0.1 to 1.0 mg/mL. Further dilutions should be made in low endotoxin medium or a buffered solution containing a carrier protein such as heat inactivated tissue culture grade HSA. It is recommended that all culture media containing supplements, such as growth factor, be sterile filtered prior to use for cell, gene, or tissue-based applications to minimize risk of contamination.	
Suggested Working Dilutions:	The optimal concentration should be determined for each specific application.	
Storage:	Lyophilized human TNF- α should be stored at 2 to 8°C, preferably desiccated. Store reconstituted human TNF- α at $\leq -20^\circ\text{C}$ (not in a frost-free freezer). Keep freeze-thaw cycles to a minimum.	
Expiration Date:	Expires one year from date of receipt when stored as instructed.	
References:	Aggarwal, B.B., W.J. Kohr, P.E. Hass, B. Moffat, S.A. Spencer, W.J. Henzel, T.S. Bringman, G.E. Nedwin, D.V. Goeddel, and R.N. Harkins (1984) Human tumor necrosis factor production, purification and characterization. <i>J. Biol. Chem.</i> 260:2345-2354. Chew, L.J., H.G. Pan, J.Y. Yu, S. Tian, W.Q. Huang, J.Y. Zhang, S. Pang, and L.Y. Li (2002) A novel secreted splice variant of vascular endothelial cell growth inhibitor. <i>FASEB J.</i> 16:742-744. Cignetti, A., E. Bryant, B. Allione, A. Vitale, R. Foa, and M.A. Cheever (1999) CD34(+) acute myeloid and lymphoid leukemic blasts can be induced to differentiate into dendritic cells. <i>Blood</i> 94(6):2048-2055. Kim, S.G., K.A. Soltysiak, Z.G. Gao, T.S. Chang, E.J. Chung, and K.A. Jacobson (2003) Tumor necrosis factor alpha-induced apoptosis in astrocytes is prevented by the activation of P2Y(6), but not P2Y(4) nucleotide receptors. <i>Biochem. Pharmacol.</i> 65(6):923-931.	

For Research Use or Non-Commercial Manufacturing of Cell Based Products for Clinical Research.

CAUTION: Not intended for direct administration into humans or animals

www.invitrogen.com

Manufactured under ISO 13485 Quality Standard

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288

For technical support or support related to CTS™ products, www.invitrogen.com/celltherapysupport

PICTS-Hu TNF- α

(Rev 07/10) DCC-10-1460

Important Licensing Information - These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or our website, www.invitrogen.com). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses.



References (cont'd):

- Kim, S.G., Z.-G. Gao, K.A. Soltysiak, T.-S. Chang, C. Brodie, and K.A. Jacobson (2003) P2Y6 nucleotide receptor activates PKC to protect 1321N1 astrocytoma cells against tumor necrosis factor-induced apoptosis. *Cell. Mol. Neurobiol.* 23(3):401-418.
- Kruszynski, M., D.J. Shealy, A.O. Leone, and G.A. Heavner (1999) Identification of TNF-alpha binding peptides from a D-amino acid hexapeptide library that specifically inhibit TNF-alpha binding to recombinant p55 receptor. *Cytokine* 11(1):37-44.
- Mazanet, M.M., K. Neote, and C.C.W. Hughs (2000) Expression of IFN-inducible T cell chemoattractant by human endothelial cells is cyclosporin A-resistant and promotes T cells adhesion: implications for cyclosporin A-resistant immune inflammation. *J. Immunol.* 164:5383-5388.
- Means, T.K., B.W. Jones, A.B. Schromm, B.A. Shurtleff, J.A. Smith, J. Keane, D.T. Golenbock, S.N. Vogel, and M.J. Fenton (2001) Differential effects of a toll-like receptor antagonist on *Mycobacterium tuberculosis*-induced macrophage responses. *J. Immunol.* 166 6):4074-4082.
- Nixon, C.S., M.J. Steffen, and J.L. Ebersole (2000) Cytokine responses to *Treponema pectinovorum* and *Treponema denticola* in human gingival fibroblasts. *Infection and Immunity* 68(9):5284-5292.
- Ottonello, L., P. Morone, P. Dapino, and F. Dallegri (1996) Monoclonal Lym-1 antibody-dependent lysis of B-lymphoblastoid tumor targets by human complement and cytokine-exposed mononuclear and neutrophilic polymorphonuclear leukocytes. *Blood* 87(12):5171-5178.
- Pati, S., C.B. Pelsler, J. Dufraigne, J.L. Bryant, M.S. Reitz Jr. and F.F. Weichold (2002) Antitumorigenic effects of HIV protease inhibitor ritonavir: inhibition of Kaposi sarcoma. *Blood* 99 (10):3771-3779.
- Qiu, D., G. Zhao, Y. Aoki, L. Shi, A. Uyei, S. Nazarian, J.C. Ng, and P.N. Kao (1999) Immunosuppressant PG490 (triptolide) inhibits T-cell interleukin-2 expression at the level of purine-box/nuclear factor of activated T-cells and NF-kappaB transcriptional activation. *J. Biol. Chem.* 274(19):13443-13450.
- Rajan, R., R. Vanderslice, S. Kapur, J. Lynch, R. Thompson, and D. Djakiew (1996) Epidermal growth factor (EGF) promotes chemomigration of a human prostate tumor cell line, and EGF immunoreactive proteins are present at sites of metastasis in the stroma of lymph nodes and medullary bone. *Prostate* 28(1):1-9.
- Rodriguez, P., M. Heyman, C. Candalh, M.A. Blaton, and C. Bouchaud (1995) Tumour necrosis factor-alpha induces morphological and functional alterations of intestinal HT29 cl.19A cell monolayers. *Cytokine* 7(5):441-448.
- Sciaky, D., W. Brazer, D.M. Center, W.W. Cruikshank, and T.J. Smith (2000) Cultured human fibroblasts express constitutive IL-16 mRNA: Cytokine induction of active IL-16 protein synthesis through a caspase-3-dependent mechanism. *J. Immunol.* 164 (7):3806-3814.
- Shields, D.C., N.G. Avgeropoulos, N.L. Banik, and W.R. Tyor (2000) Acute multiple sclerosis characterized by extensive mononuclear phagocyte infiltration. *Neurochemical Research* 25(11):1517-1520.
- Siegel, S.A., D.J. Shealy, M.T. Nakada, J. Le, D.S. Woulfe, L. Probert, G. Kollias, J. Ghraeyb, J. Vilcek, and P.E. Daddona (1995) The mouse/human chimeric monoclonal antibody cA2 neutralizes TNF in vitro and protects transgenic mice from cachexia and TNF lethality in vivo. *Cytokine* 7(1):15-25.
- Wang, H.S., H.J. Cao, V.D. Winn, L.J. Rezanka, Y. Frobert, C.H. Evans, D. Sciaky, D.A. Young, and T.J. Smith (1996) Leukoregulin induction of prostaglandin-endoperoxide H synthase-2 in human orbital fibroblasts. An in vitro model for connective tissue inflammation. *J. Biol. Chem.* 271(37):22718-22728.
- Yu, Y., M. Hagihara, K. Ando, B. Gansuud, H. Matsuzawa, T. Tsuchiya, Y. Ueda, H. Inoue, T. Hotta, and S. Kato (2001) Enhancement of human cord blood CD34(+) cell-derived NK cell cytotoxicity by dendritic cells. *J. Immunol.* 166(3):1590-1600.

Explanation of symbols

Symbol	Description	Symbol	Description
	Catalogue Number		Batch code
	Research Use Only		In vitro diagnostic medical device
	Use by		Temperature limitation
	Manufacturer		European Community authorised representative
	Without, does not contain		With, contains
	Protect from light		Consult accompanying documents
	Directs the user to consult instructions for use (IFU), accompanying the product.		

For Research Use or Non-Commercial Manufacturing of Cell Based Products for Clinical Research.

CAUTION: Not intended for direct administration into humans or animals

www.invitrogen.com

Manufactured under ISO 13485 Quality Standard

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288

For technical support or support related to CTS™ products, www.invitrogen.com/celltherapysupport

PICTS-Hu TNF- α

(Rev 07/10) DCC-10-1460

Important Licensing Information - These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or our website, www.invitrogen.com). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses.

invitrogen™