



Instruction Manual

PureLink™ miRNA Isolation Kit

For rapid, efficient purification of small RNA molecules

Catalog no. K1570-01

Version C

18 July 2005

25-0753

Table of Contents

Table of Contents.....	iii
Experienced Users Procedure	v
Kit Contents and Storage	vii
Accessory Products	viii
Introduction.....	1
Overview	1
Experimental Overview	4
Methods	5
Preparing Lysates.....	5
Purification Procedure.....	12
Analyzing RNA	15
Expected Results.....	16
Troubleshooting	17
Appendix	19
One-Column Protocol.....	19
Technical Service	23
Purchaser Notification.....	25
References.....	26

Experienced Users Procedure

Introduction

This quick reference sheet is included for experienced users of the PureLink™ miRNA Isolation Kit. If you are a first time user, follow the detailed protocol in this manual.

Step	Action
Sample Preparation	<p>Mammalian Cells (up to 1x 10⁶ cells)</p> <ol style="list-style-type: none">1. Harvest cells and add 300 µl Binding Buffer (L3) to cells.2. Mix well by vortexing and add 300 µl 70% ethanol. Mix well. <p>Mammalian Tissue (up to 5 mg)</p> <ol style="list-style-type: none">1. Add 300 µl Binding Buffer (L3) and homogenize tissue using tissue homogenizer.2. Centrifuge at 12,000 x g for 1 minute to remove any particulate materials. Transfer supernatant to a clean tube and add 300 µl 70% ethanol. Mix well by vortexing. <p>Plant Tissue (up to 100 mg)</p> <ol style="list-style-type: none">1. Grind tissue to a powder in liquid nitrogen and add 300 µl Binding Buffer (L3). Mix well by vortexing.2. Centrifuge the lysate at 12,000 x g for 2 minutes to remove any particulate materials. Transfer supernatant to a clean tube and add 300 µl 70% ethanol. Mix well by vortexing. <p>Yeast Cells (up to 2x 10⁷ cells)</p> <ol style="list-style-type: none">1. Harvest 1 ml fresh, log phase yeast cells (OD₆₆₀ ~ 1.0) by centrifugation and resuspend cell pellet in 70 µl TE Buffer.2. Add 30 units zymolase (lyticase) enzyme to the lysate and incubate for 30 minutes at 30 °C.3. Add 300 µl Binding Buffer (L3) and 210 µl of 96-100% ethanol. Mix well by vortexing. <p>Bacterial Cells (up to 5x 10⁶ cells)</p> <ol style="list-style-type: none">1. Harvest up to 5 x 10⁶ <i>E. coli</i> cells by centrifugation and resuspend cell pellet in 90 µl TE with 1 mg lysozyme. Add 1 µl 5% SDS and mix well.2. Incubate at room temperature for 5 minutes.3. Add 300 µl Binding Buffer (L3) and add 210 µl 96-100% ethanol. Mix well by vortexing. <p>Note: To process larger sample amount, see page 10 for a protocol using the TRIzol® Reagent.</p>

Continued on next page

Experienced Users Procedure, Continued

Step	Action
Purification Procedure	<p>The purification procedure is designed for use with a microcentrifuge capable of centrifuging $>10,000 \times g$ and can be completed in 10-15 minutes</p> <ol style="list-style-type: none">1. Add prepared sample from previous page to a Spin Cartridge in a collection tube.2. Centrifuge the Spin Cartridge at $12,000 \times g$ for 1 minute. Total RNA is bound to the cartridge. Keep the flow through.3. Add 96-100% ethanol to the flow through to a final concentration of 70%. Mix well by vortexing.4. Transfer 700 μl sample to a second Spin Cartridge in collection tube. Centrifuge the Spin Cartridge at $12,000 \times g$ for 1 minute. Small RNA molecules bind to the Spin Cartridge. Discard the flow through.5. Repeat Step 4 for the remaining sample. Place the Spin Cartridge in the collection tube.6. Wash the Spin Cartridge with 500 μl Wash Buffer (W5) with ethanol (page 13). Centrifuge at $12,000 \times g$ for 1 minute. Discard the flow through.7. Repeat wash step with 500 μl Wash Buffer (W5) once.8. Discard the collection tube and place the Spin Cartridge in a Wash Tube supplied with the kit.9. Centrifuge the Spin Cartridge at maximum speed for 1 minute to remove any residual Wash Buffer (W5).10. Place the Spin Cartridge in a clean 1.7-ml Recovery Tube supplied with the kit.11. Elute the RNA with 50-100 μl sterile RNase-free water supplied with the kit (add water to the center of the cartridge).12. Incubate at room temperature for 1 minute.13. Centrifuge the Spin Cartridge at maximum speed for 1 minute to elute RNA. <i>The Recovery Tube contains purified small RNA molecules. Discard the Spin Cartridge.</i>14. Store purified RNA at -80°C or use RNA for the desired downstream application.

Kit Contents and Storage

Shipping and Storage

All components of the PureLink™ miRNA Isolation Kit are shipped at room temperature. Upon receipt, store all components at room temperature.

Contents

The components included in the PureLink™ miRNA Isolation Kit are listed below.

Sufficient reagents are provided in the kit to perform 25 reactions.

Component	Amount
Binding Buffer (L3)	45 ml
Wash Buffer (W5)	10 ml
Buffer (W4) for single column purification	3 ml
Sterile, RNase-free Water	5 ml
Spin Cartridges with Collection Tubes	50
Wash Tubes (2.0 ml)	50
Recovery Tubes (1.7 ml)	50

Product Qualification

The PureLink™ miRNA Isolation Kit is functionally qualified as described below.

Yeast tRNA and 10 bp DNA Ladder are purified using the kit as described in this manual. The purified small RNA is analyzed by gel electrophoresis.

Agarose gel electrophoresis must show the co-purification of yeast tRNA with 10 bp DNA Ladder fragments less than 40 bp and DNA fragments greater than 40 bp must not be detectable.

In addition, each kit component is free of ribonuclease contamination and is lot qualified for optimal performance.

Accessory Products

Additional Products

The following products are also available from Invitrogen. For more details on these products, visit our Web site at www.invitrogen.com or contact Technical Service (page 19).

Product	Quantity	Catalog no.
RNase AWAY®	250 ml	10328-011
UltraPure™ DEPC-treated Water	1 L	750023
UltraPure™ DNase/RNase-Free Distilled Water	500 ml	10977-015
Quant-iT™ RNA Assay Kit	1000 assays	Q-33140
Novex® TBE-Urea Gels, 15% (1.0 mm, 10 well)	1 box	EC6885BOX
Novex® TBE-Urea Gels, 10% (1.0 mm, 10 well)	1 box	EC6875BOX
Novex® TBE-Urea Sample Buffer (2X)	10 ml	LC6876
Novex® TBE Running Buffer (5X)	1 L	LC6675
Invitrogen Homogenizer (clarifier; differs from tissue homogenizer. See page 8)	50	12183-026
10 bp DNA Ladder	50 µg	10821-015
TrackIt™ 10 bp DNA Ladder	20 applications	10488-019
PureLink™ Micro-to-Midi Total RNA Purification System	50 reactions	12183-018
UltraPure™ Glycogen	100 µl	10814-010
TRIzol®	100 ml	15596-026

Introduction

Overview

Introduction

The PureLink™ miRNA Isolation Kit provides a rapid and efficient method to purify small RNA molecules from biological samples for functional analysis.

The kit is specifically designed to isolate high-quality small RNA molecules including tRNA, 5S rRNA, 5.8S rRNA and regulatory RNA molecules such as microRNA (miRNA) and short interfering RNA (siRNA). Using the kit to isolate small RNA molecules results in efficient removal of large RNA molecules such as 28S rRNA, 18S rRNA, and mRNA that can inhibit expression analysis of small RNA molecules. See below for an overview of the purification protocol.

The purified RNA is suitable for northern blot analysis and microarray analysis.

For information on small RNA molecules, see next page.

System Overview

The PureLink™ miRNA Isolation Kit is based on the selective binding of small RNA molecules to silica-based membrane in the presence of 70% ethanol.

The lysate is prepared from mammalian cells and tissue, plant tissues, *E. coli* cells and yeast cells using Binding Buffer (L3) containing guanidine isothiocyanate, a chaotrope capable of protecting the RNA from endogenous RNases. Ethanol is added to the lysate to a final concentration of 35% and the lysate is processed through the first Spin Cartridge. Large RNA molecules bind to the silica-based membrane in the cartridge while small RNA molecules are recovered in the flow through fraction.

Ethanol is added to the flow through fraction to a final concentration of 70% and the sample is loaded onto a second Spin Cartridge. The small RNA molecules bind to the silica-based membrane in the cartridge and impurities are removed by thorough washing with Wash Buffer. The RNA is then eluted in sterile, RNase free water.

Continued on next page

Overview, Continued

Small RNA Molecules

Small regulatory RNA molecules play an important role in regulation of gene expression in various organisms by binding to the target mRNAs through sequence complementation (Liu, 2004). Small RNA molecules include tRNA, 5S rRNA, and 5.8S rRNA and regulatory RNA molecules such as siRNA and miRNA.

siRNAs are ~21-23 nucleotide double-stranded RNA molecules involved in post-transcriptional gene silencing using the RNAi (RNA interference) pathway (Elbashir *et al.*, 2002). miRNAs are ~21-22 nucleotide double stranded RNA molecules that play important roles in the regulation of translation and degradation of target mRNAs through base pairing to partially complementary sites in the untranslated regions of the message (Lim, 2003).

Isolation and purification of small regulatory RNA molecules allows elucidation of biologically significant pathways for regulation of gene expression and requires enrichment of small RNA molecules from biological samples for functional analysis.

The standard protocols for isolating total RNA and mRNA are not optimized for isolation of small RNA molecules and result in the loss of substantial amounts of small RNA. Co-purification of larger RNA molecules such as 28S rRNA, 18S rRNA, and mRNA with small RNA inhibits expression analysis of small RNA.

The PureLink™ miRNA Isolation Kit is specifically designed for purification of small RNA with minimal contamination from large RNA molecules.

Continued on next page

Overview, Continued

Advantages

Using PureLink™ miRNA Isolation Kit to isolate small RNA molecules provides the following advantages:

- Rapid and efficient purification of high-quality small RNA molecules using spin column-based centrifugation
 - Specifically designed to purify small RNA molecules including miRNA and siRNA from a variety of samples in less than 15 minutes
 - Minimal contamination from large RNA molecules and genomic DNA
 - Reliable performance of the purified small RNA in downstream applications
-

Specifications


Binding Capacity:	~1 mg nucleic acid
Starting Material:	Varies (see page 6)
Column Reservoir Capacity:	700 μ l
Wash Tube Capacity:	2.0 ml
Recovery Tube Capacity:	1.7 ml
Centrifuge Compatibility:	Capable of centrifuging at $>10,000 \times g$
RNA Yield:	Varies (see page 16)
RNA Recovery:	$>80\%$

Experimental Overview

Introduction

The flow chart for purifying small RNA molecules using the PureLink™ miRNA Isolation Kit is shown below.

Prepare lysate using
Binding Buffer (L3)



Add ethanol to lysate



Load sample onto the
first spin column



Add ethanol to the flow
through



Load sample onto the
second spin column



Wash column with
Wash Buffer (W5)
twice



Elute miRNA with
water into Elution
Tube

Methods

Preparing Lysates

Introduction

Instructions for preparing lysates from mammalian cells and tissues, plant tissues, yeast, and bacteria are described below.

To obtain high-quality RNA, follow the guidelines recommended below.



The PureLink™ miRNA Isolation Kit buffers contain guanidine isothiocyanate. Always wear a laboratory coat, disposable gloves, and eye protection when handling buffers.

Do not add bleach or acidic solutions directly to solutions containing guanidine isothiocyanate or sample preparation waste as it forms reactive compounds and toxic gases when mixed with bleach or acids.

General Handling of RNA

Observe the following guidelines to prevent RNase contamination:

- Use disposable, individually wrapped, sterile plastic ware
 - Use only sterile, new pipette tips and microcentrifuge tubes
 - Wear latex gloves while handling reagents and RNA samples to prevent RNase contamination from the surface of the skin
 - Always use proper microbiological aseptic techniques when working with RNA
 - Use RNase AWAY® Reagent (page viii) to remove RNase contamination from surfaces
-

Continued on next page

Preparing Lysates, Continued

Sample Amount and Quality

There are different protocols for preparing lysates depending on the starting material (sample). Based on your sample, choose an appropriate lysate preparation protocol from the table below.

To obtain high yield of small RNA molecules and minimize any degradation, collect the sample and proceed immediately to sample preparation or freeze the sample in liquid nitrogen immediately after collection.

To obtain the best results, use the appropriate protocol based on your sample and the recommended sample amount for purification as described in the table below.

If you wish to start with less amount of sample, use the appropriate protocol based on your sample without changing the volume of reagent used. Note that if you start with less amount of sample, the RNA yield may be lower.

Sample	Amount	Page no.
Mammalian cells	up to 1×10^6 cells	7
Mammalian tissues	up to 5 mg	8
Plant tissues	up to 100 mg	8
Yeast cells	2×10^7 cells	9
<i>E. coli</i> cells	5×10^6 cells	9
Dicing reaction	50 μ l	11
Larger sample amount	--	10



Note

To minimize RNA degradation and to obtain optimal RNA yield, it is important to perform complete disruption of tissue in lysis buffer **quickly** and not to exceed the recommended starting amounts for various samples listed in the above table.

Using increased amount of starting material results in co-purification of large molecules such as 18S/28S rRNA, mRNA, and possibly gDNA that did not bind to the column during the first column screening, is retained as flow through, and co-purifies with small RNA molecules.

If you wish to process large amount of starting material, see page 10 for recommended conditions.

Continued on next page

Preparing Lysates, Continued

Materials Needed

- 96-100% ethanol
 - 70% ethanol
 - Sample (see previous page)
 - Binding Buffer (L3) supplied with the kit
 - TE buffer (10 mM Tris-HCl, 1 mM EDTA, pH 8.0)
 - Zymolase (lyticase) enzyme for yeast lysate
 - Lysozyme and 5% SDS for bacterial cell lysate
 - Invitrogen Homogenizer (clarifier; for viscous or debris-containing samples)
 - General tissue homogenizer for homogenizing mammalian and plant tissues
 - TRIzol® Reagent and chloroform for processing large sample amount
 - Sterile, DNase-free microcentrifuge tubes
 - Water bath or heat block set at 30°C for yeast lysate
-

Mammalian Cells Lysate

Procedure to prepare lysate from mammalian cells is described below.

1. For adherent cells (up to 1×10^6 cells), remove the growth medium from the culture plate.
For suspension cells (up to 1×10^6 cells), harvest the cells and centrifuge the cells at $250 \times g$ for 5 minutes to pellet cells. Remove the growth medium.
 2. Resuspend cells from Step 1 in 300 μ l Binding Buffer (L3) supplied with the kit. Mix well by vortexing or pipetting up and down until the cells appear lysed.
 3. Add 300 μ l 70% ethanol to the cell lysate. Mix well by vortexing.
 4. Proceed to the **Binding Step** (page 13).
-

Continued on next page

Preparing Lysates, Continued

Mammalian Tissue Lysate

Procedure to prepare lysate from mammalian tissues is described below.

1. Place ~5 mg of minced mammalian tissue into a sterile microcentrifuge tube.
2. Add 300 μ l Binding Buffer (L3) supplied with the kit. Ensure the tissue is completely immersed in the buffer.
3. Homogenize the tissue using a tissue homogenizer.
4. Centrifuge the lysate at 12,000 \times g for 5 minutes at room temperature to remove any particulate material.
5. Transfer the supernatant to a new, sterile microcentrifuge tube.

Note: If the lysate is viscous or contains cell debris, clarify the lysate using the Homogenizer available from Invitrogen (page viii).

6. Add 300 μ l 70% ethanol to the lysate. Mix well by vortexing.
 7. Proceed to the **Binding Step** (page 13).
-

Plant Tissue Lysate

Procedure to prepare lysate from plant is described below.

1. Freeze the plant tissue in liquid nitrogen and grind the tissue to a powder.
For soft tissue, cut the tissue into small pieces.
2. Add 300 μ l Binding Buffer (L3) supplied with the kit to the tissue. For ground powder, vortex the powder until the powder is completely resuspended. For the soft tissue, homogenize using a tissue homogenizer.
3. Centrifuge the lysate at 12,000 \times g for 2 minutes at room temperature to remove any particulate material.
4. Transfer the supernatant to another sterile microcentrifuge tube.

Note: If the lysate is viscous or contains cell debris, clarify the lysate using a Homogenizer available from Invitrogen (page viii).

5. Add 300 μ l 70% ethanol to the lysate. Mix well by vortexing.
 6. Proceed to the **Binding Step** (page 13).
-

Continued on next page

Preparing Lysates, Continued

Yeast Lysate

Procedure to prepare lysate from yeast cells is described below.

1. Harvest 1 ml fresh, log-phase yeast cells ($OD_{660} = 1.04$) by centrifugation. If you are using a frozen cell pellet, proceed to Step 2.
 2. Resuspend the cell pellet in 70 μ l TE buffer (10 mM Tris-HCl, 1 mM EDTA, pH 8.0).
 3. Add 30 units zymolase (lyticase) enzyme and incubate at 30°C for 30 minutes to lyse the cells.
 4. Centrifuge at 12,000 \times g for 5 minutes at room temperature to remove any particulate material.
 5. Resuspend the pellet in 300 μ l Binding Buffer (L3) supplied with the kit.
 6. Add 210 μ l 96-100% ethanol to the lysate. Mix well by vortexing.
 7. Proceed to **Binding Step** (page 13).
-

E. coli Lysate

Procedure to prepare *E. coli* cell lysate is described below.

1. Prepare 90 μ l TE buffer (10 mM Tris-HCl, 1 mM EDTA, pH 8.0) containing 1 mg lysozyme enzyme. Store on ice until use.
 2. Harvest up to 5×10^6 *E. coli* cells by centrifugation. If you are using a frozen cell pellet, proceed to Step 3.
 3. Resuspend the cell pellet in 90 μ l TE buffer containing 1 mg lysozyme from Step 1.
 4. Add 1 μ l 5% SDS to the lysate and mix well by vortexing. Incubate at room temperature for 5 minutes.
 5. Add 300 μ l Binding Buffer (L3) supplied with the kit to the lysate. Mix well by vortexing.
 6. Add 210 μ l 96-100% ethanol to the lysate. Mix well by vortexing.
 7. Proceed to the **Binding Step** (page 13).
-

Continued on next page

Preparing Lysates, Continued

Large Sample Amount

If you wish to process large sample amount as compared to the amount listed on page 6, use TRIzol[®] Reagent to prepare lysate as described below.

Lysate Preparation

Tissues

Homogenize tissue samples in 1 ml TRIzol[®] Reagent per 50-100 mg tissue using a tissue homogenizer.

Adherent Cells

Lyse cells directly in a culture dish by adding 1 ml TRIzol[®] Reagent to a culture dish and passing the cell lysate several times through a pipette. The amount of TRIzol[®] Reagent required is based on the culture dish area (1 ml per 10 cm²) and not on the number of cells present.

Suspension Cells

Harvest cells and pellet cells by centrifugation. Use 1 ml of the TRIzol[®] Reagent per 5-10 × 10⁶ animal, plant or yeast cells, or per 1 × 10⁷ bacterial cells and lyse cells by repetitive pipetting.

Phase Separation

1. Incubate the lysate with TRIzol[®] Reagent at room temperature for 5 minutes to allow complete dissociation of nucleoprotein complexes.
2. Add 0.2 ml chloroform per 1 ml TRIzol[®] Reagent and shake the tube vigorously by hand for 15 seconds. Avoid vortexing the sample.
3. Incubate at room temperature for 2-3 minutes.
4. Centrifuge the sample at 12,000 × g for 15 minutes at 4°C. After centrifugation, the mixture separates into a lower red phenol-chloroform phase, an interphase, and a colorless upper aqueous phase. The volume of aqueous phase is ~600 µl.
5. Transfer ~ 400 µl of the colorless, upper phase containing RNA to a fresh tube.
6. Add 96-100% ethanol into 400 µl of the aqueous phase to obtain a final concentration to 35% (e.g. add 215 µl of 96-100% ethanol) and mix well by vortexing.
7. Proceed to the **Binding Step** (page 13).

Continued on next page

Preparing Lysates, Continued

Dicing Reaction

If you are purifying d-siRNA produced in the dicing reaction using the BLOCK-iT™ Dicer RNAi Kit (available from Invitrogen) or an equivalent kit, use the protocol below to prepare samples prior to purification.

To 300 µl dicing reaction, add 300 µl Binding Buffer (L3) and 300 µl 96-100% ethanol to obtain a final volume of 900 µl. Mix well by vortexing. Proceed to the **Binding Step**, page 13.

Note: You can also use 300 µl isopropanol instead of ethanol as suggested in the BLOCK-iT™ Dicer RNAi Kit manual for preparing the dicing reaction for purifying d-siRNA. For more details on the BLOCK-iT™ Dicer RNAi Kit, visit www.invitrogen.com or contact Technical Service (page 23).

Purification Procedure

Introduction

The purification procedure is designed for purifying small RNA molecules using spin column based centrifugation in a total time of **10-15 minutes**.

Materials Needed

- 96-100% ethanol
- Microcentrifuge capable of centrifuging $>10,000 \times g$

Components supplied with the kit

- Wash Buffer (W5)
 - Sterile, RNase-free Water (pH $>$ 7.0)
 - Spin Cartridge in Collection Tubes
 - Wash Tubes and Recovery Tubes
-



Follow the recommendations below to obtain the best results:

- Perform all centrifugation steps at room temperature
 - Be sure to add ethanol to a final concentration of 70% and transfer sample to a second column to isolate small RNA molecules
 - Perform the recommended wash steps to obtain high-quality RNA
 - Pipet water in the center of the Spin Cartridge and perform a 1 minute incubation
-



Note

The purification procedure described in this section is designed for purifying small RNA molecules from tissues, cells, and dicing reaction using the two-column protocol.

If you are purifying small RNA molecules from a dicing reaction, an *optional* one-column protocol followed by ethanol precipitation of RNA is described on page 19.

If you wish to purify total RNA from the same sample, the total RNA is bound to the first spin column and can be washed and eluted using the reagents and protocol described in the PureLink™ Micro-to-Midi Total RNA Purification Kit (see page viii for ordering information). The manual can be downloaded from www.invitrogen.com or by contacting Technical Service (page 23).

Continued on next page

Purification Procedure, Continued

Elution Buffer For two-column purification protocol, use sterile, RNase-free **water** supplied with the kit to elute small RNA molecules (page 14).

For one-column purification protocol, use **Buffer (W4)** for single column purification supplied with the kit to elute small RNA molecules (page 21).

Before Starting Add 40 ml 96-100% ethanol to 10 ml Wash Buffer (W5) included with the kit. Store the Wash Buffer (W5) with ethanol at room temperature.

- Binding Step**
1. Remove a Spin Cartridge in a Collection Tube from the package.
 2. Add the lysate with Binding Buffer (L3) and ethanol (total volume is 600 μ l) prepared as described on pages 5-11 to the Spin Cartridge.
 3. Centrifuge the Spin Cartridge 12,000 \times g for 1 minute at room temperature.
 4. **Keep the flow through.** Add 96-100% ethanol to the flow through to obtain a final concentration of 70% ethanol and mix well by vortexing. **Note:** The lysate already contains 35% ethanol added during sample preparation.
For example: To 600 μ l lysate, add 700 μ l 96-100% ethanol and mix well by vortexing.
Note: To isolate total RNA from the same sample, place the Spin Cartridge into a new Wash Tube and process the Spin Cartridge for isolating total RNA as described in the PureLink™ Micro-to-Midi Total RNA Purification Kit (page viii).
 5. Remove **another** Spin Cartridge in a Collection Tube from the package and transfer 700 μ l sample from Step 4 to the Spin Cartridge.
-

Continued on next page

Purification Procedure, Continued

Binding Step, continued

6. Centrifuge the Spin Cartridge at $12,000 \times g$ for 1 minute at room temperature.
 7. Transfer remaining sample from Step 4 to the Spin Cartridge from Step 6 and centrifuge at $12,000 \times g$ for 1 minute at room temperature.
 8. Discard the flow through and place the Spin Cartridge into the collection tube.
 9. Proceed to the **Washing Step**, below.
-

Washing Step

1. Add 500 μ l Wash Buffer (W5) with ethanol (page 13) to the Spin Cartridge.
 2. Centrifuge $12,000 \times g$ for 1 minute at room temperature.
 3. **Repeat** the wash step with 500 μ l Wash Buffer (W5) with ethanol.
 4. Discard the flow through and place the Spin Cartridge into the Wash Tube supplied with the kit.
 5. Centrifuge the Spin Cartridge at maximum speed for 2-3 minutes at room temperature to remove any residual Wash Buffer. Discard the Wash Tube.
 6. Proceed to the **Elution Step**, below.
-

Elution Step

1. Place the Spin Cartridge in a clean 1.7-ml Recovery Tube supplied with the kit.
 2. Add 50-100 μ l of Sterile, RNase-free water (pH >7.0) to the center of the Spin Cartridge.
 3. Incubate at room temperature for 1 minute.
 4. Centrifuge the Spin Cartridge at maximum speed for 1 minute at room temperature.
 5. The Recovery Tube contains purified small RNA molecules. Remove and discard the cartridge.
Based on the volume of elution buffer used for elution, the recovery of the elution volume will vary and is usually 90% of the elution buffer volume used.
 6. Store small RNA at -80°C or use small RNA for the desired downstream application. You may determine the quality and quantity of the purified small RNA molecules as described on the next page.
-

Continued on next page

Analyzing RNA

RNA Yield

The quantity of the purified small RNA molecules is easily quantitated using UV absorbance at 260 nm or Quant-iT™ RNA Assay Kit.

UV Absorbance

1. Dilute an aliquot of the small sample in 10 mM Tris-HCl, pH 7.0. Mix well. Transfer to a cuvette (1-cm path length).

Note: The RNA must be in a neutral pH buffer to accurately measure the UV absorbance.

2. Determine the OD₂₆₀ of the solution using a spectrophotometer blanked against 10 mM Tris-HCl, pH 7.0.

Calculate the amount of total RNA using the following formula:

$$\text{Total RNA } (\mu\text{g}) = \text{OD}_{260} \times 40 \mu\text{g} / (1 \text{ OD}_{260} \times 1 \text{ ml}) \times \text{dilution factor} \times \text{total sample volume (ml)}$$

Quant-iT™ RNA Assay Kits

The Quant-iT™ RNA Assay Kit (page viii for ordering information) provides a rapid, sensitive, and specific method for RNA quantitation with minimal interference from DNA, protein, or other common contaminants that affect UV absorbance readings.

The kit contains a state-of-the-art quantitation reagent and pre-diluted standards for standard curve. The assay is performed in a microtiter plate format and is designed for reading in standard fluorescent microplate readers.

RNA Quality

Typically, small RNA molecules isolated using the PureLink™ miRNA Isolation Kit have an OD_{260/280} of >1.8 when samples are diluted in Tris-HCl (pH 7.5). An OD_{260/280} of >1.8 indicates that RNA is reasonably clean of proteins and other UV chromophores (heme, chlorophyll, etc.) that could either interfere with downstream applications or negatively affect the stability of the stored RNA.


To visualize small RNA molecules, perform denaturing gel electrophoresis using 10% or 15% TBE acrylamide gels (page viii for ordering information). An example of denaturing TBE gel analysis is shown on the next page.

Expected Results

Results

An example of small RNA molecules isolated from various samples is shown in the figure below.

Small RNA molecules from various samples were isolated using the PureLink™ miRNA Isolation Kit as described in this manual. Samples (10 µl eluate) were analyzed on a Novex® TBE-Urea 15% Gel and RNA bands were visualized with ethidium bromide staining after electrophoresis.



Expected Yields

The yield of small RNA molecules obtained from various samples using the PureLink™ miRNA Isolation Kit is listed below. The RNA quantitation was performed with the Quant-iT™ RNA Assay Kit (page viii).

Material	Amount	Yield (µg)
HeLa cells	1x10 ⁶	1.29
293F cells	1x10 ⁶	1.95
Mouse liver	5 mg	0.71
Rat Spleen	5 mg	1.39
Spinach	60 mg	1.64
Yeast (<i>S. cerevisiae</i>)	1x10 ⁶	6.21
Bacteria (<i>E. coli</i>)	2x10 ⁶	0.55

Troubleshooting

Problem	Cause	Solution
Low RNA yield	Incomplete lysis or too much sample has clogged the filter	Decrease the amount of starting material used. For tissues, cut the tissue into smaller pieces and ensure the tissue is completely immersed in the Binding Buffer (L3) to obtain optimal lysis. Decrease the sample volume used, if cartridge is clogged or load the sample on 2 spin cartridges.
	Low RNA content	Various tissues have different small RNA content and some tissue may not contain any small RNA at detectable levels.
	Flow through from first cartridge discarded	Do not discard the flow through from the first cartridge. The flow through contains small RNA molecules.
	Incorrect binding conditions	For efficient binding of small RNA molecules, always add ethanol to the flow through to a final concentration of 70%.
	Ethanol not added to Wash Buffer (W5)	Be sure to add 96–100% ethanol to Wash Buffer (W5) as described on page 13.
	Incorrect elution conditions	Add water to the center of the cartridge and perform incubation for 1 minute with water before centrifugation. Be sure to use RNase-free water for elution using the two-column protocol and Buffer (W4) for single-column protocol.
	RNA quantitation performed with water	Be sure the RNA quantitation using UV absorbance is performed with 10 mM Tris-HCl, pH 7.0.(page 15) to accurately measure the UV absorbance.

Continued on next page

Troubleshooting, Continued

Problem	Cause	Solution
RNA degraded	RNA contaminated with RNase	Follow the guidelines on page 5 to prevent RNase contamination.
	Poor quality starting materials	Always use fresh samples or samples frozen at -80°C. For lysis, process the sample quickly to avoid degradation.
Genomic DNA contamination	Large sample amount loaded on the first column	Be sure to use the protocol with TRIzol® Reagent when processing large sample amount as described on page 10. Perform DNase I digestion on the RNA sample after elution to remove genomic DNA contamination.
Total RNA contamination	Processed large amount of starting material without using TRIzol® Reagent	Processing large amount of starting material using the standard protocols without TRIzol® Reagent will result in co-purification of large molecules such as 18S/28S rRNA and mRNA that did not bind to the first column. If you wish to process large amount of starting material, use TRIzol® Reagent to prepare samples as described on page 10.

Appendix

One-Column Protocol

Introduction

An optional one-column purification protocol is described below to purify siRNA from dicing reaction only.



Important

The one-column purification protocol is only applicable to samples obtained from dicing reaction.

The one-column protocol is not recommended with samples prepared from tissues or cells as the amount of sample used for purification may clog the spin cartridges, and the volume and amount of buffers used for wash and elution steps are not sufficient to produce high-quality RNA.

Experimental Outline

1. Add Binding Buffer (L3) and ethanol to the dicing reaction to adjust the conditions for binding.
 2. Load the sample onto the Spin Cartridge to bind small RNA molecules.
 3. Wash the Spin Cartridge with Wash Buffer (W5) containing ethanol to remove impurities.
 4. Elute the siRNA with diluted Buffer (W4) containing ethanol and RNase-free water for single column purification.
 5. Perform ethanol precipitation of siRNA from the eluate.
-

Materials Needed

- 96-100% ethanol (cold)
- 20 µg/µl Glycogen (page viii)
- 70% ethanol (cold)

Components supplied with the kit

- Binding Buffer (L3)
 - Wash Buffer (W5)
 - Buffer (W4) for single column purification
 - Spin Cartridge with Collection Tube
 - Wash Tubes and Recovery Tubes
-

Continued on next page

One-Column Protocol, Continued

Before Starting

Dilute Wash Buffer (W4) and (W5) as described below before starting the One-Column purification procedure.

- **Wash Buffer (W4):** Add 1.5 ml of RNase-free water and 1.5 ml of 96-100% ethanol to the 3 ml Wash Buffer (W4) included with the kit. The diluted buffer (W4) will be used as Elution buffer in the One-Column protocol.
 - **Wash Buffer (W5):** Add 40 ml 96-100% ethanol to 10 ml Wash Buffer (W5) included with the kit.
 - Store the prepared Wash Buffer (W4) and (W5) with ethanol at room temperature.
-

Starting Material

The **One-Column Procedure** is designed for use with **up to 8 µg of RNA in 50 µl** dicing reaction containing as starting material.

- If you wish to process **<8 µg RNA in <50 µl dicing reaction**, increase the volume of your starting material to 50 µl.
- If you wish to process **up to 8 µg in >50 µl dicing reaction**, adjust the buffer volumes of the procedure below accordingly. When loading the reaction mix with the increased volume to the column (Step 4 of the procedure, next page), perform multiple loading of the sample mix onto the Spin Cartridge (Step 4 and 5 of the procedure, next page). The reservoir capacity of the Spin Cartridge is ~700 µl.

Note: Do not load more than 8 µg of RNA from a dicing reaction to the column.

Continued on next page

One-Column Protocol, Continued

One-Column Procedure

Follow the protocol below to purify small RNA from a dicing reaction.

1. Add 150 μ l Binding Buffer (L3) to 50 μ l dicing reaction volume and mix well to obtain a final volume of 200 μ l.
2. Add 96-100% ethanol to the sample to obtain a final ethanol concentration of 75%.
For example: To 200 μ l of sample from Step 1, add 600 μ l 96-100% ethanol to obtain a final volume of 800 μ l.
3. Mix well.
4. Load sample mix onto the Spin Cartridge in a Collection Tube.
Note: The reservoir capacity of the Spin Cartridge is ~700 μ l. If the volume of your sample mix, after performing Step 3, exceeds 700 μ l, perform multiple loading of the sample mix onto the Spin Cartridge (Step 4 and 5).
5. Centrifuge the Spin Cartridge at 12,000 \times g for 1 minute at room temperature. Discard the flow through.
6. Wash the Spin Cartridge with 500 μ l Wash Buffer (W5) with ethanol. Centrifuge the Spin Cartridge at 12,000 \times g for 1 minute at room temperature.
7. Repeat wash step with Wash Buffer (W5) once.
8. Discard flow collection tube. Place the spin cartridge in the Wash Tube supplied with the kit.
9. Centrifuge the Spin Cartridge at 12,000 \times g for 1 minute to remove any residual Wash Buffer (W5). Discard the Wash Tube.
10. Place the Spin Cartridge in the Recovery Tube supplied with the kit.
11. Add 100 μ l Buffer (W4) for single-column purification supplied with the kit.
12. Incubate at room temperature for 1 minute.
13. Centrifuge the Spin Cartridge at 12,000 \times g for 1 minute at room temperature. *The eluate contains siRNA.* Discard the cartridge.
14. Proceed to ethanol precipitation of siRNA, from the eluate, next page.

Continued on next page

One-Column Protocol, Continued

Ethanol Precipitation

Perform ethanol precipitation of siRNA from the eluate, as described below.

1. Add 200 μl cold 96-100% ethanol and 1 μl glycogen solution (20 $\mu\text{g}/\mu\text{l}$) to the eluate from Step 12, previous page.
 2. Mix well and incubate at -20°C for 15 minutes.
 3. Centrifuge the tube in a microcentrifuge at maximum speed for 15 minutes at 4°C .
 4. Discard the supernatant carefully without disturbing the pellet and wash the pellet by resuspending the pellet in 0.5 ml cold 70% ethanol.
 5. Centrifuge the tube in a microcentrifuge at maximum speed for 10 minutes at 4°C .
 6. Discard the supernatant carefully without disturbing the pellet.
 7. Air-dry the pellet for ~ 5 minutes.
 8. Resuspend the dried pellet in 50 μl RNase-free water supplied with the kit. The pellet contains siRNA.
 9. Store the siRNA at -80°C or use siRNA for the desired downstream application.
-

Technical Service

World Wide Web



Visit the Invitrogen Web Resource using your World Wide Web browser. At the site, you can:

- Get the scoop on our hot new products and special product offers
- View and download vector maps and sequences
- Download manuals in Adobe® Acrobat® (PDF) format
- Explore our catalog with full color graphics
- Obtain citations for Invitrogen products
- Request catalog and product literature

Once connected to the Internet, launch your Web browser (Internet Explorer 5.0 or newer or Netscape 4.0 or newer), then enter the following location (or URL):

<http://www.invitrogen.com>

...and the program will connect directly. Click on underlined text or outlined graphics to explore. Don't forget to put a bookmark at our site for easy reference!

Contact Us

For more information or technical assistance, please call, write, fax, or email. Additional international offices are listed on our web page (www.invitrogen.com).

Corporate Headquarters:

Invitrogen Corporation
1600 Faraday Avenue
Carlsbad, CA 92008 USA
Tel: 1 760 603 7200
Tel (Toll Free): 1 800 955 6288
Fax: 1 760 602 6500
E-mail: tech_service@invitrogen.com

European Headquarters:

Invitrogen Ltd
Inchinnan Business Park
3 Fountain Drive
Paisley PA4 9RF, UK
Tel: +44 (0) 141 814 6100
Tech Fax: +44 (0) 141 814 6117
E-mail: eurotech@invitrogen.com

MSDS Requests

To request an MSDS, visit our Web site at www.invitrogen.com. On the home page, go to 'Technical Resources' and select 'Search Technical Resources'. Use the product name or catalog no. and search all technical resources or check the box and search just MSDSs.

Continued on next page

Technical Service, Continued

Limited Warranty

Invitrogen is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, please contact our Technical Service Representatives. Invitrogen warrants that all of its products will perform according to the specifications stated on the certificate of analysis. The company will replace, free of charge, any product that does not meet those specifications. This warranty limits Invitrogen Corporation's liability only to the cost of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. Invitrogen reserves the right to select the method(s) used to analyze a product unless Invitrogen agrees to a specified method in writing prior to acceptance of the order. Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore Invitrogen makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, please report it to our Technical Service Representatives. **Invitrogen assumes no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.**

Purchaser Notification

Limited Use Label License No. 171: RNAi Purification

The purchase of this product conveys to the buyer the non-transferable right to use the purchased product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of patents claimed in this product and owned by Invitrogen based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Invitrogen Corporation, 1600 Faraday Avenue, Carlsbad, California 92008. Phone (760) 603-7200. Fax (760) 602-6500.

References

- Elbashir, S. M., Harborth, J., Weber, K., and Tuschl, T. (2002). Analysis of Gene Function in Somatic Mammalian Cells Using Small Interfering RNAs. *Methods* 26, 199-213.
- Lim, L. P., Glasner, M. E., Yekta, S., Burge, C. B., Bartel, D. P. (2003). Vertebrate microRNA Genes. *Science* 299, 1540.
- Liu, C. G., Calin, G. A., Meloon, B., Gamliel, N., Sevignani, C., Ferracin, M., Dumitru, C. D., Shimizu, M., Zupo, S., Dono, M., Alder, H., Bullrich, F., Negrini, M., and Croce, C. M. (2004). An Oligonucleotide Microchip for Genome-wide microRNA Profiling in Human and Mouse Tissues. *Proc. Natl. Acad. Sci. USA* 101, 9740-9744.

©2004-2005 Invitrogen Corporation. All rights reserved.

For research use only. Not intended for any animal or human therapeutic or diagnostic use. RNase AWAY® and TRIzol® are registered trademarks of Molecular Bio-Products, Inc.



Corporate Headquarters:

Invitrogen Corporation
1600 Faraday Avenue
Carlsbad, California 92008
Tel: 1 760 603 7200
Tel (Toll Free): 1 800 955 6288
Fax: 1 760 602 6500
Email: tech_service@invitrogen.com

European Headquarters:

Invitrogen Ltd
Inchinnan Business Park
3 Fountain Drive
Paisley PA4 9RF, UK
Tel (Free Phone Orders): 0800 269 210
Tel (General Enquiries): +44 (0) 141 814 6100
Fax: +44 (0) 141 814 6260
Email: eurotech@invitrogen.com

International Offices:

Argentina 5411 4556 0844
Australia 1 800 331 627
Brazil 55 11 5051 7422
Canada 800 263 6236
China 10 6849 2578
Hong Kong 2407 8450
Japan 03 3663 7974
New Zealand 0800 600 200
Singapore 65 686 186 38
Taiwan 2 2651 6156

For other countries see our website

www.invitrogen.com

