Citrate Buffer, pH 6.0

For Heat Induced Epitope Retrieval (HIER)

<u>CAT. NO.</u>	<u>Form</u>	<u>Quantity</u>
00-5000	20X Concentrate	100 mL

INTENDED USE

For In Vitro Diagnostic Use. CAUTION: Not for human or animal therapeutic use. Uses other than the labeled intended use may be a violation of local law.

BACKGROUND

Citrate Buffer solution is used in Heat Induced Epitope Retrieval (HIER) methods. Certain antibodies require HIER of formalin-fixed paraffin-embedded tissue⁽¹⁾ and cytological preparations⁽²⁾ to 'retrieve' immunoreactivity that has been compromised by processing (see Figure 1). Our laboratory results show enhanced immunostaining when Citrate Buffer, pH 6.0 is used with most antibodies. HIER is used to reverse the loss of antigenicity that occurs with some epitopes in formalin-fixed, paraffin-embedded tissues. We have found that some antibodies require HIER, and this procedure also enhances immunostaining with many antibodies. However, immunoreactivity with some antibodies is decreased with HIER—always refer to the primary antibody product insert. During HIER, tissue slides are immersed in Citrate Buffer solution and boiled (>95 °C) by hot plate or other appropriate heat source. Refer to the *Procedure For Use With HIER* section in this product insert for the complete method.

STORAGE AND SHELF LIFE

Store at 2-8°C. Do not freeze

STABILITY

There are no obvious signs to indicate instability. It has been quality controlled to assure consistent and reliable performance. Do not use after the expiration date stamped on container. With proper storage there is no significant loss of performance. If this buffer is stored under any conditions other than those specified, it must be validated by the user.

REAGENT PREPARATION

Dilute Citrate Buffer solution 1:20 in deionized or distilled water before using.

REAGENTS AND MATERIALS REQUIRED BUT NOT PROVIDED

- Hot Plate
- Slide holder
- 1 Liter Glass Beaker (Pyrex)
- PBS (Invitrogen® Cat. No. 00-3000)
- HistoGrip™ (Invitrogen® Cat. No. 00-8050)
- Deionized or distilled water

HIER PROCEDURE FOR FORMALIN FIXED, PARAFFIN EMBEDDED TISSUE SECTIONS

1. Adhere tissue to slide using HistoGrip™ (Invitrogen® Cat. No. 00-8050) silane or poly-L-lysine.
2. Deparaffinize in xylene and rehydrate tissue in graded alcohols.
3. Rinse in deionized water.
4. Rinse in PBS.
5. Block endogenous peroxidase with 0.5% hydrogen peroxide in for 10 minutes.
6. Rinse in PBS.
7. Put the slides in a slide rack and place in a 1L glass beaker (Pyrex) containing 500 mL of working solution of citrate buffer.
8. Place beaker on hot plate. Heat the solution until it boils and keep it boiling for 15 minutes. Note: alternative heat sources such as a pressure cooker can be used with an appropriate procedure.
9. After heating, remove beaker with slides from the hot plate and allow it to cool for 25 minutes.
10. Rinse slides with PBS.
11. If required, do avidin/biotin blocking. Note: HIER may enhance endogenous avidin/biotin activity.
12. If required, do general protein blocking step.
13. Start immunostaining protocol.

www.invitrogen.com

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288

E-mail: techsupport@invitrogen.com

PI005000

(Rev 12/09) DCC-09-2058

© 2007 Invitrogen Corporation. All right reserved. These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or www.invitrogen.com). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses.

TROUBLESHOOTING

Possible causes for negative staining on positive slides:

1. Steps in the staining protocol were performed in incorrect sequence.
2. Primary or secondary antibody incubation steps were omitted.
3. Labile antigens were destroyed.
4. Specimen was improperly fixed and/or processed.
5. Specimen dehydrated during staining.

Possible causes for weak staining on all slides:

1. Specimen retained excess liquid after rinsing steps.
2. Incubation times were insufficient.
3. Substrate prepared improperly.
4. Deparaffinization was incomplete (staining may be accompanied by high background).

Possible causes for high background staining:

1. Endogenous peroxidase activity was incompletely blocked.
2. Deparaffinization was incomplete.
3. Excessive application of tissue adhesive.
4. Inadequate rinsing of slides.
5. Over-development of substrate.
6. Dehydration of specimen during staining.

REFERENCES

1. Norton A.J., et al. *J Pathol* 173 (4): 371-391, 1994.
2. Suthipintawong C, et al. *Diagn Cytopathol* 15: 167-174, 1996.

RELATED BUFFERS

	Concentration	Size	Invitrogen® Cat. No.	Recommended Antibody Use*
Citrate Buffer, pH 4.0	20X	100 mL	00-5100	Most cytokeratins, human immunoglobulins
Citrate Buffer, pH 7.1	20X	100 mL	00-5200	HER2, ER, PR, Ki-67, Bcl-2, Bcl-X _L , PCNA, Cyclin D1, p53, CD61, EGFr
EDTA	20X	100 mL	00-5500	CD3, CD4, CD8

TRADEMARKS

CAS-Block™, Clearmount™, HistoGrip™, Histomount™, and Peroxo-Block™ are trademarks of Zymed Laboratories, Inc. Zymed® and Histostain® are registered trademarks of Zymed Laboratories, Inc.

Authorized Representative for IVDD 98/79/EC: Invitrogen Ltd.
 Inchinnan Business Park
 3 Fountain Drive
 Paisley
 PA4 9RF
 UK

www.invitrogen.com

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288

E-mail: techsupport@invitrogen.com

PI005000

(Rev 12/09) DCC-09-2058

© 2007 Invitrogen Corporation. All right reserved. These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or www.invitrogen.com). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses.

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