

## Amplite™ Fluorimetric Aldehyde Quantitation Kit

### \*Blue Fluorescence\*

Ordering Information	Storage Conditions	Instrument Platform
Product Number: 10052 (200 assays)	Keep at -20 °C Avoid exposure to light	Fluorescence microplate readers

### Introduction

The formation, reactivity and toxicity of aldehydes originating from the peroxidation of lipids of cellular membranes have received great attention in recent years. Rapid and accurate measurement of aldehydes is an important task for biological research, chemical research, food industry and environmental pollution surveillance. There are a few reagents or assay kits available for quantifying the number of aldehydes. Most of the existing aldehyde test methods are based on separations either by the tedious and expensive HPLC-MS or GC-MS.

Both Amplite™ Colorimetric Aldehyde Quantitation Kit (10051) and Amplite™ Fluorimetric Aldehyde Quantitation kit (10052) are used for quantifying aldehydes at higher pH. Kit 10052 uses a proprietary fluorogenic dye that generates a strongly fluorescent product upon reacting with an aldehyde. Kit 10052 is much more sensitive than Kit 10051. This fluorimetric kit provides a sensitive mix-and-read method to detect as little as 0.3 nanomole of aldehyde in a 100 µL assay volume (3 µM). The assay can be performed in a convenient 96-well or 384-well microtiter-plate format and easily adapted to automation without a separation step. Its signal can be read by a fluorescence microplate reader at Ex/Em = 365/435 nm.

### Kit Key Features

<b>Broad Application:</b>	Used for quantifying aldehydes in a variety of applications, such as enzyme reactions.
<b>Sensitive:</b>	Detect as little as 0.3 nanomole of aldehyde in a 100 µL assay volume.
<b>Continuous:</b>	Easily adapted to automation without a separation step.
<b>Convenient:</b>	Formulated to have minimal hands-on time.
<b>Non-Radioactive:</b>	No special requirements for waste treatment.

### Kit Components

Components	Amount
Component A: AldeLight™ Blue	1 vial
Component B: Assay Buffer	1 bottle (30 mL)
Component C: Reaction Buffer	1 vial (6 mL)
Component D: Aldehyde Standard	1 vial
Component E: DMSO	1 vial (100 µL)

### Assay Protocol for One 96-Well Plate

#### Brief Summary


**Prepare enzyme reaction solution (50 µL) → Add AldeLight™ Blue reaction mixture (50 µL)**  
**→ Incubate at RT for 15 to 30 minutes → Add 25 µL of Reaction Buffer**  
**→ Monitor fluorescence increase at Ex/Em = 365/435 nm**

*Note: Thaw all the kit components to room temperature before starting the experiment.*



## Data Analysis

The fluorescence in blank wells (0  $\mu\text{M}$  Aldehyde Standard and AldeLight™ Blue reaction mixture only) is used as a control, and subtracted from the values of those wells with the aldehyde reactions. An aldehyde standard curve is shown in Figure 1.



**Figure 1.** Aldehyde dose response was measured in a solid black 96-well plate with Amplite™ Fluorimetric Aldehyde Quantitation Kit using a Gemini fluorescence microplate reader (Molecular Devices). As low as 3  $\mu\text{M}$  of aldehyde can be detected with 15 minutes incubation ( $n=3$ ). *Note: The fluorescence background increases with time, thus it is important to subtract the fluorescence intensity value of the blank wells for each data point.*

## References

1. Trevor M. Kitson. (1985) High concentrations of aldehydes slow the reaction of cytoplasmic aldehyde dehydrogenase with thiol-group modifiers *Biochem. J.* 228, 765.
2. Crabb DW, Matsumoto M, Chang D, You M (2004). Overview of the role of alcohol dehydrogenase and aldehyde dehydrogenase and their variants in the genesis of alcohol-related pathology. *The Proceedings of the Nutrition Society* 63 (1): 49.
3. Steinmetz CG, Xie P, Weiner H, Hurley TD (1997). Structure of mitochondrial aldehyde dehydrogenase: the genetic component of ethanol aversion. *Structure* 5 (5): 701.
4. O'Donnell JM, Kudej RK, LaNoue KF, Vatner SF, Lewandowski ED. (2004) Limited transfer of cytosolic NADH into mitochondria at high cardiac workload. *Am J Physiol Heart Circ Physiol*, 286, H2237.
5. Zurek G, and karst U (2000). 2,4-Dinitro-3,5,6-trideuterophenylhydrazones for the quantitation of aldehydes and ketones in air samples by liquid chromatography-mass spectrometry. *J of chromatography A*, 869, 251.
6. Ou Z, Ogamo A, Guo L, Konda Y, Harigaya Y, and Nakagawa Y. (1995). Identification and quantitation of choline glycerophospholipids that contain aldehyde residues by fluometric high-performance liquid chromatography. *Analytical biochemistry* 227, 289.

**Warning: This kit is only sold to end users. Neither resale nor transfer to a third party is allowed without written permission from AAT Bioquest®. Chemical analysis of the kit components is strictly prohibited. Please call us at 408-733-1055 or e-mail us at [info@aatbio.com](mailto:info@aatbio.com) if you have any questions.**