

Amplite™ Fluorimetric Hydrogen Peroxide Assay Kit

Near Infrared Fluorescence

Ordering Information	Storage Conditions	Instrument Platform
Product Number: 11502 (500 assays)	Keep in freezer Avoid exposure to light	Fluorescence microplate readers

Introduction

Hydrogen peroxide (H₂O₂) is a reactive oxygen metabolic by-product that serves as a key regulator for a number of oxidative stress-related states. It is involved in a number of biological events that have been linked to asthma, atherosclerosis, diabetic vasculopathy, osteoporosis, a number of neurodegenerative diseases and Down's syndrome. Perhaps the most intriguing aspect of H₂O₂ biology is the recent report that antibodies have the capacity to convert molecular oxygen into hydrogen peroxide to contribute to the normal recognition and destruction processes of the immune system. Measurement of this reactive species will help to determine how oxidative stress modulates a variety of intracellular pathways.

The Amplite™ Fluorimetric Hydrogen Peroxide Assay Kit uses our unique Amplite™ IR peroxidase substrate to quantify hydrogen peroxide in solutions and cell extracts. Amplite™ IR generates the fluorescence that is pH-independent from pH 4 to 10. It is a superior alternative to ADHP (Amplex Red™) for the detections that require low pH where ADHP has reduced fluorescence. In addition, Amplite™ IR generates a product that has maximum absorption at 647 nm with maximum emission at 670 nm. This near infrared fluorescence minimizes the assay background that is often caused by the autofluorescence of biological samples. It can also be used to detect a variety of oxidase activities through enzyme-coupled reactions.

This Amplite™ Fluorimetric Hydrogen Peroxide Assay Kit is an optimized “mix and read” assay that is compatible with HTS liquid handling instruments. It provides a sensitive, one-step fluorometric assay to detect as little as 3 picomoles of H₂O₂ in a 100 µL assay volume (30 nM). The assay can be performed in a convenient 96-well or 384-well microtiter-plate format and easily adapted to automation without a separation step. Its signal can be easily read by either a fluorescence microplate reader at Ex/Em = ~640/680 nm or an absorbance microplate reader at ~650 nm. Due to its long emission wavelength, this kit has low interference from biological samples.

Kit Key Features

Broad Application:	Can be used for quantifying hydrogen peroxide in solutions, in cell extracts and in live cells, and for detecting a variety of oxidase activities through enzyme-coupled reactions.
Sensitive:	Detect as low as 30 picomoles of H ₂ O ₂ in solution.
Continuous:	Easily adapted to automation without a separation step.
Convenient:	Formulated to have minimal hands-on time. No wash is required.
Non-Radioactive:	No special requirements for waste treatment.

Kit Components

Components	Amount
Component A: Amplite™ IR Peroxidase Substrate	1 vial
Component B: H ₂ O ₂	1 vial (3% stabilized solution, 200 µL)
Component C: Assay Buffer	1 bottle (100 mL)
Component D: Horseradish Peroxidase	1 vial (20 units)
Component E: DMSO	1 vial (0.5 mL)

Assay Protocol for One 96-Well Plate

Brief Summary

Prepare H₂O₂ reaction mixture (50 µL) → Add H₂O₂ standards or test samples (50 µL) → Incubate at room temperature for 0-30 minutes → Monitor fluorescence intensity at Ex/Em = 640/680 nm

Note: Thaw all the kit components at room temperature before starting the experiment.

1. Prepare stock solutions:

- 1.1 **100X Amplite™ IR peroxidase substrate stock solution:** Add 250 µL of DMSO (Component E) into the vial of Amplite™ IR Peroxidase Substrate (Component A). The stock solution should be used promptly; any remaining solution should be aliquoted and refrozen at -20 °C.
Note: Avoid repeated freeze-thaw cycles, and protect from light.
- 1.2 **20 U/mL peroxidase stock solution:** Add 1 mL of Assay Buffer (Component C) into the vial of Horseradish Peroxidase (Component D).
Note: The unused HRP solution should be divided into single use aliquots and stored at -20 °C.
- 1.3 **20 mM H₂O₂ stock solution:** Add 22.7 µL of 3% H₂O₂ (0.88 M, Component B) into 977 µL of Assay Buffer (Component C).
Note: The diluted H₂O₂ stock solution is not stable. The unused portion should be discarded.

2. Prepare H₂O₂ reaction mixture:

Prepare the H₂O₂ reaction mixture according to the following table and keep from light.

Table 1. H₂O₂ Reaction mixture for one 96-well plate (2X)

Components	Volume
100X Amplite™ IR Peroxidase Substrate Stock Solution (from Step 1.1)	50 µL
20 U/mL Peroxidase Stock Solution (from Step 1.2)	200 µL
Assay Buffer (Component C)	4.75 mL
Total volume	5 mL

3. Prepare serial dilutions of H₂O₂ standard (0 to 10 µM):

Warning 1: Amplite™ IR Peroxidase Substrate (Component A) is unstable in the presence of thiols such as DTT and β mercaptoethanol. If the final concentration of the thiols is higher than 10 µM, it would significantly decrease the assay dynamic range.

Warning 2: NADH and glutathione (reduced form of GSH) may interfere with the assay.

- 3.1 Add 1 µL of 20 mM H₂O₂ stock solution (from Step 1.3) into 1999 µL of Assay Buffer (Component C) to get 10 µM H₂O₂ standard.
- 3.2 Take 200 µL of 10 µM H₂O₂ standard to perform 1:3 serial dilutions to get 3, 1, 0.3, 0.1, 0.03, 0.01 and 0 µM serial dilutions of H₂O₂ standard.
- 3.3 Add serial dilutions of H₂O₂ standard and H₂O₂-containing test samples into a solid black 96-well microplate as described in Tables 2 and 3.

Data Analysis

The fluorescence in blank wells (with the assay buffer only) is used as a control, and is subtracted from the values for those wells with the H₂O₂ reactions. A H₂O₂ standard curve is shown in Figure 1.

Note: The fluorescence background increases with time, thus it is important to subtract the fluorescence intensity value of the blank wells for each data point.




Figure 1. H₂O₂ dose response was measured in a 96-well black plate with the Amplitude™ Fluorimetric Hydrogen Peroxide Assay Kit using a Gemini fluorescence microplate reader (Molecular Devices). As low as 0.03 μM H₂O₂ can be detected with 1 minute incubation (n=3). The insert shows the low levels of H₂O₂ detection.

References

1. Hoffmann, O. M., Becker, D., and Weber, J. R. (2007) *J Cereb Blood Flow Metab.*
2. Funk, R. S., and Krise, J. P. (2007) *Mol Pharm* 4, 154-9.
3. Krebs, B., Wiebelitz, A., Balitzki-Korte, B., Vassallo, N., Paluch, S., Mitteregger, G., Onodera, T., Kretzschmar, H. A., and Herms, J. (2007) *J Neurochem* 100, 358-67.
4. Yang, Y., Xu, S., An, L., and Chen, N. (2007) *J Plant Physiol.*
5. Lee, J. E., Kim, H., Jang, H., Cho, E. J., and Youn, H. D. (2007) *J Neurochem.*
6. Goth, L. (2006) *Redox Rep* 11, 281-2.
7. Laloi, C., Stachowiak, M., Pers-Kamczyc, E., Warzych, E., Murgia, I., and Apel, K. (2007) *Proc Natl Acad Sci U S A* 104, 672-7.
8. Watts, R. J., Finn, D. D., Cutler, L. M., Schmidt, J. T., and Teel, A. L. (2006) *J Contam Hydrol.*

Warning: This kit is only sold to end users. Neither resale nor transfer to a third party is allowed without written permission from AAT Bioquest. Chemical analysis of kit components is strictly prohibited. Please call us at 408-733-1055 or e-mail us at info@aatbio.com if you have any questions.