

# Cell Meter™ Phosphatidylserine Apoptosis Assay Kit

## \*Green Fluorescence Optimized for Flow Cytometry\*

Ordering Information	Storage Conditions	Instrument Platform
Product Number: 22831 (100 assays)	Keep at 4 °C and avoid exposure to light	Flow Cytometer

### Introduction

Our Cell Meter™ assay kits are a set of tools for monitoring cellular functions. There are a variety of parameters that can be used for monitoring cell viability. This particular kit is designed to monitor cell apoptosis through measuring the translocation of phosphatidylserine (PS). In apoptosis, PS is transferred to the outer leaflet of the plasma membrane. The appearance of phosphatidylserine on the cell surface is a universal indicator of the initial/intermediate stages of cell apoptosis and can be detected before morphological changes can be observed.

This kit uses our proprietary green fluorescent Apopxin™ PS sensor that specifically binds PS with affinity much higher than Annexin V ( $K_d < 10$  nM). The PS sensor used in this kit has green fluorescence upon binding to membrane PS. The stain has the spectral properties almost identical to those of FITC with much higher photostability, making it convenient to be used with the common fluorescence instruments equipped with the light sources and filters for FITC, the most common fluorophore. This kit is optimized for flow cytometry applications while Kit 22791 is optimized to be used with a fluorescence microplate reader or a fluorescence microscope. This kit is a convenient alternative to Annexin V-FITC conjugate. The kit provides all the essential components with an optimized protocol for flow cytometric applications.

### Kit Key Features

<b>Non-Radioactive:</b>	No special requirements for waste treatment.
<b>Convenient:</b>	All essential assay components are included.
<b>Optimized Performance:</b>	Provide optimal conditions for detecting the translocation of phosphatidylserine.
<b>Enhanced Value:</b>	Less expensive than the sum of individual components.

### Kit Components

Components	Amount
Component A: Apopxin™ Green (100X stock solution)	1 vial (200 µL/vial)
Component B: Assay Buffer	50 mL
Component C: 100X Propidium Iodide	1 vial (200 µL)

### Assay Protocol

#### Brief Summary


**Prepare cells with test compounds (200 µL/sample) → Add Apopxin™ Green assay solution  
→ Incubate at room temperature for 30-60 minutes → Analyze cells with  
a flow cytometer using FL1 channel (Ex/Em = 490/525 nm)**

1. Treat cells with test compounds for a desired period of time (4-6 hours for Jurkat cells treated with camptothecin) to induce apoptosis.  
*Note: Apopxin™ binding flow cytometric analysis on adherent cells is not routinely tested since specific membrane damage may occur during cell detachment or harvesting. However, methods for utilizing Annexin V for flow cytometry on adherent cell types have been previously reported by Casiola-Rosen et al. and van Engeland et al (see Refs 1 and 2).*
2. Centrifuge the cells to get  $1-5 \times 10^5$  cells/tube.
3. Resuspend cells in 200 µL of Assay Buffer (Component B).

- Add 2  $\mu$ L of Apopxin™ Green (Component A) into the cells.  
*Optional: Add 2  $\mu$ L of 100X Propidium Iodide (Component C) into the cells for necrosis cells.*
- Incubate at room temperature for 10 to 20 minutes, protected from light.
- Optional: add 200 to 300  $\mu$ L of Assay Buffer (Component B) to increase volume before analyzing the cells with a flow cytometer (See Step 7).*
- Monitor the fluorescence intensity of Apopxin™ Green using the FL1 channel (Ex/Em = 490/525 nm), and measure the cell viability with propidium iodide using the FL2 channel.

## Data Analysis

In live non-apoptotic cells, Apopxin™ Green detects innate apoptosis in non-induced cells, which is typically 2-6% of all cells. In apoptotic cells, Apopxin™ Green binds to phosphatidylserine, which is located on the outer leaflet of the cell membrane, resulted in increased staining intensity.



**Figure 1.** The detection of binding activity of Apopxin™ Green and phosphatidylserine in Jurkat cells. Jurkat cells were treated without (A, and C-Blue line) or with 20  $\mu$ M camptothecin (B and C-Red line) in a 37 °C, 5% CO<sub>2</sub> incubator for 4-5 hours, and then dye loaded with Apopxin™ Green and Propidium Iodide for 30 minutes. The fluorescence intensity of Apopxin™ Green was measured with a FACSCalibur (Becton Dickinson, San Jose, CA) flow cytometer using the FL1 channel and the fluorescence intensity of Propidium Iodide was measured using the FL2 channel.

## References

- van Engeland M, Ramaekers FCS, Schutte B, Reutelingsperger CPM: A novel assay to measure loss of plasma membrane asymmetry during apoptosis of adherent cells in culture. *Cytometry* 24:131–139, 1996.
- L Casciola-Rosen, A Rosen, M Petri, and M Schlissel. Surface blebs on apoptotic cells are sites of enhanced procoagulant activity: implications for coagulation events and antigenic spread in systemic lupus erythematosus. *Proc Natl Acad Sci U S A*. 1996 February 20; 93(4): 1624–1629.
- Hanshaw RG, Lakshmi C, Lambert TN, Johnson JR, Smith BD. (2005) Fluorescent detection of apoptotic cells by using zinc coordination complexes with a selective affinity for membrane surfaces enriched with phosphatidylserine. *Chembiochem*, 6, 2214.
- Koulov AV, Stucker KA, Lakshmi C, Robinson JP, Smith BD. (2003) Detection of apoptotic cells using a synthetic fluorescent sensor for membrane surfaces that contain phosphatidylserine. *Cell Death Differ*, 10, 1357.
- Hall MP, Burson KK, Huestis WH. (1998) Interactions of a vesicular stomatitis virus G protein fragment with phosphatidylserine: NMR and fluorescence studies. *BiochimBiophys Acta*, 1415, 101.
- Saurel O, Cezanne L, Milon A, Tocanne JF, Demange P. (1998) Influence of annexin V on the structure and dynamics of phosphatidylcholine/phosphatidylserine bilayers: a fluorescence and NMR study. *Biochemistry*, 37, 1403.

**Warning: This kit is only sold to end users. Neither resale nor transfer to a third party is allowed without written permission from AAT Bioquest. Chemical analysis of the kit components is strictly prohibited. Please call us at 408-733-1055 or e-mail us at [info@aatbio.com](mailto:info@aatbio.com) if you have any questions.**