Catalog Number: 190186, 191418, 194819, 196055, 510722, 510723, 510823

# Dimethyl Sulfoxide

### Structure:



Molecular Formula: C<sub>2</sub>H<sub>6</sub>SO (for non-heavy labels)

Molecular Weight: 78.13 (non-heavy labels)

CAS #: 67-68-5

Synonyms: DMSO; Methyl sulfoxide; Methyl sulphoxide; Sulfinylbismethane

Physical Description: Clear colorless liquid to solid. The melting point is approximately 18°C so the product may appear anywhere from the liquid phase to the solid phase when in the pure form. By heating the solid form to approximately 30°C, the product can be melted without harming it's stability. DMSO supercools easily; melts at room temperature slowly.

**Solubility:** Soluble in water, methanol, acetone, ether, benzene, chloroform. To prepare sterile solutions use a a teflon or nylon membrane to sterile-filter the DMSO - do not use a cellulose acetate membrane.

**Stability:** A thermally stable compound. DMSO is stable up to 100° C in alkaline, acidic and neutral conditions. It is stable in neutral or alkaline conditions at temperatures approaching its boiling point of 189° C. DMSO can be heated to 150° C for 24 hours with less than 0.1% loss in purity. <sup>5</sup>

DMSO is very hygroscopic. 1

Density: Approximately 1.10 g/ml

Autoprotolysis constant: Approximately 33 at 25° C.7

Dielectric constant: 45 1

Description: A dipolar, aprotic solvent. Has been shown to accelerate

strand renaturation (1-10% concentration) and is believed to give the nucleic acid thermal stability against depurination.<sup>3,8</sup>

## Typical Uses:

- Used to enhance dermal absorption of many chemicals.
- A solvent for many organic and inorganic compounds including fats, carbohydrates, dyes, resins, and polymers.
- Used in antifreeze or hydraulic fluids.
- As a cryopreservative for cell cultures.<sup>5</sup>
- Used in the oxidation of thiols and disulfides to sulfonic acids.<sup>4</sup>
- Used as a PCR cosolvent to help improve yields, especially in long PCR.

# Plastic Compatibility:

# Incompatible Moderately compatible with: Compatible with: with:

- polysulfon e
- flexible and rigid PVC tubing
- polycarbon ate
- polystyrene
- Halar ECTFE
   (ethylene-chlorotrifl
   uoroethylene
   copolymer)
- Tefzel ETFE (ethylene-tetrafluoro ethylene)
- Low-density polyethylene (LDPE)
- High-density polyethylene (HDPE)
- Polypropylene
- Polypropylene copolymer (PPCO)
- Nylon
- Teflon ETFE (ethylene-tetrafluo roethylene)

Caution: Rapidly absorbed through skin and mucous membranes.

### Availability:

| Catalog Number | Description         | Size    |
|----------------|---------------------|---------|
| 190186         | Dimethyl Sulfoxide  | 100 ml  |
|                |                     | 500 ml  |
|                |                     | 1 liter |
| 194819         | Dimethyl Sulfoxide, | 50 ml   |
|                | molecular biology   | 100 ml  |
|                | reagent             | 250 ml  |

| 191418 | Dimethyl Sulfoxide, ACS | 100 ml  |
|--------|-------------------------|---------|
|        | Reagent Grade           | 500 m1  |
|        |                         | 1 liter |
| 196055 | Dimethyl Sulfoxide,     | 25 ml   |
|        | cell culture reagent    | 100 ml  |
|        |                         | 500 ml  |
|        |                         | 1 liter |

# Heavy Labels Available:

| 510722 | CD <sub>3</sub> SOCD <sub>3</sub><br>Purity: 99.9% D atom   | 1 g<br>5 g<br>10 g<br>25 g<br>50 g |
|--------|---|------------------------------------|
| 510723 | Dimethyl Sulfoxide - D <sub>6</sub> Contains +1% TMS CAS # 2206-27-1 CD <sub>3</sub> SOCD <sub>3</sub> Purity: 99.9% D atom Density = 1.19 gm/ml MW = 84.17 | 1 g<br>5 g<br>10 g<br>25 g<br>50 g |
| 510823 | Dimethyl Sulfoxide - D <sub>6</sub><br>CAS # 2206-27-1<br>CD <sub>3</sub> SOCD <sub>3</sub><br>Purity: 99.96% D atom<br>Density = 1.19 gm/ml<br>MW = 84.17  | 0.5 ml<br>5 ml                     |

#### References:

- 1. Merck Index, 13th Ed., No. 3285.
- 2. Bretherick's Handbook of Reactive Chemical Hazards, 4th ed., p. 299-303.
- 3. Cheng, S., et al., *Proc. Natl. Acad. Sci. USA*, v. 91, 5695-5699 (1994).
- 4. Lowe, O.G., J. Org. Chem., v. 41, 2061 (1976).
- 5. Martindale, The Extra Pharmacopoeia, 29th ed., p. 1426 (1989).
- 6. Nalgene Reference/Chemical Resistance Chart (Nalgene/Nunc Life Science Products Company catalog; www.nalgenenunc.com)
- 7. Rondinini, S., et al., *Pure and Applied Chem.*, v. **59**, 1693-1702 (1987).

8. Winship, P.R., et al., "An improved method for directly sequencing PCR amplified material using dimethyl sulphoxide." *Nucleic Acids Res.*, v. 17, 1266 (1989).