Catalog Number: 103168, 103169, 103170, 103171, 152487, 194724, 194725,

194841, 194842, 806576

Tween®

Structure (Tween® 80):19

where sum of w, x, y and z = 20

Description: Tweens® are a series of nonionic surfactants derived from sorbitan esters. They are soluble or dispersible in water but differ widely in organic and oil solubilities. Used as oil-in-water emulsifiers in pharmaceuticals, cosmetics, cleaning compounds, etc. (Not sold by MP for human use).

	Polysorbate 20	Tween® 40	Tween® 60	Tween® 80	Tween® 85
CAS #	9005-64-5	9005-66-7	9005-67-8	9005-65-6	9005-70-3
Synonym	Polyoxyethyle ne- sorbitan monolaurate; Polysorbate 20	Polyoxyethyle ne- sorbitan monopalmitate ; Polysorbate 40	Polyoxyethyle ne- sorbitan monostearate; Polysorbate 60	Polyoxyethyle ne- sorbitan monooleate; Polysorbate 80	Polyoxyethyle ne- sorbitan trioleate; Polysorbate 85
Content	~55% lauric acid (remainder primarily myristic, palmitic and stearic acid)	~90% palmitic acid	~55% stearic acid (remainder primarily palmitic acid)	~75% oleic acid	~70% oleic acid (remainder primarily elaidic, linolenic and palmitic acids)
Specific Gravity	~1.1 g/ml	~1.08 g/ml	~1.1 g/ml	~1.08 g/ml	~1.0 g/ml

Viscosity @ 25°C	400 cps	600 cps		400 cps	300 cps
Acid No.	0-2.0	0-2.0	0-2.0	0-2.0	0-2.0
Saponificat ion No.	40-50	43-49	45-55	45-55	83-93
Hydroxyl No.	96-108	89-105	81-96	65-80	39-52
Water	2. 5-3. 0%	2. 5-3. 0%	2. 5-3. 0%	2. 5-3. 0%	4. 8-5. 2%
Critical Micelle Concen. (CMC)	0.06 mM (60 mg/L or 6 × 10 ⁻³ mole/L)			0.012 mM (13-15 mg/L)	
Molecular Weight	1228		1311.7	1310	

Physical Description: Yellow, oily liquid

Typical Usage:

Polysorbate 20 is typically used as an emulsifying agent for the preparation of stable oil—in—water emulsions, particularly in pharmaceutical applications. ¹⁷ Polysorbate 20 has been used in pre-extraction of membranes to remove peripheral proteins (used at 2% for extraction of membrane—bound proteins). ^{1,10,12} Several resources may be helpful in determining usage concentrations. ^{8,9,14} Polysorbate 20 has been used as a blocking agent on nitrocellulose at a typical concentration of 0.05%.

Tween® 40 can be used as a non-ionic detergent in cell lysis, nuclei isolation and cell fractionation. ^{2,13,15}

Tween® 80 has been widely used in biochemical applications including: solubilizing proteins, isolating nuclei from cells in culture^{4,7,15}, selective protein extraction^{4,15}, growing of tubercule bacilli⁶, and emulsifying and dispersing substances in medicinal and food products. It has little or no activity as an anti-bacterial agent.³ It has been shown to have an adverse effect on the antibacterial effect of methyl paraben and related compounds.¹⁶

Critical Micelle Concentration (CMC): Detergents with high CMC values are generally easy to remove by dilution; detergents with low CMC values

are advantageous for separations on the basis of molecular weight. As a general rule, detergents should be used at their CMC and at a detergent-to-protein weight ratio of approximately ten. ¹⁹⁻²⁰

Solubility: Soluble/miscible in water to give a clear yellow solution; miscible with alcohol, dioxane, and ethyl acetate; practically insoluble in liquid paraffin and fixed oils (such as mineral oil). ¹⁷ Autoclaving of solutions is not recommended. Sterile filtering is suggested with a 0.22 micron filter. Tween® may need to be warmed to about 40° C and alternated with portions of hot distilled water while being poured through the filter.

Availability:

Catalog Number	Description	Size
103168	Polysorbate 20	100 ml 500 ml 1 liter
194724	Polysorbate 20, cell culture reagent	100 ml 500 ml 4 liter
806576	Polysorbate 20	100 g
194841	Polysorbate 20, molecular biology reagent	50 ml 100 ml
103169	Tween® 40	100 ml 500 ml
103171	Tween® 60	100 ml 500 ml 1 liter
103170	Tween® 80	100 ml 500 ml 1 liter 4 liters 5 liters
194725	Tween® 80, cell culture reagent	100 ml 500 ml 4 liters
194842	Tween® 80, molecular biology reagent	50 ml 100 ml
152487	Tween® 85	100 ml 500 ml

4 liters

References:

- 1. Black, P.N., et al., *J. Biol. Chem.*, v. 262, 1412-1419 (1987).
- 2. Culvenor, J.G., et al., "Characteristics of plasma membrane isolated from a mouse T lymphoma line: comparison after nitrogen cavitation, shearing, detergent treatment, and microvesiculation." J. Cell. Biochem., v. 20, 127-138 (1982).
- 3. Dawson, et al. (eds.), *Data for Biochemical Research, 3rd Ed.*, Oxford Press, p. 289 (1989).
- 4. Deshpande, R.G., et al., "Isolation of a contact-dependent haemolysin from Mycobacterium tuberculosis." *J. Med. Microbiol.*, v. 46, 233-238 (1997).
- 5. Donbrow, M., et al., *J. Pharmaceutical Sciences*, v. 67, 1676-1681 (1978).
- 6. Dubos, R. J. and Davis, B. D., J. Exp. Med., v. 83, 409 (1946).
- 7. Fisher, H.W. and Harris, H., Proc. R. Soc. B, v. 156, 521 (1962).
- 8. Harris and Angal (eds.), *Protein Purification Methods: A Practical Approach*, IRL Press (1989).
- 9. Harris and Angal (eds.), *Protein Purification Applications: A Practical Approach*, IRL Press (1989).
- 10. Helenius, A. and Simons, K., *Biochem. Biophys. Acta*, v. **455**, 796 (1976).
- 11. Helenius, A., et al., Methods in Enzymology, v. 56, 734-749 (1979).
- 12. Karasawa, et al., *J. Biol. Chem.*, v. 274, 8655-8661 (1999).
- 13. Link, A. J. (ed.), 2-D Proteome Analysis Protocols, v. 112, Methods in Molecular Biology, Humana Press: Totowa, NJ (1999).
- 14. Neugebauer, J. M., *Methods in Enzymology*, v. 182, 239 (1992).
- 15. Storrie, B and Madden, E.A., "Isolation of subcellular organells." *Methods Enzymol.*, v. 182, 203-225 (1990).
- 16. Disinfection, Sterilization & Preservation, 4th Ed., Block, S.S. (Lea & Febiger Pub.) Chapter 4 (1991).
- 17. Martindale: The Extra Pharmacopoeia, 30th Ed., Pharmaceutical Press, p. 1030, 1247 (1993).
- 18. Sys. Analysis of Surface Active Agents, 2nd Ed., p. 533.
- 19. Womack, M. D., et al., Biochim. Biophys. Acta, v. 733, 210 (1983)
- 20. Matson, R. S. and Goheen, S. C., LC-GC, v. 4(7), 624 (1986)
- 21. Merck Index, 12th Ed., No. 7742.

Tween® is a registered trademark of Uniqema, a business unit of ICI

Americas, Inc.