



CELLULAR IMAGING & ANALYSIS

# CLIP-Cell™ Starter Kit

Instruction Manual

NEB #E9200S  
10 reactions

 NEW ENGLAND  
**BioLabs**® Inc.  
*enabling technologies in the life sciences*





## Table of Contents:

---

Introduction .....	2
Construction and Expression of CLIP-tag Fusion Proteins .....	4
Use of the CLIP <sub>f</sub> Control Plasmid pCLIP <sub>f</sub> -H2B .....	7
Labeling CLIP-tag Fusion Proteins .....	8
Use of CLIP-Cell™ Block to Block Reactivity of CLIP <sub>f</sub> Fusion Proteins .....	13
Troubleshooting .....	15
Appendix: Sequence of CLIP-tag region of pCLIP <sub>f</sub> .....	18
Ordering Information .....	19

## Kit Components:

---

*All kit components should be stored at –20°C except where noted.*

pCLIP<sub>f</sub> Vector ..... 20 µg  
Supplied in TE buffer (10 mM Tris-HCl, pH 8.0, 1 mM EDTA), 0.5 µg/µl

pCLIP<sub>f</sub>-H2B Control Plasmid ..... 20 µg  
Supplied in TE buffer (10 mM Tris-HCl, pH 8.0, 1 mM EDTA), 0.5 µg/µl

CLIP-Cell™ 505 ..... 10 nmol  
Supplied dry. Prepare a stock solution in DMSO as described on page 10.

CLIP-Cell™ TMR-Star ..... 6 nmol  
Supplied dry. Prepare a stock solution in DMSO as described on page 10.

CLIP-Cell™ Block ..... 20 nmol  
Supplied dry. Prepare a stock solution in DMSO as described on page 13.

*Note: For long-term storage, all kit components should be stored at –20°C. Plasmid solutions can be stored at 4°C for up to one week. Undissolved dye and blocking substrates can be stored at 4°C for up to 4 weeks protected from light and moisture. With proper storage at –20°C the substrates should be stable for at least two years dry or 3 months dissolved in DMSO.*

## Required Materials Not Included

Mammalian Cell Lines

DNA Transfection Reagents

Standard Tissue Culture Media and Plasticware

DMSO


Hoechst 33342 for Nuclear Staining (optional)

## Introduction:

The SNAP-tag® and CLIP-tag™ are novel tools for the specific, covalent attachment of virtually any molecule to a protein of interest, providing simplicity and extraordinary versatility to the imaging of proteins in live and fixed cells, and to the study of proteins *in vitro*. The creation of a single gene construct yields a tagged fusion protein capable of forming a covalent linkage to a variety of functional groups, including fluorophores, biotin, or beads. This system provides a powerful and unique tool to study the role of proteins in a variety of highly dynamic processes, including protein trafficking, turnover and complex formation.

The CLIP-tag is a 20 kDa mutant of the human DNA repair protein O<sup>6</sup>-alkylguanine-DNA alkyltransferase (hAGT) that reacts specifically and rapidly with benzylcytosine (BC) derivatives, leading to covalent labeling of the CLIP-tag with a synthetic probe (Figure 1). The CLIP-tag has a number of features that make it ideal for a variety of protein labeling applications. The rate of the reaction of the CLIP-tag with these derivatives is largely independent of the nature of the synthetic probe attached to BC, permitting the labeling of CLIP fusion proteins with a wide variety of functional groups. Many of these CLIP-tag substrates are cell-permeable, allowing live-cell imaging of protein expression and localization (Figure 2). The ability to turn on the signal at will, together with the availability of a cell-permeable nonfluorescent blocking agent (CLIP-Cell™ Block), allows time-resolved pulse-chase analysis of protein trafficking. Finally, the availability of orthogonal protein labeling systems from NEB permits simultaneous labeling of multiple proteins in a single cell (SNAP-tag, another hAGT variant that reacts exclusively with O<sup>6</sup>-benzylguanine substrates, and the ACP/MCP tags, small protein tags which can be enzymatically labeled on the cell surface with Coenzyme A derivatives).

Figure 1:




CLIP-tag Reaction.


The CLIP-Cell Starter Kit contains a mammalian expression plasmid (pCLIP<sub>i</sub>) encoding the CLIP-tag flanked by restriction sites for cloning a gene of interest, and two cell-permeable fluorescent CLIP-tag substrates. A positive control plasmid (pCLIP<sub>i</sub>-H2B), encoding a CLIP-tagged protein (histone H2B) with a well-characterized nuclear localization, is also included. Lastly, a negative

control “blocking agent” (CLIP-Cell Block) is included that interacts with the CLIP-tag, but is not fluorescent. There are two steps to using this system: subcloning and expression of the protein of interest as a CLIP<sub>f</sub> fusion, and labeling of the fusion with the CLIP-tag substrate of choice.

Figure 2. Live cell imaging of CLIP<sub>f</sub> fusion proteins.



*Live CHO-K1 cells transiently transfected with pCLIP<sub>f</sub>-NK1R. Cells were labeled with CLIP-Cell™ 505 (green) for 30 minutes and counterstained with Hoechst 33342 (blue).*




*Live COS-7 cells expressing CLIP<sub>f</sub>-histone H2B (CLIP<sub>f</sub>-H2B) and SNAP<sub>f</sub>-adrenergic beta-2 surface receptor (SNAP<sub>f</sub>-ADRB2) were labeled with CLIP-Cell™ TMR-Star (red) and SNAP-Surface® Alexa Fluor® 488 (green).*

# Construction and Expression of CLIP-tag Fusion Proteins:

The mammalian expression plasmid pCLIP<sub>f</sub> is intended for the cloning and stable or transient expression of CLIP-tag protein fusions in mammalian cells. This plasmid encodes the CLIP<sub>f</sub> gene, which is expressed under control of the CMV promoter. The expression vector has an IRES (internal ribosome entry site) and a neomycin resistance gene downstream of CLIP<sub>f</sub> for the efficient selection of stable transfectants. Codon usage of the gene is optimized for expression in mammalian cells. pCLIP<sub>f</sub> contains two multiple cloning sites to allow cloning of the fusion partner as a fusion to the N- or C-terminus of the CLIP-tag.

pCLIP<sub>f</sub> contains an improved version of CLIP-tag, termed CLIP<sub>f</sub>. CLIP<sub>f</sub> displays faster kinetics in *in vitro* labeling and fast, specific and efficient labeling in live and fixed cell applications, thereby rendering it a desired research tool for analysis of protein dynamics.

Figure 3.



pCLIP<sub>f</sub> plasmid map.

## Detailed Description of pCLIP<sub>f</sub>

The sequence of the cloning region can be found in the Appendix. The complete plasmid sequence can be downloaded at [www.neb.com](http://www.neb.com). This plasmid encodes the gene CLIP<sub>f</sub>, which is a mutant form of the human gene for O<sup>6</sup>-alkylguanine-DNA-alkyltransferase (hAGT). The codon usage of the gene is optimized for expression in mammalian cells. In the plasmid sequence, the CLIP<sub>f</sub> gene is encoded from bp 969 to 1514.

This plasmid is intended for the cloning and stable or transient expression of CLIP-tag protein fusions in mammalian cells. It is particularly suitable for the efficient production of stable cell lines expressing CLIP-tag gene fusions. The plasmid contains the CMV promoter followed by the genes for CLIP<sub>f</sub> and neomycin resistance separated by the IRES of the encephalomyocarditis virus (ECMV), which permits the translation of two open reading frames from one messenger RNA. After selection of stable mammalian cells for neomycin resistance, nearly all surviving colonies should therefore stably express the CLIP<sub>f</sub> fusion protein. An intron is also included because this is believed to improve expression levels. Unless your expression experiments require a pure population of cells, you can simply use the pool of resistant cells. Otherwise cell clones can be isolated and characterized using standard procedures. The plasmid also contains the beta-lactamase (Ampicillin resistance) gene for maintenance in bacteria. The gene of interest can be cloned upstream or downstream of the CLIP<sub>f</sub> coding sequence, as a fusion to the N- or C-terminus of the CLIP<sub>f</sub>. pCLIP<sub>f</sub> can also be used as an expression control plasmid, expressing the CLIP-tag alone, in which case the CLIP-tag protein is distributed throughout the cell. The CLIP<sub>f</sub> gene can be isolated from the plasmid using PCR or direct cloning in order to subclone it into a different vector system of choice.

## Cloning of CLIP-tag Fusions in pCLIP<sub>f</sub> Vector

### *Cloning by PCR*

1. To subclone the gene of interest into pCLIP<sub>f</sub> fused to the N-terminus of CLIP<sub>f</sub>, use the available restriction sites: *NheI*, *EcoRV* (blunt), *AscI*, *SwaI* (blunt), *BsrGI*, *AgeI* or *EcoRI*, which are located upstream of the CLIP-tag.
2. To subclone the gene of interest into pCLIP<sub>f</sub> fused to the C-terminus of CLIP<sub>f</sub>, use the available restriction sites downstream of the CLIP-tag: *SbfI*, *BamHI*, *PmeI* (blunt), *XhoI*, *PacI* or *NotI*.

*Note: When fusing the gene of interest to the C-terminus of CLIP<sub>f</sub>, note that there is an in-frame stop codon between the *PacI* and *NotI* sites, so *SbfI*, *BamHI*, *PmeI* (blunt), *XhoI* or *PacI* must be used as the 5' cloning site for your insert.*

*PmeI and XhoI cannot be used together for cloning because they share a cytosine as part of their recognition sequences.*

### *Primer Design and Cloning Considerations*

1. Design your PCR primers to include a sufficient overlap with the sequence of the gene you want to amplify, adding 5–6 bases on the 5' side of the introduced restriction site to ensure efficient cleavage prior to cloning.
2. A stop codon can be included at the C-terminus of the fusion (in front of the downstream cloning site) in order to terminate translation at this position.
3. For fusions upstream of CLIP<sub>+</sub>, ensure that a start codon is included. The addition of a Kozak sequence (e.g. GCCRCCATG, where the start codon is underlined) will increase the translation efficiency.
4. In general, any linker peptide between the proteins should be kept short to avoid degradation by proteases. If required, specific protease cleavage sites can be introduced into the linker peptide.
5. Care should be taken to design the cloning so that the fusion partners in the resulting construct are in frame.
6. Perform the PCR reaction and subsequent cloning steps according to established molecular biology protocols.
7. The ligated vector should be transformed into bacteria and the resulting plasmid isolated via a standard miniprep procedure.
8. After subcloning the gene of interest into pCLIP<sub>+</sub> as a fusion with the CLIP<sub>+</sub> gene, the resulting plasmid can be used for stable or transient expression of the CLIP<sub>+</sub> fusion proteins in a suitable cell line.

### *Direct Cloning*

1. Direct cloning can also be used to make fusions with the CLIP-tag. This is only possible if the fusion partner has compatible sites adjacent to the gene of interest.
2. Care should be taken to design the cloning strategy so that the fusion partners in the resulting construct are in frame.

*Note: When fusing the gene of interest to the C-terminus of CLIP<sub>+</sub>, note that there is a stop codon between the PaeI and NotI sites, so SbfI, BamHI, PmeI (blunt), XhoI or PaeI must be used as the 5' cloning site for your insert.*

*PmeI and XhoI cannot be used together for cloning because they share a cytosine as part of their recognition sequences.*

## Expression of CLIP-tag Fusions

### *Transient Expression*

Expression of the fusion protein cloned in pCLIP<sub>f</sub> can be achieved by transiently transfecting cells in culture with standard transfection protocols. The appropriate reagent and transfection time to permit adequate expression must be empirically determined, using guidelines provided by the manufacturer of the transfection reagent as a starting point. We recommend using pCLIP<sub>f</sub> or pCLIP<sub>f</sub>-H2B as expression control plasmids. Figure 2B shows that the Histone H2B-CLIP<sub>f</sub> fusion protein (from pCLIP<sub>f</sub>-H2B) gives a nuclear localized signal when labeled with CLIP-Cell substrates. If the empty pCLIP<sub>f</sub> plasmid is used as a control vector for transfection, a uniform distribution of the CLIP-tag between nucleus and cytoplasm should be seen. Both pCLIP<sub>f</sub> and the localization control plasmids have performed well in stable and transient transfection of CHO-K1, COS-7, U-2 OS and NIH 3T3 cells. Note that the intensity of the fluorescence may vary depending on cell line and labeling substrate used. We recommend using TransPass D2 (NEB #M2554) in combination with TransPass V (NEB #M2561) or Roche's FuGENE® 6 transfection reagent for both transient and stable transfections.

### *Stable Expression*

pCLIP<sub>f</sub> and pCLIP<sub>f</sub>-H2B can be transfected as described above for transient transfection or by other standard transfection methods. Twenty four to 48 hours after transfection begin selecting mammalian cultures in 600–1,200 µg/ml G418 (geneticin) depending on the cell line. It is recommended that you establish a kill curve for each cell line to determine optimal selection conditions. After 8–12 days of continuous selection, stable colonies will become visible. It is possible to use pools of stable cell populations for initial cell labeling to test for the presence of CLIP-tag expression. In addition, monoclonal cell lines can be isolated and characterized if desired.

## Use of the CLIP<sub>f</sub> Control Plasmid pCLIP<sub>f</sub>-H2B:


This control plasmid contains the gene encoding the Histone H2B protein cloned upstream of the CLIP<sub>f</sub> coding sequence in pCLIP<sub>f</sub>, as a fusion to the N-terminus of the CLIP<sub>f</sub>. Histone H2B is a member of the core histones that package DNA in the nucleus. The H2B-CLIP<sub>f</sub> fusion protein gives nuclear fluorescence when labeled with CLIP-Cell substrates (Figure 2B). The full sequence and map for pCLIP<sub>f</sub>-H2B can be downloaded at [www.neb.com](http://www.neb.com). We strongly recommend carrying out parallel expression and labeling experiments with this plasmid as a positive control for your experiment, using the procedure described above.

## Labeling CLIP-tag Fusion Proteins:

The kit includes two cell-permeable fluorescent CLIP-tag substrates, CLIP-Cell 505 and CLIP-Cell TMR-Star. Both can be used to label CLIP<sub>f</sub> fusion proteins inside living or fixed cells, on cell surfaces, or in solution.

CLIP-Cell 505 (BC-505) is a photostable green fluorescent substrate that is based on the Dyomics dye DY-505 and is suitable for standard fluorescein filter sets. It has an excitation maximum at 504 nm and an emission maximum at 532 nm (Figure 4). This kit contains 10 nmol of CLIP-Cell 505 substrate, sufficient to make 2 ml of a 5  $\mu$ M CLIP<sub>f</sub> fusion protein labeling solution.


Figure 4.



(A) Structure of CLIP-Cell 505 (MW 586.6 g/mol) (B) Excitation (dotted line) and emission (solid line) spectra of CLIP-Cell 505 after coupling to CLIP-tag in buffer at pH 7.5.

CLIP-Cell TMR-Star (BC-TMR-Star) is a red fluorescent substrate that is based on tetramethylrhodamine and is suitable for standard rhodamine filter sets. It has an excitation maximum at 554 nm and an emission maximum at 580 nm (Figure 5). This kit contains 6 nmol of CLIP-Cell TMR-Star substrate, sufficient to make 2 ml of a 3  $\mu$ M CLIP<sub>f</sub> fusion protein labeling solution.

Figure 5.



(A) Structure of CLIP-Cell TMR-Star (MW 642.7 g/mol) (B) Excitation (dotted line) and emission spectra of CLIP-Cell TMR-Star after coupling to CLIP-tag in buffer at pH 7.5.

## Instructions for Cellular Labeling

### *Preparation of Labeling Stock Solution*

Dissolve one vial of CLIP-tag substrate in 10  $\mu$ l of fresh DMSO to yield a labeling stock solution of 1 mM CLIP-Cell 505 or 0.6 mM CLIP-Cell TMR-Star. Mix by vortexing for 10 minutes until all the CLIP-tag substrate is dissolved. Store this stock solution in the dark at 4°C, or for extended storage at -20°C. Different stock concentrations can be made, depending on your requirements. The substrates are soluble up to at least 10 mM.

### *Protocol for Intracellular Labeling Reaction*

1. Dilute the labeling stock solution 1:200 in medium to yield a labeling medium of 5  $\mu$ M CLIP-Cell-505 or 3  $\mu$ M CLIP-Cell TMR-Star. Mix dye with medium thoroughly by pipetting up and down 10 times (necessary for reducing backgrounds). For best performance, add the CLIP-tag substrate to complete medium, including serum (0.5% BSA can be used for experiments carried out in serum-free media). Do not prepare more medium with CLIP-tag substrate than you will consume within one hour.

NUMBER OF WELLS IN PLATE	RECOMMENDED VOLUME FOR CELL LABELING
6	1 ml
12	500 $\mu$ l
24	250 $\mu$ l
48	100 $\mu$ l
96	50 $\mu$ l

These recommendations are for culturing cells in polystyrene plates. For confocal imaging, we recommend using chambered coverglass such as Lab-Tek II Chambered Coverglass which is available in a 1, 2, 4 or 8 well format from Nunc ([www.nuncbrand.com](http://www.nuncbrand.com)).

2. Replace the medium on the cells expressing a CLIP<sub>f</sub> fusion protein with the CLIP-tag labeling medium and incubate at 37°C, 5% CO<sub>2</sub> for 30 minutes.
3. Wash the cells three times with tissue culture medium containing serum and incubate in fresh medium for 30 minutes. Replace the medium one more time to remove unreacted CLIP-tag substrate that has diffused out of the cells.
4. Image the cells using an appropriate filter set. CLIP<sub>f</sub> fusion proteins labeled with CLIP-Cell 505 should have an excitation maximum at 504 nm and an emission maximum at 532 nm, and can be imaged with standard fluorescein filter sets. CLIP<sub>f</sub> fusion proteins labeled with CLIP-Cell TMR-Star should have an excitation maximum at 554 nm and an emission maximum at 580 nm, and can be imaged with standard rhodamine filter sets.

5. We recommend routinely labeling one well of non-transfected or mock-transfected cells as a negative control.

## Notes for Cellular Labeling

### *Blocking Unreacted CLIP-tag with CLIP-Cell Block*

In many cases the labeling of a non-transfected cell sample or a mock-transfected cell sample will be completely sufficient as a control. In some cases, however, it may be desirable to block the CLIP-tag activity in a cell sample expressing the CLIP<sub>f</sub> fusion protein to generate a control. This can be achieved using the included nonfluorescent CLIP-tag substrate, CLIP-Cell Block (bromomethylcytosine, BTC). CLIP-Cell Block may also be used in pulse-chase experiments to block the CLIP-tag reactivity during the chase between two pulse-labeling steps. A protocol for blocking with CLIP-Cell Block can be found on page 13.

### *Optimizing Labeling*

Optimal substrate concentrations and reaction times range from 1–20  $\mu$ M and 5–30 minutes, respectively, depending on experimental conditions and expression levels of the CLIP<sub>f</sub> fusion protein. Best results are usually obtained at concentrations between 1 and 5  $\mu$ M substrate and 30 minutes reaction time. Increasing substrate concentration and reaction time usually results in a higher background without necessarily increasing the signal to background ratio.

### *Stability of Signal*

The turnover rates of the CLIP<sub>f</sub> fusion protein under investigation may vary widely depending on the fusion partner. We have seen half-life values ranging from less than one hour to more than 12 hours. Where protein turnover is rapid, we recommend analyzing the cells under the microscope immediately after the labeling reaction or, if the application allows it, fixing the cells directly after labeling. As an alternative to visualize proteins with fast turnover rates, CLIP<sub>f</sub> fusion proteins can be labeled at lower temperatures (4 or 16°C). Labeling times may need to be optimized.

### *Fixation of Cells*

After labeling the CLIP<sub>f</sub> fusion proteins, the cells can be fixed with standard fixation methods such as para-formaldehyde, ethanol, methanol, methanol/acetone etc., without loss of signal. We are not aware of any incompatibility of the CLIP-tag label with any fixation method.

### *Counterstaining*

Cells can be counterstained with any live-cell dye that is compatible with the fluorescent properties of the CLIP-tag substrate for simultaneous microscopic detection. We routinely add 5  $\mu$ M Hoechst 33342 to the medium that is used prior to the final 30 minutes incubation (Step 3 above) as a DNA counterstain for nuclear visualization. Counterstaining of cells is also possible after fixation and permeabilization.

## Immunocytochemistry

Antibody labeling can be performed after CLIP-tag labeling and fixation of the cells according to standard protocols without loss of the CLIP-tag signal. The fixation conditions should be selected based on experience with the protein of interest. For example some fixation methods destroy epitopes of certain proteins and therefore do not allow antibody staining afterwards.

## Instructions for Labeling of Proteins *in vitro*

1. Dissolve the vial of CLIP-Cell 505 (10 nmol) in 10  $\mu$ l of fresh DMSO or the vial of CLIP-Cell TMR-Star (6 nmol) in 6  $\mu$ l of fresh DMSO to yield a labeling stock solution of 1 mM CLIP-tag substrate. Mix by vortexing for 10 minutes until all the CLIP-tag substrate is dissolved. Dilute this 1 mM stock solution 1:4 in fresh DMSO to yield a 250  $\mu$ M stock for labeling proteins *in vitro*.
2. Set up the reactions, in order, as follows:

COMPONENT	VOLUME	FINAL CONCENTRATION
Deionized Water	32 $\mu$ l	
5X CLIP-tag Reaction Buffer	10 $\mu$ l	1X
50 mM DTT	1 $\mu$ l	1 mM
50 $\mu$ M CLIP-tag Purified Protein	5 $\mu$ l	5 $\mu$ M
250 $\mu$ M CLIP-tag Substrate	2 $\mu$ l	10 $\mu$ M
<b>Total Volume</b>	50 $\mu$ l	

3. Incubate in the dark for 30 minutes at 37°C.
4. Run sample on an SDS-PAGE gel and detect using a fluorescent gel scanner or store samples at -20°C or -80°C in the dark.

### Removal of Unreacted Substrate (optional)

After the labeling reaction you may wish to separate the nonreacted substrate from the labeled CLIP<sub>f</sub> fusion protein. You can use gel filtration or dialysis. Please refer to the vendor's instructions for the separation tools you are using.

### Note for Labeling *in vitro*

We recommend the routine addition of 1 mM DTT to all buffers used for handling, labeling and storage of the CLIP-tag. The stability of the CLIP-tag is improved in the presence of reducing agents; however it can also be labeled in their absence (e.g. for a redox-sensitive protein) if handling at temperatures above 4°C is minimized. CLIP<sub>f</sub> fusion proteins can be purified before labeling, but the labeling reaction also works in non-purified protein solutions (including cell lysates).

## Use of CLIP-Cell Block to Block Reactivity of CLIP<sub>f</sub> Fusion Proteins:


CLIP-Cell Block (bromothienylcytosine, BTC) is a non-fluorescent compound that blocks the reactivity of the CLIP-tag in live cells (Figure 6). It can be used to generate inactive controls in live and fixed cell labeling experiments performed with CLIP<sub>f</sub> fusion proteins. CLIP-Cell Block is highly membrane permeable and once in the cell reacts with the CLIP-tag irreversibly inactivating it for subsequent labeling steps.

### Instructions for Use with CLIP-Cell Substrates

In many cases the labeling of a non-transfected cell sample or a mock-transfected cell sample will be completely sufficient as a negative control for cell labeling. In some cases, however, it may be desirable to block the CLIP-tag activity in a cell sample expressing the CLIP<sub>f</sub> fusion protein to generate a control. This is done by a pre-incubation of the cells with CLIP-Cell Block, followed by the incubation with the labeling solution. CLIP-Cell Block may also be used in pulse-chase experiments to block the CLIP-tag reactivity during the chase between two pulse-labeling steps.

*Note that CLIP-Cell Block is a potent blocker of the CLIP-tag! Always take care to avoid carryover of CLIP-Cell Block to samples that you do not wish to block.*

Figure 6.



Structure of CLIP-Cell Block (MW 286.15 g/mol)

### Preparation of Stock Solution

Dissolve one tube of CLIP-Cell Block (20 nmol) in 10  $\mu$ l of fresh DMSO to give a stock solution of 2 mM. Mix by vortexing for 10 minutes, until all the CLIP-Cell Block is dissolved. Store this stock solution in the dark at 4°C or for extended storage at -20°C. We recommend using a final concentration of 10  $\mu$ M, which is a 1:200 dilution of this stock solution.

### *Blocking CLIP-tag Activity with CLIP-Cell Block*

The following steps describe the use of CLIP-Cell Block in a typical control labeling experiment:

1. Prepare two cell samples suitable for labeling, each expressing the CLIP<sub>f</sub> fusion protein of interest.
2. Mix an appropriate amount of medium with CLIP-Cell Block stock solution in a ratio of 1:200 to give a blocking medium of 10  $\mu$ M CLIP-Cell Block. For best performance, add the dissolved CLIP-Cell Block to complete medium, including serum. Do not prepare more medium with CLIP-Cell Block than you will consume within one hour.
3. Mix an appropriate amount of medium with DMSO in a ratio of 1:200, to give a final concentration of 0.5% v/v DMSO.
4. Replace the medium on one sample of cells with the blocking medium. These are your blocked cells. Replace the medium on the other sample of cells with the medium containing DMSO. These are your test cells. Incubate both cell samples at 37°C, 5% CO<sub>2</sub> for 30 minutes.
5. Remove CLIP-Cell Block or DMSO-containing medium by washing both samples of cells twice with complete medium.
6. Label both cell samples with the CLIP-Cell substrate using the protocol on page 10.
7. Inspect both samples under the fluorescence microscope. The blocked cells should show no fluorescence, whereas the test cells should show fluorescence localized to where the CLIP<sub>f</sub> fusion protein is present in the cell.

*Note that there is a constant turnover and resynthesis of proteins in the cell. After having blocked all existing CLIP<sub>f</sub> fusion proteins within the cell, new CLIP<sub>f</sub> fusion protein molecules may be synthesized in the meantime and may get labeled during incubation with a fluorescent CLIP-tag substrate. This will give the impression that the blocking was ineffective. In order to minimize these effects of protein synthesis and protein transport, cells may have to be treated with cycloheximide and incubation with the fluorescent CLIP-tag substrate may have to be performed at 4°C.*

## Troubleshooting:

### Cloning of the Gene of Interest

If subcloning of the gene of interest into the pCLIP<sub>f</sub> vector does not work, reconfirm all the cloning steps (primer design, choice of restriction site, etc.). If all steps are confirmed as being correct, then try the cloning using different restriction sites. Be sure to include a positive and negative control for the ligation reaction. Alternatively try to subclone the CLIP<sub>f</sub> gene into an expression vector already containing your gene of interest.

### Expression

In general we have not experienced problems expressing CLIP<sub>f</sub> protein fusions. However if your fusion protein does not appear to be expressed, try expressing the H2B-CLIP<sub>f</sub> or CLIP<sub>f</sub>-NK1R protein fusions as positive controls using cells transiently transfected with pCLIP<sub>f</sub>-H2B (included) or pCLIP<sub>f</sub>-NK1R (NEB #N9216). Labeling of such cells with a fluorescent CLIP-Cell substrate should show strong nuclear or ER, Golgi and surface-localized fluorescence, respectively. The empty pCLIP<sub>f</sub> plasmid can also be used as a control (cytosolic and nuclear fluorescence). Note that the intensity of this fluorescence may vary depending on cell line and substrate used. Expression of localization controls but not your fusion protein can be due to a variety of causes. It is possible that this fusion protein may be toxic for your cell line. It is difficult to troubleshoot such instances, but the use of a different expression plasmid or cell line or tagging the opposite end of the protein may help. Signs of host cell toxicity could include slow proliferation or apoptosis. Counterstaining live cells with Hoechst 33342 or fixed cells with DAPI can be used to determine whether nuclei are healthy if toxicity is suspected.

### Problems with Cellular Labeling

#### *No Labeling*

If no labeling is seen, the most likely explanation is that the fusion protein is not expressed. Verify your transfection method to confirm that the cells contain the fusion gene of interest. If this is confirmed, check for expression of the CLIP<sub>f</sub> fusion protein. If no antibody against the fusion partner is available, Anti-SNAP-tag Antibody (NEB #P9310) can be used. Alternatively, CLIP-Vista Green (NEB #S9235) can be used to confirm presence of CLIP<sub>f</sub> fusion in cell extracts following SDS-PAGE, without the need for Western blotting.

#### *Weak Labeling*

Weak labeling may be caused by insufficient exposure of the fusion protein to the substrate. Try increasing the concentration of CLIP-tag substrate and/or the incubation time within the range of 1–20  $\mu$ M and 5–30 minutes, respectively. Alternatively the protein may be poorly expressed and/or turn over rapidly. If the protein has limited stability in the cell, it may help to analyze the samples immediately after labeling.

### *High Background*

Background fluorescence may be controlled by reducing the concentration of CLIP-tag substrate used, and by shortening the incubation time. The presence of fetal calf serum or BSA during the labeling incubation should reduce non-specific binding of substrate to surfaces.

### *Signal Strongly Reduced after Short Time*

If the fluorescence signal decreases rapidly, it may be due to instability of the fusion protein. The signal may be stabilized by fixing the cells. Alternatively try switching the CLIP-tag from the N- to the C-terminus or vice versa. Photobleaching is generally not a problem as both CLIP-Cell 505 and CLIP-Cell TMR-Star are very photostable. However, if you experience problems with photobleaching, addition of a commercially available anti-fade reagent may be helpful.

## Problems with Labeling *in vitro*

### *Solubility*

If solubility problems occur with your CLIP<sub>f</sub> fusion protein, we recommend testing a range of pH (pH 5.0–pH 10.0). The salt concentration may also need to be optimized for your particular fusion protein (50–250 mM).

### *Loss of Protein due to Aggregation or Sticking to Tube*

If stickiness of the fusion protein is a problem we recommend adding Tween 20 at a final concentration of 0.05% to 0.1%. The CLIP-tag activity is not affected by this concentration of Tween 20.

### *Incomplete Labeling*

If exhaustive labeling of a protein sample is not achieved using the recommended conditions, try the following protocol modifications: Increase the incubation time to two hours total at 25°C or to 24 hours at 4°C; or halve the volume of protein solution labeled. Both approaches may be combined. If you still have poor labeling results, we recommend checking the activity of the CLIP-tag using CLIP-Vista Green (NEB #S9235).

If the CLIP<sub>f</sub> fusion has been stored in the absence of DTT or other reducing agent, or has been stored at 4°C for a prolonged period, its activity may be compromised. Include 1 mM DTT in all solutions of the CLIP<sub>f</sub> fusion protein, and store the fusion protein at –20°C.

Using less than the recommended amount of substrate stock solution (1%) can significantly slow down the reaction rate.

### *Loss of Activity of Protein of Interest*

If your fusion protein is particularly sensitive to degradation or to loss of activity, you can try reducing the labeling time or decreasing the labeling temperature. If you label at 4°C we recommend overnight incubation.

## Appendix:

### Sequence of CLIP-tag Region of pCLIP<sub>f</sub> Vector

Unique restriction sites in the regions flanking the CLIP<sub>f</sub> gene are displayed above the coding strand. The complete sequence of pCLIP<sub>f</sub> and pCLIP<sub>f</sub>-H2B can be downloaded at [www.neb.com](http://www.neb.com).

```

      NheI   EcoRV AscI           SwaI   BsrGI   AgeI   EcoRI
...GCTAGC GATATCGGCG CGCCAGCATT TAAATCTGTA CAGACCGGTG AATTC
      CGATCG CTATAGCCGC GCGGTCGTAA ATTTAGACAT GTCTGGCCAC TTAAG...
```

```

      SbfI   BamHI           PmeI   XhoI           PacI           NotI
...CCTGCA GGCGGATCCG CGTTTAACT CGAGGTTAAT TAATGAGCGG CCGC
      GGACGT CCGCCTAGGC GCAAATTTGA GCTCCAATTA ATTACTCGCC GGCG...
```

### Note

NEB 10-beta Competent *E. coli* (High Efficiency) (NEB #C3019) is recommended for propagating and subcloning of the vector and control plasmid. If using the *Clal* site for subcloning, *dam-/dcm-* Competent *E. coli* (NEB #C2925) is recommended. The restriction endonuclease *Clal* (NEB #R0197) is methylation sensitive.

## Ordering Information

PRODUCT	NEB #	SIZE
CLIP-Cell™ Starter Kit	E9200S	10 reactions
COMPANION PRODUCTS		
CLIP-Vista Green	S9235S	50 nmol
CLIP-Biotin	S9221S	50 nmol
pCLIP <sub>f</sub> -NK1R Control Plasmid	N9216S	20 µg
pCLIP <sub>f</sub> -Cox8A Control Plasmid	N9217S	20 µg
CLIP-Cell™ 360	S9215S	50 nmol
CLIP-Cell™ 430	S9216S	50 nmol
CLIP-Cell™ 505	S9217S	50 nmol
CLIP-Cell™ Fluorescein	S9218S	50 nmol
CLIP-Cell™ TMR-Star	S9219S	30 nmol
CLIP-Cell™ Block	S9220S	100 nmol
CLIP-tag™ Purified Protein	P9311S	100 µg
Anti-SNAP-tag® Antibody (Polyclonal)	P9310S	100 µl

**For Research Use Only.** No other use is authorized, including without limitation Commercial Use. Commercial Use means any and all uses of the product and derivatives by a party for monetary or other consideration including, but not limited to: (1) product manufacture; (2) providing of a service, information or data; (3) resale of the product, derivatives or components; (4) use of the product, derivatives or components by therapeutic, diagnostic or prophylactic clinical or commercial purposes. Information regarding any other license may be obtained by contacting Business Development and Licensing, New England Biolabs, Inc., 240 County Road, Ipswich, MA 01938. Tel: 978-927-5054, Fax: 978-921-1350.

**Notice To Buyer/User:** The Buyer/User has a non-exclusive license to use this system or any component thereof for **RESEARCH AND DEVELOPMENT PURPOSES ONLY**. Commercial use of this system or any components thereof requires a license from New England Biolabs, Inc., 240 County Road, Ipswich, MA 01938. For detailed information: see [www.neb.com/cia/legal](http://www.neb.com/cia/legal).

FuGENE® is a registered trademark of Roche.







## USA

New England Biolabs, Inc.  
240 County Road  
Ipswich, MA 01938-2723  
Telephone: (978) 927-5054  
Toll Free: (USA Orders) 1-800-632-5227  
Toll Free: (USA Tech) 1-800-632-7799  
Fax: (978) 921-1350  
e-mail: [info@neb.com](mailto:info@neb.com)  
[www.neb.com](http://www.neb.com)

## Canada

New England Biolabs, Ltd.  
Telephone: (905) 837-2234  
Toll Free: 1-800-387-1095  
Fax: (905) 837-2994  
Fax Toll Free: 1-800-563-3789  
e-mail: [info@ca.neb.com](mailto:info@ca.neb.com)

## China, People's Republic

New England Biolabs (Beijing), Ltd.  
Telephone: 010-82378265/82378266  
Fax: 010-82378262  
e-mail: [info@neb-china.com](mailto:info@neb-china.com)

## France

New England Biolabs France  
Free Call: 0800/100 632  
Free Fax: 0800/100 610  
e-mail: [info@fr.neb.com](mailto:info@fr.neb.com)

## Germany

New England Biolabs GmbH  
Telephone: +49/(0)69/305 23140  
Free Call: 0800/246 5227 (Germany)  
Fax +49/(0)69/305 23149  
Free Fax: 0800/246 5229 (Germany)  
e-mail: [info@de.neb.com](mailto:info@de.neb.com)

## Japan

New England Biolabs Japan, Inc.  
Telephone: +81 (0)3 5669 6191  
Fax +81 (0)3 5669 6192  
e-mail: [info@neb-japan.com](mailto:info@neb-japan.com)

## United Kingdom

New England Biolabs (UK), Ltd.  
Telephone: (01462) 420616  
Call Free: 0800 318486  
Fax: (01462) 421057  
Fax Free: 0800 435682  
e-mail: [info@uk.neb.com](mailto:info@uk.neb.com)