

# NotI



1-800-632-7799  
info@neb.com  
www.neb.com



R0189S 052120814081

## R0189S



**500 units**      **10,000 U/ml**      **Lot: 0521208**

**RECOMBINANT**    **Store at -20°C**    **Exp: 8/14**

### Recognition Site:

5' . . . G C G G C C G C . . . 3'  
3' . . . C G C C G G C G . . . 5'

**Source:** An *E. coli* strain that carries the cloned NotI gene from *Nocardia otitidis-caviarum* (ATCC 14630)

**Also Available In  
High Fidelity (HF™) Format**

# NotI



1-800-632-7799  
info@neb.com  
www.neb.com



R0189S 052120814081

## R0189S



**500 units**      **10,000 U/ml**      **Lot: 0521208**

**RECOMBINANT**    **Store at -20°C**    **Exp: 8/14**

### Recognition Site:

5' . . . G C G G C C G C . . . 3'  
3' . . . C G C C G G C G . . . 5'

**Source:** An *E. coli* strain that carries the cloned NotI gene from *Nocardia otitidis-caviarum* (ATCC 14630)

**Also Available In  
High Fidelity (HF™) Format**

Supplied in: 200 mM NaCl, 20 mM Tris-HCl (pH 7.4), 0.1 mM EDTA, 1 mM DTT, 0.15% Triton X-100, 200 µg/ml BSA and 50% glycerol.

**Reagents Supplied with Enzyme:**  
10X NEBuffer 3, 100X BSA.

**Reaction Conditions:** 1X NEBuffer 3, supplemented with 100 µg/ml BSA. Incubate at 37°C.

**1X NEBuffer 3:**  
100 mM NaCl  
50 mM Tris-HCl  
10 mM MgCl<sub>2</sub>  
1 mM DTT  
pH 7.9 @ 25°C

**Unit Definition:** One unit is defined as the amount of enzyme required to digest 1 µg of pBC4 DNA in 1 hour at 37°C in a total reaction volume of 50 µl.

**Diluent Compatibility:** Diluent Buffer C  
250 mM NaCl, 10 mM Tris-HCl, 0.1 mM EDTA, 1 mM DTT, 0.15% Triton X-100, 200 µg/ml BSA and 50% glycerol (pH 7.4 @ 25°C).

Supplied in: 200 mM NaCl, 20 mM Tris-HCl (pH 7.4), 0.1 mM EDTA, 1 mM DTT, 0.15% Triton X-100, 200 µg/ml BSA and 50% glycerol.

**Reagents Supplied with Enzyme:**  
10X NEBuffer 3, 100X BSA.

**Reaction Conditions:** 1X NEBuffer 3, supplemented with 100 µg/ml BSA. Incubate at 37°C.

**1X NEBuffer 3:**  
100 mM NaCl  
50 mM Tris-HCl  
10 mM MgCl<sub>2</sub>  
1 mM DTT  
pH 7.9 @ 25°C

**Unit Definition:** One unit is defined as the amount of enzyme required to digest 1 µg of pBC4 DNA in 1 hour at 37°C in a total reaction volume of 50 µl.

**Diluent Compatibility:** Diluent Buffer C  
250 mM NaCl, 10 mM Tris-HCl, 0.1 mM EDTA, 1 mM DTT, 0.15% Triton X-100, 200 µg/ml BSA and 50% glycerol (pH 7.4 @ 25°C).

### Quality Control Assays

**Ligation:** After 10-fold overdigestion with NotI, > 95% of the DNA fragments can be ligated with T4 DNA Ligase (at a 5' termini concentration of 1–2 µM) at 16°C. Of these ligated fragments, > 95% can be recut.

**16-Hour Incubation:** A 50 µl reaction containing 1 µg of DNA and 200 units of enzyme incubated for 16 hours resulted in the same pattern of DNA bands as a reaction incubated for 1 hour with 1 unit of enzyme.

**Exonuclease Activity:** Incubation of 200 units of enzyme with 1 µg sonicated <sup>3</sup>H DNA (10<sup>5</sup> cpm/µg) for 4 hours at 37°C in 50 µl reaction buffer released < 0.1% radioactivity

**Endonuclease Activity:** Incubation of 200 units of enzyme with 1 µg φX174 RF I DNA for 4 hours at 37°C in 50 µl reaction buffer resulted in < 5% conversion to RF II.

### Quality Control Assays

**Ligation:** After 10-fold overdigestion with NotI, > 95% of the DNA fragments can be ligated with T4 DNA Ligase (at a 5' termini concentration of 1–2 µM) at 16°C. Of these ligated fragments, > 95% can be recut.

**16-Hour Incubation:** A 50 µl reaction containing 1 µg of DNA and 200 units of enzyme incubated for 16 hours resulted in the same pattern of DNA bands as a reaction incubated for 1 hour with 1 unit of enzyme.

**Exonuclease Activity:** Incubation of 200 units of enzyme with 1 µg sonicated <sup>3</sup>H DNA (10<sup>5</sup> cpm/µg) for 4 hours at 37°C in 50 µl reaction buffer released < 0.1% radioactivity

**Endonuclease Activity:** Incubation of 200 units of enzyme with 1 µg φX174 RF I DNA for 4 hours at 37°C in 50 µl reaction buffer resulted in < 5% conversion to RF II.

### Enzyme Properties

**Activity in NEBuffers:**  
NEBuffer 1    0%  
NEBuffer 2    50%  
NEBuffer 3    100%  
NEBuffer 4    25%

When using a buffer other than the optimal (supplied) NEBuffer, it may be necessary to add more enzyme to achieve complete digestion.

**Survival in a Reaction:** A minimum of 0.25 unit is required to digest 1 µg of substrate DNA in 16 hours.

**Heat Inactivation:** 50 units of enzyme were inactivated by incubation at 65°C for 20 minutes.

**Note:** Cleavage of mammalian genomic DNA is blocked by CpG methylation.

Supercoiled plasmids may require up to 5-fold more NotI for complete digestion than linear DNAs.

(see other side)

CERTIFICATE OF ANALYSIS

### Enzyme Properties

**Activity in NEBuffers:**  
NEBuffer 1    0%  
NEBuffer 2    50%  
NEBuffer 3    100%  
NEBuffer 4    25%

When using a buffer other than the optimal (supplied) NEBuffer, it may be necessary to add more enzyme to achieve complete digestion.

**Survival in a Reaction:** A minimum of 0.25 unit is required to digest 1 µg of substrate DNA in 16 hours.

**Heat Inactivation:** 50 units of enzyme were inactivated by incubation at 65°C for 20 minutes.

**Note:** Cleavage of mammalian genomic DNA is blocked by CpG methylation.

Supercoiled plasmids may require up to 5-fold more NotI for complete digestion than linear DNAs.

(see other side)

CERTIFICATE OF ANALYSIS

**Companion Products:**

NotI-HF™  
#R3189S 500 units  
#R3189L 2,500 units  
#R3189M 2,500 units

NotI-HF™ RE-Mix™  
#R5189S 25 reactions

 = Time-Saver™ Qualified (See [www.neb.com](http://www.neb.com) for details).

U.S. Patent No. 5,371,006

Page 2 (R0189)

**Companion Products:**

NotI-HF™  
#R3189S 500 units  
#R3189L 2,500 units  
#R3189M 2,500 units

NotI-HF™ RE-Mix™  
#R5189S 25 reactions

 = Time-Saver™ Qualified (See [www.neb.com](http://www.neb.com) for details).

U.S. Patent No. 5,371,006

Page 2 (R0189)