

ELISA for Mouse IgM

Product Code: 3845-1AD-6

CONTENTS:

Vial 1 (red top)

Monoclonal anti-IgM antibody (150 µl)

Concentration: 0.5 mg/ml

Vial 2 (yellow top)

ALP-conjugated anti-IgM antibody (150 µl)

Vial 3

Lyophilised mouse IgM standard

To ensure total recovery of stated quantity, vials have been overfilled.

STORAGE:


Shipped at ambient temperature. On arrival box 1 should be stored refrigerated at 4-8°C and box 2 should be stored frozen at -20°C.

General

Intended use: For quantitative determination of native mouse IgM in serum and plasma.

Reagents: Monoclonal anti-IgM antibody is supplied in sterile-filtered (0.2 μm) PBS with sodium azide (0.02%). ALP-conjugated anti-IgM antibody is supplied in 0.1 M Tris-buffer with 0.15% Kathon CG.

Recommended standard dilution: 0.1-500 ng/ml



Guidelines for Mouse IgM ELISA

- Day 1** 1. Coat a high protein binding ELISA plate with monoclonal anti-IgM antibody, diluted to 1 µg/ml in PBS, pH 7.4, by adding 100 µl/well. Incubate overnight at 4-8°C.
- Day 2** 2. Wash twice with PBS (200 µl/well).
3. Block plate by adding 200 µl/well of PBS with 0.05% Tween 20 (PBS-Tween) containing 0.1% BSA (incubation buffer*). Incubate for 1 hour at room temperature.
4. Wash five times with PBS-Tween.
5. Prepare mouse IgM standard by reconstituting contents of vial 3 in 500 µl PBS to make up a stock solution of 50 µg/ml. The stock solution should be used immediately or stored in aliquots at -20°C for future use. We recommend the aliquots not be refrozen after initial use. For the test, prepare dilutions of the stock using the standard range as a guideline.
6. Add 100 µl/well of samples or standards diluted in incubation buffer and incubate for 2 hours at room temperature.
7. Wash as in step 4.
8. Add 100 µl/well of anti-IgM-ALP diluted 1:500 in incubation buffer. Incubate for 1 hour at room temperature.
9. Wash as in step 4.
10. Add 100 µl/well of appropriate substrate solution e.g. p-nitrophenyl-phosphate (pNPP).
11. Measure the optical density (405 nm for pNPP) in an ELISA reader after suitable developing time.

* The same buffer is used for blocking and for dilution.

NOTE; for research use only.

MABTECH shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



MABTECH AB

Box 1233

SE-131 28 Nacka Strand

Sweden

Tel: +46 8 716 27 00

Fax: +46 8 716 27 01

E-mail: mabtech@mabtech.com

www.mabtech.com

Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the following standards:



MABTECH Inc

M.E.B. 220

3814 West Street

Cincinnati, OH 45227

USA

Tel: +1 513 871 4500

Fax: +1 513 871 7353

E-mail: mabtech.usa@mabtech.com

MABTECH AUSTRALIA Pty Ltd

44 Gresswell Road

Macleod, VIC 3085

Australia

Tel: +61 3 9459 9630

Fax: +61 3 9455 0084

E-mail: mabtech.au@mabtech.com

**MABTECH AB Büro Deutschland
Germany**

Tel: +49 40 4135 7935

Fax: +49 40 4135 7945

E-mail: mabtech.de@mabtech.com

MABTECH AB Bureau de liaison France

BP 255, 1300 route des Crêtes

06905 Sophia Antipolis

France

Tel: +33 (0)4 92 38 80 70

Fax: +33 (0)4 92 38 80 71

E-mail: mabtech.fr@mabtech.com

2012-10-25