Catalog Number: 151972

# Protease, Strain V8

**Molecular Weight:** 27,000<sup>7</sup> or 27,700 by sedimentation equilibrium.

**CAS** # 9001-92-7

Synonym: Endoproteinase Glu-C

**E.C.** 3.4.21.19

Source: Staphylococcus aureus strain V8

**Preparation:** The enzyme was isolated, by the method of Drapeau et al.<sup>4</sup>, from the culture filtrate of S. aureus V8 strain. The preparation is homogenous on ultracentrifugation and disc electrophoresis.

**Extinction Coefficient:**  $E^{1\%}_{280} = 4.26^{10}$ 

**Optimum pH:** 4.0 and 7.8 with hemoglobin substrate.<sup>4</sup>

**Inhibitors:** Diisopropyl fluorophosphate (DFP) and monovalent anions such as F<sup>-</sup>, Cl<sup>-</sup>, Br<sup>-</sup>, CH<sub>3</sub>COO<sup>-</sup> and NO<sub>3</sub><sup>-</sup>. <sup>10</sup>

Physical Description: Off-white lyophilized solid

**Activity:** Approximately 500 units/mg

Unit Definition: One unit will hydrolyze case in to change  $A_{280}$  0.001 per minute, pH 7.8, at 37°C.

**Solubility:** Soluble in aqueous buffers (such as 10 mM Tris, pH 7.8 at 1 mg/ml). Solutions can be stored refrigerated for up to 1 month or aliquoted and stored at -20°C for approximately 6 months.

**Description:** Protease V8 is used for selective cleavage of proteins for amino acid sequence determination<sup>2</sup> or peptide mapping.<sup>3,9</sup> The *Staphylococcus* strain V8 protease specifically cleaves peptide bonds on the carboxyl (COOH) terminal side of aspartic and glutamic acid residues.<sup>4</sup>

## **Assay:**

Method: Enzyme activity is determined by the casein digestion assay described by Drapeau, et al.4

### Reagents:

- 1. 1% Casein in 0.05 M Tris·PO<sub>4</sub> buffer, pH 7.8: Dissolve 1 gram Hammersten grade casein in 50 ml 0.01 N NaOH with gentle heating and stirring. Add 40 ml reagent grade water and 5.0 ml 1.0 M Tris. Adjust pH to 7.8 with H<sub>3</sub>PO<sub>4</sub> and q.s. to 100 ml.
- 2. 10% Trichloroacetic acid (TCA)

## Enzyme:

Dissolve at 1 mg/ml in reagent grade water.

#### Procedure:

Equilibrate a series of tubes with 5.0 ml of 1% casein at  $37^{\circ}$ C for 5 minutes. At zero time add 10 ul or 20 ul of enzyme. Mix. Include a reagent blank. Exactly ten minutes after adding sample, stop reaction by adding 5.0 ml TCA. Mix. Allow tubes to stand ten minutes and then filter. Read  $A_{280}$  of the filtrate.

## Calculation:

Units/mg = 
$$\frac{A_{280} \text{ (test)} - A_{280} \text{ (reagent blank)}}{10 \text{ minutes x mg of enzyme in reaction}}$$

Click Here for a list of protease inhibitors offered by MP Biomedicals and general protease inhibitor information.

### **References:**

- 1. Drapeau, G.R., Meth. Enzymol., v. 45, 469 (1976).
- 2. Drapeau, G.R., Meth. Enzymol., v. 47, 189 (1977).
- 3. Cleveland, D.W., et al., "Peptide mapping by limited proteolysis in sodium dodecyl sulfate and analysis by gel electrophoresis." *J. Biol. Chem.*, **v. 252**, 1102 (1977).
- 4. Drapeau, G.A., Boily, Y. and Houmard, J., "Purification and properties of an extracellular protease of *Staphylococcus aureus*." *J. Biol. Chem.*, v. 247, 6720 (1972).
- 5. Kunitz, M., J. Gen. Physiol., v. 30, 291 (1947).
- 6. Drapeau, G., "Protease from *Staphylococcus aureus*." in *Methods in Enzymology*, XLVB, (Lorand, L., ed.), 469 (1974).
- 7. Drapeau, G., "The primary structure of Staphylococcal protease." Can. J. Biochem., v. 56, 534 (1978).
- 8. Dugas, H., Gaudet, F. and Leduc, P., "Structural studies of Staphylococcal protease. III. Binding of anions to the spin-labeled enzyme." *Can. J. Biochem.*, **v. 56**, 7 (1977).
- 9. Hall, T., et al., "Messenger RNA in G1 protein of french bean seed: Cell-free translation and product characterization." *Proc. Natl. Acad. Sci. USA*, **v. 75**, 3196 (1978).
- 10. Houmard, J., "Kinetic investigation of the Staphylococcal protease-catalyzed hydrolysis of synthetic substrates." *Eur. J. Biochem.*, **v. 68**, 621 (1976).
- 11. Houmard, J. and Drapeau, G., "Staphylococcal protease: A proteolytic enzyme specific for glutamoyl bonds." *Proc. Natl. Acad. Sci. USA*, **v. 69**, 3506 (1972).