



PRODUCT INFORMATION

Thermo Scientific PowerCut Dicer

#F-602S 60 U

Lot _ Expiry Date _

Ordering information

Component	#F-602S	#F-602L
PowerCut Dicer, 1 U/ μ L	60 U (60 μ L)	300 U (300 μ L)
5X PowerCut Dicer Reaction Buffer	500 μ L	500 μ L

Store at -20°C

www.thermoscientific.com/onebio

Description

Thermo Scientific PowerCut Dicer is a recombinant endoribonuclease originating from *Giardia intestinalis*. It cleaves dsRNA to small interfering RNA (siRNA) with a 100% efficiency. The produced siRNA fragments have a length of 25–27 nucleotides, and are capable of triggering RNA interference when transfected into cells.

Definition of Activity Unit

One unit is defined as the amount of enzyme that is needed to completely cleave 1 μ g of 192 bp double-stranded RNA substrate to siRNA in 16 hours at 37°C.

Source

PowerCut Dicer enzyme is purified from an *E.coli* strain expressing the cloned PowerCut Dicer gene from *Giardia intestinalis*.

Storage buffer: 50 mM Tris-HCl pH 8 (25°C), 50 mM NaCl, 5 mM MgCl₂, 0.2 mg/ml BSA, 50% glycerol.

Reaction buffer: 5X PowerCut Dicer Reaction Buffer.

Rev.2

Guidelines for using PowerCut™ Dicer Reaction protocol

Make sure to avoid contaminating RNases when setting up the reactions.

Table 1. Pipetting instructions.

Component	20 μ L rxn	Final concentration
H ₂ O (RNase free)	Add to 20 μ L	
5X PowerCut Dicer Reaction Buffer	4 μ L	1X
dsRNA	X μ L	4 μ g
PowerCut Dicer (1 U/ μ L)	4 μ L	4 U

Incubate at 37°C for 16 h, preferably in a heat block.

Purify the siRNA using standard methods.

Agarose gel analysis of the reactions



Figure 1. A dicing reaction was performed using dsRNA from 193 bp furin, 812 bp eGFP and phi6 as substrates. The reactions were assembled and incubated according to the standard protocol, and unpurified aliquots were run on a 2% agarose gel. No uncut substrate is seen in the reactions demonstrating 100% cleaving efficiency. M denotes a 27 bp RNA marker.

Exonuclease contamination assay: Incubation of 1 U of PowerCut Dicer (4 h, 37°C, 50 µL) with 1 µg of sonicated [³H]-DNA (2 × 10⁵ cpm/µg) in the assay buffer released <0.5% of radioactivity.

Endonuclease contamination assay: Incubation of 1 U of PowerCut Dicer (4 h, 37°C, 50 µL) with 1 µg of ΦX174 RFI DNA in the assay buffer resulted in <10% conversion to RFII form.

Ribonuclease contamination assay: Incubation of short ssRNA in the presence or absence of 1 U of PowerCut Dicer (1 h, 37°C, 20 µL) in the assay buffer produced similar results.

Technical support

EMEA: ts.molbio.eu@thermofisher.com

Americas & APAC: ts.molbio@thermofisher.com

DISCLAIMER

The use of this product in some applications may require a license from third parties. Purchase of this product does not provide any such license. Further information on purchasing licenses to practice a specific technology must be obtained from the appropriate intellectual property owners of such technology.

PRODUCT USE LIMITATION

This product has been developed and is sold exclusively for research purposes and in *in vitro* use only. This product has not been tested for use in diagnostics or drug development, nor are they suitable for administration to humans or animals. Please refer to www.thermoscientific.com/onebio for Material Safety Data Sheet of the product.

© 2012 Thermo Fisher Scientific, Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.