

Annexin-binding Buffer for Flow Cytometry

Catalog no. V13246

Table 1. Contents and storage information.

Material	Amount	Concentration	Storage*	Stability
50 mM HEPES, 700 mM NaCl, 12.5 mM CaCl ₂ , pH 7.4	50 mL	5X	<ul style="list-style-type: none"> • 2–6°C • DO NOT FREEZE 	When stored as directed this product is stable for 1 year from the date of receipt.

Number of assays: Sufficient material is supplied for 500 assays based on the protocol below.

Introduction

Apoptosis is a carefully regulated process of cell death that occurs as a normal part of development. Inappropriately regulated apoptosis is implicated in disease states, such as Alzheimer's disease and cancer. Apoptosis is distinguished from necrosis, or accidental cell death, by characteristic morphological and biochemical changes, including compaction and fragmentation of the nuclear chromatin, shrinkage of the cytoplasm, and loss of membrane asymmetry.^{1–5}

In normal live cells, phosphatidyl serine (PS) is located on the cytoplasmic surface of the cell membrane. However, in apoptotic cells, PS is translocated from the inner to the outer leaflet of the plasma membrane, thus exposing PS to the external cellular environment.⁶ In leukocyte apoptosis, PS on the outer surface of the cell marks the cell for recognition and phagocytosis by macrophages.^{7,8} The human anticoagulant, annexin V, is a 35–36 kDa Ca²⁺-dependent phospholipid-binding protein that has a high affinity for PS.⁹ The binding of annexin V to PS requires the presence of Ca⁺⁺. Under defined salt and calcium conditions, annexin V labeled with a fluorophore or biotin can identify apoptotic cells by binding to PS exposed on the outer leaflet.¹⁰ Refer to our website at www.invitrogen.com for a wide selection of products for apoptosis research.




Figure 1. Jurkat cells (T-cell leukemia, human) treated with 10 μ M camptothecin for four hours (panel B) or untreated control (panel A). Alexa Fluor® 488 annexin V and Propidium iodide dead cell stain were used with the annexin-binding buffer in this example. Cells were stained and analyzed by flow cytometry using 488 nm excitation on the Attune™ Acoustic Cytometer, 530/30 and 575/24 bandpass filters and collected by means of a standard 100 μ L/minute collection rate. Note that the camptothecin treated cells have a higher percentage of apoptotic cells (panel B) than the basal level of apoptosis seen in the control cells (panel A). A =apoptotic cells, V = viable cells, N = necrotic cells.

Before Starting

Materials Required but Not Provided

Deionized water

Preparing 1X Buffer

Dilute 5X Annexin-binding Buffer 1:4 with deionized water before use. Store the diluted buffer at 2–8°C.

The final 1X concentration of Annexin-binding buffer is 100 mM HEPES, 140 mM NaCl, 25 mM CaCl_2 , pH 7.4.

Experimental Protocol

This assay is optimized using Jurkat cells treated with camptothecin to induce apoptosis. Some modifications may be required for use with other cell types.

This assay is optimized using traditional hydrodynamic focusing flow cytometers as well as using the Attune™ Acoustic Cytometer. Any collection rate (25 μ L/min through 1,000 μ L/min) may be used with the Attune™ Acoustic Cytometer.

1. Induce apoptosis in cells using the desired method. Prepare a negative control by incubating

cells in the absence of inducing agent.

2. Harvest the cells after the incubation period and wash in cold phosphate-buffered saline (PBS).
3. Prepare 1X annexin-binding buffer. For example, for ~10 assays, add 1 mL 5X annexin-binding buffer to 4 mL deionized water.
4. Re-centrifuge the washed cells (from step 2), discard the supernatant, and resuspend the cells in 1X annexin-binding buffer.
5. Determine the cell density and dilute in 1X annexin-binding buffer to $\sim 1 \times 10^6$ cells/mL, preparing a sufficient volume to have 100 μ L per assay.
6. Add appropriate amount of annexin V conjugate and dead cell dye to each 100 μ L of cell suspension.
7. Incubate the cells at room temperature for 15 minutes.
8. After the incubation period, add 400 μ L 1X annexin-binding buffer, mix gently and keep the samples on ice.
9. As soon as possible, analyze the stained cells by flow cytometry.

References

1. Immunol Cell Biol 76, 1 (1998); 2. Cytometry 27, 1 (1997); 3. J Pharmacol Toxicol Methods 37, 215 (1997); 4. FASEB J 9, 1277 (1995); 5. Am J Pathol 146, 3 (1995); 6. Cytometry 31, 1 (1998); 7. J Immunol 148, 2207 (1992); 8. J Immunol 151, 4274 (1993); 9. J Biol Chem 265, 4923 (1990); 10. Blood 84, 1415 (1994).

Product List Current prices may be obtained from our website or from our Customer Service Department.

Cat. no.	Product Name	Unit Size
V13241	Alexa Fluor® 488 annexin V/Dead Cell Apoptosis Kit with Alexa® Fluor 488 annexin V and PI *for flow cytometry* *50 assays*	1 kit
V13245	Alexa Fluor® 488 annexin V/Dead Cell Apoptosis Kit with Alexa® Fluor 488 annexin V and PI *for flow cytometry* *250 assays*	1 kit
Related Products		
V13240	Single Channel Annexin V/ Dead Cell Apoptosis Kit *Alexa Fluor® 488 annexin V/SYTOX® Green* *50 assays* *for flow cytometry* . .	1 kit
V13246	Annexin-binding buffer *5X concentrate* *for flow cytometry*	50 mL
V35112	PE Annexin V/ Dead Cell Apoptosis Kit *with SYTOX® Green* *50 assays* *for flow cytometry*	1 kit
V35113	APC Annexin V/Dead Cell Apoptosis Kit *with APC annexin V and SYTOX® Green* *50 assays* *for flow cytometry*	1 kit
V35114	Metabolic Activity/Annexin V/Dead Cell Apoptosis Kit *with C12 resazurin, APC annexin V, and SYTOX® Green* *50 assays* *for flow cytometry*	1 kit
V35116	Mitochondrial Membrane Potential/Annexin V Apoptosis Kit *Alexa Fluor® 488 annexin V/MitoTracker® Red CMXRos* *50 assays* *for flow cytometry*	1 kit
V35136	Violet Annexin V/Dead Cell Apoptosis Kit *Pacific Blue™ annexin V/SYTOX® AADvanced™* *for flow cytometry* *50 assays*	1 kit
V23200	Vybrant® Apoptosis Assay Kit #6 *biotin-X annexin V/Alexa Fluor® 350 streptavidin/propidium iodide* *50 assays*	1 kit
A13201	annexin V, Alexa Fluor® 488 conjugate *100 assays*	100 assays
A23204	annexin V, Alexa Fluor® 647 conjugate *100 assays*	100 assays
A35110	annexin V, allophycocyanin conjugate (APC annexin V) *50 assays*	50 assays
A35111	annexin V, R-phycoerythrin conjugate (R-PE annexin V) *50 assays*	50 assays
A35122	annexin V, Pacific Blue™ conjugate *for flow cytometry* *100 assays*	100 assays
A35137	Violet Ratiometric Membrane Asymmetry Probe/Dead Cell Apoptosis Kit *for flow cytometry* *100 assays*	1 kit

Contact Information

Molecular Probes, Inc.

29851 Willow Creek Road
Eugene, OR 97402
Phone: (541) 465-8300
Fax: (541) 335-0504

Customer Service:

6:00 am to 4:30 pm (Pacific Time)
Phone: (541) 335-0338
Fax: (541) 335-0305
probesorder@invitrogen.com

Toll-Free Ordering for USA:

Order Phone: (800) 438-2209
Order Fax: (800) 438-0228

Technical Service:

8:00 am to 4:00 pm (Pacific Time)
Phone: (541) 335-0353
Toll-Free (800) 438-2209
Fax: (541) 335-0238
probestech@invitrogen.com

Invitrogen European Headquarters

Invitrogen, Ltd.
3 Fountain Drive
Inchinnan Business Park
Paisley PA4 9RF, UK
Phone: +44 (0) 141 814 6100
Fax: +44 (0) 141 814 6260
Email: euroinfo@invitrogen.com
Technical Services: eurotech@invitrogen.com

**For country-specific contact information,
visit www.invitrogen.com.**

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Service Department in Eugene, Oregon.

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

Limited Use Label License No. 223: Labeling and Detection Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2010, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice. For Research Use Only. Not intended for any animal or human therapeutic or diagnostic use.