

Cell Therapy Systems Recombinant Human Interleukin 4 CTSTM (IL-4)

PRODUCT ANALYSIS SHEET

Catalog Number:	CTP0041	CTP0043	
Quantity:	100 µg	1 mg	
Lot Number:	See product label		
Molecular Weight:	14.9 kDa		
Purity:	>95% as determined by SDS-PAGE analysis.		
Amino Acid Sequence:	HKCDITLQEI IKTLNSLTEQ KTLCTELTVT DIFAASKNTT EKETFCRAAT VLRQFYSHHE KDTRCLGATA QQFHRHKQLI RFLKRLDRNL WGLAGLNSCP VKEANQSTLE NFLERLKTIM REKYSKCSS		
Biological Activity:	$ED_{50} = 0.05$ to 0.4 ng/mL (Specific Activity: 2.0 x 10 ⁷ to 2.5 x 10 ⁶ units/mg). The biological activity is determined by measuring the dose-dependent proliferation of human TF-1 cells. A concentration range of 0.1 to 10.0 ng/mL is effective for most <i>in vitro</i> applications.		
Formulation:	Lyophilized, carrier free.		
Sterility:	Filtered through a 0.22 micron sterile filter.		
Endotoxin:	<0.1 ng/µg		
Production:	Recombinant human IL-4 is produced in E. coli and purified via sequential chromatography.		
Reconstitution Recommendation:	We recommend that the vial be briefly centrifuged prior to opening to bring the contents to the bottom. Lyophilized human IL-4 may be reconstituted in sterile deionized H ₂ O to 0.1-1.0 mg/mL to regain full activity. These stock solutions should be apportioned into working aliquots and stored at $\leq -20^{\circ}$ C. Further dilution should be made in low endotoxin medium or buffered solution with FCS or tissue culture grade BSA. It is recommended that all culture media containing supplements, such as growth factor, be sterile filtered prior to use for cell, gene, or tissue-based applications to minimize risk of contamination.		
Suggested Working Dilutions:	The optimal concentration should be determined for each specific application.		
Storage:	Lyophilized human IL-4 should be stored at 2 to 8°C, preferably desiccated. Store reconstituted human IL-4 at ≤ -20 °C (not in a frost-free freezer). Keep freeze-thaw cycles to a minimum.		
Expiration Date:	Expires one year from date of receipt when stored as instructed.		

For Research Use or Non-Commercial Manufacturing of Cell Based Products for Clinical Research. CAUTION: Not intended for direct administration into humans or animals www.invitrogen.com

Manufactured under ISO 13485 Quality Standard

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288 For technical support or support related to CTSTM products, <u>www.invitrogen.com/celltherapysupport</u>

PICTS-Hu IL-4

(Rev 07/10) DCC-10-1460

Important Licensing Information - These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or our website, www.invitrogen.com). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses.



References:

Paul, W.E. (1991) Interleukin 4: A prototypical immunoregulatory lymphokine. Blood 77:1859-1870.

Chapoval, A.I., K. Tamada, and L.P. Chen (2000) In vitro growth inhibition of a broad spectrum of tumor cell lines by activated human dendritic cells. Blood 95(7):2346-2351.

Francisco, J.A., K.L. Donaldson, D. Chace, C.B. Siegall, and A.F. Wahl (2000) Agonistic properties and in vivo antitumor activity of the anti-CD40 antibody SGN-14. Cancer Reasearch 60(12):3225-3231.

Gubina, E., X. Luo, E. Kwon, K. Sakamoto, Y.F. Shi, and R.A. Mufson (2001) ßc cytokine receptor-induced stimulation of cAMP response element binding protein phosphorylation requires protein kinase C in myeloid cells: A novel cytokine signal transduction cascade. J. Immunol. 167(8):4303-4310.

Kahlert, H., E. Grage-Griebenow, H.T. Stuwe, O. Cromwell, and H. Fiebig (2000) T cell reactivity with allergoids: Influence of the type of APC. J. Immunol. 165(4):1807-1815.

Kim, D.-K., T.V. Lee, A. Catilleja, B.W. Anderson, G.E. Peoples, A.P. Kudelka, J.L. Murray, T. Sittisomwong, J.T. Wharton, J.-W. Kim, and C.G. Ioannides (1999) Folate binding protein peptide 191-199 presented on dendritic cells can stimulate CTL from ovarian and breast cancer patients. Anticancer Research 19:2907-2916.

Parada, N.A., D.M. Center, H. Kornfeld, W.L. Rodriguez, J. Cook, M. Vallen, and W.W. Cruikshank (1998) Synergistic activation of CD4+ T cells by IL-16 and IL-2. J. Immunol. 160(5):2115-2120.

Piccinini, G., A. Foli, G. Comolli, J. Lisziewicz, and F. Lori (2002) Complementary antiviral efficacy of hydroxyurea and protease inhibitors in human immunodeficiency virus-infected dendritic cells and lymphocytes. J. Virol.76 (5):2274-2278.

Rao, R.M., D.O. Haskard, and R.C. Landis (2002) Enhanced recruitment of TH2 and CLA-Negative lymphocytes by the S128R polymorphism of E-Selectin(1). J. Immunol. 169 (10):5860-5865.

Voburka, Z., V. Vetvicka, J. Vetvickova, and M. Fusek (2002) Cytokines affect procathepsin D-stimulated proliferation of breast cancer cells. Anticancer Res. 22(2A):913-919.

Wagers, A.J., C.M. Waters, L.M. Stoolman, and G.S. Kansas (1998) Interleukin 12 and interleukin 4 control T cell adhesion to endothelial selectins through opposite effects on alpha1,3-fucosyltransferase VII gene expression. J. Exp. Med. 188(12):2225-2231.

Wang, H.S., H.J. Cao, V.D. Winn, L.J. Rezanka, Y. Frobert, C.H. Evans, D. Sciaky, D.A. Young, and T.J. Smith (1996) Leukoregulin induction of prostaglandin-endoperoxide H synthase-2 in human orbital fibroblasts. An in vitro model for connective tissue inflammation. J. Biol. Chem. 271(37):22718-22728.

Yu, Y., M. Hagihara, K. Ando, B. Gansuvd, H. Matsuzawa, T. Tsuchiya, Y. Ueda, H. Inoue, T. Hotta, and S. Kato (2001) Enhancement of human cord blood CD34(+) cell-derived NK cell cytotoxicity by dendritic cells. J. Immunol. 166(3):1590-1600.

Explanation of symbols					
Symbol	Description	Symbol	Description		
REF	Catalogue Number	LOT	Batch code		
RUO	Research Use Only	IVD	In vitro diagnostic medical device		
X	Use by	ł	Temperature limitation		
***	Manufacturer	EC REP	European Community authorised representative		
[-]	Without, does not contain	[+]	With, contains		
from Light	Protect from light	Â	Consult accompanying documents		
[]i	Directs the user to consult instructions for use (IFU), accompanying the product.				

For Research Use or Non-Commercial Manufacturing of Cell Based Products for Clinical Research. CAUTION: Not intended for direct administration into humans or animals <u>www.invitrogen.com</u>

Manufactured under ISO 13485 Quality Standard

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288 For technical support or support related to CTS™ products, <u>www.invitrogen.com/celltherapysupport</u>

PICTS-Hu IL-4

(Rev 07/10) DCC-10-1460

Important Licensing Information - These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or our website, <u>www.invitrogen.com</u>). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses.

