Promega

Technical Bulletin

CellTiter 96[®] AQueous One Solution Cell Proliferation Assay

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MTS/MTT-Based Cell Viability Assay Systems

Product	Size	Cat.#
CellTiter 96 [®] AQueous Non-Radioactive Cell Proliferation Assay	1,000 assays	G5421
	5,000 assays	G5430
	50,000 assays	G5440
CellTiter 96 [®] AQueous MTS Reagent Powder*	250mg	G1112
	1g	G1111
CellTiter 96 [®] Non-Radioactive Cell Proliferation Assay	1,000 assays	G4000
	5,000 assays	G4100

For Laboratory Use. *PMS is not supplied with MTS Reagent Powder and must be obtained separately.

Luminescent-Based Cell Viability Assay System

Product	Size	Cat.#
CellTiter-Glo [®] Luminescent Cell Viability Assay	10ml	G7570
	10 × 10ml	G7571
	100ml	G7572
	10 × 100ml	G7573

Resazurin-Based Cell Viability Assay System

Product	Size	Cat.#
CellTiter-Blue [®] Cell Viability Assay	20ml	G8080
	100ml	G8081
	10 × 100ml	G8082

Fluorescent-Based Cell Viability Assay

Product	Size	Cat.#
CellTiter-Fluor [™] Cell Viability Assay	10ml	G6080
	5 × 10ml	G6081
	2 × 50ml	G6082

For Laboratory Use.

Cytotoxicity Assay Systems (LDH)

Product	Size	Cat.#
CytoTox-ONE™ Homogeneous Membrane Integrity Assay	200–800 assays	G7890
	1,000–4,000 assays	G7891
CytoTox 96® Non-Radioactive Cytotoxicity Assay*	1,000 assays	G1780
CytoTox-Glo™ Cytotoxicity Assay*	10ml	G9290
	5 × 10ml	G9291
	2 × 50ml	G9292

*For Laboratory Use.

Apoptosis Assay Systems

Product	Size	Cat.#
Apo-ONE® Homogeneous Caspase-3/7 Assay	1ml	G7792
	10ml	G7790
	100ml	G7791
Caspase-Glo® 2 Assay*	10ml	G0940
	50ml	G0941
Caspase-Glo® 6 Assay*	10ml	G0970
	50ml	G0971
Caspase-Glo® 3/7 Assay*	2.5ml	G8090
	10ml	G8091
	100ml	G8092
Caspase-Glo® 8 Assay*	2.5ml	G8200
	10ml	G8201
	100ml	G8202
Caspase-Glo® 9 Assay*	2.5ml	G8210
	10ml	G8211
	100ml	G8212
CaspACE™ Assay System, Colorimetric*	50 assays	G7351
	100 assays	G7220
DeadEnd™ Fluorometric TUNEL System	60 reactions	G3250
DeadEnd™ Colorimetric TUNEL System	40 reactions	G7130
	20 reactions	G7360

*For Laboratory Use.

Apoptosis Reagents

Product	Size	Cat.#
CaspACE™ FITC-VAD-FMK In Situ Marker	50µl	G7461
	125µl	G7462
Anti-ACTIVE® Caspase-3 pAb	50µl	G7481
Anti-Cytochrome c mAb	100µg	G7421
Anti-pS ⁴⁷³ Akt pAb	40µl	G7441
Anti-PARP p85 Fragment pAb	50µl	G7341
Caspase Inhibitor Z-VAD-FMK	125µl	G7232
	50µl	G7231
Caspase Inhibitor, Ac-DEVD-CHO	100µl	G5961

Viability and Cytotoxicity Assay

Product	Size	Cat.#
MultiTox-Fluor Multiplex Cytotoxicity Assay (live/dead cell protease activity determination)	10ml	G9200
	5 × 10ml	G9201
	2 × 50ml	G9202
CytoTox-Fluor™ Cytotoxicity Assay (dead cell protease activity determination)	10ml	G9260
	5 × 10ml	G9261
	2 × 50ml	G9262
MultiTox-Glo Multiplex Cytotoxicity Assay (live/dead cell protease activity determination)	10ml	G9270
	5 × 10ml	G9271
	2 × 50ml	G9272

For Laboratory Use.

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1. Barltrop, J.A. *et al.* (1991) 5-(3-carboxymethoxyphenyl)-2-(4,5-dimethylthiazolyl)-3-(4-sulfophenyl)tetrazolium, inner salt (MTS) and related analogs of 3-(4,5-dimethylthiazolyl)-2,5-diphenyltetrazolium bromide (MTT) reducing to purple water-soluble formazans as cell-viability indicators. *Bioorg. Med. Chem. Lett.* **1**, 611-4.
2. Berridge, M.V. and Tan, A.S. (1993) Characterization of the cellular reduction of 3-(4,5-dimethylthiazol-2-yl)-2,5-diphenyltetrazolium bromide (MTT): Subcellular localization, substrate dependence, and involvement of mitochondrial electron transport in MTT reduction. *Arch. Biochem. Biophys.* **303**, 474-82.
3. Cory, A.H. *et al.* (1991) Use of an aqueous soluble tetrazolium/formazan assay for cell growth assays in culture. *Cancer Commun.* **3**, 207-12.
4. Riss, T.L. and Moravec, R.A. (1992) Comparison of MTT, XTT, and a novel tetrazolium compound for MTS for in vitro proliferation and chemosensitivity assays. *Mol. Biol. Cell (Suppl.)* **3**, 184a.
5. Mosmann, T. (1983) Rapid colorimetric assay for cellular growth and survival: Application to proliferation and cytotoxicity assays. *J. Immunol. Methods* **65**, 55-63.
6. Bernabei, P.A. *et al.* (1989) In vitro chemosensitivity testing of leukemic cells: Development of a semiautomated colorimetric assay. *Hematol. Oncol.* **7**, 243-53.
7. *CellTiter 96® Non-Radioactive Cell Proliferation Assay Technical Bulletin #TB112*, Promega Corporation.
8. Chen, C.-H. *et al.* (1990) MTT colorimetric assay detects mitogen responses of spleen but not blood lymphocytes. *Int. Arch. Allergy Appl. Immunol.* **93**, 249-55.

(a) The MTS tetrazolium compound is the subject of U.S. Pat. No. 5,185,450 assigned to the University of South Florida and is licensed exclusively to Promega Corporation.

© 1996-2007 Promega Corporation. All Rights Reserved.

Anti-ACTIVE, Apo-ONE, Caspase-Glo, CellTiter 96, CellTiter-Blue, CellTiter-Glo and CytoTox 96 are registered trademarks of Promega Corporation. CaspACE, CellTiter-Fluor, CytoTox-Fluor, CytoTox-Glo, CytoTox-ONE and DeadEnd are trademarks of Promega Corporation.

Products may be covered by pending or issued patents or may have certain limitations. Please visit our Web site for more information.

NEN is a registered trademark of NEN Life Science Products, Inc.

All prices and specifications are subject to change without prior notice.

Product claims are subject to change. Please contact Promega Technical Services or access the Promega online catalog for the most up-to-date information on Promega products.