

# TnT® Quick Coupled Transcription/Translation Systems

INSTRUCTIONS FOR USE OF PRODUCTS L1170, L1171, L2080 AND L2081.

# Quick PROTOCOL

## Transcription/Translation Procedure with Plasmid DNA

### Before You Begin

Upon removal from storage at  $-70^{\circ}\text{C}$ , rapidly thaw the TnT® Quick Master Mix by hand and place on ice. Thaw all other components at room temperature and store on ice.

### Preparation of Template

The template should be free of ethanol, calcium, RNase and salt. DNA from the Wizard® Plus Minipreps DNA Purification System, the Wizard® PCR Preps System or the standard alkaline lysate method (Sambrook *et al.*) will work with TnT® reactions.

### Translation Procedure

1. Assemble the reaction components, appropriate for the label being used, in a 0.5ml microcentrifuge tube. Gently mix by pipetting or stirring with a pipette tip and, if necessary, centrifuge briefly.


Components	Standard Reaction	Transcend™ Reaction	FluoroTect™ Reaction
TnT® T7 Quick Master Mix	40µl	40µl	40µl
Methionine, 1mM	—	1µl	1µl
[ <sup>35</sup> S]methionine (1,000Ci/mmol at 10mCi/ml)	2µl	—	—
Plasmid DNA Template (0.5µg)	2µl	2µl	2µl
Transcend™ Biotin-Lysyl-tRNA	—	1–2µl	—
FluoroTect™ Green <sub>Lys</sub> tRNA	—	—	1–2µl
Nuclease-Free Water to a final volume of	<u>50µl</u>	<u>50µl</u>	<u>50µl</u>

2. Incubate the reaction at  $30^{\circ}\text{C}$  for 60–90 minutes.
3. Analyze the results. Procedures for incorporation assays, gel analysis of translation products and an assay for luciferase production in the control reactions are provided in the *TnT® Quick Coupled Transcription/Translation Systems Technical Manual #TM045*.

\*See additional protocol information in Technical Manual #TM045, available online at: [www.promega.com](http://www.promega.com)




Keep all components on ice.




Assemble reaction components. Gently mix. Return unused components to  $-70^{\circ}\text{C}$ .



Centrifuge briefly if necessary.



Incubate at  $30^{\circ}\text{C}$  for 60–90 minutes.



Analyze.

2933MA04\_0A

### ORDERING/TECHNICAL INFORMATION:

[www.promega.com](http://www.promega.com) • Phone 608-274-4330 or 800-356-9526 • Fax 608-277-2601

© 2000–2010 Promega Corporation. All Rights Reserved.



# Promega

Printed in USA. Revised 9/10.  
Part #9FB030

# TnT® Quick Coupled Transcription/Translation Systems

INSTRUCTIONS FOR USE OF PRODUCTS L1170, L1171, L2080 AND L2081.

# Quick PROTOCOL

## Transcription/Translation Procedure with PCR-Generated DNA

### Before You Begin

Upon removal from storage at  $-70^{\circ}\text{C}$ , rapidly thaw the TnT® Quick Master Mix by hand and place on ice. Thaw all other components at room temperature and store on ice.

### Template Considerations

PCR products (5–7 $\mu\text{l}$ ) can be used directly from the amplification reaction.

**Note:** For PCR-generated templates, include 1 $\mu\text{l}$  of the T7 TnT® PCR Enhancer in each 50 $\mu\text{l}$  reaction.

### Translation Procedure

1. Assemble the reaction components, appropriate for the label being used, in a 0.5ml microcentrifuge tube. Gently mix by pipetting or stirring with pipette tip and, if necessary, centrifuge briefly.

Components	Standard Reaction	Transcend™ Reaction	FluoroTect™ Reaction
TnT® T7 Quick Master Mix	40 $\mu\text{l}$	40 $\mu\text{l}$	40 $\mu\text{l}$
Methionine, 1mM	–	1 $\mu\text{l}$	1 $\mu\text{l}$
[ <sup>35</sup> S]methionine (1,000Ci/mmol at 10mCi/ml)	2 $\mu\text{l}$	–	–
PCR-Generated DNA			
Template (0.5 $\mu\text{g}$ )	2.4–5 $\mu\text{l}$	2.5–5 $\mu\text{l}$	2.5–5 $\mu\text{l}$
T7 TnT® PCR Enhance	1 $\mu\text{l}$	1 $\mu\text{l}$	1 $\mu\text{l}$
Transcend™ Biotin-Lysyl-tRNA	–	1–2 $\mu\text{l}$	–
FluoroTect™ Green <sub>Lys</sub> tRNA	–	–	1–2 $\mu\text{l}$
Nuclease-Free Water to a final volume of	<u>50<math>\mu\text{l}</math></u>	<u>50<math>\mu\text{l}</math></u>	<u>50<math>\mu\text{l}</math></u>

2. Incubate the reaction at  $30^{\circ}\text{C}$  for 60–90 minutes.
3. Analyze the results. Procedures for incorporation assays, gel analysis of translation products and an assay for luciferase production in the control reactions are provided in the *TnT® Quick Coupled Transcription/Translation Systems Technical Manual #TM045*.

\*See additional protocol information in Technical Manual #TM045, available online at: [www.promega.com](http://www.promega.com)



Keep all components on ice.




Assemble reaction components. Gently mix. Return unused components to  $-70^{\circ}\text{C}$ .



Centrifuge briefly if necessary.



Incubate at  $30^{\circ}\text{C}$  for 60–90 minutes.



Analyze.

2933MA04\_0A

### ORDERING/TECHNICAL INFORMATION:

[www.promega.com](http://www.promega.com) • Phone 608-274-4330 or 800-356-9526 • Fax 608-277-2601



**Promega**

Printed in USA. Revised 9/10.  
Part #9FB030

© 2000–2010 Promega Corporation. All Rights Reserved.