Glu-C, Sequencing Grade:

 Part No.
 Size

 V165A
 50μg (5 × 10μg)

Description: Glu-C, Sequencing Grade (Cat.# V1651), is a serine protease that specifically cleaves at the C-terminus of either aspartic or glutamic acid residues (1–4). In ammonium bicarbonate and ammonium acetate, the enzyme specificity is higher at the glutamic residues. In phosphate buffers, cleavage occurs at the aspartic and glutamic residues (1,3,4). Glu-C activity is optimal in the pH range of 4.0–9.0. This sequencing grade enzyme can be used alone or in combination with other proteases to produce protein digests for peptide mapping applications or protein identification by peptide mass fingerprinting or MS/MS spectral matching. Glu-C, Sequencing Grade, is suitable for in-solution digestion reactions.

Biological Source: S. aureus V8.

Form: Lyophilized in the presence of trehalose (1.65mg/vial).

Molecular Weight: 27kDa.

Storage Conditions: Store at 4°C. See the Product Information Label for expiration date.

Usage Note: Reconstitute with filtered double-distilled water. Reconstituted enzyme solution can be stored at 4° C for 1 week or at -20° C for 8 weeks. Reconstituted enzyme can be frozen and thawed 5 times without a loss of activity.

Quality Control Assays

This lot passes the following Quality Control specifications:

Activity: β -insulin substrate is digested at a ratio of 1:10 protease: substrate for 1 hour at 37°C. Reaction products are analyzed by HPLC. After 1 hour, all substrate is consumed and the expected Glu-C-dependent peaks are observed.

Specificity: No prominent nonspecific peaks are observed by reverse-phase HPLC analysis after 18 hours of incubation of β -insulin with Glu-C compared to the 1-hour digestion.

Usage Information on Back

References

- 1. Drapeau, G., Boily, Y. and Houmard, J. (1972) J. Biol. Chem. 247, 6720-6.
- 2. Drapeau, G. (1978) *Can. J. Biochem.* **56**, 534–44.
- 3. Drapeau, G.R. (1977) Meth. Enzymol. 47, 189-91.
- 4. Drapeau, G.R. (1976) Meth. Enzymol. 45, 469-75.
- 5. Houmard, J. (1976) Eur. J. Biochem. 68, 621-7.

Stevens

J. Stevens, Quality Assurance

Part# 9PIV165 Revised 3/13





Promega Corporation

2800 Woods Hollow Road	
Madison, WI 53711-5399	USA
Telephone	608-274-4330
Toll Free	800-356-9526
Fax	608-277-2516
Internet	www.promega.com

PRODUCT USE LIMITATIONS, WARRANTY, DISCLAIMER

Promega manufactures products for a number of intended uses. Please refer to the product label for the intended use statements for specific products. Promega products contain chemicals which may be harmful if misused. Due care should be exercised with all Promega products to prevent direct human contact.

Each Promega product is shipped with documentation stating specifications and other technical information. Promega products are warranted to meet or exceed the stated specifications. Promega's sole obligation and the customer's sole remedy is limited to replacement of products free of charge in the event products fail to perform as warranted. Promega makes no other warranty of any kind whatsoever, and SPECIFICALLY DIS-CLAIMS AND EXCLUDES ALL OTHER WARRANTIES OF ANY KIND OR NATURE WHATSOEVER, DIRECTLY OR INDIRECTLY, EXPRESS OR IMPLIED, INCLUDING, WITHOUT LIMITATION, AS TO THE SUITABILITY, PRO-DUCTIVITY, DURABILITY, FITNESS FOR A PARTICU-LAR PURPOSE OR USE, MERCHANTABILITY, CONDI-TION, OR ANY OTHER MATTER WITH RESPECT TO PROMEGA PRODUCTS. In no event shall Promega be liable for claims for any other damages, whether direct, incidental, foreseeable, consequential, or special (including but not limited to loss of use, revenue or profit), whether based upon warranty, contract, tort (including negligence) or strict liability arising in connection with the sale or the failure of Promega products to perform in accordance with the stated specifications.

© 2010, 2012, 2013 Promega Corporation. All Rights Reserved.

ProteaseMAX is a trademark of Promega Corporation

Products may be covered by pending or issued patents or may have certain limitations. Please visit our Web site for more information.

All specifications are subject to change without prior notice.

Product claims are subject to change. Please contact Promega Technical Services or access the Promega online catalog for the most up-to-date information on Promega products.

Part# 9PIV165 Printed in USA. Revised 3/13

Signed by:



Usage Information

1. Preparation of Protein

In general, proteins require efficient solubilization, denaturation, disulphide bond reduction and alkylation for optimal digestion. The following optional steps are provided as a guideline to facilitate protease digestion with this product.

- Solubilization/Denaturation: Dissolve protein in compatible buffer (Table 1). Proteins that are difficult to dissolve or require denaturation for efficient digestion can be solubilized in a minimal volume of denaturant, such as 6–8M urea or 6M guanidine HCI. Incubate at 20–37°C for up to 1 hour. Some proteins may benefit from heating the sample to 60°C for one hour, or 95°C for 15–20 minutes.
- 2. **Disulphide Reduction:** Add DTT (or β -mercaptoethanol) to a final concentration of 5mM. Heat to 50–60°C for 20 minutes.
- Alkylation: Cool to room temperature. Add iodoacetamide to a final concentration of 15mM. Incubate in the dark at room temperature for 15 minutes. Adjust the reaction volume with buffer (pH 7.5) to reduce the final component concentrations for optimal digestion (Table 1).
- Digestion: Add Glu-C, Sequencing Grade, to a final proteinase:protein ratio of 1:200 to 1:20 (w/w). Vortex to mix and centrifuge. Incubate sample at 37°C for 2–18 hours.

Table 1. Denaturant Compatibility and Reaction Conditions. The denaturing components listed were titrated into Glu-C reactions. Percent change of enzyme activity in the presence of the denaturants was determined by comparison to a reaction containing no denaturant. Denaturants were tested using a β -casein activity assay.

Component	Concentration	% Activity			
Urea	0.4M	150			
	0.8M	150			
	1.4M	3			
Guanidine HCI	0.18M	100			
	0.33M	57			
SDS	0.004%	90			
	0.009%	90			
	0.016%	75			
Acetonitrile	10%	100			
	20%	200			
	25%	250			
	30%	74			
Additional Information					
Recommended Enzyme:Protein Ratio 1:20–200		1:20–200			
		<0700			
keconninended remperature ≤37°C		S31 ⁻ U			
pH Range	Range 4.0-9.0				
Freeze-Thaw Cycles	eze-Thaw Cycles 5				

Inhibitors

2. Related Products

Product	Size	Conc.	Cat.#
Asp-N, Sequencing Grade	2µg		V1621
Arg-C, Sequencing Grade	10µg		V1881
Chymotrypsin, Sequencing Grade	25µg		V1061
	100µg (4 × 25µg)		V1062
Elastase	5mg		V1891
Endo H	10,000u	500u/µl	V4871
	50,000u	500u/µl	V4875
Endoproteinase Lys-C, Sequencing Grade	5µg		V1071
Fetuin	500µg	10mg/ml	V4961
Immobilized Trypsin	2ml		V9012
	4ml (2 × 2ml)		V9013
Pepsin	250mg		V1959
PNGase F	500u	10u/µl	V4831
ProteaseMAX [™] Surfactant, Trypsin Enhancer	1mg		V2071
	5×1 mg		V2072
Protein Deglycosylation Mix	20 reactions		V4931
rLys-C, Mass Spec Grade	15µg		V1671
Sequencing Grade Modified Trypsin	100µg (5 × 20µg)		V5111
Sequencing Grade Modified Trypsin, Frozen	100µg (5 × 20µg)		V5113
Thermolysin	25mg		V4001
Trypsin Gold, Mass Spectrometry Grade	100µg		V5280
Trypsin/Lys-C Mix, Mass Spec Grade	20µg		V5071
	100µg		V5072
	100µg (5 × 20µg)		V5073

Diisopropyl fluorophosphate (DFP) Monovalent anions such as F^- , CI^- , Br^- , CH_3COO^- , and NO_3^- (1,5)