



Promega

Technical Bulletin

pGEM[®]-7Zf(+) Vector

INSTRUCTIONS FOR USE OF PRODUCT P2251.



www.promega.com

PRINTED IN USA.
Revised 12/06



AF9TB048 1206TB048

Part# TB048

pGEM[®]-7Zf(+) Vector

All technical literature is available on the Internet at: www.promega.com/tbs/
 Please visit the web site to verify that you are using the most current version of this
 Technical Bulletin. Please contact Promega Technical Services if you have questions on use
 of this system. E-mail: techserv@promega.com

I. Description.....	1
II. Product Components and Storage Conditions	2
III. pGEM [®] -7Zf(+) Vector Multiple Cloning Region and Circle Map.....	2
IV. pGEM [®] -7Zf(+) Vector Restriction Sites	4
V. Related Products	6
VI. Reference	7

I. Description

The pGEM[®]-7Zf(+) Vector is a derivative of the pGEM[®]-3Zf(+) Vector and contains the origin of replication of the filamentous phage f1. The plasmid serves as a standard cloning vector, as a template for in vitro transcription and as a template for the production of circular ssDNA. The plasmid contains SP6 and T7 RNA polymerase promoters flanking a multiple cloning region within the α -peptide coding region of β -galactosidase (1). Insertional inactivation of the α -peptide allows recombinant clones to be directly identified by color screening on indicator plates. The multiple cloning region is unique and includes restriction sites for ApaI, AatII, SphI, XbaI, XhoI, EcoRI, KpnI, SmaI, Csp45I, ClaI, HindIII, BamHI, SacI, BstXI and NsiI. This arrangement is designed specifically for generation of unidirectional deletions with the Erase-a-Base[®] System. The polylinker contains restriction enzyme sites that produce 5' overhangs or blunt ends (sensitive to Exonuclease III) flanked on both sides by blocks of restriction sites that generate 3' overhangs (resistant to Exonuclease III).

For induction of ssDNA, bacterial cells containing pGEM[®]-7Zf(+) Vector recombinants are infected with an appropriate helper phage. The plasmid then enters the f1 replication mode, and the resulting ssDNA is exported from the cell as an encapsidated virus-like particle. The sequence of the ssDNA rescued upon infection with helper phage is complementary to the sequence shown in Figure 1. The exported ssDNA can be used for mutagenesis in vitro or can be sequenced using the T7 Promoter Primer (Cat.# Q5021) or pUC/M13 Forward Primer (Cat.# Q5391, Q5601).

The sequences of Promega vectors are available online at:
www.promega.com/vectors/ and from the GenBank[®] database.

II. Product Components and Storage Conditions

Product	Size	Cat.#
pGEM®-7Zf(+) Vector	20µg	P2251

The pGEM®-7Zf(+) Vector is provided with a glycerol stock of bacterial strain JM109. The JM109 cells do not contain the vector and are not competent.

Storage Conditions: Store the pGEM®-7Zf(+) Vector at -20°C and the glycerol stock of JM109 cells at -70°C.

III. pGEM®-7Zf(+) Vector Multiple Cloning Region and Circle Map




Figure 1. pGEM®-7Zf(+) Vector promoter and multiple cloning region sequence. The sequence shown corresponds to RNA synthesized by T7 RNA polymerase and is complementary to RNA synthesized by SP6 RNA polymerase. The strand shown is complementary to the ssDNA strand produced by this vector.




Figure 2. pGEM®-7Zf(+) Vector circle map and sequence reference points. The pGEM®-7Zf(+) and pGEM®-7Zf(-) Vectors are identical except for the orientation of the f1 origin. Use the T7 Promoter Primer or pUC/M13 Forward Primer to sequence ssDNA produced by the pGEM®-7Zf(+) Vector.

pGEM®-7Zf(+) Vector sequence reference points:

T7 RNA polymerase transcription initiation site	1
SP6 RNA polymerase transcription initiation site	123
T7 RNA polymerase promoter (-17 to +3)	2981-3
SP6 RNA polymerase promoter (-17 to +3)	121-140
multiple cloning region	10-110
binding site of pUC/M13 Reverse Sequencing Primer	158-174
<i>lacZ</i> start codon	162
<i>lac</i> operon sequences	2818-2978; 148-377
<i>lac</i> operator	182-198
β -lactamase (Amp^r) coding region	1319-2179
phage f1 region	2362-2817
binding site of pUC/M13 Forward Sequencing Primer	2938-2954

Specialized applications of the pGEM®-7Zf(+) Vector:

- Used with the Erase-a-Base® System
- ssDNA production
- Blue/white screening for recombinants
- Transcription in vitro from dual-opposed promoters (For protocol information, please request the *Riboprobe® in vitro Transcription Systems Technical Manual*, #TM016.)
- Translation in vitro (For protocol information, please request the *TnT® Quick Coupled Transcription/Translation System Technical Manual*, #TM045.)

Note: All Promega technical literature is available on the Internet at:
www.promega.com

IV. pGEM®-7Zf(+) Vector Restriction Sites

The following restriction enzyme tables were constructed using DNASTAR® sequence analysis software. Please note that we have not verified this information by restriction digestion with each enzyme listed. The location given specifies the 3' end of the cut DNA (the base to the left of the cut site). For more information on the cut sites of these enzymes, or if you identify a discrepancy, please contact your local Promega Branch or Distributor. In the U.S., contact Promega Technical Services at 800-356-9526. The vector sequence is available in the GenBank® database (GenBank®/EMBL Accession Number X65310) and on the Internet at: www.promega.com/vectors/

Table 1. Restriction Enzymes That Cut the pGEM®-7Zf(+) Vector Between 1 and 5 Times.

Enzyme	# of Sites	Location	Enzyme	# of Sites	Location
AatII	1	20	BstXI	1	100
AccIII	1	81	Cfr10I	2	1472, 2687
Acc65I	1	49	ClaI	1	67
Acyl	2	17, 1929	Csp45I	1	61
AflIII	2	96, 499	DdeI	4	774, 1183, 1349, 1889
Alw26I	2	1453, 2229	DraI	3	1258, 1277, 1969
Alw44I	2	813, 2059	DraIII	1	2586
AlwNI	1	915	DrdI	2	607, 2541
Apal	1	14	EaeI	3	338, 1780, 2967
AspHI	4	91, 817, 1978, 2063	EarI	3	383, 2187, 2875
AvaI	2	37, 54	EclHKI	1	1392
AvaII	2	1530, 1752	EcoICRI	1	89
BamHI	1	78	EcoRI	1	43
BanI	4	49, 243, 1340, 2623	FokI	5	116, 1358, 1539, 1826, 2913
BanII	3	14, 91, 2661	FspI	2	1614, 2837
BbuI	1	26	HaeII	4	377, 747, 2737, 2745
BglI	2	1512, 2830	HgaI	4	610, 1188, 1918, 2803
BsaI	1	1453	HindIII	1	72
BsaAI	1	2586	Hsp92I	2	17, 1929
BsaHI	2	17, 1929	KpnI	1	53
BsaJI	5	53, 54, 238, 659, 2933	MaeI	5	32, 994, 1247, 1582, 2737
BsaOI	5	415, 839, 1762, 1911, 2858	MluI	1	96
Bsp120I	1	10	MspAII	5	323, 841, 1086, 2027, 2887
BspHI	2	1219, 2227	NaeI	1	2689
BssSI	2	672, 2056			
BstOI	5	239, 527, 648, 661, 2934			

Table 1. Restriction Enzymes That Cut the pGEM®-7Zf(+) Vector Between 1 and 5 Times (continued).

Enzyme	# of Sites	Location	Enzyme	# of Sites	Location
NciI	5	55, 56, 879, 1575, 1926	ScaI	1	1872
NgoMIV	1	2687	SinI	2	1530, 1752
NsiI	1	109	SmaI	1	56
NspI	2	26, 503	SphI	1	26
PaeR7I	1	37	SspI	2	2196, 2378
Ppu10I	1	105	TfiI	2	334, 474
PspAI	1	54	VspI	3	270, 329, 1564
PvuI	2	1762, 2858	XbaI	1	31
PvuII	2	323, 2887	XhoI	1	37
RsaI	2	51, 1872	XmaI	1	54
SacI	1	91	XmnI	1	1991

Table 2. Restriction Enzymes That Do Not Cut the pGEM®-7Zf(+) Vector.

AccI	BsmI	Eco8II	PacI	SpeI
AccB7I	BspMI	EcoNI	PflMI	SpII
AflII	BsrGI	EcoRV	PinAI	SrfI
AgeI	BssHIII	EheI	PmeI	Sse8387I
AscI	Bst1107I	FseI	PmlI	StuI
AvrII	Bst98I	HincII	PpuMI	StyI
BalI	BstEII	HindII	PshAI	SwaI
BbeI	BstZI	HpaI	Psp5II	Tth111I
BbrPI	Bsu36I	I-PpoI	PstI	XcmI
BbsI	CspI	KasI	RsrII	
BclI	DraII	NarI	SacII	
BglII	DsaI	NcoI	SalI	
BlpI	EagI	NdeI	SfiI	
Bpu1102I	Eco47III	NheI	SgfI	
BsaBI	Eco52I	NotI	SgrAI	
BsaMI	Eco72I	NruI	SnaBI	

Table 3. Restriction Enzymes That Cut the pGEM®-7Zf(+) Vector 6 or More Times.

AcI	CfoI	HphI	MspI	SfaNI
AluI	DpnI	Hsp92II	NdeII	TaqI
BbvI	DpnII	MaeII	NlaIII	Tru9I
Bsp1286I	Fnu4HI	MaeIII	NlaIV	XhoII
BsrI	HaeIII	MboI	PleI	
BsrSI	HhaI	MboII	Sau3AI	
Bst7II	HinfI	MnII	Sau96I	
BstUI	HpaII	MseI	ScrFI	

Note: The enzymes listed in boldface type are available from Promega.

V. Related Products

pGEM® Vectors

Product	Size	Cat.#
pGEM®-3Z Vector	20µg	P2151
pGEM®-4Z Vector	20µg	P2161
pGEM®-3Zf(+) Vector	20µg	P2271
pGEM®-3Zf(-) Vector	20µg	P2261
pGEM®-5Zf(+) Vector	20µg	P2241
pGEM®-5Zf(-) Vector	20µg	P2351
pGEM®-7Zf(-) Vector	20µg	P2371
pGEM®-9Zf(-) Vector	20µg	P2391
pGEM®-11Zf(+) Vector	20µg	P2411
pGEM®-11Zf(-) Vector	20µg	P2421
pGEM®-13Zf(+) Vector	20µg	P2541

All pGEM® Vectors are provided with a glycerol stock of bacterial strain JM109. The JM109 cells do not contain vector and are not competent.

Other Vectors

Product	Size	Cat.#
pSP64 Poly(A) Vector	20µg	P1241
pSP72 Vector	20µg	P2191
pSP73 Vector	20µg	P2221

Sequencing Primers

Product	Size	Cat.#
SP6 Promoter Primer	2µg	Q5011
T7 Promoter Primer	2µg	Q5021
pUC/M13 Primer, Reverse (17mer)	2µg	Q5401
pUC/M13 Primer, Forward (17mer)	2µg	Q5391
pUC/M13 Primer, Forward (24mer)	2µg	Q5601
pUC/M13 Primer, Reverse (22mer)	2µg	Q5421

Related Systems

Product	Cat.#
Erase-a-Base® System	E5750
Riboprobe® System – SP6*	P1420
Riboprobe® System – T7*	P1440
TNT® T7 Quick Coupled Transcription/Translation System*	L1170
TNT® SP6 Quick Coupled Transcription/Translation System*	L2080

*For Laboratory Use.

VI. Reference

1. Yanish-Perron, C. *et al.* (1985) Improved M13 phage cloning vectors and host strains: Nucleotide sequences of the M13mp18 and pUC19 vectors. *Gene* **33**, 103-19.

© 1988–2006 Promega Corporation. All Rights Reserved.

Erase-a-Base, pGEM, Riboprobe and TNT are registered trademarks of Promega Corporation.

DNASTAR is a registered trademark of DNASTAR, Inc. GenBank is a registered trademark of the U.S. Department of Health and Human Services.

Products may be covered by pending or issued patents or may have certain limitations. Please visit our Web site for more information.

All prices and specifications are subject to change without prior notice.

Product claims are subject to change. Please contact Promega Technical Services or access the Promega online catalog for the most up-to-date information on Promega products.