



Promega

Technical Bulletin

pGEM[®]-3Zf(-) Vector

INSTRUCTIONS FOR USE OF PRODUCT P2261.



www.promega.com

PRINTED IN USA.
Revised 12/06



AF9TB045 1206TB045

Part# TB045

pGEM[®]-3Zf(-) Vector

All technical literature is available on the Internet at: www.promega.com/tbs/
Please visit the web site to verify that you are using the most current version of this Technical Bulletin. Please contact Promega Technical Services if you have questions on use of this system. E-mail: techserv@promega.com

I. Description.....	1
II. Product Components and Storage Conditions	2
III. pGEM [®] -3Zf(-) Vector Multiple Cloning Region and Circle Map	2
IV. pGEM [®] -3Zf(-) Vector Restriction Sites.....	4
V. Related Products	6
VI. Reference	7

I. Description

The pGEM[®]-3Zf(-) Vector is a derivative of the pGEM[®]-3Z Vector and contains the origin of replication of the filamentous phage f1. The plasmid serves as a standard cloning vector, as a template for in vitro transcription and as a template for the production of circular ssDNA.

The pGEM[®]-3Zf(-) Vector contains SP6 and T7 RNA polymerase promoters flanking the multiple cloning region within the α -peptide coding region of β -galactosidase (1). Insertional inactivation of the α -peptide allows recombinant clones to be directly identified by color screening on indicator plates. For induction of ssDNA, bacterial cells containing pGEM[®]-3Zf(-) recombinants are infected with an appropriate helper phage. The plasmid then enters the f1 replication mode, and the resulting ssDNA is exported from the cell as an encapsidated virus-like particle. The sequence of the ssDNA rescued upon infection with helper phage is identical to the sequence shown in Figure 1. The exported ssDNA can be used for mutagenesis in vitro or can be sequenced using the SP6 Promoter Primer (Cat.# Q5011) or pUC/M13 Reverse Primer (Cat.# Q5401, Q5421).

The sequences of Promega vectors are available online at:
www.promega.com/vectors/ and from the GenBank[®] database.

II. Product Components and Storage Conditions

Product	Size	Cat.#
pGEM [®] -3Zf(-) Vector	20µg	P2261

The pGEM[®]-3Zf(-) Vector is provided with a glycerol stock of bacterial strain JM109. The JM109 cells do not contain the vector and are not competent.

Storage Conditions: Store the pGEM[®]-3Zf(-) Vector at -20°C and the glycerol stock of JM109 cells at -70°C.

III. pGEM[®]-3Zf(-) Vector Multiple Cloning Region and Circle Map




Figure 1. pGEM[®]-3Zf(-) Vector promoter and multiple cloning region sequence.

The sequence shown corresponds to RNA synthesized by T7 RNA polymerase and is complementary to RNA synthesized by SP6 RNA polymerase. The strand shown is the same as the ssDNA strand produced by this vector.




Figure 2. pGEM[®]-3Zf(-) Vector circle map and sequence reference points. The pGEM[®]-3Zf(-) and pGEM[®]-3Zf(+) Vectors are identical except for the orientation of the f1 origin. Use the SP6 or pUC/M13 Reverse prime to sequence ssDNA produced by the pGEM[®]-3Zf(-) Vector.

pGEM[®]-3Zf(-) Vector sequence reference points:

T7 RNA polymerase transcription initiation site	1
multiple cloning region	5-61
SP6 RNA polymerase promoter (-17 to +3)	67-86
SP6 RNA polymerase transcription initiation site	69
<i>lac</i> operon sequences	94-323; 3018-3178
binding site of pUC/M13 Reverse Sequencing Primer	104-120
<i>lacZ</i> start codon	108
<i>lac</i> operator	128-144
β -lactamase (Amp ^r) coding region	1265-2125
phage f1 region	2562-3017
binding site of pUC/M13 Forward Sequencing Primer	3138-3154
T7 RNA polymerase promoter (-17 to +3)	3181-3

Specialized applications of the pGEM[®]-3Zf(-) Vector:

- ssDNA production.
- Blue/white screening for recombinants.
- Transcription in vitro from dual-opposed promoters (For protocol information, please request the *Riboprobe[®] in vitro Transcription Systems Technical Manual*, #TM016.)
- Translation in vitro (For protocol information, please request the *TNT[®] Quick Coupled Transcription/Translation System Technical Manual*, #TM045.)

Note: All Promega technical literature is available on the Internet at:
www.promega.com

IV. pGEM[®]-3Zf(-) Vector Restriction Sites

The following restriction enzyme tables were constructed using DNASTAR[®] sequence analysis software. Please note that we have not verified this information by restriction digestion with each enzyme listed. The location given specifies the 3' end of the cut DNA (the base to the left of the cut site). For more information on the cut sites of these enzymes, or if you identify a discrepancy, please contact your local Promega Branch or Distributor. In the U.S., contact Promega Technical Services at 800-356-9526. Vector sequences are available in the GenBank[®] database (GenBank[®]/EMBL Accession Number X65307) and on the Internet at: www.promega.com/vectors/

Table 1. Restriction Enzymes That Cut the pGEM[®]-3Zf(-) Vector Between 1 and 5 Times.

Enzyme	# of Sites	Location	Enzyme	# of Sites	Location
AatII	1	2260	Cfr10I	2	1418, 2687
AccI	1	39	DraI	3	1204, 1223, 1915
Acc65I	1	17	DraII	1	2314
AcyI	2	1875, 2257	DraIII	1	2795
AflIII	1	445	DrdI	3	553, 2422, 2839
Alw26I	4	1399, 2175, 2328, 2370	EaeI	3	284, 1726, 3167
Alw44I	3	759, 2005, 2502	EarI	3	329, 2133, 3075
AlwNI	1	861	EclHKI	1	1338
AspHI	5	15, 763, 1924, 2009, 2506	EcoICRI	1	13
AvaI	1	21	EcoRI	1	5
AvaII	2	1476, 1698	FokI	5	1304, 1485, 1772, 2415, 3113
BamHI	1	26	FspI	2	1560, 3037
BanI	4	17, 189, 1286, 2751	HaeII	4	323, 693, 2637, 2645
BanII	2	15, 2721	HgaI	5	556, 1134, 1864, 2422, 2570
BbuI	1	54	HincII	1	40
BglI	2	1458, 3030	HindII	1	40
Bsa I	1	1399	HindIII	1	56
BsaAI	1	2792	Hsp92I	2	1875, 2257
BsaHI	2	1875, 2257	KpnI	1	21
BsaJI	5	21, 22, 184, 605, 3133	MaeI	5	33, 940, 1193, 1528, 2639
BsaOI	5	361, 785, 1708, 1857, 3058	NaeI	1	2689
BspHI	3	1165, 2173, 2278	NdeI	1	2509
BspMI	1	51	NgoMIV	1	2687
BssSI	3	618, 2002, 2309	NspI	3	54, 449, 2366
BstOI	5	185, 473, 594, 607, 3134	PspAI	1	21
			PstI	1	48
			PvuI	2	1708, 3058

Table 1. Restriction Enzymes That Cut the pGEM®-3Zf(-) Vector Between 1 and 5 Times (continued).

Enzyme	# of Sites	Location	Enzyme	# of Sites	Location
PvuII	2	269, 3087	SspI	2	2142, 3000
RsaI	3	19, 1818, 2494	TaqI	5	9, 39, 545, 1989, 2757
SacI	1	15	TfiI	2	280, 420
SalI	1	38	VspI	3	216, 275, 1510
Scal	1	1818	XbaI	1	32
SinI	2	1476, 1698	XmaI	1	21
SmaI	1	23	XmnI	1	1937
SphI	1	54			
Sse8387I	1	48			

Table 2. Restriction Enzymes That Do Not Cut the pGEM®-3Zf(-) Vector.

AccII	BlpI	Bsu36I	FseI	PfIMI	SnaBI
AccB7I	Bpu1102I	Clal	HpaI	PinAI	SpeI
AflII	BsaBI	CspI	I-PpoI	PmeI	SplI
AgeI	BsaMI	Csp45I	KasI	PmlI	SrfI
ApaI	BsmI	DsaI	MluI	Ppu10I	StuI
AscI	Bsp120I	EagI	NarI	PpuMI	StyI
AvrII	BsrGI	Eco47III	NcoI	PshAI	Swal
BalI	BssHII	Eco52I	NheI	Psp5II	Tth111I
BbeI	Bst1107I	Eco72I	NotI	RsrII	XcmI
BbrPI	Bst98I	Eco81I	NruI	SacII	XhoI
BbsI	BstEII	EcoNI	NsiI	SfiI	
BclI	BstXI	EcoRV	PacI	SgfI	
BglII	BstZI	EheI	Paer7I	SgrAI	

Table 3. Restriction Enzymes That Cut the pGEM®-3Zf(-) Vector 6 or More Times.

AcI	CfoI	HpaII	MseI	Sau3AI
AluI	DdeI	HphI	MspI	Sau96I
BbvI	DpnI	Hsp92II	MspAII	ScrFI
Bsp1286I	DpnII	MaeII	NciI	SfaNI
BsrI	Fnu4HI	MaeIII	NdeII	Tru9I
BsrSI	HaeIII	MboI	NlaIII	XhoII
Bst7II	HhaI	MboII	NlaIV	
BstUI	HinfI	MnlI	PleI	

Note: The enzymes listed in boldface type are available from Promega.

V. Related Products

pGEM® Vectors

Product	Size	Cat.#
pGEM®-3Z Vector	20µg	P2151
pGEM®-4Z Vector	20µg	P2161
pGEM®-3Zf(+) Vector	20µg	P2271
pGEM®-5Zf(+) Vector	20µg	P2241
pGEM®-5Zf(-) Vector	20µg	P2351
pGEM®-7Zf(+) Vector	20µg	P2251
pGEM®-7Zf(-) Vector	20µg	P2371
pGEM®-9Zf(-) Vector	20µg	P2391
pGEM®-11Zf(+) Vector	20µg	P2411
pGEM®-11Zf(-) Vector	20µg	P2421
pGEM®-13Zf(+) Vector	20µg	P2541

All pGEM®-Zf Vectors are provided with a glycerol stock of bacterial strain JM109. The JM109 cells do not contain the vector and are not competent.

Other Vectors

Product	Size	Cat.#
pSP64 Poly(A) Vector	20µg	P1241
pSP72 Vector	20µg	P2191
pSP73 Vector	20µg	P2221

Sequencing Primers

Product	Size	Cat.#
SP6 Promoter Primer	2µg	Q5011
T7 Promoter Primer	2µg	Q5021
pUC/M13 Primer, Reverse (17mer)	2µg	Q5401
pUC/M13 Primer, Forward (17mer)	2µg	Q5391
pUC/M13 Primer, Forward (24mer)	2µg	Q5601
pUC/M13 Primer, Reverse (22mer)	2µg	Q5421

Riboprobe® in vitro Transcription Systems

Product	Cat.#
Riboprobe® System – SP6	P1420
Riboprobe® System – T3	P1430
Riboprobe® System – T7	P1440

For Laboratory Use.

TNT® Quick Coupled Transcription/Translation Systems

Product	Cat.#
TNT® T7 Quick Coupled Transcription/Translation System	L1170
TNT® T7 Quick Coupled Transcription/Translation System, Trial Size	L1171
TNT® SP6 Quick Coupled Transcription/Translation System	L2080
TNT® SP6 Quick Coupled Transcription/Translation System, Trial Size	L2081

For Laboratory Use.

VI. Reference

1. Yanisch-Perron, C., Vieira, J. and Messing, J. (1985) Improved M13 phage cloning vectors and host strains: Nucleotide sequences of the M13mp18 and pUC19 vectors. *Gene* **33**, 103-19.

© 1988-2006 Promega Corporation. All Rights Reserved.

pGEM, Riboprobe and TNT are registered trademarks of Promega Corporation.

DNASTAR is a registered trademark of DNASTAR, Inc. GenBank is a registered trademark of U.S. Department of Health and Human Services.

Products may be covered by pending or issued patents or may have certain limitations. Please visit our Web site for more information.

All prices and specifications are subject to change without prior notice.

Product claims are subject to change. Please contact Promega Technical Services or access the Promega online catalog for the most up-to-date information on Promega products.