



**Promega**

## Technical Bulletin

---

# pCI and pSI Mammalian Expression Vectors

INSTRUCTIONS FOR USE OF PRODUCTS E1721 AND E1731.



[www.promega.com](http://www.promega.com)

# pCI and pSI Mammalian Expression Vectors

All technical literature is available on the Internet at: [www.promega.com/tbs/](http://www.promega.com/tbs/)  
Please visit the web site to verify that you are using the most current version of this  
Technical Bulletin. Please contact Promega Technical Services if you have questions on use  
of this system. E-mail: [techserv@promega.com](mailto:techserv@promega.com).

1. Description.....	1
2. Product Components and Storage Conditions .....	2
3. pCI Mammalian Expression Vector Multiple Cloning Sequence and Map .....	2
4. pSI Mammalian Expression Vector Multiple Cloning Sequence and Map .....	4
5. Vector Components.....	5
A. Enhancer/Promoter Regions.....	5
B. Chimeric Intron.....	5
C. T7 Promoter.....	5
D. Multiple Cloning Region.....	6
E. SV40 Late Polyadenylation Signal .....	6
F. $\phi$ 1 Origin of Replication and Plasmid Replicon .....	6
6. Related Products .....	7
7. References .....	7
8. Appendix .....	9
A. pCI Mammalian Expression Vector Restriction Enzyme Sites.....	9
B. pSI Mammalian Expression Vector Restriction Enzyme Sites .....	11

## 1. Description

The pCI<sup>(a)</sup> and pSI Mammalian Expression Vectors are designed to promote constitutive expression of cloned DNA inserts in mammalian cells. The major difference between the two vectors is the enhancer/promoter region controlling the expression of the inserted gene. The pSI Expression Vector contains the simian virus 40 (SV40) enhancer and early promoter region, whereas the pCI Expression Vector contains the human cytomegalovirus (CMV) major immediate-early gene enhancer/promoter region.

The pSI and pCI Vectors can be used for both transient and stable expression of genes. For stable expression, the pSI or pCI Vectors must be co-transfected with an expression vector containing a selectable gene for mammalian cells.


## 2. Product Components and Storage Conditions

Product	Size	Cat. #
pCI Mammalian Expression Vector	20µg	E1731
pSI Mammalian Expression Vector	20µg	E1721


pCI and pSI Mammalian Expression Vectors are supplied frozen in TE buffer (pH 8.0).

**Storage Conditions:** Store the vectors at -20°C.

## 3. pCI Mammalian Expression Vector Multiple Cloning Sequence and Map



**Figure 1. Multiple cloning region sequence and T7 promoter of the pCI Mammalian Expression Vector.** The sequence shown corresponds to RNA synthesized by T7 RNA polymerase. The strands shown are the same as the ssDNA strands produced by the pCI Vector.




**pCI Mammalian Expression Vector sequence reference points:**

Cytomegalovirus immediate-early enhancer/promoter region	1-742
Chimeric intron	857-989
T7-EEV sequencing primer binding site	1020-1041
T7 RNA Polymerase Promoter (-17 to +2)	1034-1052
T7 promoter transcription start site	1051
Multiple cloning region	1052-1104
SV40 late polyadenylation signal	1111-1332
Phage f1 region	1422-1877
$\beta$ -lactamase ( $Amp^r$ ) coding region	2314-3174
<i>ColEI</i> -derived origin of replication	3936




Use the T7 EEV Promoter Primer (Cat.# Q6700) to sequence ssDNA produced by the pCI and pSI Vectors. **Do Not** use the T7 Promoter Primer (Cat.# Q5021) to sequence the pCI or pSI Vectors. There is a sequence difference between the primer and the promoter sequences.

#### 4. pSI Mammalian Expression Vector Multiple Cloning Sequence and Map



0771MA09\_4B

**Figure 2. Multiple cloning region sequence and T7 promoter of the pSI Mammalian Expression Vector.** The sequence shown corresponds to RNA synthesized by T7 RNA polymerase. The strands shown are the same as the ssDNA strands produced by the pSI Vector.



0864VA06\_4B

#### pSI Mammalian Expression Vector sequence reference points:

SV40 enhancer and early promoter region	1-419
Chimeric intron	483-615
T7 EEV sequencing primer binding site	646-667
T7 RNA Polymerase Promoter (-17 to +2)	660-678
T7 promoter transcription start site	677
Multiple cloning region	678-730
SV40 late polyadenylation signal	737-958
Phage f1 region	1048-1503
β-lactamase (Amp <sup>r</sup> ) coding region	1940-2800



Use T7 EEV Promoter Primer (Cat.# Q6700) to sequence ssDNA produced by the pCI and pSI Vectors. **Do Not** use the T7 Promoter Primer (Cat.# Q5021) to sequence the pCI or pSI Vectors. There is a sequence difference between the primer and the promoter sequences.

## 5. Vector Components

### 5.A. Enhancer/Promoter Regions

The SV40 enhancer/promoter region present in the pSI Vector and the CMV enhancer/promoter region present in the pCI Vector allow strong, constitutive expression in many cell types. The promiscuous nature of the CMV enhancer/promoter has been demonstrated in transgenic mice, where expression of the chloramphenicol acetyltransferase (CAT) gene under the regulation of the CMV enhancer/promoter was observed in 24 of the 28 tissues examined (1). The pSI Vector contains the SV40 origin of replication, which allows transient, episomal replication in cells expressing the SV40 large T antigen such as COS-1 or COS-7 cells (2).

### 5.B. Chimeric Intron

Downstream of the enhancer/promoter region is a chimeric intron composed of the 5'-donor site from the first intron of the human  $\beta$ -globin gene and the branch and 3'-acceptor site from the intron that is between the leader and the body of an immunoglobulin gene heavy chain variable region (3). The sequences of the donor and acceptor sites, along with the branchpoint site, have been changed to match the consensus sequences for splicing (4).

Transfection studies have demonstrated that the presence of an intron flanking the cDNA insert frequently increases the level of gene expression (5-8). The increase in expression level due to the presence of the intron depends on the particular cDNA insert. For example, in transient transfections of 293s cells, we have found that the presence of the chimeric intron in the pSI Vector results in a 21-fold increase in expression of the CAT gene. In contrast, the chimeric intron increases the gene expression level from the luciferase cDNA by only 3-fold.


The intron is located 5' to the cDNA insert in order to prevent utilization of possible cryptic 5'-donor splice sites within the cDNA sequence (9). In transgenic experiments, the presence of an intron is necessary to promote a high level of expression for virtually all cDNA inserts (10-12).

### 5.C. T7 Promoter

A T7 promoter is located downstream of the intron (i.e., immediately upstream of the multiple cloning region). This promoter can be used to synthesize RNA transcripts in vitro using T7 RNA Polymerase (Cat.# P2075).

#### 5.D. Multiple Cloning Region

The multiple cloning region is immediately downstream from the T7 promoter. The unique restriction sites available in the multiple cloning region of the pCI and pSI Vectors are identical, **except that the pSI Vector does not contain a unique *Kpn I* site within this region.** The sites in the multiple cloning region are compatible with subcloning cDNAs that have been prepared with the Universal RiboClone® cDNA Synthesis Systems (Cat.# C4360).

 **Note:** There are no ATG sequences in either the multiple cloning region or between the transcription start site and the multiple cloning region. Thus, an ATG for the initiation of translation must be present in the inserted DNA.

#### 5.E. SV40 Late Polyadenylation Signal

Polyadenylation signals cause the termination of transcription by RNA polymerase II and signal the addition of approximately 200-250 adenosine residues to the 3'-end of the RNA transcript (13). Polyadenylation has been shown to enhance RNA stability and translation (14,15). The late SV40 polyadenylation signal is extremely efficient and has been shown to increase the steady-state level of RNA approximately 5-fold more than the early SV40 polyadenylation signal (16).

#### 5.F. f1 Origin of Replication and Plasmid Replicon

The backbones for the pCI and pSI Vectors were derived from the pGEM®-3Zf(+) Vector. As a result, these vectors are high-copy plasmids and contain the origin of replication of the filamentous phage f1. For generation of single-stranded DNA (ssDNA) from the f1 origin, bacteria transformed with the pSI or pCI Vectors carrying the DNA insert of interest are infected with an appropriate helper phage. The plasmid then enters the f1 replication mode, and the resulting ssDNA is exported from the cell as an encapsidated virus particle. The ssDNA molecule exported has the sequence of the strand shown for the multiple cloning region (Figures 1 and 2). For further information, please contact your local Promega Branch Office or Distributor. In the U.S., contact Technical Services at 1-800-356-9526.

The "poison" sequence present in pBR322 that has been shown to inhibit replication of SV40 origin-containing vectors in COS cells has been deleted in the pCI and pSI Vectors (17). This results in more efficient expression of the cloned cDNAs in COS cells and other cells that have been transformed with the SV40 large T antigen.

## 6. Related Products

Product	Size	Cat. #
Flexi® Cloning System Entry/Transfer	5 entry, 20 transfer reactions	C8640
Flexi® Cloning System Transfer	100 transfer reactions	C8820
Carboxy Flexi® System, Transfer	50 transfer reactions	C9320
pCI-neo Mammalian Expression Vector	20µg	E1841
pAdVAntage™ Vector	20µg	E1711
T7 EEV Promoter Primer	2µg	Q6700
ProFection® Mammalian Transfection System—Calcium Phosphate	40 transfections	E1200
TransFast™ Transfection Reagent	1.2mg	E2431
Tfx™-50 Transfection Reagent	2.1mg	E1811
Tfx™-20 Transfection Reagent	4.8mg	E2391

Product	Size	Cat.#
Universal RiboClone® cDNA Synthesis System	1 system	C4360
T7 RNA Polymerase*	1,000 units	P2075

\*For Laboratory Use.

## 7. References

- Schmidt, E.V. *et al.* (1990) The cytomegalovirus enhancer: A pan-active control element in transgenic mice. *Mol. Cell. Biol.* **10**, 4406-11.
- Gluzman, Y. (1981) SV40-transformed simian cells support the replication of early SV40 mutants. *Cell* **23**, 175-82.
- Bothwell, A.L. *et al.* (1981) Heavy chain variable region contribution to the NPb family of antibodies: Somatic mutation evident in a gamma 2a variable region. *Cell* **24**, 625-37.
- Senapathy, P., Shapiro, M.B. and Harris, N.L. (1990) Splice junctions, branch point sites, and exons: Sequence statistics, identification, and applications to genome project. *Meth. Enzymol.* **183**, 252-78.
- Gross, M.K., Kainz, M.S. and Merrill, G.F. (1987) Introns are inconsequential to efficient formation of cellular thymidine kinase mRNA in mouse L cells. *Mol. Cell. Biol.* **7**, 4576-81.
- Buchman, A.R. and Berg, P. (1988) Comparison of intron-dependent and intron-independent gene expression. *Mol. Cell. Biol.* **8**, 4395-405.



## 7. References (continued)

7. Evans, M.J. and Scarpulla, R.C. (1989) Introns in the 3'-untranslated region can inhibit chimeric CAT and  $\beta$ -galactosidase gene expression. *Gene* **84**, 135-42.
8. Huang, M.T. and Gorman, C.M. (1990) Intervening sequences increase efficiency of RNA 3' processing and accumulation of cytoplasmic RNA. *Nucl. Acids Res.* **18**, 937-47.
9. Huang, M.T. and Gorman, C.M. (1990) The simian virus 40 small-t intron, present in many common expression vectors, leads to aberrant splicing. *Mol. Cell. Biol.* **10**, 1805-10.
10. Brinster, R.L. *et al.* (1988) Introns increase transcriptional efficiency in transgenic mice. *Proc. Natl. Acad. Sci. USA* **85**, 836-40.
11. Choi, T. *et al.* (1991) A generic intron increases gene expression in transgenic mice. *Mol. Cell. Biol.* **11**, 3070-4.
12. Palmiter, R.D. *et al.* (1991) Heterologous introns can enhance expression of transgenes in mice. *Proc. Natl. Acad. Sci. USA* **88**, 478-82.
13. Proudfoot, N. (1991) Poly(A) signals. *Cell* **64**, 671-4.
14. Bernstein, P. and Ross, J. (1989) Poly(A), poly(A) binding protein and the regulation of mRNA stability. *Trends Biochem. Sci.* **14**, 373-7.
15. Jackson, R.J. and Standart, N. (1990) Do the poly(A) tail and 3' untranslated region control mRNA translation? *Cell* **62**, 15-24.
16. Carswell, S. and Alwine, J.C. (1989) Efficiency of utilization of the simian virus 40 late polyadenylation site: Effects of upstream sequences. *Mol. Cell. Biol.* **9**, 4248-58.
17. Lusky, M. and Botchan, M. (1981) Inhibition of SV40 replication in simian cells by specific pBR322 DNA sequences. *Nature* **293**, 79-81.

## 8. Appendix

### 8.A. pCI Mammalian Expression Vector Restriction Enzyme Sites

The following restriction enzyme tables were constructed using DNASTAR® sequence analysis software. Please note that we have not verified this information by restriction digestion with each enzyme listed. The location given specifies the 3' end of the cut DNA (the base to the left of the cut site). For more information on the cut sites of these enzymes or if you identify a discrepancy, please contact your local Promega Branch or Distributor. In the U.S., contact Promega Technical Services at 800-356-9526. Vector sequences are available in the GenBank® database (GenBank®/EMBL Accession Number U47119) and online at: [www.promega.com/vectors/](http://www.promega.com/vectors/)

**Table 1. Restriction Enzymes That Cut the pCI Vector Between 1 and 5 Times.**

<b>Enzyme</b>	<b># of Sites</b>	<b>Location</b>	<b>Enzyme</b>	<b># of Sites</b>	<b>Location</b>
<b>AatII</b>	5	278, 331, 414, 600, 2182	<b>Bst98I</b>	2	820, 1017
<b>AccI</b>	1	1088	<b>BstOI</b>	5	243, 436, 3830, 3843, 3964
<b>Acc65I</b>	1	1075	<b>BstZI</b>	1	1098
<b>AflIII</b>	2	820, 1017	<b>Cfr10I</b>	2	1547, 3016
<b>AflIII</b>	1	1069	<b>ClaI</b>	1	1336
<b>Alw44I</b>	3	1932, 2429, 3675	<b>DraI</b>	4	1302, 2523, 3215, 3234
<b>AlwNI</b>	1	3580	<b>DraII</b>	1	2121
<b>AspHI</b>	5	721, 1936, 2433, 2518, 3679	<b>DraIII</b>	1	1655
<b>AvaI</b>	2	1058, 1092	<b>DrdI</b>	4	809, 1699, 2018, 3887
<b>AvaII</b>	2	2737, 2959	<b>DsaI</b>	1	513
<b>BalI</b>	2	10, 64	<b>EaeI</b>	4	8, 62, 1098, 2708
<b>BamHI</b>	1	1343	<b>EagI</b>	1	1098
<b>BanI</b>	5	618, 943, 1075, 1611, 3148	<b>EarI</b>	2	1360, 2302
<b>BanII</b>	2	721, 1581	<b>EclHKI</b>	1	3101
<b>BbsI</b>	1	928	<b>Eco52I</b>	1	1098
<b>BglIII</b>	1	4001	<b>EcoICRI</b>	1	719
<b>BsaI</b>	2	882, 3035	<b>EcoI</b>	1	1063
<b>BsaOI</b>	5	1101, 1382, 2583, 2732, 3655	<b>FokI</b>	5	950, 2019, 2662, 2949, 3130
<b>BsaAI</b>	2	493, 1652	<b>FspI</b>	2	1401, 2878
<b>BsaBI</b>	1	1342	<b>HaeII</b>	3	1497, 1505, 3749
<b>BsaJI</b>	3	513, 1092, 3829	<b>HincII</b>	3	669, 1089, 1241
<b>BsaMI</b>	2	1162, 1255	<b>HindII</b>	3	669, 1089, 1241
<b>BsmI</b>	2	1162, 1255	<b>HindIII</b>	1	748
<b>BspHI</b>	3	2156, 2261, 3269	<b>HpaI</b>	1	1241
<b>BspI</b>	1	844	<b>KpnI</b>	1	1079
<b>BsrGI</b>	1	96	<b>MluI</b>	1	1069
<b>BssSI</b>	3	2125, 2432, 3816			

**Note:** The enzymes listed in boldface type are available from Promega.

**Table 1. Restriction Enzymes That Cut the pCI Vector Between 1 and 5 Times (continued).**

Enzyme	# of Sites	Location	Enzyme	# of Sites	Location
MspA1	4	1999, 2465, 3406, 3651	SalI	1	1087
NaeI	1	1549	ScaI	2	1030, 2620
NcoI	1	513	SinI	2	2737, 2959
NdeI	2	387, 1927	SmaI	1	1094
NgoMIV	1	1547	SnaBI	1	493
NheI	1	1052	SpeI	1	152
NotI	1	1098	SspI	4	5, 52, 1860, 2296
NspI	1	2076	StyI	1	513
PaeR7I	1	1058	VspI	2	160, 2926
PspA1	1	1092	XbaI	1	1081
PstI	1	830	XhoI	1	1058
PvuI	2	1382, 2732	XmaI	1	1092
SacI	1	721	XmnI	1	2501

**Table 2. Restriction Enzymes That Do Not Cut the pCI Vector.**

AccB7I	Bsp120I	EcoRV	PmlI	SpII
AccIII	BssHIII	EheI	Ppu10I	SrfI
AgeI	Bst1107 I	FseI	PpuMI	Sse8387 I
ApaI	BstEII	I-PpoI	PshAI	StuI
AscI	BstXI	KasI	Psp5II	SwaI
AvrII	Bsu36I	NarI	PvuII	TfiI
BbeI	CspI	NruI	RsrII	Tth111 I
BbrPI	Csp45I	NsiI	SacII	XcmI
BbuI	Eco47III	Pacl	SfiI	
BclI	Eco72I	PfiMI	SgfI	
BlpI	Eco81I	PinAI	SgrAI	
Bpu1102 I	EcoNI	PmeI	SphI	

**Table 3. Restriction Enzymes That Cut the pCI Vector 6 or More Times.**

Acil	Bst7II	HinfI	MnlI	Sau96I
AcyI	BstUI	HpaII	MseI	ScrFI
AluI	CfoI	HphI	MspI	SfaNI
Alw26I	DdeI	Hsp92I	NciI	TaqI
BbvI	DpnI	Hsp92II	NdeII	Tru9I
BglI	DpnII	MaeI	NlaIII	XhoII
BsaHI	Fnu4HI	MaeII	NlaIV	
Bsp1286 I	HaeIII	MaeIII	PleI	
BsrI	Hgal	MboI	RsaI	
BsrSI	HhaI	MboII	Sau3A I	

**Note:** The enzymes listed in boldface type are available from Promega.

## 8.B. pSI Mammalian Expression Vector Restriction Enzyme Sites

The following restriction enzyme tables were constructed using DNASTAR® sequence analysis software. Please note that we have not verified this information by restriction digestion with each enzyme listed. The location given specifies the 3' end of the cut DNA (the base to the left of the cut site). For more information on the cut sites of these enzymes or if you identify a discrepancy, please contact your local Promega Branch or Distributor. In the U.S., contact Promega Technical Services at 800-356-9526. Vector sequences are available in the GenBank® database (GenBank®/EMBL Accession Number U47121) and online at: [www.promega.com/vectors/](http://www.promega.com/vectors/)

**Table 4. Restriction Enzymes That Cut the pSI Vector Between 1 and 5 Times.**

<b>Enzyme</b>	<b># of Sites</b>	<b>Location</b>	<b>Enzyme</b>	<b># of Sites</b>	<b>Location</b>
<b>AatII</b>	1	1808	BspHI	3	1782, 1887, 2895
<b>AccI</b>	1	714	BspMI	1	470
<b>Acc65I</b>	2	48, 701	BssSI	3	1751, 2058, 3442
AcyI	2	1805, 2187	<b>Bst98I</b>	2	446, 643
AflII	2	446, 643	<b>BstZI</b>	1	724
AflIII	1	695	Cfr10I	2	1173, 2642
<b>Alw44I</b>	3	1558, 2055, 3301	<b>Clal</b>	1	962
AlwNI	1	3206	<b>DraI</b>	4	928, 2149, 2841, 2860
AspHI	4	1562, 2059, 2144, 3305	DraII	1	1747
<b>AvaI</b>	2	684, 718	DraIII	1	1281
<b>AvaII</b>	2	2363, 2585	DrdI	4	435, 1325, 1644, 3513
AvrII	1	398	DsaI	2	9, 305
<b>BamHI</b>	1	969	EaeI	2	724, 2334
<b>BanI</b>	5	48, 569, 701, 1237, 2774	EagI	1	724
<b>BanII</b>	1	1207	Ear I	2	986, 1928
BbsI	1	554	<b>EclHKI</b>	1	2727
<b>BbuI</b>	2	146, 218	<b>Eco52I</b>	1	724
<b>BglI</b>	3	351, 1037, 2609	<b>EcoRI</b>	1	689
<b>BglII</b>	1	3627	FspI	2	1027, 2504
BsaI	2	508, 2661	<b>HaeII</b>	3	1123, 1131, 3375
BsaOI	5	727, 1008, 2209, 2358, 3281	HgaI	5	1056, 1637, 2195, 2925, 3503
BsaAI	1	1278	<b>HincII</b>	2	715, 867
BsaBI	1	968	HindII	2	715, 867
BsaHI	2	1805, 2187	<b>HindIII</b>	1	414
<b>BsaMI</b>	2	788, 881	<b>HpaI</b>	1	867
BsmI	2	788, 881	<b>Hsp92I</b>	2	1805, 2187
<b>Bsp1286I</b>	5	1207, 1562, 2059, 2144, 3305	<b>KpnI</b>	2	52, 705
			<b>MluI</b>	1	695

**Note:** The enzymes listed in boldface type are available from Promega.

**Table 4. Restriction Enzymes That Cut the pSI Vector Between 1 and 5 Times (continued).**

Enzyme	# of Sites	Location	Enzyme	# of Sites	Location
<b>MspAII</b>	5	74, 625, 2091, 3032, 3277	<b>RsaI</b>	5	50, 656, 703, 1570, 2246
<b>NaeI</b>	1	1175	<b>SaII</b>	1	713
<b>NcoI</b>	2	9, 305	<b>ScaI</b>	2	656, 2246
<b>NdeI</b>	1	1553	<b>SfiI</b>	1	351
<b>NgoMIV</b>	1	1173	<b>SinI</b>	2	2363, 2585
<b>NheI</b>	1	678	<b>SmaI</b>	1	720
<b>NotI</b>	1	724	<b>SphI</b>	2	146, 218
<b>NsiI</b>	2	148, 220	<b>SspI</b>	2	1486, 1922
<b>NspI</b>	3	146, 218, 1702	<b>StuI</b>	1	397
<b>PaeR7I</b>	1	684	<b>StyI</b>	3	9, 305, 398
<b>Ppu10I</b>	2	144, 216,	<b>TfiI</b>	1	420
<b>PspAI</b>	1	718	<b>VspI</b>	1	2552
<b>PstI</b>	1	456	<b>XbaI</b>	1	707
<b>PvuI</b>	2	1008, 2358	<b>XhoI</b>	1	684
<b>PvuII</b>	1	74	<b>XmaI</b>	1	718
			<b>XmnI</b>	1	2127

**Table 5. Restriction Enzymes That Do Not Cut the pSI Vector.**

<b>Acc B7I</b>	Bsp120I	Eco81I	PfiMI	SgrAI
<b>AccIII</b>	BsrGI	<b>EcoICRI</b>	PinAI	<b>SnaBI</b>
<b>AgeI</b>	<b>BssHIII</b>	EcoNI	PmeI	<b>SpeI</b>
<b>ApaI</b>	Bst1107 I	<b>EcoRV</b>	PmlI	SplI
AscI	<b>BstEII</b>	EheI	PpuMI	SrfI
<b>BalI</b>	<b>BstXI</b>	FseI	PshAI	Sse8387I
BbeI	<b>Bsu36I</b>	<b>I-PpoI</b>	Psp5II	Swal
BbrPI	<b>CspI</b>	KasI	RsrII	<b>Tth111I</b>
<b>BclI</b>	<b>Csp45I</b>	<b>NarI</b>	<b>SacI</b>	XcmI
BlpI	<b>Eco47III</b>	<b>NruI</b>	<b>SacII</b>	
Bpu1102 I	Eco72I	Pacl	<b>SgfI</b>	

**Note:** The enzymes listed in boldface type are available from Promega.

**Table 6. Restriction Enzymes That Cut the pSI Vector 6 or More Times.**

AcI	<b>CfoI</b>	<b>HpaII</b>	MseI	ScrFI
<b>AluI</b>	<b>DdeI</b>	HphI	<b>MspI</b>	SfaNI
<b>Alw26I</b>	<b>DpnI</b>	<b>Hsp92II</b>	<b>NciI</b>	<b>TaqI</b>
BbvI	DpnII	MaeI	<b>NdeII</b>	<b>Tru9I</b>
BsaJI	Fnu4HI	MaeII	NlaIII	<b>XhoII</b>
BsrI	FokI	MaeIII	NlaIV	
<b>Bsr SI</b>	<b>HaeIII</b>	MboI	PleI	
Bst71I	<b>HhaI</b>	<b>MboII</b>	<b>Sau3AI</b>	
BstOI	<b>HinfI</b>	Mnl I	Sau96I	

**Note:** The enzymes listed in boldface type are available from Promega.

<sup>(4)</sup> The CMV promoter and its use are covered under U.S. Pat. Nos. 5,168,062 and 5,385,839 owned by the University of Iowa Research Foundation, Iowa City, Iowa, and licensed FOR RESEARCH USE ONLY. Research Use includes contract research for which monetary or other consideration may be received. Other commercial users must obtain a license to these patents directly from the University of Iowa Research Foundation.

© 1994–2006 Promega Corporation. All Rights Reserved.

Flexi, Transfectam, pGEM, ProFection and RiboClone are registered trademarks of Promega Corporation. pAdvAntage, Tfx and TransFast are trademarks of Promega Corporation.

DNASTAR is a registered trademark of DNASTAR, Inc. GenBank is a registered trademark of U.S. Dept. of Health and Human Services.

Products may be covered by pending or issued patents or may have certain limitations. Please visit our Web site for more information.

All prices and specifications are subject to change without prior notice.

Product claims are subject to change. Please contact Promega Technical Services or access the Promega online catalog for the most up-to-date information on Promega products.