

## Platinum® PCR SuperMix

<b>Cat. nos.</b>	<b>Size</b>	<b>Store at -30°C to -10°C</b>
11306-016	100 reactions	
11306-081	5000 reactions	
<b>Pub. Part no.</b> 11306.pps	MAN0000950	<b>Rev. Date</b> 10 Apr 2012

### Description

Platinum® PCR SuperMix provides qualified reagents for the amplification of nucleic acid templates by polymerase chain reaction (PCR). The mixture contains anti-*Taq* DNA polymerase antibody, Mg<sup>++</sup>, dNTPs, and recombinant *Taq* DNA polymerase at concentrations sufficient to allow amplification during PCR. Platinum® PCR SuperMix is supplied at 1.1X concentration to allow approximately 10% of the final reaction volume to be used for the addition of primer and template solutions. Reagents sufficient for 100 or 5000 amplification reactions of 50 µL each are provided.

Anti-*Taq* DNA polymerase antibody inhibits polymerase activity providing an automatic “hot start” (Chou, 1992; Sharkey, 1994) and permits ambient temperature set-up. Antibody-mediated hot starts improve PCR specificity and yield (Westfall, 1997). Due to specific binding of the antibody, Platinum® PCR SuperMix is present in an inactive form and is reactivated after a denaturation step in PCR cycling at 94°C.

Platinum® PCR SuperMix may be stored at either -30°C to -10°C or 2°C to 8°C. Storage at 2°C to 8°C avoids the necessity of thawing the mix before assembling the PCR. No detectable reduction of PCR performance or enzyme activity is observed after storage of Platinum® PCR SuperMix for 12 months at 2°C to 8°C. Repeated freeze-thaw cycles do not reduce performance or activity.

Component	100-rxn size	5000-rxn size
Platinum® PCR SuperMix	4 × 1.125 mL	4 × 56.25 mL

**Product Use:** For research use only.

Not intended for any animal or human therapeutic or diagnostic use.

## Components

22 U/mL complexed recombinant *Taq* DNA polymerase with Platinum® *Taq* Antibody, 22 mM Tris-HCl (pH 8.4), 55 mM KCl, 1.65 mM MgCl<sub>2</sub>, 220 μM dGTP, 220 μM dATP, 220 μM dTTP, 220 μM dCTP, and stabilizers.

## Guidelines and Recommendations

Since PCR is a powerful technique capable of amplifying trace amounts of DNA, all appropriate precautions should be taken to avoid cross-contamination. Ideally, amplification reactions should be assembled in a DNA-free environment. Use of aerosol-resistant barrier tips is recommended. Take care to avoid contamination of PCR SuperMix with the primers or template DNA used in individual reactions. PCR products should be analyzed in an area separate from the reaction assembly area.

A standard 50-μL reaction uses 45 μL of Platinum® PCR SuperMix and 5 μL of primer and template solutions. For the primer sets used in the development of Platinum® PCR SuperMix, no decrease in product yield was observed if the amount of template and primer solution added is between a fraction of a microliter and 20 μL. Lower yield occurs as the Mg<sup>++</sup> concentration drops to a suboptimal level. If the final Mg<sup>++</sup> concentration is adjusted to 1.5 mM, the volume of primer and template solution that can be added to 45 μL of Platinum® PCR SuperMix can exceed 50 μL.

## Protocol

The following protocol is suggested as a starting point and guideline when using Platinum® PCR SuperMix. We recommend assembling reactions on ice from pre-chilled components. This protocol is for a reaction size of approximately 50 µL. The reaction size may be adjusted as desired.

**Note:** For multiple reactions with common components, prepare a master mix of the components common to all reactions to reduce pipetting errors.

1. Set up the reaction tubes/plates on ice.
2. Add the following components in any order to each reaction vessel.
  - 45-µL Platinum® PCR SuperMix
  - Primers (200 nM final concentration per primer is recommended)\*
  - Template DNA solution\*

\*Total volume of primer and template solution can be 0.5–20 µL.

3. Cap reaction vessels and load in thermal cycler at 94°C.
4. Incubate tubes in a thermal cycler at 94°C for 30 seconds–2 min to completely denature the template and activate the enzyme.
5. Perform 25–35 cycles of PCR amplification as follows:
  - Denature: 94°C for 15–30 seconds
  - Anneal: 55°C for 15–30 seconds
  - Extend: 72°C for 1 minute per kb

## References

- Chou, Q., Russel, M., Birch, D., Raymond, J., Bloch, W. (1992) *Nucl. Acids Res.*, 20, 1717.
- Sharkey, D.J., Scalice, E.R., Christy, K.G., Atwood, S.M., Daiss, J.L. (1994) *BioTechnology*, 12, 506.
- Westfall, B., Sitaraman, K., Solus, J., Hughes, J., Rashtchian, A. (1997) *Focus*®, 19.2, 46.

## Certificate of Analysis

The Certificate of Analysis (CofA) provides detailed quality control information for each product. The CofA is available on our website at [www.lifetechnologies.com/cofa](http://www.lifetechnologies.com/cofa), and it is searchable by product lot number, which is printed on each box.

## Limited Use Label License No. 358: Research Use Only

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com) or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

## Limited Product Warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at [www.lifetechnologies.com/termsandconditions](http://www.lifetechnologies.com/termsandconditions). If you have any questions, please contact Life Technologies at [www.lifetechnologies.com/support](http://www.lifetechnologies.com/support).

©2012 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

**For support visit** [www.lifetechnologies.com/support](http://www.lifetechnologies.com/support) or email [techsupport@lifetech.com](mailto:techsupport@lifetech.com)

[www.lifetechnologies.com](http://www.lifetechnologies.com)

*life*  
technologies™