

# Protein Expression Sample Preparation Kit

Insert PN 4405447 Rev. A  
Printed in USA



For Research Use Only.  
Not for use in diagnostic procedures.

## Product

Product	Part number
Protein Expression Sample Preparation Kit	4405443

For safety and biohazard guidelines, refer to the "Safety" section in the *Protein Expression Sample Preparation Kit Protocol*. Read the MSDS and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

## Reagent characteristics

The Protein Expression Sample Preparation Kit contains the reagents required for the preparation of target and control cell lysates for use with TaqMan® Protein Expression Assay experiments. The **Cell Lysis Reagent** is a buffered detergent solution, which is combined with user-supplied protease and phosphatase inhibitors, and used to lyse adherent or nonadherent cultured cells under nondenaturing conditions. The **Cell Resuspension Buffer** is used to resuspend the cells before adding the lysis reagent and to dilute the lysates in preparation for the assay-binding step.

For more information on the Sample Preparation Kit or the TaqMan® Protein Expression Assays chemistry, see the *Protein Expression Sample Preparation Kit Protocol*.

## Components and storage conditions

Part	Part number	Quantity	Storage
<b>Cell Lysis Reagent, 2X</b>	4405374	3 × 1 mL	2 to 8 °C
<b>Cell Resuspension Buffer</b>	4405377	1 × 25 mL	

When the components are stored at the recommended temperatures, this product maintains performance at least until the expiration date on the tube label.

## Preparation for use

Before use, gently swirl or invert the **Cell Lysis Reagent** and **Cell Resuspension Buffer** to mix the reagents.

## Usage guidelines

- See the *Protein Expression Sample Preparation Kit Protocol* for instructions on performing the cell lysis and preparing the lysates for the assay binding step.
- Prepare the **Cell Lysis Buffer** with protease and phosphatase inhibitors immediately prior to use. Discard any unused solution after 2 hours.
- After each use, store the protease and phosphatase inhibitor stock solutions at –15 to–25 °C.

**For Research Use Only. Not for use in diagnostic procedures.**

NOTICE TO PURCHASER: PLEASE REFER TO THE PROTEIN EXPRESSION SAMPLE PREPARATION KIT PROTOCOL FOR LIMITED LABEL LICENSE OR DISCLAIMER INFORMATION.

© 2009, Life Technologies Corporation. All rights reserved. Trademarks of Life Technologies Corporation and its affiliated companies: AB® (Design), Applied Biosystems®. TaqMan is a registered trademark of Roche Molecular Systems, Inc.

04/2009

