



AcTEV™ Protease

Cat. No. 12575-015

Size: 1,000 units

Cat. No. 12575-023

Size: 10,000 units

Description

AcTEV™ Protease is an enhanced form of Tobacco Etch Virus (TEV) protease that is highly site-specific, active, and more stable than native TEV protease (1). AcTEV™ Protease recognizes the seven-amino-acid sequence Glu-Asn-Leu-Tyr-Phe-Gln-Gly (2-5) and cleaves between Gln and Gly with high specificity. The protease is used to cleave affinity tags from fusion proteins. The optimal temperature for cleavage is 30°C; however, the enzyme is active over wide ranges of temperature (see table on page 3) and pH (pH 6.0-8.5). Following digestion, AcTEV™ Protease is easily removed from the cleavage reaction by affinity chromatography using the polyhistidine tag at the N-terminus of the protease. AcTEV™ Protease is purified from *E. coli* by affinity chromatography using the polyhistidine tag.

Components

Item	Composition	Quantity	
		1000 U	10,000 U
AcTEV™ Protease (10U/μl)	AcTEV™ Protease in: 50 mM Tris-HCl, pH 7.5 1 mM EDTA 5 mM DTT 50% (v/v) glycerol 0.1% (w/v) Triton® X-100	100 μl	1.0 ml
20X TEV Buffer	1 M Tris-HCl, pH 8.0 10 mM EDTA	1.0 ml	8.0 ml
0.1 M DTT	0.1 M DTT in deionized water	500 μl	2.0 ml

Store AcTEV™ Protease at -20°C or at -80°C for long-term storage.

Store 20X TEV Buffer at 4°C or -20°C and store 0.1 M DTT at -20°C.

Part No. 12575.pps

Rev. Date: 07/25/03

This product is distributed for laboratory research only. CAUTION: Not for diagnostic use. The safety and efficacy of this product in diagnostic or other clinical uses has not been established.

For technical questions about this product, call the Invitrogen Tech-Line™ U.S.A. 800 955 6288

Unit Definition

One unit of AcTEV™ cleaves ≥85% of 3 µg control substrate in 1 h at 30°C.

Unit Assay Conditions

The AcTEV™ Protease assay is performed in 1X TEV Buffer (50 mM Tris-HCl, pH 8.0, 0.5 mM EDTA) and 1 mM DTT with 1 unit enzyme and 3 µg control substrate (6) at 30°C for 1 hour in a total volume of 30 µl.

Quality Control

In a quality control assay, this product must demonstrate functional absence of any non-specific protease activity.

Recommended Conditions for Cleavage of a Fusion Protein

Example of a time course experiment with 10 units AcTEV™ Protease is shown below. If the protein of interest is heat-labile, incubate at 4°C with longer incubation times and/or more enzyme (see table on next page).

1. Add the following to a microcentrifuge tube:

Fusion Protein	20 µg
20X TEV Buffer	7.5 µl
0.1 M DTT	1.5 µl
AcTEV™ Protease, (10 units)	1.0 µl
Water	to 150 µl

2. Incubate at 30°C. Remove 30 µl aliquots at 1, 2, 4, and 6 hours.

3. Add 30 µl 2X SDS sample buffer (125 mM Tris-HCl, pH 6.8; 4% SDS; 1.4 M β-mercaptoethanol; 20% (v/v) glycerol; 0.01% bromophenol blue). Keep samples at -20°C until experiment is complete.

4. Analyze 40 µl of sample by SDS-PAGE using a suitable gel.

The percent protein cleavage is determined by analyzing the amount of cleaved products formed and amount of uncleaved protein remaining after digestion.

After evaluating the initial results, you may optimize the cleavage reaction for your specific protein by optimizing the amount of AcTEV™ Protease, incubation temperature, or reaction time.

Varying Parameters for Cleavage

The percent of 3 µg control substrate hydrolyzed by one unit of AcTEV™ Protease at various temperatures was examined (see table below). More cleaved protein is formed with AcTEV™ by increasing the incubation time. If time is critical, add more AcTEV™ Protease to increase hydrolysis.

Percentage Substrate Hydrolyzed				
Time	4°C	16°C	21°C	30°C
0.5 h	36	62	72	85
1 h	58	85	99	99
2 h	77	99	99	99
3 h	88	99	99	99

Cloning of TEV Cleavage Site into Vectors

You may clone your gene of interest into expression vectors available from Invitrogen that contain an affinity tag and a TEV cleavage site such as Champion™ pET151/D-TOPO (cat. no. K151-01) for *E. coli* expression and pFastBac™ HT (cat. no. 10584-027) for baculovirus expression.

To introduce the TEV cleavage site (Glu-Asn-Leu-Tyr-Phe-Gln-Gly) into an expression vector of choice, design a primer with nucleotides encoding the recognition site amino acids. You may use the free OligoPerfect™ Designer software available on our web site at www.invitrogen.com for primer design. The software conveniently adds a TEV cleavage site into your primer with the click of a button.

Removal of AcTEV™ Protease after Cleavage

The AcTEV™ Protease contains a polyhistidine tag at the N-terminus. After cleavage of the fusion protein, remove AcTEV™ Protease from the cleavage reaction by affinity chromatography on a nickel chelating resin such as ProBond™ Resin (cat. no. K801-01). Dilute cleavage reaction in the binding buffer for ProBond™ and perform binding and elution as described in the ProBond™ Purification manual available at www.invitrogen.com. **The cleaved native protein will be in the flow-through fractions.**

References

1. Nayak, S. *et al.* (2003) *Focus*[®] in press.
2. Dougherty, W.G., *et al.* (1988) *EMBO* 7, 1281.
3. Carrington, J.C. and Dougherty, W.G. (1988) *Proc. Natl. Acad. Sci. USA* 85, 3391.
4. Dougherty, W.G., *et al.* (1989) *Virology* 172, 302.
5. Dougherty, W.G., and Parks, T.D. (1989) *Virology* 172, 145.
6. Van Hoy, M., *et al.* (1993) *Cell* 72, 587.

Limited Use Label Licenses 135 and 163:

135: This product is licensed under U.S. Patent No. 5,532,142. **163:** The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patent(s) based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Invitrogen Corporation, 1600 Faraday Avenue, Carlsbad, California 92008. Phone (760) 603-7200. Fax (760) 602-6500.

©2003 Invitrogen Corporation. All rights reserved.

Triton[®] is a registered trademark of Rohm & Haas Co.