

Lipofectamine® 2000 CD Reagent

Part no. 12566014.pps

MAN0001090

Rev. Date 14 July 2011

Cat. no. 12566-014

Size 1.0 mL

Store at 4°C (do not freeze)

Description

Lipofectamine® 2000 CD is a proprietary *animal origin-free* formulation for transfecting nucleic acids into eukaryotic cells. A Drug Master File (DMF) has been submitted to the FDA. Contact Life Technologies Technical Service or your local Sales Representative for permission to cross-reference the DMF. Using Lipofectamine® 2000 CD for transfection provides the following advantages:

- Highest transfection efficiency in many cell types and formats. Refer to the Cell Lines database at www.invitrogen.com for a list of cell types transfected.
- The animal origin-free formulation ensures that mammalian cell culture and bioproduction processes are free of animal-derived materials.
- DNA-Lipofectamine® 2000 CD complexes can be added directly to cells in culture medium.
- It is not necessary to remove complexes or change/add medium after transfection, but complexes may be removed after 4–6 hours.

Important Guidelines for Transfection

- Culture cells in serum-free media that are free of animal-derived components. Test serum-free media for compatibility with Lipofectamine® 2000 CD since some serum-free formulations (e.g. CD 293, 293 SFM II, CD Hybridoma) may inhibit cationic lipid-mediated transfection.
- For consistent animal origin-free transfection, use OptiPro™ SFM (Cat. no. 12309-019) to dilute DNA and Lipofectamine® 2000 CD before complexing.
- Transfect cells at high cell density: For adherent cells, we recommend 90–95% confluence at the time of transfection for high efficiency, high expression levels, and to minimize cytotoxicity. For suspension cultures, transfect cells at a density of $0.8\text{--}1.1 \times 10^6$ cells/mL. Optimization may be necessary. Maintain the same seeding conditions between experiments.
- *Do not* add antibiotics to media during transfection since this causes cell death.
- *Do not* add Pluronic® or charged media extracts (e.g. dextran sulfate) to media during transfection, as these reagents can inhibit transfection.

Intended Use: For research use only.

Not intended for any animal or human therapeutic or diagnostic use.

Transfecting Adherent Cells

Use the following procedure to transfect mammalian cells in a *24-well format*. For other formats, see **Scaling Up or Down Transfections** (page 3). Use the recommended Lipofectamine® 2000 CD amount as a starting point, and then optimize conditions for your cell line and serum-free medium, as needed. All amounts and volumes are given on a per-well basis.

1. One day before transfection, plate $0.5\text{--}2 \times 10^5$ cells in 500 μL of growth medium without antibiotics so that they will be 90–95% confluent at the time of transfection.
2. *For each transfection sample*, prepare complexes as follows:
 - a. Dilute DNA in 50 μL of OptiPro™ SFM. Mix gently.
 - b. Mix Lipofectamine® 2000 CD gently before use, then dilute the appropriate amount in 50 μL of OptiPro™ SFM. Incubate for 5 minutes at room temperature.

Note: Combine diluted Lipofectamine® 2000 CD with diluted DNA within 30 minutes.
 - c. After 5 minute incubation, combine the diluted DNA with diluted Lipofectamine® 2000 CD (total volume = 100 μL). Mix gently and incubate for 20 minutes at room temperature (solution may appear cloudy).

Note: Complexes are stable for 6 hours at room temperature.
3. Add the 100 μL of complexes to a well containing cells and medium. Mix gently by rocking the plate back and forth.
4. Incubate cells at 37°C in a CO₂ incubator for 18–48 hours prior to testing for transgene expression. It is not necessary to change the medium, but medium may be replaced after 4–6 hours.
5. *For stable cell lines:* Passage cells at a 1:10 (or higher dilution) into fresh growth medium 24 hours after transfection. Add selective medium the following day.

Scaling Up or Down Transfections

To transfect cells in different tissue culture formats, vary the amounts of Lipofectamine® 2000 CD, DNA, cells, and medium used in proportion to the relative surface area, as shown in the following table. With automated, high-throughput systems, we recommend a complexing volume of 50 µL for transfections in 96-well plates.

Culture vessel	Surf. area per well (cm ²)	Relative surf. area vs. 24-well	Vol. of plating medium	DNA (µg) in media vol. (µL)	Lipofectamine® 2000 CD (µL) in media vol. (µL)
96-well	0.3	0.2	100 µL	0.2 µg in 25 µL	0.5 µL in 25 µL
24-well	2	1	500 µL	0.8 µg in 50 µL	2.0 µL in 50 µL
12-well	4	2	1 mL	1.6 µg in 100 µL	4.0 µL in 100 µL
35-mm	10	5	2 mL	4.0 µg in 250 µL	10 µL in 250 µL
6-well	10	5	2 mL	4.0 µg in 250 µL	10 µL in 250 µL
60-mm	20	10	5 mL	8.0 µg in 0.5 mL	20 µL in 0.5 mL
10-cm	60	30	15 mL	24 µg in 1.5 mL	60 µL in 1.5 mL

Note: Surface areas are determined from actual measurements of tissue culture vessels, and may vary depending on the manufacturer.

Transfecting Cells in Suspension

Use the following procedure to transfect suspension cells at any scale. Before beginning, make sure that cells are healthy and > 90% viable.

1. On the day of transfection, prepare a single-cell suspension from stock cells at less than 3×10^6 cells/mL. Seed cells at a density of $0.8\text{--}1.1 \times 10^6$ cells/mL in growth media without antibiotics.
2. *For each transfection sample*, prepare complexes as described in step 2 of **Transfecting Adherent Cells** on page 2, using the following reagent amounts and volumes *for every mL of cells transfected*:
 - Dilute 0.5–1.5 µg DNA in 34 µL of OptiPro™ SFM
 - Dilute 1–10 µL of Lipofectamine® 2000 CD in 34 µL of OptiPro™ SFM
3. Add the complexes to the flask containing cells and media.
4. Incubate the cells at 37°C in a CO₂ incubator on an orbital shaker rotating at 125 rpm for 24–96 hours prior to testing for transgene expression.

Optimizing Transfection

To obtain the highest transfection efficiency and low non-specific effects, optimize transfection conditions by varying cell density and Lipofectamine® 2000 CD concentrations.

Certificate of Analysis

The Certificate of Analysis provides detailed quality control information for each product. Certificates of Analysis are available on our website. Go to www.lifetechnologies.com/support and search for the Certificate of Analysis by product lot number, which is printed on the box.

Limited Use Label License

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

©2011 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

Pluronic is a registered trademark of BASF Corporation.

For support visit www.lifetechnologies.com/support
or email techsupport@lifetech.com

www.lifetechnologies.com

