

**BacMam Histone H3 [pSer10] Cellular Assay**

Catalog Number: A12898

Literature Part Number: A12898PIS (MAN0003207)

Literature Lot Number: V1

Revision date: 13 August 2010

**FAST FACTS**


For first-time BacMam users, we recommend using cells like U-20S (ATCC, HTB-96) which can be transduced exceptionally well.

BacMam Enhancer Solution (PV5835) is **not** required for this assay.

**Optimal Virus %:** We recommend that you perform a titration of the BacMam Histone H3 Reagent to determine the optimal percentage of virus for your transduction in your cell background of interest or when you receive a new lot of virus. Select the lowest percentage of BacMam Reagent that yields the largest assay window (response ratio). See example below.

| Kit Components                                     | Amount    | Storage | Handling   |
|--|-----------|---------|--|
| BacMam Histone H3 Reagent                          | 2 × 25 mL | 4°C     | <ul style="list-style-type: none"> <li>Do not freeze</li> <li>Minimize exposure to ambient light</li> <li>Use sterile technique</li> <li>Aliquot to minimize handling, if necessary</li> </ul> |
| LanthaScreen® Tb-anti-Histone H3 [pSer10] Antibody | 10 µg     | -20°C   | Aliquot if necessary to avoid multiple freeze/thaw cycles  |
| 6X LanthaScreen® Cellular Assay Lysis Buffer       | 6 mL      | 4°C     | Supplement with inhibitors and antibody  |

**Titration of BacMam Histone H3 Reagent in U-20S cells** (Detection of Histone H3 phosphorylation at Ser19 induced by calyculin A)



| Additional Materials Required, but not provided   | Source   | Part #  |
|---|--|---------|
| Positive Control Agonist<br><b>We recommend:</b> Calyculin A (prepare stock solution at 1,000 µM in DMSO) | EMD Biosciences  | 208851  |
| Cell Line of Interest   | Various  | Various |
| DMSO  | Fluka  | 41647   |
| Protease Inhibitor  | Sigma  | P8340   |
| Phosphatase Inhibitor   | Sigma  | P0044   |
| White tissue culture-treated, 384-well assay plates   | Corning  | 3570    |
| Fluorescence plate reader with top-read and TR-FRET capability  | <a href="http://www.invitrogen.com/instrumentsetup">www.invitrogen.com/instrumentsetup</a> for details |         |
| <b>Optional:</b> Clear-bottom, tissue-culture treated, 384-well plates                                    | Corning  | 3712    |

**Detailed Protocols and Additional Assay Performance Data Available**

Please visit [www.invitrogen.com](http://www.invitrogen.com) and search for A12898 to download the full protocol and validation packet for this assay. Protocols and validation packets are located under the "How to Use" tab on the product page. Validation Packets include assay performance under variable experimental conditions.

**Technical Support**

For additional assistance in running this BacMam-enabled Cellular Assay, please contact our technical support team at [drugdiscoverytech@lifetech.com](mailto:drugdiscoverytech@lifetech.com) or 760-603-7200, extension 40266.

## Quick Reference Protocol for U-2 OS Cells—BacMam Transduction and Agonist Assays

This protocol is designed for experienced users using U-2 OS cells. Conditions may need to be optimized for different cell types. See the online protocol for more detailed information.

|  |   | Cell-Free Control Wells  | Unstimulated Control Wells   | Stimulated Control Wells  | Test Compound Wells                               |
|--|---|--|--|---|---|
| BacMam Transduction                          | <b>Step 1</b><br>Grow, Harvest Cells and Transduce  | <ul style="list-style-type: none"> <li>Grow cells in Growth Medium* to 70–90% confluency (~0.75 to 1.0×10<sup>5</sup> cells/cm<sup>2</sup>).</li> <li>Harvest, wash cells and resuspend in Assay Medium** at 3.5 × 10<sup>5</sup> cells/mL.</li> <li>Perform five 3-fold serial dilutions of BacMam reagent in Assay Medium.</li> <li>Add 0.4 mL undiluted or serial diluted BacMam reagent to 1 mL cells to generate a virus titration range of 0.1% to 30% (v/v) final concentration.</li> </ul> |  |   |   |
|  | <b>Step 2</b><br>Plate Cells/virus mixture  | Add 20 µL/well Assay Media only  | Add 20 µL cells and BacMam mixture per well (~5000 cells/well), and pulse spin the plate |   |   |
|  | <b>Step 3</b><br>Incubate Cells   | Incubate the plate at 37°C/5% CO <sub>2</sub> for 20–24 hours  |  |   |   |
| LanthaScreen® Assay                          | <b>Step 4</b><br>Add media, control agonist (calyculin A), or test compounds  | Add 10 µL/well of 0.3% DMSO in Assay Media   |  | Add 10 µL/well of 3X Agonist (for calyculin, 0.3 µM) in Assay Media | Add 10 µL/well of 3X Test Compound in Assay Media |
|  | <b>Step 5</b><br>Stimulate Cells  | Incubate the plate at 37°C/5% CO <sub>2</sub> for 1 hour   |  |   |   |
|  | <b>Step 6</b><br>Prepare Complete 6X Lysis Buffer   | To 1 mL of 6X Lysis Buffer, add 30 µL 100X protease inhibitor, 30 µL 100X phosphatase inhibitor, and Tb-anti-Histone H3 [pSer10] Antibody to 12 nM   |  |   |   |
|  | <b>Step 7</b><br>Add Lysis Buffer (including Tb-Ab)   | Add 6 µL/ well of Complete 6X Lysis Buffer to each well  |  |   |   |
|  | <b>Step 8</b><br>Cell Lysis/ Assay Equilibration  | Incubate plate for ~3 hours at room temperature in the dark  |  |   |   |
| <b>Step 9</b><br>Read Plate and Analyze Data | See <b>LanthaScreen® Detection</b> in the online protocol—Excitation filter: 337 nm (30 nm bandwidth); Emission filters: 490 nm (10 nm bandwidth) and 520 nm (25 nm bandwidth)*** |  |  |   |   |

\* **Growth Media:** McCoy's 5A Media (Invitrogen 16600) with 10% dFBS, 10 mM HEPES, 0.1 mM NEAA, 1 mM Sodium Pyruvate, and 100 U/mL Penicillin/ 100 µg/mL Streptomycin

\*\* **Assay Media:** Opti-MEM® I (Invitrogen 11058) with 0.5% cdFBS (or dialyzed FBS), 0.1 mM NEAA, 1 mM Sodium Pyruvate, and 100 U/mL Penicillin/ 100 µg/mL Streptomycin

\*\*\* **Instrument Specific Set-up Guides:** describe instrument specific filters, dichroic mirrors and reading parameters, available at [www.invitrogen.com/instrumentsetup](http://www.invitrogen.com/instrumentsetup).

**We would like your feedback**

Please let us know how you like this Product Information Sheet.

Visit [www.invitrogen.com/feedbacksurvey](http://www.invitrogen.com/feedbacksurvey) to take a quick 5-minute survey about your experience and get a free gift upon completion.\*

**All customers who complete this survey will receive a free gift.**

\*Available while supplies last. Offer void where prohibited by federal, state, or local laws or regulation or agency/institutional policy.

This product is subject to one or more limited use label licenses. Please refer to [www.invitrogen.com](http://www.invitrogen.com) for the corresponding limited use label licenses.

The performance of this product is guaranteed for six months from the date of purchase if stored and handled properly.

The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

©2010 Life Technologies. All rights reserved.