

Version No. 1

Representative Data for BacMam Histone H3 [pSer10] Cellular Assays

BacMam Histone H3 [pSer10] Cellular Assays

Overview

Introduction to BacMam. LanthaScreen® Cellular Assays can now be used together with a gene delivery method known as BacMam for the transient expression of the desired GFP-substrate fusion in many cell lines of interest. The BacMam delivery system utilizes modified baculovirus as a vehicle to deliver and express target genes in mammalian cells. This technology has several advantages over commonly used transient methods for heterologous gene expression, including: 1) high transduction efficiency across a broad range of cell types, including primary and stem cells; 2) little-to-no observable cytopathic effects; 3) reproducible and titratable target gene expression; 4) virus is non-replicating in mammalian cells, therefore it is designated as Biosafety Level 1 (BSL1); and 5) compatible with simultaneous delivery of multiple genes. Thus, combination of the LanthaScreen® Cellular Assay technology and BacMam-mediated gene delivery provides a fast, convenient, and robust method for interrogating specific signal transduction events in a cell background of choice. Please see Kost, et al for more information pertaining to BacMam gene expression in cells.



Schematic Illustration of the workflow of BacMam-enabled Cellular Assay Using LanthaScreen[®] Technology. Cells are treated with BacMam virus encoding a GFP-fusion protein and plated in a 384-well assay plate. 24-48 hours post-transduction, the cells are stimulated to induce the specific post-translational modifications (such as phosphorylation as shown) of the GFP-substrate fusion. Cells are then lysed in the presence of terbium-labeled anti-modification specific antibody and TR-FRET from terbium to GFP can be detected.

BacMam-enabled Cellular Histone H3 Assay Using LanthaScreen[®] **Technology.** The eukaryotic nucleosome, composed of histones H2A, H2B, H3, and H4, regulates the structure of chromatin and consequently modulates gene transcription profiles in a concerted manner. Nucleosome function is directly regulated by a multitude of posttranslational modifications on amino-terminal tails of core histones, including acetylation, phosphorylation, methylation and ubiquitination. The combination of baculovirus-mediated gene delivery (BacMam) with LanthaScreen® Cellular Assay technology enables a platform for the analysis of specific posttranslational modifications of histones. BacMam provides a convenient genetic delivery tool for a GFP-H3 fusion protein in the cell line of interest. This kit describes an HTS-compatible cellular immunoassay measuring phosphorylation of GFP-H3 to set the transmitted of the set of the transmitted of GFP-H3 to set the transmitted of the transmitted of



Version No. 1

Representative Data for BacMam Histone H3 [pSer10] Cellular Assays

Validation Table of Contents

Validation experiments of BacMam Histone H3 [pSer10] Cellular Assay were performed in a 384-well format using U-2 OS as a host cell line. Following treatment with Calyculin A, the phosphorylation status of GFP-Histone H3 was then measured using LanthaScreen® Tb-anti-H3 [pSer10] antibody as a detection reagent. Included are representative results using this assay, including an example of the assay optimization process, which can be applied to your cell line of choice.

- 1. Representative assay optimization and validation for Histone H3 pSer10 assay
 - a. Comparison of emission ratio and response ratio graphs
 - b. Optimization of BacMam Histone H3 transduction
 - c. Optimization of cell number per well
 - d. Optimization of stimulation time
 - e. Testing effect of carrier flask confluence
 - f. Lysis buffer equilibration time
 - g. Aurora Kinase Inhibitor Testing
 - h. Cell-cycle dependent Histone H3 Ser10 phosphorylation
 - i. Localization of GFP-Histone H3
- 2. Assay Portability
 - a. Alternate cell backgrounds tested: HeLa, NIH-3T3



Version No. 1

Page 3 of 6

Representative Data for BacMam Histone H3 [pSer10] Cellular Assays

Figure 1a — Comparison of emission ratio and response ratio graphs



[Calyculin A] (nM)

Measurement of Histone H3 Phosphorylation using the BacMam Histone H3 [pSer10] Cellular Assay. U-2 OS cells were harvested, resuspended in assay medium (Opti-MEM I + 0.1% charcoal-stripped FBS) and mixed with BacMam Histone H3 Reagent to achieve 10% virus (v/v) concentration. The mix was applied to wells of a 384-well plate, yielding a cell density of ~5,000 cells per well. Following a 24 hour incubation, cells were treated for 1 hour with a serial dilution of Calyculin A. Cells were lysed by adding fully-supplemented 6X lysis buffer containing Tb-anti-Histone H3 [pSer10] Antibody and equilibrated for 3 hours at room temperature before TR-FRET was measured on a BMG Labtech Pherastar plate reader. Raw TR-FRET Emission Ratios (520nm/490nm) and Response Ratios (measure of assay window) normalized to untreated controls are shown.

Figure 1b — Optimization of BacMam Histone H3 transduction



Measurement of Histone H3 Phosphorylation using the BacMam Histone H3 [pSer10] Cellular Assay. U-2 OS cells were harvested, resuspended in assay medium and mixed with BacMam Histone H3 Reagent to achieve the indicated % virus (v/v) concentrations. The mixes were applied to designated wells of a 384well plate, yielding a cell density of ~5,000 cells per well. Following a 24 hour incubation, cells were treated for 1 hour with a serial dilution of Calyculin A. Cells were lysed by adding fully-supplemented 6X lysis buffer containing Tb-anti-Histone H3 [pSer10] Antibody and equilibrated for 3 hours at room temperature before TR-FRET was measured on a BMG Labtech Pherastar plate reader.

Figure 1c – Optimization of cell number



Measurement of Histone H3 Phosphorylation using the BacMam Histone H3 [pSer10] Cellular Assay. U-2 OS cells were harvested, resuspended in assay medium and mixed with BacMam Histone H3 Reagent to achieve 10% virus (v/v) concentration. The mix was applied to wells of a 384-well plate, yielding varying cell densities per well. Following a 24 hour incubation, cells were treated for 1 hour with a serial dilution of Calyculin A. Cells were lysed by adding fullysupplemented 6X lysis buffer containing Tb-anti-Histone H3 [pSer10] Antibody and equilibrated for 3 hours at room temperature before TR-FRET was measured on a BMG Labtech Pherastar plate reader.

Tel: 800 955 6288 x40266 • E-mail: drugdiscoverytech@invitrogen.com • Web: www.invitrogen.com 19July2010



Version No. 1

Representative Data for BacMam Histone H3 [pSer10] Cellular Assays





[Calyculin A] (nM)

Measurement of Histone H3 Phosphorylation using the BacMam Histone H3 [pSer10] Cellular Assay. U-2 OS cells were harvested, resuspended in assay medium and mixed with BacMam Histone H3 Reagent to achieve 10% virus (v/v) concentration. The mix was applied to wells of a 384-well plate, yielding a cell density of ~5,000 cells per well. Following a 24 hour incubation, cells were treated for the indicated times with a serial dilution of Calyculin A. Cells were lysed by adding fully-supplemented 6X lysis buffer containing Tb-anti-Histone H3 [pSer10] Antibody and equilibrated for 3 hours at room temperature before TR-FRET was measured on a BMG Labtech Pherastar plate reader.

Figure 1e — Testing effect of carrier flask confluence



[Calyculin A] (nM)

Measurement of Histone H3 Phosphorylation using the BacMam Histone H3 [pSer10] Cellular Assay. Two T225 flasks of U-2 OS cells at varying confluence were harvested, resuspended in assay medium and mixed with BacMam Histone H3 Reagent to achieve 10% virus (v/v) concentration. The mix was applied to wells of a 384-well plate, yielding a cell density of ~5,000 cells per well. Following a 24 hour incubation, cells were treated for 1 hour with a serial dilution of Calyculin A. Cells were lysed by adding fullysupplemented 6X lysis buffer containing Tb-anti-Histone H3 [pSer10] Antibody and equilibrated for 3 hours at room temperature before TR-FRET was measured on a BMG Labtech Pherastar plate reader.



BacMam Histone H3 Lysis Buffer Equilibration Time



[Calyculin A] (nM)

Measurement of Histone H3 Phosphorylation using the BacMam Histone H3 [pSer10] Cellular Assay. U-2 OS cells were harvested, resuspended in assay medium and mixed with BacMam Histone H3 Reagent to achieve 10% virus (v/v) concentration. The mix was applied to wells of a 384-well plate, yielding a cell density of ~5,000 cells per well. Following a 24 hour incubation, cells were treated for 1 hour with a serial dilution of Calyculin A. Cells were lysed by adding fully-supplemented 6X lysis buffer containing Tb-anti-Histone H3 [pSer10] Antibody and equilibrated for 3 hours at room temperature before TR-FRET was measured on a BMG Labtech Pherastar plate reader.





Aurora Kinase Inhibitor VX-680 Profile. U-2 OS cells were transduced and plated onto a 384-well assay plate as described in Figure 1. Following a 24 hour incubation, cells were treated for 1 hour with a serial dilution of VX-680 followed by one hour incubation with 80 nM of Calyculin A. Cells were lysed by adding fully-supplemented 6X lysis buffer containing Tb-anti-Histone H3 [pSer10] Antibody and equilibrated for 3 hours at room temperature before TR-FRET was measured on a BMG Labtech Pherastar plate reader.



Version No. 1

Representative Data for BacMam Histone H3 [pSer10] Cellular Assays





Cell-cycle arrest at M phase leading to increased Ser10 phosphorylation of Histone H3. U-2 OS cells were harvested, resuspended in assay medium and mixed with BacMam Histone H3 Reagent to achieve 10% virus (v/v) concentration. The mix was applied to wells of a 384-well plate, yielding a cell density of ~5,000 cells per well. Following a 24 hour incubation, cells were treated for 3 hours with serial dilutions of nocodazole. Cells were lysed by adding fully-supplemented 6X lysis buffer containing Tb-anti-Histone H3 [acLys9] antibody or [pSer10] antibody and equilibrated for 3 hours at room temperature before TR-FRET was measured on a BMG Labtech Pherastar plate reader.

Figure 1i — Localization of GFP-Histone H3



Colocalization of GFP-Histone H3 with DNA in the

nucleus. U-2 OS cells were transduced with 10% (v/v) BacMam GFP-Histone H3 reagent for 20 hours prior to Hoechst staining and imaging.

Figure 2a — Alternate Cell Backgrounds

Measuring Histone H3 Phosphorylation in HeLa Cells



[Calyculin A] (nM)

Measurement of Histone H3 Phosphorylation using the BacMam Histone H3 [pSer10] Cellular Assay. HeLa cells were harvested, resuspended in assay medium and mixed with BacMam Histone H3 Reagent to achieve 30% virus (v/v) concentration. The mix was then applied to wells of a 384-well plate, yielding a cell density of ~10,000 cells per well. Following a 24 hour incubation, cells were treated for 1 hour with a serial dilution of Calyculin A. Cells were lysed by adding fully-supplemented 6X lysis buffer containing Tb-anti-Histone H3 [pSer10] Antibody and equilibrated for 3 hours at room temperature before TR-FRET was measured on a BMG Labtech Pherastar plate reader.

Measuring Histone H3 Phosphorylation in NIH-3T3 Cells



Measurement of Histone H3 Phosphorylation using the BacMam Histone H3 [pSer10] Cellular Assay. NIH-3T3 cells were seeded and simultaneously transduced overnight in a T75 flask using 25% (v/v) BacMam Histone H3 Reagent in combination with 0.5 μ M Trichostatin A. After 24 hours transduction, cells were harvested, resuspended in assay medium, and then applied to wells of a 384-well plate at a cell density of ~10,000 cells per well. Following an additional 24 hour incubation, cells were treated for 1 hour with a serial dilution of Calyculin A. Cells were lysed by adding fully-supplemented 6X lysis buffer containing Tb-anti-Histone H3 [pSer10] Antibody and equilibrated for 3 hours at room temperature before TR-FRET was measured on a BMG Labtech Pherastar plate reader.



Version No. 1

Page 6 of 6

Representative Data for BacMam Histone H3 [pSer10] Cellular Assays

General References:

- Carlson, C. B.; Mashock, M.J.; Bi, K., BacMam-enabled LanthaScreen cellular assays for PI3K/Akt pathway compound profiling in disease-relevant cell backgrounds. *J. Biomol. Screen* 2010, *15*, 327-334.
- Robers, M. B.; Horton, R. A.; Bercher, M. R.; Vogel, K. W.; Machleidt, T., High-throughput cellular assays for regulated post-translational modifications. *Anal. Biochem.* 2008, 372, 189–197.
- Huwiler, K.G.; Machleidt, T.; Chase, L.; Hanson, B.; Robers, M. B., Characterization of serotonin 5hydroxytryptamine-1A receptor activation using a phosphorextracellular-signal regulated kinase 2 sensor. *Anal. Biochem.* 2009, 393, 95-104.
- Kost, T. A.; Condreay, J. P.; Ames, R. S.; Rees, S.; Romanos, M. A., Implementation of BacMam virus gene delivery technology in a drug discovery setting. *Drug Disc. Today* 2007, *12*, 396-403.