

## Human Epidermal Keratinocytes, neonatal (HEKn)

Cat. no. C-001-5C

Quantity: ≥500,000 viable cells/vial

### Product Description

HEKn are human epidermal keratinocytes isolated from neonatal foreskin. Each vial of this product contains ≥5 x 10<sup>5</sup> viable cells that have been cryopreserved at the end of the primary culture stage in a medium containing 10% DMSO. An independent laboratory tests the cells for the presence of mycoplasma, Hepatitis B, Hepatitis C, and HIV-1 viruses. These agents were not detected. In our laboratory, each lot of cells is performance tested by culturing the cells through multiple passages in EpiLife<sup>®</sup> Medium supplemented with Human Keratinocyte Growth Supplement (HKGS) in the absence of antibiotics and antimycotics. During this culture period, no contamination by bacteria, yeast, or other fungi was detected. Upon thawing, the cells are guaranteed to be ≥70% viable (trypan blue) and to have a potential of ≥30 population doublings when handled according to the directions provided in this document. For recommended precautions to be taken when handling human cells, please read the caution statement.

### Intended Use

Cryopreserved HEKn are intended for use by researchers investigating the molecular and biochemical bases of various normal and disease processes. **This product is for research use only. Not for use in animals, humans, or diagnostic procedures.**

### Storage and Stability

Cryopreserved HEKn should arrive frozen on dry ice. If the cells are not to be used immediately, the user should prepare a space for storage of the vial in the vapor phase of a liquid nitrogen freezer. While wearing protective eyewear, gloves, and a laboratory coat, remove the vial from its shipping container and place immediately in the liquid nitrogen freezer. Although the viability of cryopreserved cells decreases with time in storage, useful cultures can usually be established even after 2 years of storage at liquid nitrogen temperatures.

### Caution

Although cryopreserved cells have been tested for the presence of various hazardous agents, diagnostic tests are not necessarily 100% accurate. In addition, human cells may harbor other known or unknown agents or organisms which could be harmful to your health or cause fatal illness. The user should treat all human cells as potential pathogens. Wear protective clothing and eyewear. Practice appropriate disposal techniques for potentially pathogenic or biohazardous materials.

Keratinocyte Culture Systems from Cascade Biologics			
	Standard system	Extended-lifespan systems	
	Undefined	Undefined	Defined
Basal Medium	Medium 154	EpiLife <sup>®</sup> Medium	EpiLife <sup>®</sup> Medium
Growth Supplement	HKGS (S-001-5) or HKGS Kit (S-001-K)	HKGS (S-001-5) or HKGS Kit (S-001-K)	EDGS* (S-012-5)
Coating Matrix Kit	N/A	N/A	Coating Matrix Kit (R-011-K)
Subculture Reagent	Trypsin/EDTA (R-001-100)	Trypsin/EDTA (R-001-100)	Trypsin/EDTA (R-001-100)
Subculture Reagent	Trypsin Neutralizer (R-002-100)	Trypsin Neutralizer (R-002-100)	Defined Trypsin Inhibitor (R-007-100)
Expected lifespan from HEKn (C-001-5C)	16-25 population doublings	50-70 population doublings	40-60 population doublings
Expected lifespan from HEKa (C-005-5C)	16-20 population doublings	35-45 population doublings	30-40 population doublings

\*For optimal cell growth, EDGS should be used in conjunction with Coating Matrix Kit and Defined Trypsin Inhibitor.

### Initiating Cultures from Cryopreserved Cells

We recommend seeding cells recovered from cryopreservation at a density of 2.5 x 10<sup>3</sup> viable cells/cm<sup>2</sup>. For example, three 75 cm<sup>2</sup> or nine 25 cm<sup>2</sup> tissue culture flasks can usually be established from one vial containing ≥5 x 10<sup>5</sup> HEKn. The procedure given below is a sample protocol for establishing cultures from the contents of one vial.

- From the table above, determine which basal medium and growth supplement you wish to use and prepare a bottle of supplemented medium according to the instructions that accompany that product.
- Remove a vial of HEKn from liquid nitrogen storage, taking care to protect hands and eyes.
- Dip the lower half of the vial into a 37°C water bath to thaw.
- When the contents of the vial have thawed, wipe the outside of the vial with disinfecting solution and move to a Class II, type A laminar flow culture hood.
- Open the vial and pipette the suspension up and down with a 1 ml pipette to disperse the cells.
- Remove 20 µl from the vial and dilute the cell suspension in 20 µl of trypan blue solution (for example, Sigma Cat. no. T8154).
- Use a hemacytometer to determine the number of viable cells per ml.
- Dilute the contents of the vial (1 ml) to a concentration of 1.25 x 10<sup>4</sup> viable cells/ml using the supplemented medium from step 1, above.

**For research use only.**

Life Technologies Corporation • 5791 Van Allen Way • Carlsbad • CA 92008 • Tel: 800.955.6288 • [www.invitrogen.com](http://www.invitrogen.com) • E-mail: tech\_support@invitrogen.com

9. Add 5 ml of cell suspension to each 25 cm<sup>2</sup> culture flask or 15 ml of cell suspension to each 75 cm<sup>2</sup> culture flask.
10. Following inoculation, swirl the medium in the flasks to distribute the cells. HEK<sub>293</sub> attach to culture surfaces quickly, and if the medium is not distributed immediately following inoculation, the cells may grow in uneven patterns.
11. Incubate the cultures in a 37°C, 5% CO<sub>2</sub>/95% air, humidified cell culture incubator. For best results, do not disturb the culture for at least 24 hours after the culture has been initiated.

## Maintenance of Stock Cultures

1. Change the culture medium to freshly supplemented medium, 24 to 36 hours after establishing a secondary culture from cryopreserved cells. For subsequent subcultures, change the medium 48 hours after establishing the subculture.
2. Change the medium every other day thereafter, until the culture is approximately 50% confluent.
3. Once the culture reaches 50% confluence, change the medium every day until the culture is approximately 80% confluent.

### Note

- To achieve the highest cell densities, the culture medium should be changed every day as the cultures approach confluence. To obtain rapidly proliferating subcultures, HEK<sub>293</sub> should be subcultured before they become more than 80% confluent. If HEK<sub>293</sub> reach confluency, the cells mitotically arrest and some of the cells leave the proliferating pool. Allowing HEK<sub>293</sub> cultures to arrest will decrease the long-term potential yield from a cryopreserved vial. The number of subcultures (passages) that can be achieved will vary with the starting cell density and the methods employed by individual investigators.
- HEK<sub>293</sub> cultures seeded at 2.5 x 10<sup>3</sup> cells/cm<sup>2</sup> from cryopreserved cells should reach 80% confluence in 5-7 days. In this culture, most of the cells should have an epithelioid morphology and be associated with each other in colonies. Some irregularly sized and shaped cells may be observed. Occasionally, small numbers of melanocytes persist in the secondary culture. Melanocytes do not readily proliferate in medium supplemented with HKGS or EDGS, and should be virtually absent in subsequent cultures.

## Subculture of HEK<sub>293</sub>

View the culture under the microscope to confirm that it is subconfluent, and that there are mitotic cells present. This protocol is designed for the subculture of one 25 cm<sup>2</sup> culture flask. If different-sized culture vessels are to be used, reagent volumes should be adjusted accordingly.

1. Assemble the appropriate supplemented medium and subculture reagents according to the table on page 1.

**Note:** We do NOT recommend warming the reagents prior to use.

2. Assemble the appropriate culture vessels, sterile pipettes, and sterile 15 ml conical tubes (not provided).
3. Remove all of the culture medium from the flask.
4. Add 3 ml of Trypsin/EDTA solution to the flask. Rock the flask to ensure that the entire surface is covered.
5. Immediately remove all 3 ml of Trypsin/EDTA solution from the flask.
6. Add 1 ml of fresh Trypsin/EDTA solution to the flask.
7. View the culture under a microscope. Incubate the flask at room temperature until the cells have become completely round, approximately 8-10 minutes.
8. Rap the flask very gently to dislodge cells from the surface of the flask.
9. Add 3 ml of Trypsin Neutralizer solution or Defined Trypsin Inhibitor solution to the flask and transfer the detached cells to a sterile 15 ml conical tube.
10. Add 3 ml additional Trypsin Neutralizer solution or Defined Trypsin Inhibitor solution to the flask and pipette the solution over the flask surface several times to remove any remaining cells. Add this solution to the 15 ml conical tube.
11. Centrifuge the cells at 180 x g for 7 minutes. Observe the cell pellet.
12. Remove the supernatant from the tube, being careful not to dislodge the cell pellet.
13. Resuspend the cell pellet in 4 ml supplemented medium. Pipette the cells up and down with a 10 ml pipette to ensure a homogeneous cell suspension.
14. Determine the concentration of cells in the suspension.
15. Dilute the cells in supplemented medium and seed new culture vessels with 2.5 x 10<sup>3</sup> cells/cm<sup>2</sup>.
16. Incubate the cultures in a 37°C, 5% CO<sub>2</sub>/95% air, humidified cell culture incubator.

### Notes

- Damage to cultured HEK<sub>293</sub> can occur during trypsinization. This damage may result from exposure of the cells to the Trypsin/EDTA solution for excessive lengths of time, trypsinization at temperatures exceeding room temperature and/or excessive mechanical agitation. Check to make sure that the temperature of trypsinization is appropriate and, if necessary, alter the incubation time of the procedure.
- Another common source of damage is centrifugation at excessive g forces. Check to make sure that the speed of the centrifuge is appropriate. One manifestation of cellular damage that may be evident after centrifugation is strings of cells (and debris) that do not pellet in the bottom of the tube. This is due to the presence of DNA from lysed cells in the solution. If this condition exists, the cell pellet may be lost upon aspiration of the supernatant containing the DNA strings. In many cases, viable cells can be rescued by pipetting the cells (and DNA) up and down in a 10 ml pipette to shear the DNA, and centrifuging the suspension again to recover the cells.

*For research use only.*

Life Technologies Corporation • 5791 Van Allen Way • Carlsbad • CA 92008 • Tel: 800.955.6288 • [www.invitrogen.com](http://www.invitrogen.com) • E-mail: tech\_support@invitrogen.com

**Limited Use Label License No. 5: Invitrogen Technology**

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. For products that are subject to multiple limited use label licenses, the terms of the most restrictive limited use label license shall control. Life Technologies Corporation will not assert a claim against the buyer of infringement of patents owned or controlled by Life Technologies Corporation which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Life Technologies is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, California 92008. Phone (760) 603-7200. Fax (760) 602-6500. Email: [outlicensing@invitrogen.com](mailto:outlicensing@invitrogen.com).

©2009 Life Technologies Corporation. All rights reserved.

For research use only. Not intended for any animal or human therapeutic or diagnostic use.

EpiLife is a Registered Trademark of Cascade Biologics, Inc.

***For research use only.***

Life Technologies Corporation • 5791 Van Allen Way • Carlsbad • CA 92008 • Tel: 800.955.6288 • [www.invitrogen.com](http://www.invitrogen.com) • E-mail: [tech\\_support@invitrogen.com](mailto:tech_support@invitrogen.com)