

TABLE OF CONTENTS

1.0	Materials Supplied	1
2.0	Materials Required, but Not Supplied	2
2.1	Optional Equipment and Materials	2
3.0	Cell Culture Conditions	3
3.1	Media Required	3
3.2	Methods Used	3
4.0	Assay Procedure	4
4.1	Controls	4
4.2	Cell Stimulation and Assay Plate Setup	5
4.3	Substrate Loading	6
4.4	Incubation	6
4.5	Detection	6
5.0	Data Analysis	8
5.1	Background Subtraction	8
5.2	Response Ratio Calculations	8
5.3	Visual Observation of Intracellular LiveBLAzer™ FRET B/G Substrate (CCF4-AM)	9
6.0	Purchaser Notification	10

1.0 MATERIALS SUPPLIED

Description:	This cell line utilizes GeneBLAzer® beta-lactamase technology. CellSensor® NFκB- <i>bla</i> Jurkat cells contain a beta-lactamase reporter gene under control of Nuclear Factor kappa B (NFκB) response element stably integrated into Jurkat cells. This cell line can be used to detect agonists/antagonists of the NFκB signaling pathway. NFκB- <i>bla</i> Jurkat T-cells have been shown to respond to Tumor Necrosis Factor alpha (TNFα).
Cat. no.:	K1167
Shipping Condition:	Dry ice
Storage Condition:	Liquid nitrogen. NOTE: Immediately upon receipt, cells must be stored in liquid N₂. Cells CANNOT be stored at -80°C, as they will quickly lose viability.
Quantity:	~5,000,000 cells (5 × 10 ⁶ cells/ml)
Cell Line Name:	NFκB- <i>bla</i> Jurkat
Application:	This cell line can be used to detect agonists/antagonists of the NFκB signaling pathway.
Growth Properties:	Suspension
Cell Phenotype:	Lymphoblast
Selection Marker:	Blasticidin (5 μg/ml)
Vector Used:	pLenti- <i>bsd</i> /NFκB- <i>bla</i> Vector
Mycoplasma Testing:	Negative
Biosafety Level:	1

2.0 MATERIALS REQUIRED, BUT NOT SUPPLIED

Media/Reagents	Recommended Source	Cat. no.
LiveBLAzer® Loading Kit containing: LiveBLAzer™ FRET B/G Substrate (CCF4-AM) substrate DMSO for Solution A Solution B Solution C	Invitrogen	K1095 K1030 Other sizes and kits are available
RPMI 1640	Invitrogen	22400-089
DMSO	Fluka	41647
Fetal bovine serum (FBS), dialyzed, tissue-culture grade (DO NOT SUBSTITUTE)	Invitrogen	26400-044
Non-essential amino acids (NEAA)	Invitrogen	11140-050
Penicillin/Streptomycin	Invitrogen	15140-122
Sodium pyruvate	Invitrogen	11360-070
Tumor Necrosis Factor alpha (TNFα)	BD Biosciences	354066
Blasticidin antibiotic	Invitrogen	R210-01

Consumables	Recommended Source	Cat. no.
Black-wall, clear-bottom, 384-well assay plates (with low fluorescence background)	Costar	3712
Compressed air	Various	—
Conical tubes, 15 ml, sterile	Various	—
Internally threaded cryogenic vials, 1.8 ml	Various	—
Tissue culture flasks	Various	—

Equipment	Recommended Source
Class II biological safety cabinet	Various
Dual wave-length bottom-reading fluorescence plate reader	Various
Filters (see Sections 4.5.1 and 5.3)	Chroma Technologies
Hemocytometer (or another cell counting method)	Various
Humidified 37°C /5% CO ₂ incubator	Various
Inverted microscope (with phase contrast capabilities)	Various
Liquid nitrogen tank	Various

2.1 Optional Equipment and Materials

- Epifluorescence or fluorescence-equipped microscope, equipped with appropriate filters
- Microplate centrifuge

3.0 CELL CULTURE CONDITIONS

3.1 Media Required

Component	Thaw Medium	Assay Medium	Growth Medium	Freezing Medium
RPMI 1640	90%	90%	90%	80%
Dialyzed FBS	10%	10%	10%	10%
NEAA	0.1 mM	0.1 mM	0.1 mM	0.1 mM
Sodium pyruvate	1 mM	1 mM	1 mM	1 mM
Penicillin	100 U/ml	100 U/ml	100 U/ml	--
Streptomycin	100 µg/ml	100 µg/ml	100 µg/ml	--
DMSO	--	--	--	10%
Blasticidin Antibiotic	--	--	5 µg/ml	--

Note: We prepare our media by adding the listed components directly to the medium bottle. Blasticidin can be added directly to the cell culture flask to reach 5 µg/mL. Similar methods are suitable.

Note: Unless otherwise stated, all media and solutions should be at least room temperature (37°C is best) before adding them to the cells.

3.2 Methods Used

Follow these methods exactly, as they have been validated specifically for optimal performance of this cell line. At first opportunity, create and store an aliquot of cells for back-up.

3.2.1 Thawing Method

- Place 14 ml of Thaw Medium (without Blasticidin) into a T75 flask.
- Place the flask in a 37°C/5% CO₂ incubator for 15 minutes to allow medium to equilibrate to the proper pH and temperature.
- Remove vial of cells to be thawed from liquid nitrogen and rapidly thaw by placing at 37°C in a water bath with gentle agitation for 1–2 minutes. Do not submerge vial in water.
- Decontaminate the vial by wiping with 70% ethanol before opening in a Class II biological safety cabinet.
- Transfer the vial contents dropwise into 10 ml of Thaw Medium in a sterile 15 ml conical tube.
- Centrifuge cells at 200 × g for 5 minutes.
- Aspirate supernatant and resuspend the cell pellet in 1 ml of fresh Thaw Medium (without Blasticidin).
- Transfer contents to the T75 tissue culture flask containing pre-equilibrated Thaw Medium and place flask in the 37°C/5% CO₂ incubator.
- At first passage**, switch to Growth Medium with Blasticidin.

3.2.2 Propagation Method

- Cells should be passaged or fed at least twice a week. Cells should be maintained between 2 × 10⁵ and 2 × 10⁶ cells/ml.
- To passage cells, transfer cell suspension to a conical tube.
- Spin down cells and resuspend in Growth Medium.

3.2.3 Freezing Method

- Harvest the cells as described in **Section 3.2.2**, Steps 2 and 3. Count the cells, then spin cells down and resuspend in 4°C Freezing Medium (Section 3.1) at 5 × 10⁶ cells/ml.
- Dispense 1.0 ml aliquots into cryogenic vials.
- Place in an insulated container for slow cooling and store overnight at –80°C.
- Transfer to liquid nitrogen the next day for storage.

4.0 ASSAY PROCEDURE

The following instructions outline the recommended procedure for determining activity of compounds as modulators (agonists or antagonists) of the NFκB signaling pathway using beta-lactamase as the readout.

4.1 Controls

Plate layouts and experimental outlines will vary; in screening mode, we recommend using at least five wells for each control: Unstimulated Control, Stimulated Control, and Cell-Free Control.

Note: Certain solvents may affect assay performance. The effect of a test compound solvent should be assessed prior to screening. The cell stimulation described below is carried out in the presence of 0.5% DMSO to simulate the effect that the test compound solvent may have on the assay. If other solvents and/or solvent concentrations are used, the following assay should be changed accordingly.

Data point	Agonist Assay	Antagonist Assay	Purpose of Data Point
Unstimulated Control Wells	1. 32 µl cells diluted in Assay Medium 2. 4 µl 5% DMSO stock solution 3. 4 µl Assay Medium	1. 32 µl cells diluted in Assay Medium 2. 4 µl 5% DMSO stock solution 3. 4 µl Assay Medium	Provides the Unstimulated Emission Ratio, which is a portion of the Response Ratio. Unstimulated Emission Ratio = (Unstimulated 460 nm/Unstimulated 530 nm). See Section 5.0 (Data Analysis).
Stimulated Control Wells	1. 32 µl cells diluted in Assay Medium 2. 4 µl 5% DMSO stock solution 3. 4 µl 10X stock solution of TNFα	1. 32 µl cells diluted in Assay Medium 2. 4 µl 5% DMSO stock solution 3. 4 µl 10X stock solution of TNFα	Provides the Stimulated Emission Ratio, which is a portion of the Response Ratio. Stimulated Emission Ratio = (Stimulated 460 nm/Stimulated 530 nm). See Section 5.0 (Data Analysis).
Cell-Free Control Wells	1. 36 µl Assay Medium 2. 4 µl 5% DMSO stock solution	1. 36 µl Assay Medium 2. 4 µl 5% DMSO stock solution	Provides the background blue and green emission values to be subtracted from both unstimulated and stimulated blue and green emission values to yield the net unstimulated and stimulated blue and green emission values.
Test Compound Wells	1. 32 µl cells diluted in Assay Medium 2. 4 µl 10X stock of test compounds in 5% DMSO 3. 4 µl Assay Medium	1. 32 µl cells diluted in Assay Medium 2. 4 µl 10X stock of test compounds in 5% DMSO 3. 4 µl 10X stock solution of TNFα	Provides the experimental data on whether or not a test compound is active as an agonist or antagonist in the assay.

4.2 Cell Stimulation and Assay Plate Setup

4.2.1 Precautions

1. Work on a dust-free, clean surface. Always handle the 384-well, black-wall, clear-bottom assay plate by the sides; do not touch the clear bottom of the assay plate.
2. If pipetting manually, it may be necessary to centrifuge the plate briefly at room temperature (for 1 min. at 14 x g) after additions to ensure all the assay components are on the bottom of the wells.

4.2.2 Plating Cell

1. Harvest cells from culture in Growth Medium as described in **Section 3.2.2**, Step 2 and resuspend in Assay Medium at a density of 6.25×10^5 cells/ml.
2. Add 36 μ l per well of the Assay Medium to the Cell-Free Control wells. Add 32 μ l (20,000 cells) per well of the cell suspension to the Test Compound wells, the Unstimulated Control wells, and Stimulated Control wells.

Note: Cells were plated on day of assay. Cells can be plated on previous day at half the concentration in Step 1.

4.2.3 Preparation of Stock Solution

1. Prepare a stock solution of 5% DMSO.
2. Prepare a 10X stock of test compounds in 5% DMSO.
3. Prepare a 10X stock of TNF α in Assay Medium. For an antagonist screen, use ~18.7 ng/ml final concentration or EC₈₀ (10X = 187 ng/ml). For an agonist screen, use ~150 ng/ml final concentration or max stim (10X = 1500 ng/ml). We recommend running a dose response curve to determine the EC₈₀ for your 10X TNF α solution.

4.2.4 Agonist Assay Plate Setup

Note: This subsection provides directions for performing an agonist assay. Directions for performing an antagonist assay can be found in **Section 4.2.5**.

1. Add 4 μ l of the stock solution of 5% DMSO to the Unstimulated Control wells, the Stimulated Control wells, and to the Cell-Free Control wells.
2. Add 4 μ l of Assay Medium to the Test Compound and Unstimulated Control wells.
3. Add 4 μ l of the 10X stock solution of TNF α to Stimulated Control wells.
4. Add 4 μ l of the 10X stock of test compounds to the Test Compound wells.
5. Incubate the agonist assay plate in a humidified 37°C/5% CO₂ incubator for 5 hours.

4.2.5 Antagonist Assay Plate Setup

Note: This subsection provides directions for performing an antagonist assay. Directions for performing an agonist assay are provided in **Section 4.2.4**.

1. Add 4 μ l of the 10X stock of test compounds to the Test Compound wells.
2. Add 4 μ l of the stock solution of 5% DMSO to the Stimulated Control wells, the Unstimulated Control wells, and the Cell-Free Control wells.
3. If desired, incubate the test compounds with the cells at 37°C/5% CO₂ before proceeding. Typically, a one-hour incubation is sufficient.
4. Add 4 μ l of the 10X stock solution of TNF α to the Test Compound wells and the Stimulated Control wells.
5. Add 4 μ l of Assay Medium to the Unstimulated Control wells.
6. Incubate the antagonist assay plate in a humidified 37°C/5% CO₂ incubator for 5 hours.

4.3 Substrate Loading

Preparation of 6X CCF2-AM or LiveBLAzer™ FRET B/G Substrate (CCF4-AM) Loading Solution and cell loading should be done in the absence of direct strong lighting. Turn off the light in the hood.

4.3.1 LiveBLAzer™ FRET B/G Substrate (CCF4-AM) Loading Protocol

This protocol is designed for loading cells with LiveBLAzer™ FRET B/G Substrate (CCF4-AM) using the following solutions. The following protocol is written using LiveBLAzer™ FRET B/G (CCF4-AM); however, CCF2-AM may be substituted.

Reagents	Storage and Handling
LiveBLAzer™ FRET B/G Substrate (CCF4-AM)	The product is supplied as a dried powder. The molecular weight of CCF2-AM is 1082 g/mol; the molecular weight of LiveBLAzer™ FRET B/G Substrate (CCF4-AM) is 1096 g/mol. Store the product at -20°C, desiccated and protected from light.
Dry DMSO	Used to dissolve LiveBLAzer™ FRET B/G Substrate (CCF4-AM) for preparation of Solution A.
Solution A	LiveBLAzer™ FRET B/G Substrate (CCF4-AM) stock solution; 1 mM LiveBLAzer™ FRET B/G Substrate in dry DMSO. Prepare a 1 mM LiveBLAzer™ FRET B/G Substrate stock solution in dry DMSO. Store the LiveBLAzer™ FRET B/G Substrate stock solution at -20°C, desiccated and protected from light. Before each use, let frozen stock solution warm to room temperature and remove desired amount of reagent. Immediately recap the vial after each use to reduce moisture uptake and return to -20°C storage. Stored under these conditions, Solution A is stable for approximately three months. Once thawed, Solution A may appear slightly yellow. This is normal.
Solution B	Store the reagent at room temperature (18–22°C) protected from direct light. Under cold lab conditions [colder than 18°C (65°F)], the solution may freeze or a white precipitate may form. In this case, warm and mix the solution (~35°C) until thawed and the precipitate dissolves. Mix thoroughly before use.
Solution C	Keep tightly closed and store in a cool, dry place. Store the reagent at room temperature (18–22°C), protected from direct light.

4.3.2 Preparation of 6X LiveBLAzer™ FRET B/G Substrate (CCF4-AM) Loading Solution

1. Add 6 µl of Solution A to 60 µl of Solution B and vortex.
2. Add 934 µl Solution C to the combined solutions from Step 1 with vortexing.

4.3.3 Cell Loading

1. Remove assay plate from incubator and allow it to equilibrate to room temperature prior to loading.
2. Add 8 µl of the 6X Loading Buffer to each well [C_i LiveBLAzer™ FRET B/G (CCF4-AM) = 1 µM].

4.4 Incubation

1. Cover the plate to protect it from light and evaporation.
2. Incubate at room temperature for 120 minutes.
Note: Handle the plate gently and do not touch the bottom.

4.5 Detection

All measurements are made at room temperature from the bottom of the wells, preferably in 384-well black-wall, clear-bottom assay plates with low fluorescence background. Before reading the plate, remove dust from the bottom with compressed air.

Note: Some plates/fluorescence plate readers experience edge effects, which may affect data. If edge effects are noticed, plate layout should be considered when setting up the assay.

4.5.1 Instrumentation, Filters, and Plates

- Almost all dual-wavelength, bottom-read fluorescence plate readers can be used to detect beta-lactamase using LiveBLAzer™ FRET B/G (CCF4-AM).
Note: If you are uncertain of whether your instrument can be used for this assay, please contact Technical Support for assistance.

- Filters for fluorescence plate reader:

For ratiometric readout using a fluorescence plate reader:

Excitation filter: 409/20 nm

Emission filter: 460/40 nm

Emission filter: 530/30 nm

4.5.2 Reading an Assay Plate

1. Set the fluorescence plate reader to bottom-read mode and for eight scans per cycle.
2. Allow the lamp in the fluorescence plate reader to warm up for at least 10 min. before making measurements.
3. Use the following filter selections:

	Scans 1–4	Scans 5–8
Purpose:	Measure fluorescence in the blue channel	Measure fluorescence in the green channel
Excitation filter:	409/20 nm	409/20 nm
Emission filter:	460/40 nm (using the gain determined during calibration)	530/30 nm (using the gain determined during calibration)

Typically, set the fluorescence plate reader to 5 reads/well.

5.0 DATA ANALYSIS

5.1 Background Subtraction

Background subtraction for both channels (460 nm and 530 nm) is essential to obtain meaningful data. This can be accomplished either automatically using software connected to the fluorescence plate reader, or manually after each assay plate has been read.

Use the assay plate layout to identify the location of the Cell-Free Control wells. These control wells are used for background subtraction.

5.1.1 Automatic background subtraction

Use the assay plate layout to designate appropriate Cell-Free Control wells and enable background subtraction for both sets of emission scans.

5.1.2 Manual background subtraction

1. Determine the average emission from the Cell-Free Control wells at both 460 nm (Average Em₄₆₀ Background) and 530 nm (Average Em₅₃₀ Background).
2. Calculate the Net Stimulated Em₄₆₀ and Em₅₃₀ values and the Net Unstimulated Em₄₆₀ and Em₅₃₀ values, as shown below:

Stimulated Cells

Net Stimulated Em₄₆₀ = Em₄₆₀ of stimulated sample – Average Em₄₆₀ Background

Net Stimulated Em₅₃₀ = Em₅₃₀ of stimulated sample – Average Em₅₃₀ Background

Unstimulated Cells

Net Unstimulated Em₄₆₀ = Em₄₆₀ of unstimulated sample – Average Em₄₆₀ Background

Net Unstimulated Em₅₃₀ = Em₅₃₀ of unstimulated sample – Average Em₅₃₀ Background

5.2 Response Ratio Calculations

1. Calculate the Stimulated Emission Ratio for each well as shown below.

$$\text{Stimulated Emission Ratio} = \frac{\text{Net Stimulated Em}_{460}}{\text{Net Stimulated Em}_{530}}$$

2. Calculate the Unstimulated Emission Ratio for each well as shown below.

$$\text{Unstimulated Emission Ratio} = \frac{\text{Net Unstimulated Em}_{460}}{\text{Net Unstimulated Em}_{530}}$$

3. Calculate the Mean Unstimulated Emission Ratio for each assay plate as shown below.

$$\text{Mean Unstimulated Emission Ratio} = \frac{\text{Sum of Unstimulated Emission Ratios}}{\text{Total Number of Wells Containing Unstimulated Cells}}$$

4. Calculate the Response Ratio for each well of interest as shown below.

$$\text{Response Ratio} = \frac{\text{Stimulated Emission Ratio (Em}_{460}/\text{Em}_{530})}{\text{Mean Unstimulated Emission Ratio (Em}_{460}/\text{Em}_{530})}$$

Note: All Response Ratios are calculated using the Mean Unstimulated Emission Ratio. The Response Ratio determined with 150 ng/ml TNFα and LiveBLAzer™ FRET B/G Substrate (CCF4-AM) was 14.3. The calculated EC₅₀ and EC₈₀ values should be ~2.37 ng/ml and ~18.7 ng/ml, respectively.

5.3 Visual Observation of Intracellular LiveBLAzer™ FRET B/G Substrate (CCF4-AM)

Note: Microscopic visualization of cells will cause photobleaching. Always read the assay plate in the fluorescence plate reader before performing microscopic visualization.

An inverted microscope equipped for epifluorescence and either a xenon or mercury excitation lamp is typically required to view the LiveBLAzer™ FRET B/G Substrate (CCF4-AM) signal in cells. To inspect the cells, a long-pass filter passing blue and green fluorescence light is needed so that you can visually identify whether the cells are fluorescing green or blue.

Recommended filter sets for observing beta-lactamase activity are described below and are available from Chroma Technologies (800-824-7662).

Chroma Set # 41031

Excitation filter: HQ405/20x (405 ±10)

Dichroic mirror: 425 DCXR

Emission filter: HQ435LP (435 long-pass)

Filter sizes vary for specific microscopes and need to be specified when the filters are ordered. For epifluorescence microscopes, a long-pass dichroic mirror is needed to separate excitation and emission light and should be matched to the excitation filter (to maximally block the excitation light around 405 nm, yet allow good transmission of the emitted light).

6.0 PURCHASER NOTIFICATION

Limited Use Label License No. 108: Lentiviral Technology

The Lentiviral Technology (based upon the lentikat™ system) is licensed from Cell Genesys, Inc., under U.S. patents and patents and applications in other countries for internal research purposes only. Use of this technology for gene therapy applications or bioprocessing other than for non-human research use requires a license from Cell Genesys (Cell Genesys, Inc. 342 Lakeside Drive, Foster City, California 94404). The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity), including non-gene therapy research and target validation applications in laboratory animals.

Limited Use Label License No. 109: Retroviral Helper Lines

Retroviral helper cell lines are licensed from Wisconsin Alumni Research Foundation, under U.S. patents and patents and applications in other countries for internal research purposes only. Use of these cell lines for Commercial Purposes requires a license from Invitrogen.

Limited Use Label License No. 150: GeneBLAzer® Technology

This product and/or its use is the subject of U.S. and foreign patents licensed to Invitrogen Corporation. The right to use this product for internal research specifically excludes the right to use this product to identify, discover, and profile compounds that act as a flavor, fragrance or taste-enhancers and modify a target identified in taste, olfaction, or pheromone detection, which compound does not require FDA approval of an NDA for claims of safety and efficacy. The right to use methods claimed in the foregoing patents with this product for research purposes can only be acquired by the use of GeneBLAzer® substrates purchased from Invitrogen Corporation or its authorized distributors.

Use of Genetically Modified Organisms (GMO)

Information for European Customers. The NFκB-*bla* Jurkat cell line(s) are genetically modified with the plasmid pLenti-*bsd*/NFκB-*bla*. As a condition of sale, use of this product must be in accordance with all applicable local legislation and guidelines including EC Directive 90/219/EEC on the contained use of genetically modified organisms.

©2006, 2010 Invitrogen Corporation. All rights reserved. Reproduction forbidden without permission.