| ٠ | • | •• | |
|---|-----|-------|------|
| Ø | IN\ | /itrc | ngen |

CellSensor[®] LEF/TCF-*bla* FreeStyle 293F Cell-based Assay

Cat. nos. K1192

Shipping: Dry Ice Storage: Liquid Nitrogen

Protocol part no. K1192.pps

Rev. date: 28 May 2008

Table of Contents

| 1. | Overview of GeneBLAzer [®] Technology | . 1 |
|----|--|---------------------|
| 2. | Materials Supplied | .1 |
| 3. | Materials Required, but Not Supplied 3.1 Optional Equipment and Materials | .2 2 |
| 4. | Cell Culture Conditions 4.1 Media Required 4.2 Growth Conditions | .3 3 |
| 5. | Assay Procedure | .4 4 5 |
| 6. | Data Analysis 6.1 Background Subtraction 6.2 Visual Observation of Intracellular Beta-lactamase Activity 6.3 Representative Data | .6 6 6 |
| 7. | Detailed Growth Methods. 7.1 Thawing Method | .7 7 7 |
| 8. | References | .7 |
| 9. | Purchaser Notification | . 8 |

1. Overview of GeneBLAzer[®] Technology

GeneBLAzer[®] Beta-lactamase Reporter Technology provides a highly accurate, sensitive and easy to use method of monitoring cellular response to drug candidates or other stimuli (1). The core of the GeneBLAzer[®] Technology is a Förster Resonance Energy Transfer (FRET) substrate that generates a ratiometric reporter response with minimal experimental noise. In addition to the dual-color (blue/green) readout of stimulated and unstimulated cells, this ratiometric method reduces the absolute and relative errors which can mask the underlying biological response of interest. Such errors include variations in cell number, transfection efficiency, substrate concentration, excitation path length, fluorescence detectors and volume changes. The GeneBLAzer[®] Beta-lactamase Reporter Technology has been proven in high-throughput screening campaigns for a range of target classes, including G-protein coupled receptors (2, 3), nuclear receptors (4-6) and kinase signaling pathways (7).

2. Materials Supplied

| Cell Line Name: | LEF/TCF- <i>bla</i> FreeStyle [™] 293F | |
|---------------------|---|--|
| Description: | CellSensor [®] LEF/TCF- <i>bla</i> FreeStyle [™] 293F cells contain a beta-lactamase reporter gene under control of the β- catenin/LEF/TCF binding elements stably integrated into FreeStyle [™] 293F cells. This cell line can be used to detect agonists/antagonists of the Wnt/β-catenin signaling pathway. LEF/TCF- <i>bla</i> FreeStyle [™] 293F cells have been shown to respond to mouse Wnt3a stimulation. | |
| Product Number: | K1192 | |
| Shipping Condition: | Dry Ice | |
| Storage Condition: | Liquid nitrogen. Immediately upon receipt, cells must be stored in liquid nitrogen. Cells cannot be stored at -80°C, as they can quickly lose viability. | |
| Quantity: | ~8,000,000 (8 × 10 ⁶ cells/ml) | |
| Application: | This cell line can be used to detect agonists/antagonists of the Wnt/ β -catenin signaling pathway. | |
| Growth Properties: | Adherent | |
| Cell Phenotype: | Epithelial | |
| Selection Marker: | Blasticidin (5 μg/ml) | |
| Vector Used: | pLenti- <i>bsd</i> /LEF/TCF- <i>bla</i> Vector | |
| Mycoplasma Testing: | Negative | |
| BioSafety Level: | 2 | |

3. Materials Required, but Not Supplied

| Media/Reagents | Recommended Source | Cat. no. |
|---|-----------------------|---|
| LiveBLAzer [™] Loading Kit LiveBLAzer [™] -FRET B/G Substrate (CCF4-AM), 5 mg DMSO for Solution A Solution B Solution C | Invitrogen | K1030 Other sizes and Loading Kits are available |
| Recovery [™] Cell Culture Freezing Medium | Invitrogen | 12648-010 |
| DMEM (high-glucose) | Invitrogen | 10569-010 |
| DMSO | Fluka | 41647 |
| Opti-MEM® Reduced Serum Medium | Invitrogen | 11058-021 |
| Fetal bovine serum, (FBS), dialyzed, tissue-culture grade (DO NOT SUBSTITUTE!) | Invitrogen | 26400-036 |
| Non-essential amino acids (NEAA) | Invitrogen | 11140-050 |
| Penicillin/Streptomycin | Invitrogen | 15140-122 |
| Phosphate-buffered saline without calcium and magnesium [PBS(-)] | Invitrogen | 14190-136 |
| HEPES (1M, pH 7.3) | Invitrogen | 15630-080 |
| Sodium pyruvate | Invitrogen | 11360-070 |
| Mouse Wnt3a | R&D Systems | 1324-WN |
| Lithium Chloride | Sigma | L4408 |
| 0.05% Trypsin/EDTA | Invitrogen | 25300-054 |
| Blasticidin antibiotic | Invitrogen | R210-01 |

| Consumables | Recommended Source | Catalog # |
|---|-----------------------|-----------|
| BIOCOAT Poly-D-Lysine CellWare [™] | BD Biosciences | 354663 |
| Compressed Air | Various | — |

| Equipment | Recommended Source |
|--|-------------------------|
| Fluorescence plate reader with bottom-read capability | Various |
| Filters if required for plate reader (see Section 5.3) | Chroma Technology Corp. |

3.1 Optional Equipment and Materials

- Epifluorescence or fluorescence-equipped microscope with appropriate filters
- Microplate centrifuge

4. Cell Culture Conditions

4.1 Media Required

| Component | Growth Medium | Thaw Medium | Assay Medium | Freezing Medium |
|---------------------------------|---------------|-------------|--------------|-----------------|
| DMEM | 90% | 90% | - | - |
| Opti- MEM® | - | - | 96.5% | - |
| Dialyzed FBS | 10% | 10% | 0.5% | - |
| NEAA | 0.1 mM | 0.1 mM | 0.1 mM | - |
| Sodium pyruvate | - | - | 1 mM | - |
| HEPES (pH 7.3) | 25 mM | 25 mM | 10 mM | - |
| Penicillin | 100 U/ml | 100 U/ml | 100 U/ml | - |
| Streptomycin | 100 µg/ml | 100 µg/ml | 100 µg/ml | - |
| Blasticidin antibiotic | 5 μg/ml | - | - | - |
| Cell Culture Freezing Medium | - | - | - | 100% |

Note: Unless otherwise stated, all media and solutions should be at least room temperature (37°C is recommended for optimal performance) before adding them to the cells.

4.2 Growth Conditions

See Section 7 for detailed growth and maintenance conditions.

- 1. Cells should be thawed in Thaw Medium without Blasticidin and grown in Growth Medium with Blasticidin. Cells should be passed or fed at least twice a week and maintained in a 37°C/5% CO₂ incubator. Cells should be maintained between 10% and 90% confluency. Do not allow cells to reach confluence.
 - *Note:* We recommend passaging cells three times after thawing before using them in the beta-lactamase assay.
- 2. Cells should be frozen at 2×10^6 cells/ml in Freezing Medium.

5. Assay Procedure

The following instructions outline the recommended procedure for monitoring the Wnt/β -catenin signaling pathway using beta-lactamase as the readout.

5.1 Quick Reference Guide

For more detailed protocol, see Section 5.2.

Plate layouts and experimental outlines will vary; in screening mode, we recommend using at least 3 wells for each condition: Unstimulated, Stimulated, and Cell-free.

Note: Certain solvents may affect assay performance. The effect of a test compound solvent should be assessed prior to screening. The cell stimulation described below is carried out in the presence of 0.5% DMSO to simulate the effect that the test compound solvent may have on the assay. If other solvents and/or solvent concentrations are used, the following assay should be changed accordingly.

| | Unstimulated Wells | Stimulated Wells | Cell-free wells |
|---|---|--|--|
| Step 1 Plating Cells | 20 μl cells in Assay Medium (10,000 cells/well). | 20 μl cells in Assay Medium (10,000 cells/well). | 20 µl Assay Medium (no cells) |
| Step 2 Incubation | Incubate the p | plate in the 37°C incubator fo | or 2 to 3 hours. |
| Step 3 Adding LiCl | 10 µl Assay Medium with 30 mM LiCl. | 10 µl Assay Medium with 30 mM LiCl. | 10 µl Assay Medium with 30 mM LiCl. |
| Step 4 Incubation | Incubate the | plate in the 37°C incubator | for 15 hours. |
| Step 4 Adding Agonist | 10 µl Assay Medium with 2% DMSO. | 10 µl 4 x mWnt3a in Assay Medium with 2% DMSO. | 10 µl Assay Medium with 2% DMSO. |
| Step 5 Incubation | Incubate the plate in the 37°C incubator for 5 hours. | | |
| Step 6 Prepare 6X Substrate Mixture | 6 µl of 1 mM LiveBLAzer [™] -FRET B/G substrate (CCF4-AM) + 60 µl of solution B, mix 934 µl of Solution C, mix | | |
| Step 7 Substrate Loading | 8 μl per well | | |
| Step 8 Substrate Incubation | 3 hrs at Room Temperature in the dark | | |
| Step 9 Detection | See Section 5.3 | | |
| Step 10 Data Analysis | See Section 6 | | |

5.2 Detailed Assay Protocol

5.2.1 Precautions

- 1. Work on a dust-free, clean surface. Always handle the poly-D-Lysine 384-well, black-wall, clear-bottom assay plate by the sides; do not touch the clear bottom of the assay plate.
- 2. If pipetting manually, it may be necessary to centrifuge the plate briefly at room temperature (for 30 seconds at $14 \times g$) after additions to ensure all the assay components are on the bottom of the wells.

5.2.2 Plating Cells

- 1. Harvest cells from culture in Growth Medium and resuspend in Assay Medium at a density of 5×10^5 cells/ml.
- 2. Add 20 µl per well of the Assay Medium to the cell free control wells. Add 20 µl per well of the cell suspension to the Unstimulated and Stimulated wells.
- 3. After plating, put the plate in 37°C incubator and incubate for 2–3 hours.

- 4. Prepare 30 mM LiCl solution in the Assay Medium.
- 5. Add 10 µl per well of the 30 mM LiCl in Assay Medium to each well.
- 6. Incubate the plate in the 37°C incubator for 15 hours.

5.2.3 Preparation of Stock Solution

- 1. Prepare Assay Medium with 2% DMSO.
- 2. Prepare a 4X mWnt3a in Assay Medium with 2% DMSO. We recommend running a dose response curve to determine the EC₈₀ for your Stimulation Solution.

5.2.4 Cell Stimulation

- 1. Add 10 µl of Assay Medium with 2% DMSO to the Unstimulated wells and to the Cell-free wells.
- 2. Add 10 µl of 4X mWnt3a to Stimulated wells.
- 3. Incubate the assay plate in a humidified $37^{\circ}C/5\%$ CO₂ incubator for 5 hours.

5.2.5 Substrate Loading and Incubation

This protocol is designed for loading cells with LiveBLAzer[™]-FRET B/G Substrate (CCF4-AM) or CCF2-AM. If alternative substrates are used please follow the loading protocol provided with the substrate.

Preparation of 6X LiveBLAzer[™]-FRET B/G Substrate (CCF4-AM) or CCF2-AM Loading Solution and cell loading should be done in the absence of direct strong lighting. Turn off the light in the hood.

- 1. Prepare Solution A: 1 mM LiveBLAzer[™]-FRET B/G Substrate (CCF4-AM, MW=1096) stock solution in dry DMSO. Store the aliquots of the stock solution at -20°C until use.
- 2. Prepare 6X LiveBLAzer[™]-FRET B/G (CCF4-AM) Substrate Mixture.

Add 6 μl of Solution A to 60 μl of Solution B and vortex.

Add 934 µl of Solution C to the combined solutions from above step with vortexing.

- 3. Remove assay plate from the incubator.
- 4. Add 8 µl of 6X Substrate Mixture from **Step 2** to each well.
- 5. Cover the plate to protect it from light and evaporation.
- 6. Incubate at room temperature for 3 hours.

5.3 Detection

All measurements are made at room temperature from the bottom of the wells, preferably in 384-well, black-wall, clear-bottom assay plates with low fluorescence background. Before reading the plate, remove dust from the bottom with compressed air.

Note: If the fluorescence plate reader with bottom read capabilities is not available, please call our tech support for options of other beta-lactamase substrates that can be read with top-reading instruments.

5.3.1 Instrumentation, Filters, and Plates

- Fluorescence plate reader with bottom reading capabilities.
- Recommended filters for fluorescence plate reader:

| Excitation filter: | 409/20 nm |
|--------------------|-----------|
| Emission filter: | 460/40 nm |
| Emission filter: | 530/30 nm |

5.3.2 Reading an Assay Plate

- 1. Set the fluorescence plate reader to bottom-read mode.
- 2. Allow the lamp in the fluorescence plate reader to warm up for at least 10 min. before making measurements.
- 3. Use the following filter selections:

| | Scan 1 | Scan 2 |
|--------------------|--|---|
| Purpose: | Measure fluorescence in the blue channel | Measure fluorescence in the green channel |
| Excitation filter: | 409/20 nm | 409/20 nm |
| Emission filter: | 460/40 nm | 530/30 nm |

6. Data Analysis

6.1 Background Subtraction

Background subtraction for both emission channels (460 nm and 530 nm) is recommended.

- 1. Use the assay plate layout to identify the location of the Cell-Free wells. These control wells are used for background subtraction.
- 2. Determine the average emission from the Cell-Free wells at both 460 nm (Average Blue Background) and 530 nm (Average Green Background).
- 3. Subtract the Average Blue background from all of the Blue emission data.
- 4. Subtract the Average Green background from all of the Green emission data.
- 5. Calculate Blue/Green Emission Ratio for each well, by dividing the background subtracted blue emission values by the background subtracted green emission values.

6.2 Visual Observation of Intracellular Beta-lactamase Activity

Note: Microscopic visualization of cells will cause photobleaching. Always read the assay plate in the fluorescence plate reader before performing microscopic visualization.

An inverted microscope equipped for epifluorescence and either a xenon or mercury excitation lamp is typically required to view the LiveBLAzerTM-FRET B/G Substrate (CCF4-AM) signal in cells. To visually inspect the cells, a long-pass filter passing blue and green fluorescence light is needed so that your eye can visually identify whether the cells are fluorescing green or blue.

Recommended filter sets for observing beta-lactamase activity are described below and are available from Chroma Technologies (800-824-7662, www.chroma.com).

Chroma Set # 41031

| Excitation filter: | $HQ405/20x (405 \pm 10)$ |
|--------------------|--------------------------|
| Dichroic mirror: | 425 DCXR |
| Emission filter: | HQ435LP (435 long-pass) |

Filter sizes vary for specific microscopes and need to be specified when the filters are ordered. For epifluorescence microscopes, a long-pass dichroic mirror is needed to separate excitation and emission light and should be matched to the excitation filter (to maximally block the excitation light around 405 nm, yet allow good transmission of the emitted light).

6.3 Representative Data



| EC ₅₀ | 15 ng/ml |
|-------------------------|-----------|
| EC ₈₀ | 46 ng/ml |
| EC ₁₀₀ | 450 ng/ml |
| Z' at EC ₁₀₀ | 0.67 |

Dose Response of LEF/TCF-*bla* **FreeStyle**[™] **293F cells to mWnt3a.** LEF/TCF-*bla* FreeStyle[™] 293F cells treated overnight with 10 mM LiCl were stimulated with indicated amount of mWnt3a for 5 hours before β-Lactamase assay was performed.

7. Detailed Growth Methods

7.1 Thawing Method

- 1. Place 29 ml of Growth Medium without Blasticidin into a T225 flask.
- 2. Place the flask in a $37^{\circ}C/5\%$ CO₂ incubator for 15 minutes to allow medium to equilibrate to the proper pH and temperature.
- 3. Remove vial of cells to be thawed from liquid nitrogen and rapidly thaw by placing at 37°C in a water bath with gentle agitation for 1-2 minutes. Do not submerge vial in water.
- 4. Decontaminate the vial by wiping with 70% ethanol before opening in a Class II biological safety cabinet.
- 5. Transfer the vial contents dropwise into 10 ml of Growth Medium without Blasticidin in a sterile 15 ml conical tube.
- 6. Centrifuge cells at $200 \times g$ for 5 minutes.
- 7. Aspirate supernatant and resuspend the cell pellet in 1 ml of fresh Growth Medium without Blasticidin.
- 8. Transfer contents to the T225 tissue culture flask containing pre-equilibrated Growth Medium without Blasticidin and place flask in the $37^{\circ}C/5\%$ CO₂ incubator.
- 9. At first passage switch to Growth Medium with Blasticidin.

7.2 Propagation Method

- 1. Cells should be passaged or fed at least twice a week. Cells should be maintained between 10% and 90% confluence. Do not allow cells to reach confluence.
- 2. To passage cells, aspirate medium, rinse once in PBS, add Trypsin/EDTA (3 ml for a T75 flask and 5 ml for a T175 flask and 8 ml for T225 flask) and swirl to coat the cells evenly. Cells usually detach after ~2–5 minutes exposure to Trypsin/EDTA. Add an equal volume of Growth Medium to inactivate Trypsin.
- 3. Verify under a microscope that cells have detached and clumps have completely dispersed.
- 4. Spin down cells and resuspend in Growth Medium.

7.3 Freezing Method

- 1. Harvest the cells as described in Section 7.2. After detachment, count the cells, then spin cells down and resuspend in 4° C Cell Culture Freezing Medium at 2×10^{6} cells/ml.
- 2. Dispense 1.0 ml aliquots into cryogenic vials.
- 3. Place in an insulated container for slow cooling and store overnight at -80°C.
- 4. Transfer to liquid nitrogen the next day for storage.

8. References

- 1. Zlokarnik, G., *et al*, Quantitation of Transcription and Clonal Selection of Single Living Cells with Beta-Lactamase as Reporter, (1998) *Science*; 279: p84-88.
- Kunapuli P., Ransom R., Murphy K., Pettibone D., Kerby J., Grimwood S., Zuck P., Hodder P., Lacson R., Hoffman I., Inglese J., Strulovici B., Development of an Intact Cell Reporter Gene Beta-lactamase Assay for G Protein-coupled Receptors, (2003) Analytical Biochem.; 314: p16-29.
- 3. Xing, H., Pollok, B., et al, A Fluorescent Reporter Assay For The Detection of Ligands Acting Through G1 Proteincoupled Receptors, (2000) J. Receptor & Signal Transduction Research; 20: p189-210.
- 4. Qureshi, S., *et al*, A One-Arm Homologous Recombination Approach for Developing Nuclear Receptor Assays in Somatic Cells, (2003) *Assay and Drug Dev. Tech*; **1**: p755-766.
- Peekhaus, N. et al, A Beta-Lactamase-Dependent Gal4-Estrogen Receptor Transactivation Assay for the Ultra-High Throughput Screening of Estrogen Receptor Agonists in a 3,456-Well Format, (2003) Assay and Drug Dev Tech; 1: p789-800.
- 6. Chin, J., et al, Miniaturization of Cell-Based Beta-Lactamase-Dependent FRET Assays to Ultra-High Throughput Formats to Identify Agonists of Human Liver X Receptors, (2003) Assay and Drug Dev. Tech.; 1: p777-787.
- Whitney M., Rockenstein E., Cantin G., Knapp T., Zlokarnik G., Sanders P., Durick K., Craig F.F., Negulescu P.A., A Genome-wide Functional Assay of Signal Transduction in Living Mammalian Cells, (1998) Nat. Biotechnol.; 16: p1329-1333.

9. Purchaser Notification

Limited Use Label License No. 5: Invitrogen Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of patents owned or controlled by Invitrogen Corporation which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Invitrogen Corporation, 1600 Faraday Avenue, Carlsbad, California 92008. Phone (760) 603-7200. Fax (760) 602-6500. Email: outlicensing@invitrogen.com.

Limited Use Label License No. 150: GeneBLAzer® Technology

This product and/or its use is the subject of one or more of U.S. Patent Nos. 5,741,657, 5,955,604, 6,291,162, and 6,472,205 and foreign equivalents licensed to Invitrogen Corporation. The right to use this product for internal research specifically excludes the right to use this product to identify, discover, and profile compounds that act as a flavor, fragrance or taste-enhancers and modify a target identified in taste, olfaction, or pheromone detection, which compound does not require FDA approval of an NDA for claims of safety and efficacy. The right to use methods claimed in the foregoing patents with this product for research purposes can only be acquired by the use of GeneBLAzer® substrates purchased from Invitrogen Corporation or its authorized distributors.

Limited Use Label License No. 51: Blasticidin and the Blasticidin Selection Marker

Blasticidin and the blasticidin resistance gene (*bsd*) are the subject of U.S. Patent No. 5,527,701 sold under patent license for research purposes only. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Invitrogen Corporation, 1600 Faraday Avenue, Carlsbad, California 92008. Phone (760) 603-7200. Fax (760) 602-6500.

Use of Genetically Modified Organisms (GMO)

Information for European Customers: The CellSensor[®] LEF/TCF-*bla* HCT-116 cell line is genetically modified with the plasmid pLenti*bsd*/LEF/TCF-*bla* Vector. As a condition of sale, use of this product must be in accordance with all applicable local legislation and guidelines including EC Directive 90/219/EEC on the contained use of genetically modified organisms.

©2008, Invitrogen Corporation. All rights reserved. Reproduction forbidden without permission.