# High Capacity cDNA Reverse Transcription Kits

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© Copyright 2006, 2007, 2010 Applied Biosystems. All rights reserved.

**For Research Use Only. Not for use in diagnostic procedures.**

Information in this document is subject to change without notice. Applied Biosystems assumes no responsibility for any errors that may appear in this document.

APPLIED BIOSYSTEMS DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL APPLIED BIOSYSTEMS BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

**DISCLAIMER OF LICENSE:**

**No license is conveyed with the purchase of this product under U.S. Patents, and corresponding patent claims outside the U.S., relating to the 5' nuclease and dsDNA-binding dye processes. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.**

**TRADEMARKS:**

AB (Design), Applied Biosystems, and GeneAmp are registered trademarks and Applera, MicroAmp, and MultiScribe are trademarks of Applied Biosystems or its subsidiaries in the U.S. and/or certain other countries.

TaqMan is a registered trademark of Roche Molecular Systems, Inc.

All other trademarks are the sole property of their respective owners.

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**Worldwide Sales and Support**

Applied Biosystems vast distribution and service network, composed of highly trained support and applications personnel, reaches 150 countries on six continents. For sales office locations and technical support, please call our local office or refer to our Web site at [www.appliedbiosystems.com](http://www.appliedbiosystems.com).

Applied Biosystems is committed to providing the world's leading technology and information for life scientists.

**Headquarters**

850 Lincoln Centre Drive Foster City,  
CA 94404 USA Phone: +1  
650.638.5800  
Toll Free (In North America): +1 800.345.5224  
Fax: +1 650.638.5884

06/2010