

## SYTOX® Red Dead Cell Stain

**Table 1.** Contents and storage information.

Material	Amount	Concentration	Storage	Stability
SYTOX® Red dead cell stain	1 mL	5 µM solution in DMSO	<ul style="list-style-type: none"> <li>• ≤-20°C *</li> <li>• Store vial upright</li> <li>• Desiccate</li> <li>• Protect from light</li> </ul>	When stored as directed, product is stable for at least 1 year.

\* Before refreezing, seal the vial tightly. The DMSO solution may be subjected to many freeze-thaw cycles without reagent degradation.


**Approximate fluorescence excitation and emission maxima:** 640/658 nm, bound to DNA

## Introduction


SYTOX® Red dead cell stain is a high-affinity nucleic acid stain that easily penetrates cells with compromised plasma membranes but will not cross uncompromised cell membranes. After brief incubation with SYTOX® Red stain, the nucleic acids of dead cells fluoresce bright red when excited with 633 or 635 nm red laser light. These properties, combined with its >500-fold fluorescence enhancement upon nucleic acid binding, make the SYTOX® Red stain a simple and quantitative single-step dead-cell indicator for use with red laser-equipped flow cytometers (Figure 1). Using 633 or 635 nm excitation, SYTOX® Red dead cell stain is distinct from other dead cell probes like 7-AAD and PI, which are excited using 488 nm. The emission of SYTOX® Red stain is limited to one channel with minimal spectral overlap. This effectively frees all the channels of the 488 nm laser line.

### Spectral Characteristics

The fluorescence excitation and emission spectra of the SYTOX® Red stain are shown in Figure 1. These spectra were obtained from samples of the dye bound to DNA. The SYTOX® Red stain exhibits a fluorescence enhancement of greater than 500-fold. The SYTOX® Red stain/DNA complex has fluorescence excitation and emission maxima of 640 nm and 658 nm, respectively.



**Figure 1.** Fluorescence excitation and emission spectra of the SYTOX® Red stain bound to DNA.



**Figure 2.** A mixture of heat-killed and untreated Jurkat cells were stained with 5 nM SYTOX® Red stain for 15 minutes. Cells were analyzed on a flow cytometer equipped with a 633 nm laser and a 660/20 nm bandpass filter. Live cells are easily distinguished from the dead cell population

## Experimental Protocol

---

The following procedure was developed using Jurkat cells (a human T-cell line) but can be adapted for any cell type. Growth medium, cell density, cell type variations, and other factors may influence staining. In initial experiments, it may be best to try a range of dye concentrations to determine the one that yields optimal staining for the given cell type and experimental conditions.

If SYTOX® Red stain is used in combination with other dyes for multicolor applications, we recommend that the other stain(s) is applied to the sample first, following all manufacturer's instructions, including washes. SYTOX® Red stain should be the last stain applied to the sample, and samples should not be washed prior to flow cytometric analysis.


### Caution

No data are available addressing the mutagenicity or toxicity of this reagent. However, SYTOX® Red stain binds to nucleic acids and should be treated as a potential mutagen and used with appropriate care. The DMSO stock solution should be handled with particular caution as DMSO is known to facilitate the entry of organic molecules into tissues. Please dispose of the reagents in compliance with all pertaining local regulations.

- 1.1 Remove the vial containing the SYTOX® Red stain from the freezer and allow the contents to equilibrate to room temperature.
- 1.2 Harvest the cell sample(s). Using an appropriate buffer, adjust the cell concentration of the sample(s) to be from  $1 \times 10^5$  to  $5 \times 10^7$  cells/mL.
- 1.3 Prepare flow cytometry tubes each containing 1 mL of cell suspension.
- 1.4 Add 1  $\mu$ L of SYTOX® Red stain (Component A) to each flow cytometry tube. The final concentration of dye will be 5 nM.
- 1.5 Incubate flow cytometry tubes for a minimum of 15 minutes at room temperature or 2–6°C, protected from light.
- 1.6 Analyze samples without washing or fixing, using 633 nm or 635 nm excitation and emission monitored at ~660 nm.

## Multicolor Staining

SYTOX® Red dead cell stain has little spectral overlap with fluorophores excited by other laser lines, and may be easily combined with other dyes. SYTOX® Red dead cell stain may be easily combined with many annexin-V conjugates to distinguish live, apoptotic, and dead cells (Figure 3).



**Figure 3.** Staining pattern of Jurkat cells (T-cell leukemia, human) treated with 10  $\mu$ M camptothecin for 4 hours (panel B) or untreated (as control, panel A). Cells were suspended in annexin-binding buffer and then treated with R-phycoerythrin (R-PE) annexin V conjugate and SYTOX® Red dead cell stain for 30 minutes, followed by flow cytometric analysis. Three cell types are identified: live (L), apoptotic (A), and dead (D). SYTOX® Red fluorescence was collected in 660/20 bandpass (633 nm red laser) and R-PE fluorescence was collected in 585/42 bandpass (488 nm blue laser).

## Product List Current prices may be obtained from our website or from our Customer Service Department.

---

Cat #	Product Name	Unit Size
S34857	SYTOX® Blue dead cell stain *for flow cytometry* *1000 assays* *1 mM solution in DMSO* .....	1 mL
S7020	SYTOX® Green nucleic acid stain *5 mM solution in DMSO* .....	250 µL
S34859	SYTOX® Red dead cell stain *for 633 or 635 nm excitation* *5 µM solution in DMSO* .....	1 mL

## Contact Information

---

### Molecular Probes, Inc.

29851 Willow Creek Road  
Eugene, OR 97402  
Phone: (541) 465-8300  
Fax: (541) 335-0504

### Customer Service:

6:00 am to 4:30 pm (Pacific Time)  
Phone: (541) 335-0338  
Fax: (541) 335-0305  
probesorder@invitrogen.com

### Toll-Free Ordering for USA:

Order Phone: (800) 438-2209  
Order Fax: (800) 438-0228

### Technical Service:

8:00 am to 4:00 pm (Pacific Time)  
Phone: (541) 335-0353  
Toll-Free (800) 438-2209  
Fax: (541) 335-0238  
probetech@invitrogen.com

### Invitrogen European Headquarters

Invitrogen, Ltd.  
3 Fountain Drive  
Inchinnan Business Park  
Paisley PA4 9RF, UK  
Phone: +44 (0) 141 814 6100  
Fax: +44 (0) 141 814 6260  
Email: euroinfo@invitrogen.com  
Technical Services: eurotech@invitrogen.com

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Service Department in Eugene, Oregon.

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

### Limited Use Label License No. 223: Labeling and Detection Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2006, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.