

GFP Counterstain (BODIPY® TR Methyl Ester)

Catalog no. C34556

Table 1. Contents and storage information.

Material	Amount	Concentration	Storage*	Stability
GFP counterstain (BODIPY® TR methyl ester)	1 mL	5 mM solution in DMSO	<ul style="list-style-type: none"> • Room temperature • Protect from light 	When stored as directed this product is stable for 1 year.
Number of assays: Sufficient material is supplied to stain ~2,000 zebrafish embryos based on the protocol below.				
Approximate fluorescence excitation/emission maxima: GFP counterstain BODIPY® TR methyl ester : 598/625 in nm.				

Introduction

Many research and biotechnological applications require detailed three- and four-dimensional visualization of embryonic cells labeled with green-fluorescent protein (GFP) within their native tissue environments. Fluorescent counterstains that label all the cells in a living embryo provide a histological context for the GFP-expressing cells in the specimen.

GFP Counterstain (BODIPY® TR methyl ester) is an excellent counterstain for GFP-expressing cells and tissues. This dye readily permeates cell membranes and localizes in endomembranous organelles such as ER, the Golgi apparatus, and mitochondria, but does not appear to localize in the plasma membrane. These localization properties make the dye an ideal vital stain that can be used to reveal the:

- Location and shapes of cell nuclei
- Shapes of cells within embryonic tissues
- Boundaries of organ-forming tissues within the whole embryo

Furthermore, GFP Counterstain (BODIPY® TR methyl ester) staining is retained after paraformaldehyde fixation, and the dye does not appear to produce any teratogenic effects on embryonic development.¹ The emission spectra of enhanced GFP (EGFP) and BODIPY® TR methyl ester are well separated, with peaks at 508 nm and 625 nm, respectively (Figure 1), allowing simultaneous dual-channel confocal imaging without significant overspill of GFP fluorescence into the BODIPY® TR methyl ester detection channel.




Figure 1. Normalized absorption (---) and fluorescence emission (----) spectra of EGFP and BODIPY® TR methyl ester.

Before Starting

Materials Required but Not Provided

- Appropriate buffer (for staining cultured cells)
- 4% formaldehyde (for staining cultured cells)

Caution

No data are available addressing the mutagenicity or toxicity of this reagent. Handle the DMSO stock solution with caution, as DMSO is known to facilitate the entry of organic molecules into tissues. Dispose of the stain in compliance with all pertaining local regulations.

Experimental Protocols

Staining Zebrafish Embryos

- 1.1 Add 5 μL of the GFP Counterstain (BODIPY® TR methyl ester) stock solution directly to 250 μL of embryo rearing medium.

The resulting staining solution contains 100 μM dye and 2% DMSO and is sufficient to simultaneously stain approximately 10 zebrafish embryos.

- 1.2 Stain embryos for 1 hour, then wash three times with plain embryo rearing medium or saline solution.

Staining Cultured Cells

- 1.1 Prepare a staining solution from the GFP Counterstain (BODIPY® TR methyl ester) stock solution.

For example, prepare a 1 μM staining solution by diluting 0.2 μL of the stock solution into 1 mL of buffer.

Note: We have used dye concentrations from 0.01 to 1.0 μM in Hanks' balanced salt solution containing 10 mM HEPES, pH 7.4, with the best results at $\sim 0.1 \mu\text{M}$. We have tested this protocol on MRC-5 human fibroblasts and bovine pulmonary aorta endothelial (BPAE) cells with generally similar results. Optimize the dye concentration for your particular cell type.

- 1.2 Incubate live cells in the staining solution for 10 minutes at 37°C.
- 1.3 Wash and resuspend the cells in balanced salt solution or medium.
- 1.4 Fix cells in 4% formaldehyde (see **Note** below).

Note: Permeabilization of cells or tissue with detergents, methanol, or acetone may remove the dye from the sample, altering the pattern or intensity of labeling.

Fluorescence Imaging

Filter sets suitable for Texas Red[®] dye are recommended for imaging samples labeled with BODIPY[®] TR methyl ester. For confocal laser-scanning microscopy (CLSM), excitation at 568 nm (argon/krypton-ion laser) or 561 nm (yellow diode laser) is recommended.

Reference

1. Personal communication, Mark Cooper,, University of Washington.

Product List

Current prices may be obtained from our website or from our Customer Service Department.

Cat. no.	Product Name	Unit Size
C34556	GFP Counterstain (BODIPY [®] TR Methyl Ester) *solution in DMSO*	1 mL

Contact Information

Molecular Probes, Inc.

29851 Willow Creek Road
Eugene, OR 97402
Phone: (541) 465-8300
Fax: (541) 335-0504

Customer Service:

6:00 am to 4:30 pm (Pacific Time)
Phone: (541) 335-0338
Fax: (541) 335-0305
probesorder@invitrogen.com

Toll-Free Ordering for USA:

Order Phone: (800) 438-2209
Order Fax: (800) 438-0228

Technical Service:

8:00 am to 4:00 pm (Pacific Time)
Phone: (541) 335-0353
Toll-Free (800) 438-2209
Fax: (541) 335-0238
probestech@invitrogen.com

Invitrogen European Headquarters

Invitrogen, Ltd.
3 Fountain Drive
Inchinnan Business Park
Paisley PA4 9RF, UK
Phone: +44 (0) 141 814 6100
Fax: +44 (0) 141 814 6260
Email: euroinfo@invitrogen.com
Technical Services: eurotech@invitrogen.com

**For country-specific contact information,
visit www.invitrogen.com.**

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Service Department in Eugene, Oregon.

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

Limited Use Label License No. 223: Labeling and Detection Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2010, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.