

Qty: 5ml

rec-Protein A-Sepharose® 4B Conjugate

Catalog No. 10-1141

Lot No.

rec-Protein A-Sepharose[®] 4B Conjugate

FORM

This product is supplied as a 50% suspension (5 ml Protein A-Sepharose 4B + 5 ml PBS) of recombinant (rec) Protein A-Sepharose B 4B in 10 mM phosphate-buffered saline, pH 7.4, containing 0.1% Proclin as a preservative. This product is prepared by covalently coupling highly purified *Staph. aureus* (Cowan strain) Protein A (produced in *E. Coli*) to Sepharose 4B.

DESCRIPTION:

Protein A and Protein G bind specifically to Fc regions of many mammalian immunoglobulins. Protein A and Protein G conjugates are commonly used as affinity adsorbents to purify immunoglobulins (antibodies) and immunoglobulin subtypes from serum, hybridoma ascites fluids, tissue culture supernatants and other biological fluids. These reagents are also commonly used to capture immune complexes generated in immunoprecipitation experiments. Because Protein A and Protein G vary widely in their abilities to bind individual immunoglobulins, consult the table below before proceeding with your application. It should be noted that the Fc binding characteristics of naturally produced Protein A and Protein G are identical to their recombinant counterparts. However, also note that Invitrogen's recombinant Protein G products are deleted for two albumin binding sites present in the wild type molecule. For certain procedures, it may be advantageous to use a light chain binding reagent (e.g. Invitrogen's Kappalock™ products) instead of Protein A or Protein G.

SPECIFICATIONS

Protein-A/Sepharose® ratio: ~2.7 mg Protein A per ml beads Binding capacity: ~22 mg Rabbit lgG per ml beads

Protein A molecular weight: 32 kDa (recombinant)

Ig Binding Sites/molecule Protein A: 4 per molecule

Albumin binding sites/molecule: none

Binding Characteristics of Immunoglobulins to Protein A and Protein G

| Species | Protein A | Protein G | Species | Protein A | Protein G | Species | Protein A | Protein G |
|-------------------|-----------|-----------|-------------|-----------|-----------|-------------|-----------|-----------|
| Cat IgG | +++ | + | Human IgA2 | + | - | Mouse IgG2b | +++ | +++ |
| Chicken IgG | = | - | Human IgD | - | - | Mouse IgG3 | ++ | +++ |
| Cow IgG | + | +++ | Human IgE | ++ | - | Mouse IgM | - | - |
| Cow IgG1 | + | +++ | Human IgG | +++ | +++ | Rat IgG | + | ++ |
| Cow IgG2 | +++ | +++ | Human IgG1 | +++ | +++ | Rat IgG1 | - | + |
| Dog IgG | +++ | + | Human IgG2 | +++ | +++ | Rat IgG2a | - | +++ |
| Goat IgG | + | ++ | Human IgG3 | - | +++ | Rat IgG2b | - | + |
| Goat IgG1 | + | +++ | Human IgG4 | +++ | +++ | Rat IgG2c | ++ | ++ |
| Goat IgG2 | +++ | +++ | Human IgM | + | - | Rabbit IgG | +++ | +++ |
| Guinea Pig IgG | +++ | + | Monkey IgG | +++ | +++ | Sheep IgG | + | ++ |
| Hamster IgG | + | ++ | Pig IgG | +++ | ++ | Sheep IgG1 | + | ++ |
| Horse IgG | + | +++ | Mouse IgG | ++ | ++ | Sheep IgG2 | +++ | +++ |
| Human IgA | + | - | Mouse IgG1 | + | ++ | | | |
| Human IgA1 | + | - | Mouse IgG2a | +++ | +++ | | | |

+++, ++, +, --: strong binding, moderate binding, weak binding, no binding, respectively.

www.invitrogen.com

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288 E-mail: techsupport@invitrogen.com

PROCEDURES:

Recommended procedures for immunoglobulin affinity purification and immunoprecipitation using our Protein A and Protein G products are available on our Web site at www.invitrogen.com/methods. Hard copy is available on request from Invitrogen Technical Services by phone at 800.955.6288 or by email at tech_support@invitrogen.com. An excellent source for information on antibodies in general and on nearly all routine laboratory procedures involving antibodies is Antibodies, A.cold-Spring Harbor Laboratory (1999).

STORAGE:

Store at 2-8°C. Do Not Freeze. Freezing this product will damage or destroy its performance.

SELECTED INVITROGEN PRODUCT REFERENCES

Immunoprecipitation

- Dikic I, et al. Nature 383(6600):547-550, 1996.
- 2. Streuli CH, et al. Cell 69:927-938, 1995.
- 3. Goi K, et al. Cancer Research 57(10):1895-1902, 1997.
- 4. Elliott JG, et al. J Biol Chem 272(21):13849-13855, 1997.
- 5. Inohara H, et al. Cancer Research 56:4530-4534, 1996.
- 6. McConnell KR, et al. J Immunol 158:2083-2089, 1997.

Affinity Puriication

- 7. Šakuai T, et al. *J Cell Biol* 136:907-918, 1997.
- 8. Lasky L, et al. Cell 69:927-938, 1992.

RELATED PRODUCTS

| Product | Cat. No. |
|--|-----------|
| Protein A-Sepharose® (5 ml) | 10-1041 |
| Protein A-Sepharose (25 ml) | 10-1042 |
| Rec-Protein A-Sepharose (25 ml) | 10-1142 |
| MonoAb rec-Protein A Column Kit | 10-1160 |
| Zysorbin™ (fixed, killed Protein A-bearing Staph. aureus for IP, 10 ml) | 10-1051-1 |
| Zysorbin [™] (fixed, killed Protein A-bearing Staph. aureus for IP, 100 ml) | 10-1051-3 |
| Rec-Protein G-Sepharose (2 ml) | 10-1241 |
| Rec-Protein G-Sepharose (5 ml) | 10-1242 |
| Zysorbin-G (fixed, killed Protein G-bearing Streptococci for IP, 10 ml) | 10-1352 |
| Kappalock™-Sepharose (binds to mammalian kappa light chains, 2 ml) | 10-1841 |

Zymed[®] is a registered trademark of Zymed Laboratories Inc. Sepharose[®] is a registered trademark of Amersham Biosciences Limited. Zysorbin[™] is a Zymed trademark. Kappalock[™] is a trademark of Aaston Inc.

For Research Use Only