

# CaptureSelect™ Antibody Affinity Resins

Pub. No. 4485776 Rev. D

- Product information ..... 1
- Specifications ..... 1
- Conditions for use ..... 2
- Instructions for use ..... 2
- Example application with CaptureSelect™ IgA affinity resin ..... 2
- Ordering information ..... 3
- Customer and technical support ..... 3
- For more information ..... 3
- Limited product warranty ..... 4

## Product information

CaptureSelect™ affinity resins can be used for the purification and isolation of proteins and/or antibodies and antibody subtypes from complex sources such as plasma, serum, and cell culture supernatants.

 **WARNING!** Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from [thermofisher.com/techresources](http://thermofisher.com/techresources).

## Specifications

All resins have a pressure limit of 3 bar.

| CaptureSelect™ Affinity Resin | Binding specificity                      | Resin and particle size           | Dynamic binding capacity per mL of resin |
|-------------------------------|--|-----------------------------------|--|
| IgA                           | Human IgA (all subclasses)               | Aldehyde-activated agarose, 70 µm | >8.0 mg IgA                              |
| IgA-CH1                       | CH1 domain of Human IgA (all subclasses) | Aldehyde-activated agarose, 70 µm | >6 mg IgA                                |
| IgG1 (Human)                  | Human IgG (subclass 1)                   | Aldehyde-activated agarose, 35 µm | >8 mg IgG1                               |
| IgG3 (Human)                  | Human IgG (subclass 3)                   | Aldehyde-activated agarose, 35 µm | >6 mg IgG3                               |
| IgG4 (Human)                  | Human IgG (subclass 4)                   | NHS-activated agarose, 90 µm      | >6 mg IgG4                               |
| IgG-CH1                       | CH1 domain of human IgG                  | Aldehyde-activated agarose, 70 µm | >15 mg IgG                               |
| FcXL                          | Human IgG (all subclasses)               | Aldehyde-activated agarose, 70 µm | >20 mg IgG                               |

| CaptureSelect™ Affinity Resin  | Binding specificity                         | Resin and particle size             | Dynamic binding capacity per mL of resin |
|--|---|-------------------------------------|--|
| IgG-Fc (Multi-species) (Human, primates, rat, mouse, guinea pig, bovine, horse, goat, sheep) | IgG Fc                                      | Aldehyde-activated agarose, 35 µm   | >15 mg IgG                               |
| IgM  | Human, mouse, and rat IgMs                  | NHS-activated agarose, 90 µm        | >2.5 mg IgM                              |
| POROS™ IgM   | Human, mouse, and rat IgMs                  | CDI-activated agarose, POROS™ 50 µm | >5.0 mg IgM                              |
| KappaXL  | Human Igs (kappa light chain)               | Aldehyde-activated agarose, 70 µm   | >20 mg IgG                               |
| LC-kappa (Murine)  | Murine Igs (kappa light chain)              | Aldehyde-activated agarose, 35 µm   | >10 mg IgG                               |
| LC-lambda (Human)  | Human Igs (lambda light chain)              | NHS-activated agarose, 90 µm        | >10 mg IgG                               |
| LC-lambda (Mouse)  | Mouse Igs (lambda light chain)              | Aldehyde-activated agarose, 70 µm   | >10 mg IgG                               |
| LC-lambda (Rat)  | Rat and guinea pig Igs (lambda light chain) | Aldehyde-activated agarose, 70 µm   | >10 mg IgG                               |
| LC-lambda (Ungulate) (Hoofed animals such as horse, sheep, cow)                              | Ungulate Igs (lambda light chain)           | Aldehyde-activated agarose, 35 µm   | >10–15 mg IgG, depending on species      |

## Conditions for use

All resins have a flow rate of 150 cm/h.

| CaptureSelect™ Affinity Resin | Equilibration/wash buffer                             | Elution buffer  |
|-------------------------------|---|---|
| IgA                           | PBS, pH 7.2–7.4 (physiological pH and ionic strength) | 0.1 M Glycine, pH 3.0   |
| IgA-CH1                       |   | Neutral pH: 20 mM Tris, 2.0 M MgCl <sub>2</sub> , pH 7.0<br>Acidic: 0.1 M Glycine, pH 3.0                               |
| IgG1 (Human)                  |   | 0.1 M Glycine, pH 3.0   |
| IgG3 (Human)                  |   | 0.1 M Glycine, pH 3.0   |
| IgG4 (Human)                  |   | 0.1 M Glycine, pH 3.0   |
| IgG-CH1                       |   | 0.1 M Glycine, pH 3.0   |
| FcXL                          |   | Neutral pH: 20 mM Tris, 1.0 M MgCl <sub>2</sub> , 40% propylene glycol, pH 5.0–7.0<br>Acidic: 0.1 M Glycine, pH 4.0–3.0 |
| IgG-Fc (Multi-species)        |   | 0.1 M Glycine, pH 3.0   |
| IgM                           |   | 0.1 M Glycine, pH 3.0   |
| POROS™ IgM                    |   | 0.1 M Glycine, pH 3.0   |
| KappaXL                       |   | 0.1 M Glycine, pH 4.0–3.0   |
| LC-kappa (Murine)             |   | 0.1 M Glycine, pH 3.0   |
| LC-lambda (Human)             |   | 0.1 M Glycine, pH 3.0–2.0   |
| LC-lambda (Mouse)             |   | 0.1 M Glycine, pH 3.0   |
| LC-lambda (Rat)               |   | 0.1 M Glycine, pH 3.0   |
| LC-lambda (Ungulate)          |   | 0.1 M Glycine, pH 3.0   |

## Instructions for use

- Pack the column.
- Equilibrate with 5 to 10 column volumes (CV) of the equilibration/wash buffer recommended in “Conditions for use” on page 2.
- Prepare and load the sample.  
The sample loading volume depends on the concentration of the target molecule and the dynamic binding capacity of the resin. See “Specifications” on page 1.
  - Dissolve, dilute, or exchange samples into the equilibration buffer. This is particularly important for large samples (greater than 25% of the column volume).
  - Centrifuge and filter samples (0.22 or 0.45 µm) before injection.
- Wash with 5 to 10 CV of the equilibration/wash buffer recommended in “Conditions for use” on page 2, or until you see a stable baseline.
- Elute with 5 to 10 CV of the elution buffer recommended in “Conditions for use” on page 2, or until you see a stable baseline.
- Re-equilibrate with 5 to 10 CV of the equilibration/wash buffer recommended in “Conditions for use” on page 2, or until you see a stable baseline.
- Re-equilibrate in equilibration/wash buffer.  
If the column will not be used immediately, store the resin in 20% ethanol at 4°C (39°F), stable for up to 1 year.

## Example application with CaptureSelect™ IgA affinity resin

Refer to [www.lifetechnologies.com/captureselect](http://www.lifetechnologies.com/captureselect) for additional examples.

Resin: IgA

Sample: Human serum (200 µL)

Sample preparation: Loaded on 400-µL CaptureSelect™ IgA affinity resin


Equilibration and wash buffer: PBS, pH 7.4

Elution buffer: 0.1 M Glycine, pH 3.0

Flow: 150 cm/h


Elution fractions neutralized with 0.1 volume of 1 M Tris, pH 8.0

Starting material, flow through, and elution fractions analysis: SYPRO™ Ruby gel-stained non-reduced 4–20% acrylamide Tris-Glycine gel



1 Cycle is:

- 10 column volumes (CV) equilibration
- Sample loading
- 10 CV wash out unbound sample
- 5 CV elution
- 10 CV re-equilibration



- 1: Molecular weight marker  
 2: Human serum  
 3: Flow through IgA affinity resin  
 4: Elution IgA affinity resin

## Ordering information

| CaptureSelect™ Affinity Resin | Cat. no.   |
|-------------------------------|--|
| IgA                           | 194288005 (5 mL)<br>194288010 (10 mL)<br>194288050 (50 mL)   |
| IgA-CH1                       | 194311005 (5 mL)<br>194311010 (10 mL)<br>194311050 (50 mL)   |
| IgG1 (Human)                  | 191303005 (5 mL)<br>191303010 (10 mL)<br>191303050 (50 mL)   |
| IgG3 (Human)                  | 191304005 (5 mL)<br>191304010 (10 mL)<br>191304050 (50 mL)   |
| IgG4 (Human)                  | 290005 (5 mL)<br>290010 (10 mL)<br>290050 (50 mL)            |
| IgG-CH1                       | 194320005 (5 mL)<br>194320010 (10 mL)<br>194320050 (50 mL)   |
| IgG-CH1 pre-packed columns    | 494320001 (five 1-mL columns)<br>494320005 (one 5-mL column) |
| FcXL                          | 194328005 (5 mL)<br>194328010 (10 mL)<br>194328050 (50 mL)   |
| FcXL pre-packed columns       | 494328001 (five 1-mL columns)                                |

| CaptureSelect™ Affinity Resin | Cat. no.   |
|-------------------------------|--|
|                               | 494328005 (one 5-mL column)  |
| IgG-Fc (Multi-species)        | 191285505 (5 mL)<br>191285510 (10 mL)<br>191285550 (50 mL)<br>1912855250 (250 mL)<br>1912855500 (500 mL) |
| IgM                           | 289005 (5 mL)<br>289010 (10 mL)<br>289050 (50 mL)  |
| POROS™ IgM                    | 195289005 (5 mL)<br>195289010 (10 mL)<br>195289050 (50 mL)   |
| KappaXL                       | 194321005 (5 mL)<br>194321010 (10 mL)<br>194321050 (50 mL)   |
| KappaXL pre-packed columns    | 494321001 (five 1-mL columns)<br>494321005 (one 5-mL column)   |
| LC-kappa (Murine)             | 191315005 (5 mL)<br>191315010 (10 mL)<br>191315050 (50 mL)   |
| LC-lambda (Human)             | 084905 (5 mL)<br>084910 (10 mL)<br>084950 (50 mL)  |
| LC-lambda (Mouse)             | 194323005 (5 mL)<br>194323010 (10 mL)<br>194323050 (50 mL)   |
| LC-lambda (Rat)               | 194324005 (5 mL)<br>194324010 (10 mL)<br>194324050 (50 mL)   |
| LC-lambda (Ungulate)          | 191314005 (5 mL)<br>191314010 (10 mL)<br>191314050 (50 mL)   |

## Customer and technical support

Visit [thermofisher.com/techresources](http://thermofisher.com/techresources) for the latest in services and support, including:

- Worldwide contact telephone numbers
- Product support, including:
  - Product FAQs
  - Software, patches, and updates
- Order and web support
- Product documentation, including:
  - User guides, manuals, and protocols
  - Certificates of Analysis
  - Safety Data Sheets (SDSs; also known as MSDSs)

**Note:** For SDSs for reagents and chemicals from other manufacturers, contact the manufacturer.

## For more information

For more information on CaptureSelect™ and POROS™ products, go to [www.lifetechnologies.com/captureselect](http://www.lifetechnologies.com/captureselect).

## Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at

[www.lifetechnologies.com/termsandconditions](http://www.lifetechnologies.com/termsandconditions). If you have any questions, please contact Life Technologies at [www.lifetechnologies.com/support](http://www.lifetechnologies.com/support).

The information in this guide is subject to change without notice.

### DISCLAIMER

TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2015 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

For support visit [thermofisher.com/techresources](http://thermofisher.com/techresources) or email [techsupport@lifetech.com](mailto:techsupport@lifetech.com)

[thermofisher.com/lifescience](http://thermofisher.com/lifescience)