

EfficientFeed™ B+ AGT™ Supplement

Description

EfficientFeed™ B+ AGT™ Supplement is a dry-format, single-part product for use in fed batch culture of multiple cell lines. EfficientFeed™ B+ AGT™ Supplement is an animal-origin free, chemically defined formulation that contains no proteins, hydrolysates, or components of unknown composition. The product provides an easy to reconstitute, pH neutral solution that can be concentrated from the suggested 52.7 g/L (1X) up to 158.1 g/L (3X). As a result, key nutrients can be delivered in smaller volumes.

Product	Catalog No.	Amount	Storage	Shelf life
EfficientFeed™ B+ AGT™ Supplement	A25030-04	1 × 1 L	2°C to 8°C; Store dark and dry	12 months
	A25030-05	1 × 10 L		
	A25030-01	1 × 100 L		

Product use

Caution: For manufacturing, processing, or repacking.

Safety information

Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

Reconstitute EfficientFeed™ B+ AGT™ Supplement

1. Measure 80% of the final volume deionized or distilled water at room temperature (15°C to 30°C).
2. Add EfficientFeed™ B+ AGT™ Supplement (52.7 grams/L) to water. Mix until completely dissolved and for a minimum of 20 minutes. Solution will be cloudy, improving with time.
3. Using a calibrated vessel, dilute to final volume with deionized or distilled water. Mix for an additional 10 minutes.
4. Measure pH and osmolality. Confirm final pH of 6.8–7.3. Final osmolality should be 345–385 mOsm/kg.
5. Sterilize immediately by membrane filtration (positive pressure recommended).

Note: Once product is filtered, use immediately or product can be stored at 2°C to 8°C for up to 12 months. Protect from light.

Flexibility in reconstitution

EfficientFeed™ B+ AGT™ Supplement can be flexibly concentrated following these recommendations:

Concentration	Dissolution time
1X (52.9 g/L)	20 min
3X (158.1 g/L)	50 min

Use

EfficientFeed™ B+ AGT™ Supplement can be added as multiple discrete additions every 1–3 days or as a continuous feed. Total feed volume can range from 10–45% (1X) to 3–15% (3X) of the original culture volume depending on the cell line and process. For best results, it is recommended to initiate feeding only after reaching mid to late exponential phase.

Feeding examples

EfficientFeed™ B+ AGT™ Supplement example feeding strategies:

Conc.	Culture days										
	0	1	2	3	4	5	6	7	8	9	10
1X (52.7 g/L)	-	-	-	-	10%	-	10%	-	10%	-	10%
3X (158.1 g/L)	-	-	-	-	3.3%	-	3.3%	-	3.3%	-	3.3%

Related products

Product	Catalog No.
CD OptiCHO™ Medium, (1X) liquid	12681
CD CHO Medium, (1X) liquid	10743
CD FortiCHO™ Medium	A11483
EfficientFeed™ A+ AGT™ Supplement	A25023
EfficientFeed™ C+ AGT™ Supplement	A25031
Trypan Blue Solution	15250
Water, Distilled	15230

Explanation of symbols and warnings

The symbols present on the product label are explained below:

				
Caution	Manufacturer	Consult instructions for use	Keep away from light	Temperature limitation
				
Use by	Lot/batch #	Cat. no.	Sterilized using aseptic processing techniques	

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

Limited Use Label License No. 517:

Internal Research and Bioproduction Use

Notice to Purchaser: The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product to (a) perform internal research for the sole benefit of the purchaser; (b) manufacture protein (or other biological material) for resale; and (c) perform research or manufacturing services conducted by the purchaser on a fee for service or contract basis for or on behalf of third parties. However, the purchaser may transfer this product, its components, or materials made using this product to a third party (including contract research/manufacturing organizations), provided that each such third party agrees in writing to use such product, components, or materials solely on behalf of the purchaser, and such third party is restricted from further transferring any such product, components, or materials to any individual or entity other than the purchaser. No additional rights are granted. By purchasing this product, the purchaser agrees not to: (1) resell the product in any form; (2) use the product as a therapeutic agent or diagnostics test component; (3) reverse engineer the product or cause the product to be reverse engineered; or (4) use the product for purposes other than what is indicated in this Limited Use Label License. Life Technologies is not aware of Intellectual Property ("IP") that would be infringed by the manufacture or sale of its products. Customers are urged to perform an IP search and analysis specific to their manufacturing processes and target biologic products. Should that analysis require information about Life Technologies products, Life Technologies will provide to the customer's IP counsel the information necessary to analyze any relevant IP in a manner that protects Life Technologies proprietary information. The purchaser is responsible for obtaining all regulatory approvals necessary for any therapeutic or diagnostic use of the protein (or biological material) manufactured using this product. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Corporate entity
Life Technologies | Carlsbad, CA 92008 USA | Toll Free in USA 1.800.955.6288

© 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

For support visit thermofisher.com/support or email techsupport@lifetech.com

thermofisher.com

ThermoFisher
SCIENTIFIC