



## Instruction Manual

### **WesternBreeze®**

Chromogenic Western Blot  
Immunodetection Kit

Catalog nos. WB7103, WB7105, WB 7107

**Version F**

June 4, 2003

*IM-1004*



# Table of Contents

Table of Contents.....	i
Important Information.....	1
General Guidelines.....	3
Preparing Solutions.....	5
WesternBreeze <sup>®</sup> Chromogenic Immunodetection Protocol...	7
Troubleshooting.....	8
Technical Service.....	13



# Important Information

---

## Introduction

The WesternBreeze® Chromogenic Immunodetection Kit contains complete, optimized, ready-to-use or ready-to-dilute reagents for sensitive immunodetection of western blots or dot blots. The kit is used to detect primary antibodies immobilized on nitrocellulose (NC) or polyvinylidene difluoride (PVDF) membranes.

---

## Contents

The components included in the WesternBreeze® Chromogenic Immunodetection Kit are listed below. Sufficient reagents are supplied to detect 20 mini-blots.

Item	Description	Amount
Blocker/Diluent (Part A)	Concentrated buffered saline solution containing detergent	80 ml
Blocker/Diluent (Part B)	Concentrated Hammersten casein solution	80 ml
Antibody Wash (16X)	Concentrated buffered saline solution containing detergent	2 x 100 ml
Chromogenic Substrate	Ready-to-use solution of BCIP/NBT substrate for alkaline phosphatase.	100 ml
Secondary Antibody Solution	Ready-to-use solution of alkaline phosphatase-conjugated, affinity purified, anti-species IgG (anti-mouse and anti-rabbit kits: goat conjugate, anti-goat kits: rabbit conjugate)	2 x 100 ml
Incubation Dishes	100 mm x 100 mm x 15 mm plastic dishes with lids (Fisher No. O8-757-11A)	2

---

*Continued on next page*

# Important Information, Continued

---

## Storage

All solutions are proprietary formulations and contain chlorobutanol as a preservative; once diluted, they should be used that day. Store all solutions at 2-8°C.

---

## Safety Issues

The WesternBreeze® Chromogenic Immunodetection Kit, contains a solution of 5-bromo-4-chloro-3-indolyl-1-phosphate (BCIP) and nitro blue tetrazolium (NBT) which is a possible carcinogen and dilute solutions of irritating chemicals. We recommend wearing gloves, safety glasses, and a lab coat when using the WesternBreeze® Chromogenic Immunodetection Kit.

---

## Product Qualification

Invitrogen qualifies the WesternBreeze® Chromogenic Immunodetection Kit as follows: 10 -100 pg human IgG in NuPAGE® LDS sample buffer is electrophoresed on a 4-12% NuPAGE® Novex Bis-Tris Gel using NuPAGE® MES SDS Running Buffer. After electrophoresis, proteins are transferred onto a nitrocellulose membrane using NuPAGE® Transfer Buffer. Immunodetection is performed as described in this manual using anti-human IgG primary antibody. The amount of antigen as shown in the table below must be detected within 40 minutes of color development.

<b>Kit</b>	<b>Amount</b>
Rabbit	10 pg
Mouse	100 pg
Goat	100 pg

---

# General Guidelines

---

## General Guidelines

To obtain the best results with WesternBreeze<sup>®</sup> Chemiluminescent Immunodetection Kit:

- Use a single, clean dish for each blot.
- Avoid touching the working surface of the membrane, even with gloves.
- Use pure water, free from alkaline phosphatase activity. Stored water should be autoclaved or ultra-filtered to remove alkaline phosphatase activity. Freshly ultra-filtered water is preferred.
- Avoid cross-contamination of system solutions especially with the alkaline phosphatase substrate solution.
- Perform all washing, blocking, and incubating steps on a rotary shaker platform rotating at 1 revolution/second.
- Work quickly when changing solutions as PVDF membranes dry quickly. If the membrane dries, re-wet the membrane with methanol and rinse with water before proceeding (see page 4).
- Add solutions to the trays slowly, at the membrane edge, to avoid bubbles forming under the membrane. Decant from the same corner of the dish to ensure complete removal of previous solutions.



### Note

Since the amount of reagents used for the detection of PVDF and nitrocellulose membranes is different, it is normal for some reagents to be in excess in the kit after detecting 20 mini-blots.

---

*Continued on next page*

# General Guidelines, Continued

---

## Materials Needed

- Blotted membranes containing applied antigen or primary antibody samples
  - Primary antibody to detect applied antigen, if appropriate
  - Purified water autoclaved, sterile, or ultra filtered to remove alkaline phosphatase activity from all solutions used in the procedure
  - Clean flasks for preparing solutions
  - Forceps for manipulating blotted membranes.
  - Orbital shaker platform
- 

## Preparing the Membrane

For western blots or western dot blots, from SDS-PAGE gels to NC or PVDF membranes, wash the membranes twice for 5 minutes with 20 ml of pure water to remove gel and transfer buffer components and some weakly bound proteins. Proceed to immunodetection (see page 7).

If you are using water-washed and dried NC membranes, proceed to immunodetection (see page 7). For water-washed and dried PVDF membranes, re-wet the membrane in methanol followed by two 20 ml water washes for 5 minutes. Proceed to immunodetection (see page 7).

For native-PAGE western blots or western dot blots, a drying step is recommended to improve protein binding to the membrane. Wash the dried NC membranes twice with 20 ml water for 5 minutes before proceeding to immunodetection (see page 7). Re-wet the dried PVDF membranes in methanol, followed by two water washes of 20 ml for 5 minutes before proceeding to immunodetection (see page 7).

---

# Preparing Solutions

---

## Nitrocellulose Membranes

For western blots from Novex<sup>®</sup> or other mini gels, prepare the solutions for nitrocellulose membranes as described in the table below for ~ 60 cm<sup>2</sup> membrane.

**Note:** The Blocking Solution is used for blocking and as a Primary Antibody Diluent for NC membranes.

Solution	For NC Membrane
Blocking Solution	Ultra filtered Water 14 ml
	Blocker/Diluent (Part A) 4 ml
	<u>Blocker/Diluent (Part B) 2 ml</u>
	Total Volume 20 ml
Primary Antibody Diluent	Dilute your primary antibody according to the manufacturers recommendations into 10 ml of NC Blocking Solution (see above). Typically, commercial primary antibody preparations are diluted 1:1000 to 1:5000 to a concentration of about 1 to 0.2 µg/ml
Antibody Wash	Ultra filtered Water 150 ml
	<u>Antibody Wash Solution (16X) 10 ml</u>
	Total Volume 160 ml

---

*Continued on next page*

## Preparing Solutions, Continued

---

### PVDF Membranes

For western blots from Novex<sup>®</sup> or other mini gels, prepare the solutions for PVDF membranes as described in the table below for ~ 60 cm<sup>2</sup> membrane.

Solution	For PVDF Membrane
Blocking Solution	Ultra filtered Water 5 ml
	Blocker/Diluent (Part A) 2 ml
	<u>Blocker/Diluent (Part B) 3 ml</u>
	Total Volume 10 ml
Primary Antibody Diluent	Ultra filtered Water 7 ml
	Blocker/Diluent (Part A) 2 ml
	<u>Blocker/Diluent (Part B) 1 ml</u>
	Total Volume 10 ml
Dilute your primary antibody into this diluent according to the manufacturer's recommendations. Typically, commercial primary antibody preparations are diluted 1:1000 to 1:5000 to a concentration of about 1 to 0.2 µg/ml.	
Antibody Wash	Ultra filtered Water 150 ml
	<u>Antibody Wash Solution (16X) 10 ml</u>
	Total Volume 160 ml

---

# WesternBreeze<sup>®</sup> Chromogenic Immunodetection Protocol

---

## Protocol for Small Membranes (60 cm<sup>2</sup>)

1. Place the membrane in 10 ml of the appropriate Blocking Solution in a covered, plastic dish provided in the kit. Incubate for 30 minutes on a rotary shaker set at 1 revolution/sec. Decant the Blocking Solution.
2. Rinse the membrane with 20 ml of water for 5 minutes, then decant. Repeat once.
3. Incubate the membrane with 10 ml of Primary Antibody Solution for 1 hour, then decant.
4. Wash the membrane for 5 minutes with 20 ml of prepared Antibody Wash, then decant. Repeat 3 times.
5. Incubate the membrane in 10 ml of Secondary Antibody Solution for 30 minutes, then decant.
6. Wash the membrane for 5 minutes with 20 ml of Antibody Wash, then decant. Repeat 3 times.
7. Rinse the membrane with 20 ml of water for 2 minutes, then decant. Repeat twice.
8. Incubate the membrane in 5 ml of Chromogenic Substrate until purple bands develop on the membrane. Development is complete in 1 to 60 minutes.
9. Rinse the membrane with 20 ml of water for 2 minutes. Repeat twice.
10. Dry the membrane on a clean piece of filter paper to open air, by a stream of slightly warm air, or under an infrared lamp.

---

## Protocol for Large Membranes (200 cm<sup>2</sup>)

To blot standard size gels (~200 cm<sup>2</sup>), scale up the required solution volumes by 3.3X. Use a tray that closely matches the dimensions of the membrane for most efficient use of the solutions.

---

# Troubleshooting

---

## Introduction

Review the information provided below to troubleshoot your experiments.

---

## High Background

- **Membrane not completely wetted**

Follow instructions for pre-wetting the membrane. Use an incubation dish which is small enough to allow thorough coverage of membrane to prevent drying out. Shake or agitate during each step.
- **Membrane is contaminated**

Use only clean, new membranes. Wear clean gloves at all times and use forceps when handling membranes.
- **Blocking solution is incorrect**

Nitrocellulose and PVDF require different blocking conditions; make sure the appropriate solution volumes were used for preparation.
- **Blocking time or washing time is too short**

Make sure that each step is performed for the specified amount of time. Strictly adhere to all wash times.
- **Solutions or incubation tray is contaminated**

Use clean glassware and purified water to prepare solutions. Replace or clean the tray thoroughly with a glassware-cleaning detergent. Rinse thoroughly with purified water. Wear clean gloves at all times.
- **Blot is overdeveloped**

Follow recommended developing time or remove blot from substrate when signal-to-noise ratio is acceptable.

---

*Continued on next page*

# Troubleshooting, Continued

---

## High Background, continued

- **Concentrated Primary antibody used**  
Follow the supplier's recommended dilution or determine the optimum concentration by dot-blotting.
- 

## Non-Specific Binding

- **Membrane contaminated by fingerprints or keratin proteins**  
Wear clean gloves at all times and use forceps when handling membranes. Always handle membranes around the edges.
  - **Primary antibody too concentrated**  
Follow the supplier's recommended dilution or determine the optimum concentration by dot-blotting.
  - **Insufficient removal of SDS/weakly bound proteins from membrane after blotting**  
Follow instructions for membrane preparation before immunodetection (see page 4).
  - **Short blocking time or long washing time**  
Make sure that each step is performed for the specified amount of time.
  - **Affinity of the primary antibody for the protein standards**  
Check with the protein standard manufacturer for homologies with primary antibody.
- 

*Continued on next page*

# Troubleshooting, Continued

---

## **Weak or No Signal**

- **Poor or incomplete transfer**  
Refer to Western Blotting instructions (IM-9051) and repeat blot. After blotting, stain membrane to measure transfer efficiency. Use positive control and/or molecular weight marker.
- **Membrane not completely wet**  
Follow instructions for pre-wetting the membrane. The incubation dish must be small enough to allow thorough coverage of membrane to prevent drying out. Shake or agitate during each step.
- **Primary antibody concentration too low**  
Follow the supplier's recommended dilution or determine the optimum concentration by dot-blotting.
- **Inactive primary antibody**  
Determine activity by performing a dot-blot.
- **Low Affinity of primary antibody to antigen**  
Obtain a higher affinity primary antibody.
- **Contaminated secondary antibody solution**  
Wear gloves at all times and keep bottles tightly capped when not in use. Use only purified water when preparing reagents.
- **Protein of interest ran off the gel**  
Match gel separation range to size of protein being transferred.
- **Poor retention of proteins**  
Match gel separation range to size of protein being transferred. Use a molecular weight marker with relevant size proteins. Larger proteins require more transfer time, smaller proteins less. Use membrane with the appropriate binding capacity (see Western Blotting instructions).

---

*Continued on next page*

## Troubleshooting, Continued

---

### **Weak or No Signal,** continued

- **Sample improperly prepared; antigenicity weakened, or destroyed**  
SDS and reducing agents may interfere with some antibody / antigen affinities.
- **Sample too dilute**  
Load a higher concentration or amount of protein onto the gel.
- **Protein weakly bound to membrane**  
Ensure that transfer buffer contains 10-20% methanol.
- **Incorrect blocking solution**  
Nitrocellulose and PVDF require different blocking conditions; make sure the appropriate solution volumes were used for preparation.
- **Insufficient substrate incubation**  
Perform each step for the specified amount of time or remove blot from substrate when signal-to-noise ratio is acceptable.
- **Substrate is contaminated**  
Wear gloves at all times and keep bottles tightly capped when not in use.
- **Blots are too old**  
Protein may have broken down over time. Use freshly prepared blots.

---

*Continued on next page*

## Troubleshooting, Continued

---

### "Spotted" Membrane

- **Poor or incomplete transfer**  
Refer to Western Blotting instructions (IM-9051) and repeat blot.
  - **Membrane pads are dirty or contaminated**  
Soak with detergent and rinse thoroughly with purified water before use. Replace pads when they become worn or discolored.
  - **Membrane not completely wetted**  
Follow instructions for pre-wetting the membrane. The incubation dish must be small enough to allow thorough coverage of membrane to prevent drying out. Shake or agitate during each step.
  - **Membrane is contaminated by fingerprints or keratin proteins**  
Wear clean gloves at all times and use forceps when handling membranes. Always handle membranes around the edges.
  - **Uneven blocking**  
The incubation dish must be small enough to allow thorough coverage of membrane. Shake or agitate during each step.
- 

### Large, Scattered Signal

- **Protein is overloaded**  
Reduce load or dilute concentration of sample.
  - **Poor or incomplete transfer**  
Refer to Western Blotting instructions (IM-9051) and repeat blot.
  - **Primary antibody is too concentrated**  
Follow the supplier's recommended dilution or determine the optimum concentration by dot-blotting.
-

# Technical Service

---

## World Wide Web



Visit the Invitrogen Web Resource using your World Wide Web browser. At the site, you can:

- Get the scoop on our hot new products and special product offers
- View and download vector maps and sequences
- Download manuals in Adobe® Acrobat® (PDF) format
- Explore our catalog with full color graphics
- Obtain citations for Invitrogen products
- Request catalog and product literature

Once connected to the Internet, launch your Web browser (Internet Explorer 5.0 or newer or Netscape 4.0 or newer), then enter the following location (or URL):

<http://www.invitrogen.com>

...and the program will connect directly. Click on underlined text or outlined graphics to explore. Don't forget to put a bookmark at our site for easy reference!

---

## Contact Us

For more information or technical assistance, call, write, fax, or email. Additional international offices are listed on our web page ([www.invitrogen.com](http://www.invitrogen.com)).

### United States Headquarters:

Invitrogen Corporation  
1600 Faraday Avenue  
Carlsbad, CA 92008 USA

Tel: 1 760 603 7200

Tel (Toll Free): 1 800 955 6288

Fax: 1 760 602 6500

E-mail: [tech\\_service@invitrogen.com](mailto:tech_service@invitrogen.com)

### European Headquarters:

Invitrogen Ltd  
Inchinnan Business Park  
3 Fountain Drive

Paisley PA4 9RF, UK

Tel: +44 (0) 141 814 6100

Tech Fax: +44 (0) 141 814 6117

E-mail: [eurotech@invitrogen.com](mailto:eurotech@invitrogen.com)

---

## MSDS Requests

To request an MSDS, visit our Web site at [www.invitrogen.com](http://www.invitrogen.com). On the home page, go to 'Technical Resources', select 'MSDS', and follow instructions on the page.

---

## Technical Service, Continued

---

### Limited Warranty

Invitrogen is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, please contact our Technical Service Representatives. Invitrogen warrants that all of its products will perform according to the specifications stated on the certificate of analysis. The company will replace, free of charge, any product that does not meet those specifications. This warranty limits Invitrogen Corporation's liability only to the cost of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. Invitrogen reserves the right to select the method(s) used to analyze a product unless Invitrogen agrees to a specified method in writing prior to acceptance of the order. Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore Invitrogen makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, please report it to our Technical Service Representatives. **Invitrogen assumes no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.**

---

©1999-2003 Invitrogen Corporation. All rights reserved.

For research use only. Not intended for any animal or human therapeutic or diagnostic use.





**Corporate Headquarters:**

*Invitrogen Corporation  
1600 Faraday Avenue  
Carlsbad, California 92008  
Tel: 1 760 603 7200  
Tel (Toll Free): 1 800 955 6288  
Fax: 1 760 603 7229  
Email: tech\_service@invitrogen.com*

**European Headquarters:**

*Invitrogen Ltd  
3 Fountain Drive  
Inchinnan Business Park  
Paisley PA4 9RF, UK  
Tel (Free Phone Orders): 0800 269 210  
Tel (General Enquiries): 0800 5345 5345  
Fax: +44 (0) 141 814 6287  
Email: eurotech@invitrogen.com*

**International Offices:**

*Argentina 5411 4556 0844  
Australia 1 800 331 627  
Austria 0800 20 1087  
Belgium 0800 14894  
Brazil 0800 11 0575  
Canada 800 263 6236  
China 10 6849 2578  
Denmark 80 30 17 40*

*France 0800 23 20 79  
Germany 0800 083 0902  
Hong Kong 2407 8450  
India 11 577 3282  
Italy 02 98 22 201  
Japan 03 3663 7974  
The Netherlands 0800 099 3310  
New Zealand 0800 600 200  
Norway 00800 5456 5456*

*Spain & Portugal 900 181 461  
Sweden 020 26 34 52  
Switzerland 0800 848 800  
Taiwan 2 2651 6156  
UK 0800 838 380  
For other countries see our Web site*

[www.invitrogen.com](http://www.invitrogen.com)

