



ELISA Kit  
Catalog #KHR9081

***Human***  
**EGFR\***  
**[pY1068]**

[www.invitrogen.com](http://www.invitrogen.com)

Invitrogen Corporation  
542 Flynn Road, Camarillo, CA 93012  
Tel: 800-955-6288

**E-mail: [techsupport@invitrogen.com](mailto:techsupport@invitrogen.com)**

\*Patent Pending



## TABLE OF CONTENTS

Introduction .....	4
Principle of the Method.....	6
Reagents Provided.....	7
Supplies Required But Not Provided .....	8
Procedural Notes/Lab Quality Control .....	8
Safety .....	10
Directions for Washing.....	10
Procedure for Extraction of Proteins from Cells.....	11
Reagent Preparation and Storage.....	13
Reconstitution and Dilution of Human EGFR [pY1068] Standard .	13
Dilution of Human EGFR [pY1068] Standard .....	14
Storage and Final Dilution of Anti-Rabbit IgG HRP (100X).....	15
Dilution of Wash Buffer .....	16
Assay Method.....	16
Typical Data .....	19
Limitations of the Procedure.....	20
Performance Characteristics .....	21
Sensitivity .....	21
Precision .....	22
Recovery.....	23
Parallelism .....	23
Linearity of Dilution .....	24
Specificity.....	25
References .....	28

## INTRODUCTION

The epidermal growth factor receptor (EGFR) belongs to the family of receptor tyrosine kinases (RTKs), which regulate cell growth, survival, proliferation and differentiation. EGFR, also known as ErbB1, is related to the other members of the EGFR family of RTKs including HER2/ErbB2/neu, HER3/ErbB3 and HER4/ErbB4. EGFR is a 170 kDa type I transmembrane glycoprotein which consists of an extracellular ligand-binding domain, a single hydrophobic transmembrane region, and an intracellular segment harboring the highly conserved, tyrosine kinase domain. Several deletions in the extra- and intracellular domain of the EGFR have been found in a number of tumors. For example, EGFRvIII is a 145 kDa protein with a deletion of exons 2-7 in EGFR mRNA.

Selective compounds have been developed that target either the extracellular ligand-binding domain of EGFR or the intracellular tyrosine kinase region, resulting in interference with the signaling pathways that modulate mitogenic and other cancer-promoting responses.

EGFR is activated by the binding of ligands such as EGF, transforming growth factor  $\alpha$  (TGF- $\alpha$ ), amphiregulin,  $\beta$ -cellulin, heparin binding EGF-like growth factor (HB-EGF) and epiregulin. The binding causes EGFR homo- and heterodimerization and rapid activation of its intrinsic tyrosine kinase followed by autophosphorylation of multiple tyrosine residues in the cytoplasmic domain including pY845, pY992, pY1068, pY1086, and pY1173. The phosphorylated tyrosine residues in the COOH-terminus serve as binding sites for cytosolic signaling proteins containing Src homology 2 (SH2) domains. Tyrosine 1068 within the cytoplasmic domain of the receptor is an autophosphorylation site that

allows binding of Grb2 and activation of the Ras-Raf-ERK1 and ERK2 signaling pathway.

The Invitrogen Human EGFR [pY1068] ELISA is designed to specifically detect and quantify the levels of phosphorylated EGFR at tyrosine 1068. It does not detect non-phosphorylated EGFR protein or the EGFR when phosphorylated at other tyrosine residues such as pY845, pY1086 and pY1173. This assay is intended to detect EGFR [pY1068] from lysates of cells and tissues. For normalization of the EGFR content of the samples, an ELISA kit for human full length EGFR (Cat. #KHR9061) is available in which detection is independent of phosphorylation status.

**For Research Use Only. CAUTION: Not for human or animal therapeutic or diagnostic use.**

**READ ENTIRE PROTOCOL BEFORE USE**

## PRINCIPLE OF THE METHOD

The Invitrogen Human EGFR [pY1068] kit is a solid phase sandwich Enzyme Linked-Immuno-Sorbent Assay (ELISA). A monoclonal antibody specific for human EGFR (regardless of phosphorylation state) has been coated onto the wells of the microtiter strips provided. Samples, including a standard containing human EGFR [pY1068], control specimens, and unknowns, are pipetted into these wells. During the first incubation, the human EGFR antigen binds to the immobilized (capture) antibody. After washing, an antibody specific for EGFR phosphorylated at tyrosine 1068 is added to the wells. During the second incubation, this antibody serves as a detection antibody by binding to the immobilized EGFR protein captured during the first incubation. After removal of excess detection antibody, a horseradish peroxidase-labeled Anti-Rabbit IgG (Anti-Rabbit IgG HRP) is added. This binds to the detection antibody to complete the four-member sandwich. After a third incubation and washing to remove all the excess Anti-Rabbit IgG HRP, a substrate solution is added, which is acted upon by the bound enzyme to produce color. The intensity of this colored product is directly proportional to the concentration of human EGFR [pY1068] present in the original specimen.

## REAGENTS PROVIDED

**Note:** Store all reagents at 2 to 8°C.

<i>Reagent</i>	<i>96 Test Kit</i>
<i>Hu EGFR [pY1068] Standard.</i> Contains 0.1% sodium azide. Refer to vial label for quantity and reconstitution volume.	2 vials
<i>Standard Diluent Buffer.</i> Contains 0.1% sodium azide; 25 mL per bottle.	1 bottle
<i>Hu EGFR Antibody Coated Wells, 12 x 8 Well Strips.</i>	1 plate
<i>Hu EGFR [pY1068] Detection Antibody.</i> Contains 0.1% sodium azide; 11 mL per bottle.	1 bottle
<i>Anti-Rabbit IgG HRP (100X).</i> Contains 3.3 mM thymol; 0.125 mL per vial.	1 vial
<i>HRP Diluent.</i> Contains 3.3 mM thymol; 25 mL per bottle.	1 bottle
<i>Wash Buffer Concentrate (25X);</i> 100 mL per bottle.	1 bottle
<i>Stabilized Chromogen, Tetramethylbenzidine (TMB);</i> 25 mL per bottle.	1 bottle
<i>Stop Solution;</i> 25 mL per bottle.	1 bottle
<i>Plate Covers, adhesive strips.</i>	3

**Disposal Note:** This kit contains materials with small quantities of sodium azide. Sodium azide reacts with lead and copper plumbing to form explosive metal azides. Upon disposal, flush drains with a large volume of water to prevent azide accumulation. Avoid ingestion and contact with eyes, skin and mucous membranes. In case of contact, rinse affected area with plenty of water. Observe all federal, state and local regulations for disposal.

### **SUPPLIES REQUIRED BUT NOT PROVIDED**

1. Microtiter plate reader capable of measurement at or near 450 nm.
2. Calibrated adjustable precision pipettes, preferably with disposable plastic tips. (A manifold multi-channel pipette is desirable for large assays.)
3. Cell extraction buffer (see Recommended Formulation, p. 11).
4. Distilled or deionized water.
5. Plate washer: automated or manual (squirt bottle, manifold dispenser, etc.).
6. Data analysis and graphing software.
7. Glass or plastic tubes for diluting and aliquoting standard.
8. Absorbent paper towels.
9. Calibrated beakers and graduated cylinders in various sizes.

### **PROCEDURAL NOTES/LAB QUALITY CONTROL**

1. When not in use, kit components should be refrigerated. All reagents should be warmed to room temperature before use.
2. **Microtiter plates should be allowed to come to room temperature before opening the foil bags.** Once the desired number of strips has been removed, immediately reseal the bag and store at 2 to 8°C to maintain plate integrity.



3. Samples should be frozen if not analyzed shortly after collection. Avoid multiple freeze-thaw cycles of frozen samples. Thaw completely and mix well prior to analysis.
4. If particulate matter is present, centrifuge or filter prior to analysis.
5. All standards, controls and samples should be run in duplicate.
6. Samples containing EGFR [pY1068] protein extracted from cells should be diluted at least 1:10 with *Standard Diluent Buffer*. This dilution is necessary to reduce the matrix effect of the cell lysis buffer.
7. When pipetting reagents, maintain a consistent order of addition from well-to-well. This ensures equal incubation times for all wells.
8. Cover or cap all reagents when not in use.
9. **Do not mix or interchange different reagent lots from various kit lots.**
10. Do not use reagents after the kit expiration date.
11. Read absorbances within 2 hours of assay completion.
12. In-house controls should be run with every assay. If control values fall outside pre-established ranges, the accuracy of the assay is suspect.
13. All residual wash liquid must be drained from the wells by efficient aspiration or by decantation followed by tapping the plate forcefully on absorbent paper. *Never* insert absorbent paper directly into the wells.
14. Because *Stabilized Chromogen* is light sensitive, avoid prolonged exposure to light. Also avoid contact between *Stabilized Chromogen* and metal, or color may develop.

## **SAFETY**

All blood components and biological materials should be handled as potentially hazardous. Follow universal precautions as established by the Centers for Disease Control and Prevention and by the Occupational Safety and Health Administration when handling and disposing of infectious agents.

## **DIRECTIONS FOR WASHING**

**Incomplete washing will adversely affect the test outcome.** All washing must be performed with *Wash Buffer Concentrate (25X)* provided.

Washing can be performed manually as follows: completely aspirate the liquid from all wells by gently lowering an aspiration tip (aspiration device) into the bottom of each well. Take care not to scratch the inside of the well.

After aspiration, fill the wells with at least 0.4 mL of diluted wash solution. Let soak for 15 to 30 seconds, then aspirate the liquid. Repeat as directed under **ASSAY METHOD**. After the washing procedure, the plate is inverted and tapped dry on absorbent tissue.

Alternatively, the wash solution may be put into a squirt bottle. If a squirt bottle is used, flood the plate with wash buffer, completely filling all wells. After the washing procedure, the plate is inverted and tapped dry on absorbent tissue.

If using an automated washer, the operating instructions for washing equipment should be carefully followed. If your automated washer allows, 30 second soak cycles should be programmed into the wash cycle.

## PROCEDURE FOR EXTRACTION OF PROTEINS FROM CELLS

### A. Recommended Formulation of Cell Extraction Buffer:

10 mM Tris, pH 7.4

100 mM NaCl

1 mM EDTA

1 mM EGTA

1 mM NaF

20 mM  $\text{Na}_4\text{P}_2\text{O}_7$

2 mM  $\text{Na}_3\text{VO}_4$

1% Triton X-100

10% glycerol

0.1% SDS

0.5% deoxycholate

1 mM PMSF (stock is 0.3 M in DMSO)

Protease inhibitor cocktail (e.g., Sigma Cat. # P-2714) (reconstituted according to manufacturer's guideline). Add 250  $\mu\text{L}$  per 5 mL Cell Extraction Buffer.

This buffer is stable for 2 to 3 weeks at 4°C or for up to 6 months when aliquoted (without protease inhibitors and PMSF added) and stored at -20°C. When stored frozen, the Cell Extraction Buffer should be thawed on ice. Important: add the protease inhibitors just before using. The stability of protease inhibitor supplemented Cell Extraction Buffer is 24 hours at 4°C. PMSF is very unstable and must be added prior to use, even if added previously.

## B. Protocol for Cell Extraction

This protocol has been successfully applied to several cell lines. Researchers should optimize the cell extraction procedures for their own applications.

1. Collect cells in PBS by centrifugation (non-adherent) or scraping from culture flasks (adherent).
2. Wash cells twice with cold PBS.
3. Remove and discard the supernatant and collect the cell pellet. (At this point the cell pellet can be frozen at  $-80^{\circ}\text{C}$  and lysed at a later date).
4. Lyse the cell pellet in Cell Extraction Buffer for 30 minutes, on ice, with vortexing at 10 minute intervals. The volume of Cell Extraction Buffer depends on the cell number in cell pellet and expression of EGFR [pY1068]. For example,  $10^8$  A431 cells grown in RPMI plus 10% FBS and treated with 20 ng/mL of EGF can be extracted in 1 mL of Extraction Buffer. Under these conditions, use of 0.1-1  $\mu\text{L}$  of the clarified cell extract diluted to a volume of 100  $\mu\text{L}$ /well in Standard Diluent Buffer (See **Assay Method**) is sufficient for the detection of EGFR [pY1068].
5. Transfer extract to microcentrifuge tubes and centrifuge at 13,000 rpm for 10 minutes at  $4^{\circ}\text{C}$ .
6. Aliquot the clear lysate to clean microfuge tubes. These samples are ready for assay. Lysates can be stored at  $-80^{\circ}\text{C}$ . Avoid multiple freeze-thaws.

## REAGENT PREPARATION AND STORAGE

### A. Reconstitution and Dilution of Human EGFR [pY1068] Standard

**Note:** This *Hu EGFR [pY1068] Standard* was prepared from full length EGFR purified from A431 cells. The purified protein was allowed to autophosphorylate in the presence of 2 mM ATP and 1x Autophosphorylation Buffer (15 mM HEPES, 6 mM MnCl<sub>2</sub>, 15 mM MgCl<sub>2</sub>). 1 unit of standard is equivalent to the amount of EGFR [pY1068] derived from 42 pg of EGFR allowed to autophosphorylate. Subsequent lots of standard will be normalized to this lot of material to allow consistency of EGFR [pY1068] quantitation.

1. Reconstitute *Hu EGFR [pY1068] Standard* with *Standard Diluent Buffer*. Refer to standard vial label for instructions. Swirl or mix gently and allow to sit for 10 minutes to ensure complete reconstitution. Label as 100 Units/mL human EGFR [pY1068]. Use standard within 1 hour of reconstitution.
2. Add 0.25 mL of *Standard Diluent Buffer* to each of 6 tubes labeled 50, 25, 12.5, 6.25, 3.12 and 1.6 Units/mL human EGFR [pY1068].
3. Make serial dilutions of the standard as described in the following dilution table. Mix thoroughly between steps.

**B. Dilution of Human EGFR [pY1068] Standard**

<b>Standard:</b>	<b>Add:</b>	<b>Into:</b>
100 Units/mL	Prepare as described in step 1	
50 Units/mL	0.25 mL of the 100 Units/mL std.	0.25 mL of the Diluent Buffer
25 Units/mL	0.25 mL of the 50 Units/mL std.	0.25 mL of the Diluent Buffer
12.5 Units/mL	0.25 mL of the 25 Units/mL std.	0.25 mL of the Diluent Buffer
6.25 Units/mL	0.25 mL of the 12.5 Units/mL std.	0.25 mL of the Diluent Buffer
3.12 Units/mL	0.25 mL of the 6.25 Units/mL std.	0.25 mL of the Diluent Buffer
1.6 Units/mL	0.25 mL of the 3.12 Units/mL std.	0.25 mL of the Diluent Buffer
0 Units/mL	0.25 mL of the Diluent Buffer	An empty tube

Discard all remaining reconstituted and diluted standards after completing assay. Return the *Standard Diluent Buffer* to the refrigerator.

### C. Storage and Final Dilution of Anti-Rabbit IgG HRP (100X)

**Please Note:** The *Anti-Rabbit IgG HRP (100X)* is in 50% glycerol. This solution is viscous. To ensure accurate dilution, allow *Anti-Rabbit IgG HRP (100X)* to reach room temperature. Gently mix. Pipette *Anti-Rabbit IgG HRP (100X)* slowly. Remove excess concentrate solution from pipette tip by gently wiping with clean absorbent paper.

1. Dilute 10  $\mu\text{L}$  of this 100x concentrated solution with 1 mL of *HRP Diluent* for each 8-well strip used in the assay. Label as Anti-Rabbit IgG HRP Working Solution.

For Example:

<u># of 8-Well Strips</u>	<u>Volume of Anti-Rabbit IgG HRP (100X)</u>	<u>Volume of Diluent</u>
2	20 $\mu\text{L}$ solution	2 mL
4	40 $\mu\text{L}$ solution	4 mL
6	60 $\mu\text{L}$ solution	6 mL
8	80 $\mu\text{L}$ solution	8 mL
10	100 $\mu\text{L}$ solution	10 mL
12	120 $\mu\text{L}$ solution	12 mL

2. Return the unused *Anti-Rabbit IgG HRP (100X)* to the refrigerator.

#### D. Dilution of Wash Buffer

Allow the *Wash Buffer Concentrate (25X)* to reach room temperature and mix to ensure that any precipitated salts have redissolved. Dilute 1 volume of the *Wash Buffer Concentrate (25X)* with 24 volumes of deionized water (e.g., 50 mL may be diluted up to 1.25 liters, 100 mL may be diluted up to 2.5 liters). Label as Working Wash Buffer.

Store both the concentrate and the Working Wash Buffer in the refrigerator. The diluted buffer should be used within 14 days.

#### ASSAY METHOD: PROCEDURE AND CALCULATIONS

**Be sure to read the *Procedural Notes/Lab Quality Control* section before carrying out the assay.**

Allow all reagents to reach room temperature before use. Gently mix all liquid reagents prior to use.

**Note:** A standard curve must be run with each assay.

1. Determine the number of 8-well strips needed for the assay. Insert these in the frame(s) for current use. (Re-bag extra strips and frame. Store these in the refrigerator for future use.)
2. Add 100  $\mu\text{L}$  of the *Standard Diluent Buffer* to zero wells. Well(s) reserved for chromogen blank should be left empty.
3. Add 100  $\mu\text{L}$  of standards to the appropriate microtiter wells. Samples prepared in Cell Extraction Buffer must be diluted 1:10 or greater in *Standard Diluent Buffer* (for example, 10  $\mu\text{L}$  sample into 90  $\mu\text{L}$  buffer). While a 1:10 sample dilution has been found to be satisfactory, higher dilutions such as 1:25 or 1:50 are found to be optimal for sample resolution. The dilution chosen should be



optimized for each investigation's system. Tap gently on side of plate to thoroughly mix. (See **REAGENT PREPARATION AND STORAGE**, Section B.)

4. Cover wells with *plate cover* and incubate for **2 hours at room temperature**. Alternatively, the plate may be incubated overnight at 4°C.
5. Thoroughly aspirate or decant solution from wells and discard the liquid. Wash wells 4 times. See **DIRECTIONS FOR WASHING**.
6. Pipette 100 µL of *Hu EGFR [pY1068] Detection Antibody* solution into each well except the chromogen blank(s). Tap gently on the side of the plate to mix.
7. Cover wells with *plate cover* and incubate for **1 hour at room temperature**.
8. Thoroughly aspirate or decant solution from wells and discard the liquid. Wash wells 4 times. See **DIRECTIONS FOR WASHING**.
9. Add 100 µL *Anti-Rabbit IgG HRP Working Solution* to each well except the chromogen blank(s). (Prepare the working dilution as described in **REAGENT PREPARATION AND STORAGE**, Section C.)
10. Cover wells with the *plate cover* and incubate for **30 minutes at room temperature**.
11. Thoroughly aspirate or decant solution from wells and discard the liquid. Wash wells 4 times. See **DIRECTIONS FOR WASHING**.
12. Add 100 µL of *Stabilized Chromogen* to each well. The liquid in the wells will begin to turn blue.

13. Incubate for **30 minutes at room temperature and in the dark.** **Please Note: Do not cover the plate with aluminum foil or metalized mylar.** The incubation time for chromogen substrate is often determined by the microtiter plate reader used. Many plate readers have the capacity to record a maximum optical density (O.D.) of 2.0. The O.D. values should be monitored and the substrate reaction stopped before the O.D. of the positive wells exceed the limits of the instrument. The O.D. values at 450 nm can only be read after the *Stop Solution* has been added to each well. If using a reader that records only to 2.0 O.D., stopping the assay after 20 to 25 minutes is suggested.
14. Add 100  $\mu\text{L}$  of *Stop Solution* to each well. Tap side of plate gently to mix. The solution in the wells should change from blue to yellow.
15. Read the absorbance of each well at 450 nm having blanked the plate reader against a chromogen blank composed of 100  $\mu\text{L}$  each of *Stabilized Chromogen* and *Stop Solution*. Read the plate within 2 hours after adding the *Stop Solution*.
16. Plot on graph paper the absorbance of the standards against the standard concentration. (Optimally, the background absorbance may be subtracted from *all* data points, including standards, unknowns and controls, prior to plotting.) Draw the best smooth curve through these points to construct the standard curve. If using curve fitting software, the four parameter algorithm provides the best curve fit.
17. Read the human EGFR [pY1068] concentrations for unknown samples and controls from the standard curve plotted in step 16. **Multiply value(s) obtained for sample(s) by dilution factor to correct for the dilution in step 3.** (Samples still producing signals

higher than the highest standard (100 Units/mL) should be further diluted in *Standard Diluent Buffer* and reanalyzed, multiplying the concentration found by the appropriate dilution factor.)

18. Values of EGFR [pY1068] should be normalized for full length EGFR content by parallel measurement with the Invitrogen EGFR Full Length ELISA Kit (Cat. # KHR9061).

### **TYPICAL DATA**

The following data were obtained for the various standards over the range of 0 to 100 Units/mL human EGFR [pY1068].

<u>EGFR [pY1068] Standard (Units/mL)</u>	<u>Optical Density (450 nm)</u>
0	0.208
1.6	0.322
3.12	0.394
6.25	0.526
12.5	0.758
25	1.171
50	1.939
100	3.160

## LIMITATIONS OF THE PROCEDURE

Do not extrapolate the standard curve beyond the 100 Units/mL standard point; the dose-response is non-linear in this region and accuracy is difficult to obtain. Dilute samples >100 Units/mL with *Standard Diluent Buffer*; reanalyze these and multiply results by the appropriate dilution factor.

The influence of various extraction buffers has not been thoroughly investigated. The rate of degradation of native EGFR [pY1068] in various matrices has not been investigated. Although EGFR [pY1068] degradation in the Cell Extraction Buffer described in this protocol has not been seen to date, the possibility of this occurrence cannot be excluded.

**For Research Use Only. CAUTION: Not for human or animal therapeutic or diagnostic use.**

## PERFORMANCE CHARACTERISTICS

### SENSITIVITY

The analytical sensitivity of this assay is <0.3 Units/mL of human EGFR [pY1068]. This was determined by adding two standard deviations to the mean O.D. obtained when the zero standard was assayed 30 times.

The sensitivity of this ELISA was compared to Western blotting using known quantities of EGFR [pY1068]. The data presented in Figure 1 show that the sensitivity of the ELISA is at least 8x greater than that of Western blotting. The bands shown in the Western blot data were developed using rabbit anti-EGFR [pY1068], an alkaline phosphatase conjugated anti-rabbit IgG followed by chemiluminescent substrate and autoradiography.

**Figure 1**

Western  
blotting  
(170 kDa)



ELISA:  
OD 450 nm

0.208	0.322	0.394	0.526	0.758	1.171	1.939	3.160
-------	-------	-------	-------	-------	-------	-------	-------

EGFR  
units/test

0	0.15	0.31	0.63	1.25	2.5	5	10
---	------	------	------	------	-----	---	----

## PRECISION

### 1. Intra-Assay Precision

Samples of known EGFR [pY1068] concentration were assayed in replicates of 16 to determine precision within an assay.

	Sample 1	Sample 2	Sample 3
Mean (Units/mL)	39.9	12.10	2.58
SD	1.40	0.54	0.20
%CV	3.51	4.49	7.58

SD = Standard Deviation

CV = Coefficient of Variation

### 2. Inter-Assay Precision

Samples were assayed 48 times in multiple assays to determine precision between assays.

	Sample 1	Sample 2	Sample 3
Mean (Units/mL)	40.55	11.95	2.49
SD	1.74	0.68	0.23
%CV	4.30	5.68	9.27

SD = Standard Deviation


CV = Coefficient of Variation

## RECOVERY

The recovery of EGFR [pY1068] added to Cell Extraction Buffer followed by 1:10 dilution in *Standard Diluent Buffer*, averaged 103%.

## PARALLELISM

Natural EGFR [pY1068] from EGF treated-A431 cell lysate was serially diluted in *Standard Diluent Buffer*. The optical density of each dilution was plotted against the EGFR [pY1068] standard curve. Parallelism was demonstrated by the figure below and indicated that the Standard accurately reflects full length EGFR [pY1068] content in samples.



## LINEARITY OF DILUTION

Extract Buffer was spiked with EGFR [pY1068] and serially diluted in *Standard Diluent Buffer* over the range of the assay. Linear regression analysis of samples versus the expected concentration yielded a correlation coefficient of 0.99.


<b>Dilution</b>	<b>Cell Lysate</b>		
	<b>Measured (Units/mL)</b>	<b>Expected (Units/mL)</b>	<b>% Expected</b>
Neat	61.68	-	-
1/2	31.11	30.8	101%
1/4	15.16	15.4	98%
1/8	8.06	7.7	104%
1/16	3.92	3.9	102%



## SPECIFICITY

This EGFR [pY1068] ELISA kit is specific for measurement of phosphorylated human EGFR at tyrosine 1068. The kit does not detect non-phosphorylated EGFR protein, as shown in Figure 2.

**Figure 2**



The specificity of this assay for EGFR phosphorylated at tyrosine 1068 was confirmed by peptide competition. Phosphorylated EGFR was serially diluted and quantitated in the assay as usual except that the detection antibody was preincubated with EGFR-derived peptides at a concentration of 1  $\mu\text{g}/\text{mL}$ . The data presented in Figure 3 show that only the peptide corresponding to the region

surrounding tyrosine 1068, containing the phospho-tyrosine, could block the ELISA signal. The same sequence containing non-phosphorylated tyrosine at position 1068 did not block the signal.

**Figure 3**





Figure 4 shows that EGFR [pY1068] phosphorylation in HeLa cells is dependent on levels of EGF stimulation. HeLa cells (approximately 50% confluent) were treated with EGF at varying concentrations (0-500 ng/mL) for 10 minutes, lysed and quantitated in parallel for EGFR content (both full length and pY1068). The amount of full length EGFR remains constant, while the level of phosphorylation at tyrosine 1068 decreases with diminishing EGF dosage.

**Figure 4**
















## REFERENCES

1. Prenzel, N., et al. (2001) The epidermal growth factor receptor family as a central element for cellular signal transduction and diversification. *Endocr. Relat. Cancer* 8:11-31.
2. Raymond, E., et al. (2000) Epidermal growth factor receptor tyrosine kinase as a target for anticancer therapy. *Drugs* 60:15-23.
3. Tang, C.K., et al. (2000) Epidermal growth factor receptor vIII enhances tumorigenicity in human breast cancer. *Cancer Res.* 60:3081–3087.
4. Biscardi, J.S., et al. (1999) c-Src-mediated phosphorylation of the epidermal growth factor receptor on Tyr845 and Tyr1101 is associated with modulation of receptor function. *J. Biol. Chem.* 274(12):8335-8343.
5. Schooler, K., and Wiley, H.S. (2000) Ratiometric assay of epidermal growth factor receptor tyrosine kinase activation. *Anal. Biochem.* 277(1):135-142.
6. Sibilio, M., et al. (2000) The EGF receptor provides an essential survival signal for SOS-dependent skin tumor development. *Cell* 102:211-220.
7. Keilhack, H., et al. (1998) Phosphotyrosine 1173 mediates binding of the protein-tyrosine phosphatase SHP-1 to the epidermal growth factor receptor and attenuation of receptor signaling. *J. Biol. Chem.* 273:24839-24846.
8. Rodrigues, G.A., et al. (2000) A novel positive feedback loop mediated by the docking protein Gab1 and phosphatidylinositol 3-kinase in epidermal growth factor receptor signaling. *Mol. Cell. Biol.* 20:1448-1459.

9. Chen, K., et al. (2001) c-Jun N-terminal kinase activation by hydrogen peroxide in endothelial cells involves SRC-dependent epidermal growth factor receptor transactivation. *J. Biol. Chem.* 276:16045-16050.

**Important Licensing Information** - These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or our website, [www.invitrogen.com](http://www.invitrogen.com)). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

#### Explanation of symbols

Symbol	Description	Symbol	Description
	Catalogue Number		Batch code
	Research Use Only		<i>In vitro</i> diagnostic medical device
	Use by		Temperature limitation
	Manufacturer		European Community authorised representative
	Without, does not contain		With, contains
	Protect from light		Consult accompanying documents
	Directs the user to consult instructions for use (IFU), accompanying the product.		

Copyright © Invitrogen Corporation. 21 May 2010

## NOTES



## Human EGFR [pY1068] Assay Summary

Incubate 100  $\mu$ L Standard or Cell Extract (>1:10)  
for 2 hours at RT



**aspirate and wash 4x**

Incubate 100  $\mu$ L of Detection Antibody  
for 1 hour at RT



**aspirate and wash 4x**

Incubate 100  $\mu$ L of HRP Anti-Rabbit Antibody  
for 30 minutes at RT



**aspirate and wash 4x**

Incubate 100  $\mu$ L of Stabilized Chromogen  
for 30 minutes at RT



Add 100  $\mu$ L of Stop Solution and read at 450 nm  
**Total time: 4 hours**

