

Recommended universal buffers for double digestion



Double digestion (cutting DNA with two restriction enzymes simultaneously) is widely used to save time. In this table, the suitable buffer type, dilution rate, and additive are shown to perform double digestion with the restriction enzymes shown below. The letters L, M, H, T, and K show the type of universal buffer recommended, which is supplied at a 10-fold (10X) concentration. When "0.5X" is recommended, dilute 10X buffer 20-fold; when "1X" is recommended, dilute 10-fold; and when "2X" is recommended, dilute 5-fold for use in double digestion. As BSA is also supplied at 10X concentration, dilute 10-fold to a final concentration of 0.1% when it is used.

Enzyme	Accl	BamHI	BglII	Clal	EcoRI	EcoRV	HincII	HindIII	KpnI	NcoI	NdeI	NotI	PstI	PvuI	SacI	Sall	SmaI	SpeI	SphI	XbaI	XhoI
Supplied buffer	10X M	10X K	10X H	10X M	10X H	10X H	10X M	10X M	10X L	10X K +BSA	10X H	10X H +BSA +Triton	10X H	10X K +BSA	10X L	10X H	10X T +BSA	10X M	10X H	10X M +BSA	10X H
Accl	–	0.5X K	1X T	1X M	1X M	0.5X K	1X M	1X M	1X M	1X M +BSA	1X T	0.5X K +BSA	1X M	0.5X K	1X M	1.5X T	1X T +BSA	1X M	0.5X K	1X M	1X M
BamHI	0.5X K	–	1X K	1X K	1X K	1X K	0.5X K	1X K	0.5X K	1X K +BSA	1X K	0.5X K +BSA	1X K	1X K	0.5X K	1.5X T	0.5X T +BSA	1X K	1X K	0.5X K	1X M
BglII	1X T	1X K	–	1X H	1X H	1X H	2X K	1X K	1X T	1X K +BSA	1X H	1X H +BSA	1X H	1X K	0.5X K	1X H	1X T +BSA	1X H	1X H	2X T	1X H
Clal	1X M	1X K	1X H	–	1X H	1X H	1X M	1X M	1X M	1X K +BSA	1X H	1X H +BSA	1X H	1X K	1X M	1X H	1X T +BSA	1X M	1X H	1X H	1X H
EcoRI	1X M	1X K	1X H	1X H	–	1X H	1X M	1X M	1X M	1X K +BSA	1X H	1X H +BSA	1X H	1X K	1X M	1X H	1X T +BSA	1X H	1X H	1X M	1X H
EcoRV	0.5X K	1X K	1X H	1X H	1X H	–	2X T	1X K	0.5X K	1X K +BSA	1X H	1X H +BSA	1X H	1X K	0.5X K	1X H	0.5X K +BSA	1X H	1X H	2X T	1X H
HincII	1X M	0.5X K	2X K	1X M	1X M	2X T	–	1X M	1X M	1X M +BSA	1X T	0.5X K +BSA	1X M	0.5X K	1X M	1.5X K	1X T +BSA	1X M	2X T	1X M	1X M
HindIII	1X M	1X K	1X K	1X M	1X M	1X K	1X M	–	1X M	1X K +BSA	1X K	0.5X K +BSA	1X M	1X K	1X M	1.5X K	1X T +BSA	1X M	1X K	1X M	1X M
KpnI	1X M	0.5X K	1X T	1X M	1X M	0.5X K	1X M	1X M	–	0.5X K +BSA	1X T	0.5X K +BSA	1X M	0.5X K	1X L	1.5X T +BSA	1X T +BSA	1X M	0.5X K	1X M	1X M
NcoI	1X M +BSA	1X K +BSA	1X K +BSA	1X K +BSA	1X K +BSA	1X K +BSA	1X M +BSA	1X K +BSA	0.5X K +BSA	–	1X K +BSA	0.5X K +BSA	1X K +BSA	1X K +BSA	0.5X K +BSA	1.5X T +BSA	1X T +BSA	1X K +BSA	1X K +BSA	1X M +BSA	1X K +BSA
NdeI	1X T	1X K	1X H	1X H	1X H	1X H	1X T	1X K	1X T	1X K +BSA	–	1X H +BSA	1X H	1X K	1X T	1X T +BSA	1X H	1X H	1X T	1X H	1X H
NotI	0.5X K +BSA	0.5X K +BSA	1X H +BSA	1X H +BSA	1X H +BSA	1X H +BSA	0.5X K +BSA	0.5X K +BSA	0.5X K +BSA	0.5X K +BSA	1X H +BSA	–	1X H +BSA	2X K +BSA	0.5X K +BSA	1X H +BSA	0.5X T +BSA	1X H +BSA	1X H +BSA	0.5X K +BSA	1X H +BSA
PstI	1X M	1X K	1X H	1X H	1X H	1X H	1X M	1X M	1X M	1X K +BSA	1X H	1X H +BSA	–	1X K	1X M	1X H	0.5X T +BSA	1X H	1X H	1X M	1X H
PvuI	0.5X K	1X K	1X K	1X K	1X K	1X K	0.5X K	1X K	0.5X K	1X K +BSA	1X K	2X K +BSA	1X K	–	0.5X K	1.5X K +BSA	1X K +BSA	1X K	1X K	0.5X K	1X K
SacI	1X M	0.5X K	0.5X K	1X M	1X M	0.5X K	1X M	1X M	1X L	0.5X K +BSA	1X T	0.5X K +BSA	1X M	0.5X K	–	1.5X T +BSA	1X T +BSA	1X M	0.5X K	1X M	1X M
Sall	1.5X T	1.5X T	1X H	1X H	1X H	1X H	1.5X K	1.5X K	1.5X T +BSA	1.5X T +BSA	1X H	1X H +BSA	1X H	1.5X K +BSA	1.5X T +BSA	–	1.5X T +BSA	1X H	1X H	1.5X T	1X H
SmaI	1X T +BSA	0.5X T +BSA	1X T +BSA	1X T +BSA	1X T +BSA	0.5X K +BSA	1X T +BSA	1X T +BSA	1X T +BSA	1X T +BSA	1X T +BSA	0.5X T +BSA	0.5X T +BSA	1X K +BSA	1X T +BSA	1.5X T +BSA	–	1X T +BSA	0.5X T +BSA	1X T +BSA	1X T +BSA
SpeI	1X M	1X K	1X H	1X M	1X H	1X H	1X M	1X M	1X M	1X K +BSA	1X H	1X H +BSA	1X H	1X K	1X M	1X H	1X T +BSA	–	1X H	1X M	1X H
SphI	0.5X K	1X K	1X H	1X H	1X H	1X H	2X T	1X K	0.5X K	1X K +BSA	1X H	1X H +BSA	1X H	1X K	0.5X K	1X H	0.5X T +BSA	1X H	–	2X T	1X H
XbaI	1X M	0.5X K	2X T	1X M	1X M	2X T	1X M	1X M	1X M	1X M +BSA	1X T	0.5X K +BSA	1X M	0.5X K	1X M	1.5X T	1X T +BSA	1X M	2X T	–	1X M
XhoI	1X M	1X K	1X H	1X H	1X H	1X H	1X M	1X M	1X M	1X K +BSA	1X H	1X H +BSA	1X H	1X K	1X M	1X H	1X T +BSA	1X H	1X H	1X M	–

Note:

1. It is confirmed that 10 units of each enzyme completely digests 1 µg of DNA at 37°C in one hour in 50 µL reaction mixture.
2. The concentration of glycerol should be less than 10% to minimize star activity.
3. DNA may not be digested completely when recognition sequences of two enzymes are close to each other or when DNA takes high-structure conformation.



www.invitrogen.com

For research use only. Not intended for any animal or human therapeutic or diagnostic use, unless otherwise stated.

© 2010 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners. These products may be covered by one or more Limited Use Label Licenses (see Invitrogen catalog or www.invitrogen.com). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. CO13950 0710

Relative activity for each reaction buffer



Restriction enzyme	Cut site	Relative activity (%)					
		Buffer L	Buffer M	Buffer H	Buffer K	Buffer T (+BSA)	Buffer
AccI	5'-GT↓(Py)(Pu)AC-3' 3'-CA(Pu)(Py)TTG-5'	20	100	<20	<20*	160	-
AfaI (Rsal)	5'-GT↓AC3' 3'-CATTG-5'	60	60	40	60	100	-
AluI	5'-AG↓CT-3' 3'-TCTAG-5'	100	100	<20	40	200	-
ApaI	5'-GGGCC↓C-3' 3'-CTCGGG-5'	100	<20	<20	<20	<20	-
BalI	5'-TGG↓CCA-3' 3'-ACCTGGT-5'	20*	20*	<20*	<20*	40*	100 (Ball buffer)
BamHI	5'-G↓GATC C-3' 3'-C CTAG↑G-5'	<20*	<20	40	100	<20*	-
BglII	5'-A↓GATC T-3' 3'-T CTAG↑A-5'	<20	20	100	100*	60*	-
BmeT120I (AvaI)	5'-C↓(Py)C G(Pu) G-3' 3'-G (Pu)G C(Py)TC-5'	<20	<20	20	100	<20	-
BshTI (AgeI)	5'-A↓AGCTT-3' 3'-TT CGATA-5'	<20	20	80*	50	50	-
Clal	5'-AT↓CGAT-3' 3'-TAGCTTA-5'	40	100	120	100	60	-
DdeI	5'-C↓ITNA G-3' 3'-G ANTC-5'	60	80	100	100	80	-
DpnI	5'-GmA↓ TC-3' 3'-C T↑mAG-5'	60	60	120	140	100*	-
DraI (AhaIII)	5'-TTT↓AAA-3' 3'-AAA↑TTT-5'	100	100	60	100	80	-
EcoRI	5'-G↓AATT C-3' 3'-C TTAAT↑G-5'	20*	100*	100	120*	80*	-
EcoRV	5'-GAT↓ATC-3' 3'-CTATTAG-5'	<20*	40*	100	120*	40*	-
EcoT22I (Ava II)	5'-ATGCA↓T-3' 3'-TTACGTA-5'	<20	20	100	140*	20*	-
HaeIII	5'-GG↓CC-3' 3'-CC↑GG-5'	60	100	100	60	100	-
HapII (HpaII)	5'-C↓ICG G-3' 3'-G GC↑G-5'	100	60	<20	<20	100	-
HhaI	5'-GCG↓C-3' 3'-CTGCG-5'	80	100	100	120	120	-
HincII (HindII)	5'-GT↓(Py)↓(Pu)AC-3' 3'-CA(Pu)↑(Py)TG-5'	20	100	20	40	100	-
HindIII	5'-A↓AGCT T-3' 3'-T TCGA↑A-5'	60*	100	<20	200	100*	-
Hinfl	5'-G↓AANT C-3' 3'-C T↑TG-5'	80	100	100	160	60	-
HpaI	5'-GTT↓AAC-3' 3'-CAAT↑TG-5'	<20	40*	20	100	80*	-
KpnI	5'-GGT↓AC↓C-3' 3'-CTCATGG-5'	100	60	<20	<20	100*	-
MboI (Sau3AI)	5'-↓GATC-3' 3'-CTAG↑T-5'	20	40	60	100	40	-
MluI	5'-A↓ICGCGT-3' 3'-TGCCTA-5'	60	60	100	100*	60	-
MseI	5'-T↓TA A-3' 3'-A AT↑T-5'	40	40	20	40	100	-
NcoI	5'-C↓CATGG-3' 3'-GGTACT↑G-5'	40*	60*	20	60*	60*	-
NdeI	5'-CA↓TATG-3' 3'-GTATTAC↑-5'	<20	40	100	10	80	-
NheI	5'-G↓CTAGC-3' 3'-CGATC↑G-5'	120*	100	<20	<20	160*	-
NotI	5'-G↓C↓GGCCGC-3' 3'-CGCCGG↑TCG-5'	<20*	<20*	20**	<20	<20*	-
NruI	5'-TCG↓CGA-3' 3'-AGCT↑GCT-5'	0*	<20*	20*	20*	<20*	100 (NruI buffer)
PstI	5'-C TGA↓G-3' 3'-G↑TACGT C-5'	<20*	60*	100	80	20*	-
PvuI	5'-CGA↓TCG-3' 3'-GCTTAGC↑-5'	<20*	20*	40*	80**	40*	-
PvuII	5'-CAG↓CTG-3' 3'-GTC↑GAC-5'	80*	100	40	<20	40*	-
SacI	5'-GAGCT↓C-3' 3'-CTCGAG↑-5'	100	60	<20	<20	80	-
SacII	5'-CCG↓CGG-3' 3'-GG↑CGCC-5'	40	20	<20	<20	100	-
Sall	5'-G↓ITCGA C-3' 3'-C AGCT↑TG-5'	<20	<20	100	20*	<20	-
Scal	5'-AGT↓ACT-3' 3'-TCA↑TGA-5'	<20*	<20*	100	60*	<20*	-
SmaI	5'-CC↓GGG-3' 3'-GGG↑CCC-5'	<20	<20	<20	<20	100	-
SpeI	5'-A↓LCTAGT-3' 3'-TGTACTA↑-5'	80*	100	80	100	80*	-
SphI	5'-GCA↓TGC-3' 3'-CTGTACG↑-5'	20*	40*	100	120	20*	-
SspI	5'-AAT↓ATT-3' 3'-TAT↑TAA-5'	<20*	60*	40	100*	80*	100 (SSpI buffer)
StuI	5'-AGG↓CCT-3' 3'-TCC↑TGA-5'	60	100	60	80	140	-
TaqI (TthHB8I)	5'-T↓JCGA-3' 3'-AGCT↑T-5'	40*	80*	60*	60*	80*	100 (TaqI buffer)
XbaI	5'-T↓CTAGA-3' 3'-AGATC↑T-5'	<20	80*	20	<20	120	-
XhoI (PaeR7I)	5'-C↓TTCGAG-3' 3'-GAGCT↑TC-5'	<20	60	100	160	60	-

*Weak star activity is detected.

**100% activity is obtained by addition of 0.01% BSA.

For more information, visit www.invitrogen.com/restrictionenzymes

Buffer compositions

10X Buffer H

500 mM Tris-HCl, pH 7.5

100 mM MgCl₂

10 mM dithiothreitol (DTT)

1000 mM NaCl

10X Buffer K

200 mM Tris-HCl, pH 8.5

100 mM MgCl₂

10 mM dithiothreitol (DTT)

1000 mM KCl

10X Buffer L

100 mM Tris-HCl, pH 7.5

100 mM MgCl₂

10 mM dithiothreitol (DTT)

10X Buffer M

100 mM Tris-HCl, pH 7.5

100 mM MgCl₂

10 mM dithiothreitol (DTT)

500 mM NaCl

10X Buffer T (BSA-free)

330 mM Tris-acetate, pH 7.9

100 mM magnesium acetate

5 mM dithiothreitol (DTT)

660 mM potassium acetate