



43-4311 2 ml
FITC-Streptavidin

FITC-STREPTAVIDIN

BIOTIN/AVIDIN REAGENTS

FORM:

FITC-streptavidin is provided in liquid form in 10 mM phosphate-buffered saline, pH 7.4, containing 1% bovine serum albumin, 40% glycerol and 0.1% sodium azide. Streptavidin (from *Streptomyces avidinii*) is conjugated with FITC under optimum condition. The conjugate is purified by molecular sieving and is completely free of unbound FITC.

APPLICATIONS:

FITC-streptavidin can be used for the detection or localization of biotinylated biomolecules. Working dilution for specific application should be determined by the investigator to obtain the best conditions.

STORAGE:

Store at 2-8°C. Avoid exposure to light.

For Research Use Only. CAUTION: Not for human or animal therapeutic or diagnostic use.

www.invitrogen.com

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288

E-mail: techsupport@invitrogen.com

PI434311

(Rev 10/09) DCC-09-1531

Important Licensing Information - These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or our website, www.invitrogen.com). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.



43-4311 2 ml
FITC-Streptavidin

FITC-STREPTAVIDIN

BIOTIN/AVIDIN REAGENTS

FORM:

FITC-streptavidin is provided in liquid form in 10 mM phosphate-buffered saline, pH 7.4, containing 1% bovine serum albumin, 40% glycerol and 0.1% sodium azide. Streptavidin (from *Streptomyces avidinii*) is conjugated with FITC under optimum condition. The conjugate is purified by molecular sieving and is completely free of unbound FITC.

APPLICATIONS:

FITC-streptavidin can be used for the detection or localization of biotinylated biomolecules. Working dilution for specific application should be determined by the investigator to obtain the best conditions.

STORAGE:

Store at 2-8°C. Avoid exposure to light.

For Research Use Only. CAUTION: Not for human or animal therapeutic or diagnostic use.

www.invitrogen.com

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288

E-mail: techsupport@invitrogen.com

PI434311

(Rev 10/09) DCC-09-1531

Important Licensing Information - These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or our website, www.invitrogen.com). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.