

5' Pyrophosphohydrolase (RppH)

in Links

info@neb.com

www.neb.com



M0356S 🗰 🖼 🔤

200 units	5,000 U/ml	Lot: 0011208
RECOMBINANT	Store at -20°C	Exp: 8/14

Description: The bacterial RNA 5´ Pyrophosphohydrolase (RppH) removes pyrophosphate from the 5´ end of triphosphorylated RNA to leave a 5´ monophosphate RNA (1). The RppH protein was also known as NudH/YgdP which can split Ap_sA to ADP and ATP (2).

Source: An *E. coli* strain containing a clone of the *E. coli* RppH gene.

Supplied in: 200 mM NaCl, 0.1 mM EDTA, 20 mM Tris-HCl (pH 7.5), 1 mM DTT, 0.01% Triton X-100 and 50% olycerol.

Applications:

- Convert RNA transcript (ppp-RNA) to 5' monophosphate RNA (p-RNA) for 5' RNA ligation or for removal by XRN-1.
- RNA 5' end structure analysis.

Reagents Supplied with Enzyme: 10X NEBuffer 2

1X NEBuffer 2: 50 mM NaCl 10 mM Tris-HCl 10 mM MgCl₂ 1 mM dithiothreitol pH 7.9 @ 25°C

Unit Definition: One unit is the amount of enzyme that converts 1 μ g 300 mer RNA transcript into a XRN-1 digestible RNA in 30 minutes at 37°C.

Quality Control Assays

RNase Assay: Incubation of a 10 µl reaction containing 5 units of RppH with 40 ng of 300 mer RNA transcript for 4 hours at 37°C resulted in less than 10% degradation of RNA as determined by denaturing PAGE analysis.

Endonuclease Assay: Incubation of a 10 µl reaction containing 5 units of RppH with 300 ng of supercoiled plasmid for 4 hours at 37°C produced less than 10% nicked or linear molecules as determined by agarose gel electrophoresis.

Exonuclease Assay: Incubation of a 50 µl reaction containing 25 units of RppH with 1 µg of a mixture of single and double-stranded [³H] *E. coli* DNA for 4 hours at 37°C released < 0.5% of the total radioactivity.

Phosphatase Activity: Incubation of a 100 μ l reaction containing 25 units of RppH with 4-nitrophenyl phosphate under standard conditions (NEBuffer 2) at 37°C for 4 hours, no phosphatase activity was detected.

Notes On Use:

Since RppH is active in the presence of $Mg^{\rm 2+},$ it can be incubated with RNA ligase or XRN-1.

References:

- 1. Deana, A. et al. (2008) *Nature*, 451, 355–358.
- 2. Bessman, M.J. et al (2001) JBC, 276, 37834.

CERTIFICATE OF ANALYSIS

-

BioLabs

1-800-632-7799 info@neb.com

www.neb.com

R R NEB 2

5 <i>′</i>	Pyrophosphohydrolase	
	(RppH)	

RNA



M0356S

200 units	5,000 U/ml	Lot: 0011208
RECOMBINANT	Store at -20°C	Exp: 8/14

Description: The bacterial RNA 5' Pyrophosphohydrolase (RppH) removes pyrophosphate from the 5' end of triphosphorylated RNA to leave a 5' monophosphate RNA (1). The RppH protein was also known as NudH/YgdP which can split Ap_sA to ADP and ATP (2).

Source: An *E. coli* strain containing a clone of the *E. coli* RppH gene.

Supplied in: 200 mM NaCl, 0.1 mM EDTA, 20 mM Tris-HCl (pH 7.5), 1 mM DTT, 0.01% Triton X-100 and 50% glycerol.

Applications:

- Convert RNA transcript (ppp-RNA) to 5' monophosphate RNA (p-RNA) for 5' RNA ligation or for removal by XRN-1.
- RNA 5' end structure analysis.

Reagents Supplied with Enzyme: 10X NEBuffer 2

1X NEBuffer 2:

50 mM NaCl 10 mM Tris-HCl 10 mM MgCl₂ 1 mM dithiothreitol pH 7.9 @ 25°C

Unit Definition: One unit is the amount of enzyme that converts 1 μ g 300 mer RNA transcript into a XRN-1 digestible RNA in 30 minutes at 37°C.

Quality Control Assays

RNase Assay: Incubation of a 10 µl reaction containing 5 units of RppH with 40 ng of 300 mer RNA transcript for 4 hours at 37°C resulted in less than 10% degradation of RNA as determined by denaturing PAGE analysis.

Endonuclease Assay: Incubation of a 10 μ l reaction containing 5 units of RppH with 300 ng of supercoiled plasmid for 4 hours at 37°C produced less than 10% nicked or linear molecules as determined by agarose gel electrophoresis.

Exonuclease Assay: Incubation of a 50 μ l reaction containing 25 units of RppH with 1 μ g of a mixture of single and double-stranded [³H] *E. coli* DNA for 4 hours at 37°C released < 0.5% of the total radioactivity.

Phosphatase Activity: Incubation of a 100 μ l reaction containing 25 units of RppH with 4-nitrophenyl phosphate under standard conditions (NEBuffer 2) at 37°C for 4 hours, no phosphatase activity was detected.

Notes On Use: Since RppH is active in the presence of Mg²⁺, it can be incubated with RNA ligase or XRN-1.

References:

- 1. Deana, A. et al. (2008) Nature, 451, 355-358.
- 2. Bessman, M.J. et al (2001) JBC, 276, 37834.