

SalI



1-800-632-7799
info@neb.com
www.neb.com



R0138S 054121014101

R0138S



2,000 units **20,000 U/ml** **Lot: 0541210**
RECOMBINANT **Store at -20°C** **Exp: 10/14**

Recognition Site:

5'...GTCGAC...3'
3'...CAGCTG...5'

Source: An *E. coli* strain that carries the cloned SalI gene from *Streptomyces albus* G (ATCC 49789)

**Also Available In
High Fidelity (HF™) Format**

SalI



1-800-632-7799
info@neb.com
www.neb.com



R0138S 054121014101

R0138S



2,000 units **20,000 U/ml** **Lot: 0541210**
RECOMBINANT **Store at -20°C** **Exp: 10/14**

Recognition Site:

5'...GTCGAC...3'
3'...CAGCTG...5'

Source: An *E. coli* strain that carries the cloned SalI gene from *Streptomyces albus* G (ATCC 49789)

Supplied in: 50 mM KCl, 10 mM Tris-HCl (pH 7.5), 0.1 mM EDTA, 1 mM dithiothreitol, 300 µg/ml BSA and 50% glycerol.

Reagents Supplied with Enzyme:
10X NEBuffer 3, 100X BSA.

Reaction Conditions: 1X NEBuffer 3, supplemented with 100 µg/ml BSA. Incubate at 37°C.

1X NEBuffer 3:
100 mM NaCl
50 mM Tris-HCl
10 mM MgCl₂
1 mM dithiothreitol
pH 7.9 @ 25°C

Unit Definition: One unit is defined as the amount of enzyme required to digest 1 µg of λ DNA (HindIII digest) in 1 hour at 37°C in a total reaction volume of 50 µl.

Diluent Compatibility: Diluent Buffer A
50 mM KCl, 10 mM Tris-HCl, 0.1 mM EDTA, 1 mM dithiothreitol, 200 µg/ml BSA and 50% glycerol (pH 7.4 @ 25°C).

Supplied in: 50 mM KCl, 10 mM Tris-HCl (pH 7.5), 0.1 mM EDTA, 1 mM dithiothreitol, 300 µg/ml BSA and 50% glycerol.

Reagents Supplied with Enzyme:
10X NEBuffer 3, 100X BSA.

Reaction Conditions: 1X NEBuffer 3, supplemented with 100 µg/ml BSA. Incubate at 37°C.

1X NEBuffer 3:
100 mM NaCl
50 mM Tris-HCl
10 mM MgCl₂
1 mM dithiothreitol
pH 7.9 @ 25°C

Unit Definition: One unit is defined as the amount of enzyme required to digest 1 µg of λ DNA (HindIII digest) in 1 hour at 37°C in a total reaction volume of 50 µl.

Diluent Compatibility: Diluent Buffer A
50 mM KCl, 10 mM Tris-HCl, 0.1 mM EDTA, 1 mM dithiothreitol, 200 µg/ml BSA and 50% glycerol (pH 7.4 @ 25°C).

Quality Control Assays

Ligation: After 50-fold overdigestion with SalI, > 95% of the DNA fragments can be ligated with T4 DNA Ligase (at a 5' termini concentration of 1–2 µM) at 16°C. Of these ligated fragments, > 95% can be recut.

16-Hour Incubation: A 50 µl reaction containing 1 µg of pBR322 DNA and 150 units of SalI incubated for 16 hours at 37°C resulted in a DNA pattern free of detectable nuclease degradation as determined by agarose gel electrophoresis.

Exonuclease Activity: Incubation of 200 units of enzyme with 1 µg sonicated ³H DNA (10⁵ cpm/µg) for 4 hours at 37°C in 50 µl reaction buffer released < 0.1% radioactivity.

Endonuclease Activity: Incubation of 40 units with 1 µg φX174 RF I DNA for 4 hours at 37°C in 50 µl reaction buffer resulted in < 20% conversion to RF II.

Quality Control Assays

Ligation: After 50-fold overdigestion with SalI, > 95% of the DNA fragments can be ligated with T4 DNA Ligase (at a 5' termini concentration of 1–2 µM) at 16°C. Of these ligated fragments, > 95% can be recut.

16-Hour Incubation: A 50 µl reaction containing 1 µg of pBR322 DNA and 150 units of SalI incubated for 16 hours at 37°C resulted in a DNA pattern free of detectable nuclease degradation as determined by agarose gel electrophoresis.

Exonuclease Activity: Incubation of 200 units of enzyme with 1 µg sonicated ³H DNA (10⁵ cpm/µg) for 4 hours at 37°C in 50 µl reaction buffer released < 0.1% radioactivity.

Endonuclease Activity: Incubation of 40 units with 1 µg φX174 RF I DNA for 4 hours at 37°C in 50 µl reaction buffer resulted in < 20% conversion to RF II.

Blue/White Screening Assay: This enzyme has been tested to determine the integrity of the DNA ends produced after digestion with an excess of enzyme. An appropriate vector is digested at a unique site within *lacZ'* gene with a 10-fold excess of enzyme, ligated, transformed and plated on XGal/IPTG/Amp plates. Successful expression of β-galactosidase is a function of how intact its gene remains after cloning, an intact gene gives rise to a blue colony, an interrupted gene (i.e. degraded DNA end) gives rise to a white colony. Enzymes must produce fewer than 3% white colonies to be Blue/White Certified.

Enzyme Properties

Activity in NEBuffers:

NEBuffer 1	0%
NEBuffer 2	0%
NEBuffer 3	100%
NEBuffer 4	0%

When using a buffer other than the optimal (supplied) NEBuffer, it may be necessary to add more enzyme to achieve complete digestion.

(see other side)

CERTIFICATE OF ANALYSIS

Blue/White Screening Assay: This enzyme has been tested to determine the integrity of the DNA ends produced after digestion with an excess of enzyme. An appropriate vector is digested at a unique site within *lacZ'* gene with a 10-fold excess of enzyme, ligated, transformed and plated on XGal/IPTG/Amp plates. Successful expression of β-galactosidase is a function of how intact its gene remains after cloning, an intact gene gives rise to a blue colony, an interrupted gene (i.e. degraded DNA end) gives rise to a white colony. Enzymes must produce fewer than 3% white colonies to be Blue/White Certified.

Enzyme Properties

Activity in NEBuffers:

NEBuffer 1	0%
NEBuffer 2	0%
NEBuffer 3	100%
NEBuffer 4	0%

When using a buffer other than the optimal (supplied) NEBuffer, it may be necessary to add more enzyme to achieve complete digestion.

(see other side)

CERTIFICATE OF ANALYSIS

Survival in a Reaction: A minimum of 0.13 unit is required to digest 1 µg of substrate DNA in 16 hours.

Heat Inactivation: 50 units of enzyme were inactivated by incubation at 65°C for 20 minutes.

Notes: Supercoiled forms of pBR322 and pUC require 10-fold overdigestion with Sall to achieve complete digestion.


Cleavage of mammalian genomic DNA is blocked by CpG methylation.

Conditions of low ionic strength, high enzyme concentration, glycerol concentration > 5% or pH > 8.0 may result in star activity.

Companion Products:

Sall-HF™	
#R3138S	2,000 units
#R3138L	10,000 units
#R3138T	2,000 units
#R3138M	10,000 units

Sall-HF™ RE-Mix™	
#R5138S	100 reactions

 = Time-Saver™ Qualified (See www.neb.com for details).

Survival in a Reaction: A minimum of 0.13 unit is required to digest 1 µg of substrate DNA in 16 hours.

Heat Inactivation: 50 units of enzyme were inactivated by incubation at 65°C for 20 minutes.

Notes: Supercoiled forms of pBR322 and pUC require 10-fold overdigestion with Sall to achieve complete digestion.


Cleavage of mammalian genomic DNA is blocked by CpG methylation.

Conditions of low ionic strength, high enzyme concentration, glycerol concentration > 5% or pH > 8.0 may result in star activity.

Companion Products:

Sall-HF™	
#R3138S	2,000 units
#R3138L	10,000 units
#R3138T	2,000 units
#R3138M	10,000 units

Sall-HF™ RE-Mix™	
#R5138S	100 reactions

 = Time-Saver™ Qualified (See www.neb.com for details).