

XbaI



1-800-632-7799
info@neb.com
www.neb.com



R0145S 041120814081

R0145S



3,000 units **20,000 U/ml** **Lot: 0411208**
RECOMBINANT **Store at -20°C** **Exp: 8/14**

Recognition Site:

5'...T[▼]CTAGA...3'
3'...AGATC[▲]T...5'

Source: An *E. coli* strain that carries the cloned XbaI gene from *Xanthomonas badrii* (ATCC 11672)

New Reaction Buffer

XbaI



1-800-632-7799
info@neb.com
www.neb.com



R0145S 041120814081

R0145S



3,000 units **20,000 U/ml** **Lot: 0411208**
RECOMBINANT **Store at -20°C** **Exp: 8/14**

Recognition Site:

5'...T[▼]CTAGA...3'
3'...AGATC[▲]T...5'

Source: An *E. coli* strain that carries the cloned XbaI gene from *Xanthomonas badrii* (ATCC 11672)

New Reaction Buffer

Supplied in: 50 mM NaCl, 10 mM Tris-HCl (pH 7.4), 0.1 mM EDTA, 1 mM DTT, 200 µg/ml BSA and 50% glycerol.

Reagents Supplied with Enzyme:
10X NEBuffer 4, 100X BSA.

Reaction Conditions: 1X NEBuffer 4, supplemented with 100 µg/ml BSA. Incubate at 37°C.

1X NEBuffer 4:
50 mM potassium acetate
20 mM Tris-acetate
10 mM magnesium acetate
1 mM DTT
pH 7.9 @ 25°C

Unit Definition: One unit is defined as the amount of enzyme required to digest 1 µg of λ DNA (dam⁻/Hind III digest) in 1 hour at 37°C in a total reaction volume of 50 µl.

Diluent Compatibility: Diluent Buffer A
50 mM KCl, 10 mM Tris-HCl, 0.1 mM EDTA, 1 mM DTT, 200 µg/ml BSA and 50% glycerol (pH 7.4 @ 25°C)

Supplied in: 50 mM NaCl, 10 mM Tris-HCl (pH 7.4), 0.1 mM EDTA, 1 mM DTT, 200 µg/ml BSA and 50% glycerol.

Reagents Supplied with Enzyme:
10X NEBuffer 4, 100X BSA.

Reaction Conditions: 1X NEBuffer 4, supplemented with 100 µg/ml BSA. Incubate at 37°C.

1X NEBuffer 4:
50 mM potassium acetate
20 mM Tris-acetate
10 mM magnesium acetate
1 mM DTT
pH 7.9 @ 25°C

Unit Definition: One unit is defined as the amount of enzyme required to digest 1 µg of λ DNA (dam⁻/Hind III digest) in 1 hour at 37°C in a total reaction volume of 50 µl.

Diluent Compatibility: Diluent Buffer A
50 mM KCl, 10 mM Tris-HCl, 0.1 mM EDTA, 1 mM DTT, 200 µg/ml BSA and 50% glycerol (pH 7.4 @ 25°C)

Quality Control Assays

Ligation: After 20-fold overdigestion with XbaI, > 95% of the DNA fragments can be ligated with T4 DNA Ligase (at a 5' termini concentration of 1–2 µM) at 16°C. Of these ligated fragments, > 95% can be recut.

16-Hour Incubation: A 50 µl reaction containing 1 µg of DNA and 200 units of enzyme incubated for 16 hours resulted in the same pattern of DNA bands as a reaction incubated for 1 hour with 1 unit of enzyme.

Exonuclease Activity: Incubation of 200 units of enzyme with 1 µg sonicated ³H DNA (10⁵ cpm/µg) for 4 hours at 37°C in 50 µl reaction buffer released < 0.1% radioactivity.

Endonuclease Activity: Incubation of 200 units of enzyme with 1 µg φX174 RF I DNA for 4 hours at 37°C in 50 µl reaction buffer resulted in < 10% conversion to RF II.

Quality Control Assays

Ligation: After 20-fold overdigestion with XbaI, > 95% of the DNA fragments can be ligated with T4 DNA Ligase (at a 5' termini concentration of 1–2 µM) at 16°C. Of these ligated fragments, > 95% can be recut.

16-Hour Incubation: A 50 µl reaction containing 1 µg of DNA and 200 units of enzyme incubated for 16 hours resulted in the same pattern of DNA bands as a reaction incubated for 1 hour with 1 unit of enzyme.

Exonuclease Activity: Incubation of 200 units of enzyme with 1 µg sonicated ³H DNA (10⁵ cpm/µg) for 4 hours at 37°C in 50 µl reaction buffer released < 0.1% radioactivity.

Endonuclease Activity: Incubation of 200 units of enzyme with 1 µg φX174 RF I DNA for 4 hours at 37°C in 50 µl reaction buffer resulted in < 10% conversion to RF II.

Blue/White Screening Assay: This enzyme has been tested to determine the integrity of the DNA ends produced after digestion with an excess of enzyme. An appropriate vector is digested at a unique site within *lacZ^α* gene with a 10-fold excess of enzyme, ligated, transformed and plated on XGal/IPTG/Amp plates. Successful expression of β-galactosidase is a function of how intact its gene remains after cloning, an intact gene gives rise to a blue colony, an interrupted gene (i.e. degraded DNA end) gives rise to a white colony. Enzymes must produce fewer than 3% white colonies to be Blue/White Certified.

Survival in a Reaction: A minimum of 0.13 unit is required to digest 1 µg of substrate DNA in 16 hours.

(see other side)

CERTIFICATE OF ANALYSIS

Blue/White Screening Assay: This enzyme has been tested to determine the integrity of the DNA ends produced after digestion with an excess of enzyme. An appropriate vector is digested at a unique site within *lacZ^α* gene with a 10-fold excess of enzyme, ligated, transformed and plated on XGal/IPTG/Amp plates. Successful expression of β-galactosidase is a function of how intact its gene remains after cloning, an intact gene gives rise to a blue colony, an interrupted gene (i.e. degraded DNA end) gives rise to a white colony. Enzymes must produce fewer than 3% white colonies to be Blue/White Certified.

Survival in a Reaction: A minimum of 0.13 unit is required to digest 1 µg of substrate DNA in 16 hours.

(see other side)

CERTIFICATE OF ANALYSIS

Enzyme Properties

Activity in NEBuffers:

NEBuffer 1	0%
NEBuffer 2	100%
NEBuffer 3	75%
NEBuffer 4	100%

When using a buffer other than the optimal (supplied) NEBuffer, it may be necessary to add more enzyme to achieve complete digestion.

Heat Inactivation: 200 units of enzyme were inactivated by incubation at 65°C for 20 minutes.

Plasmid Cleavage: Number of units required to cleave 1 µg of supercoiled plasmid DNA in one hour: pUC 19 = 2 units.

Note: Blocked by overlapping *dam* methylation.

Companion Products:

XbaI RE-Mix™

#R5145S 150 reactions

dam⁻/*dcm*⁻ Competent *E. coli*

#C2925H 20 transformation reactions

#C2925I 24 transformation reactions

 = Time-Saver™ Qualified (See www.neb.com for details).

U.S. Patent No. 4,983,542

Enzyme Properties

Activity in NEBuffers:

NEBuffer 1	0%
NEBuffer 2	100%
NEBuffer 3	75%
NEBuffer 4	100%

When using a buffer other than the optimal (supplied) NEBuffer, it may be necessary to add more enzyme to achieve complete digestion.

Heat Inactivation: 200 units of enzyme were inactivated by incubation at 65°C for 20 minutes.

Plasmid Cleavage: Number of units required to cleave 1 µg of supercoiled plasmid DNA in one hour: pUC 19 = 2 units.

Note: Blocked by overlapping *dam* methylation.

Companion Products:

XbaI RE-Mix™

#R5145S 150 reactions

dam⁻/*dcm*⁻ Competent *E. coli*

#C2925H 20 transformation reactions

#C2925I 24 transformation reactions

 = Time-Saver™ Qualified (See www.neb.com for details).

U.S. Patent No. 4,983,542