

Monoclonal antibodies to Human Interleukin-12/-23 (p40)

Antibody:	MT86/221 The reagent consists of a combination of two different monoclonal antibodies.
Product code:	3450-3-250
Quantity:	250 µg
Immunogen:	Recombinant human IL-23
Isotype:	Mouse IgG1
Specificity:	Native and recombinant p40 subunit of human IL-12 and IL-23. The antibodies recognise p40 monomer, p40 homodimer, IL-12 p70 heterodimer (p40+p35) and IL-23 heterodimer (p40+p19). The antibodies cross-react with native p40 from rhesus and cynomolgus macaques and recombinant p40 from rhesus macaques. Please contact Mabtech for information on reactivity with other monkey species.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing.
Applications:	For quantification of IL-12/-23 (p40), IL-12 (p70) and IL-23 in solution e.g. cell culture supernatants using ELISA and for enumeration of IL-12/-23 (p40), IL-12 (p70) and IL-23 producing cells using ELISpot. MT86/221 are recommended as coating mAbs for analysis of: - IL-12/-23 (p40) in combination with mAb MT618 (product code 3450-6) - IL-12 (p70) in combination with mAb IL12-III (product code 3450-8) - IL-23 in combination with mAb MT155 (product code 3457-6)
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



MABTECH AB
Box 1233
SE-131 28 Nacka Strand
Sweden
Tel +46 8 716 27 00
Fax +46 8 716 27 01
mabtech@mabtech.com
www.mabtech.com

Mabtech, Inc.
M.E.B. Suite 220
3814 West Street
Mariemont, OH 45227
USA
Toll Free: +1 866 ELI-SPOT
Tel: +1 513 871 4500
Fax: +1 513 871 7353
mabtech.usa@mabtech.com

MABTECH AUSTRALIA
Pty Ltd
resolvingIMAGES
Unit 22, 196 Settlement Rd
Thomastown Victoria 3074
Australia
Tel: +61 3 9466 4007
Fax: +61 3 9466 4003
mabtech.au@mabtech.com

Mabtech AB
Büro Deutschland
GERMANY
Tel: +49 40 4135 7935
Fax: +49 40 4135 7945
mabtech.de@mabtech.com

Mabtech AB
Bureau de liaison France
BP 255, 1300 route des
Crêtes
06905 Sophia Antipolis
FRANCE
Tel: +33 4 92 38 80 70
Fax: +33 4 92 38 80 71
mabtech.fr@mabtech.com