

# Assay Buffer

- Product code:** 3652-J2
- Contents:** Protein-containing buffer, ready-to-use, with 0.15% Kathon CG as preservative.
- Applications:** To prevent interference by heterophilic antibodies in serum and plasma samples in ELISA. The Assay buffer is recommended for use with MIP-1 $\beta$ , apoA1, apoB, apoE and apoH ELISA tests. Please note that the Assay buffer is optimized for use with antibodies from Mabtech.
- Instructions:** The Assay buffer should be used for dilution of samples, standard and detection antibody. All samples should be diluted at least 2x in Assay buffer. Do not use undiluted serum/plasma samples in the assay. Sample dilution can be made either in separate tubes or directly in the plate. If sample dilution is made directly in the plate, Assay buffer should be added prior to addition of the sample.
- Important!** Do not use the Assay buffer for dilution of Streptavidin conjugates. We recommend the use of an incubation buffer e.g. PBS with 0.05% Tween 20 containing 0.1% BSA for dilution of conjugates.
- Storage:** Shipped at ambient temperature. On arrival store product at 4-8°C. Adjust to room temperature before use.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



MABTECH AB  
Box 1233  
SE-131 28 Nacka Strand  
Sweden  
Tel +46 8 716 27 00  
Fax +46 8 716 27 01  
mabtech@mabtech.com  
www.mabtech.com

Mabtech, Inc.  
M.E.B. Suite 220  
3814 West Street  
Mariemont, OH 45227  
USA  
Toll Free: +1 866 ELI-SPOT  
Tel: +1 513 871 4500  
Fax: +1 513 871 7353  
mabtech.usa@mabtech.com

MABTECH AUSTRALIA  
Pty Ltd  
resolvingIMAGES  
Unit 22  
Thomastown Victoria 3074  
Australia  
Tel: +61 3 9466 4007  
Fax: +61 3 9466 4003  
mabtech.au@mabtech.com

Mabtech AB  
Büro Deutschland  
GERMANY  
Tel: +49 40 4135 7935  
Fax: +49 40 4135 7945  
mabtech.de@mabtech.com

Mabtech AB  
Bureau de liaison France  
BP 255, 1300 route des  
Crêtes  
06905 Sophia Antipolis  
FRANCE  
Tel: +33 4 92 38 80 70  
Fax: +33 4 92 38 80 71  
mabtech.fr@mabtech.com